

Table of contents
-------------------

\\CRC

protocols

studies

v5motion

[localizer](#)  
[AAHead\\_Scout\\_32ch-head-coil](#)  
[b1map\\_sag\\_p2](#)  
[gre\\_field\\_mapping\\_2mm\\_AABrain](#)  
[cmrr\\_mbep2d\\_p3\\_mb2\\_1.6iso\\_AABrain](#)  
[cmrr\\_mbep2d\\_p3\\_mb2\\_1.6iso\\_reversePE\\_AABrain](#)  
[t1\\_mp2rage\\_sag\\_p3\\_0.75mm\\_AAbasis](#)  
[mp2rage\\_siemens\\_ipat3\\_0.69mm\\_AAbasis](#)  
[cmrr\\_mbep2d\\_p3\\_mb2\\_1.6iso\\_Neuro3D\\_Visual\\_AABrain](#)  
[cmrr\\_mbep2d\\_p3\\_mb2\\_1.6iso\\_Neuro3D\\_Audio\\_AABrain](#)  
[cmrr\\_mbep2d\\_p3\\_mb2\\_1.6iso\\_Neuro3D\\_Visual\\_v2\\_AABrain](#)  
[cmrr\\_mbep2d\\_p3\\_mb2\\_1.6iso\\_Neuro3D\\_Audio\\_v2\\_AABrain](#)  
[cmrr\\_mbep2d\\_p3\\_mb2\\_1iso\\_AABrain](#)  
[cmrr\\_mbep2d\\_p4\\_mb2\\_750um\\_AABrain](#)  
[cmrr\\_mbep2d\\_p3\\_mb2\\_1.6iso\\_AABrain](#)  
[cmrr\\_mbep2d\\_p3\\_mb2\\_1.6iso\\_reversePE\\_AABrain](#)  
[gre\\_field\\_mapping\\_2mm\\_eufind](#)  
[cmrr\\_mbep2d\\_p4\\_mb2\\_750um\\_AABrain](#)  
[cmrr\\_mbep2d\\_p4\\_mb2\\_750um\\_reversePE\\_AABrain](#)  
[cmrr\\_mbep2d\\_p4\\_mb2\\_750um\\_GRAPPA-GRE\\_FA75](#)  
[cmrr\\_mbep2d\\_p4\\_mb2\\_750um\\_AABrain](#)  
[cmrr\\_mbep2d\\_p4\\_mb2\\_750um\\_GRAPPA\\_GRE\\_FA90\\_B1148](#)  
[cmrr\\_mbep2d\\_p4\\_mb2\\_750um\\_GRAPPA-GRE\\_FA75\\_Band\\_T1](#)  
[cmrr\\_mbep2d\\_p4\\_mb2\\_750um\\_GRAPPA-GRE\\_FA75\\_Band\\_T1](#)  
[cmrr\\_mbep2d\\_p4\\_mb2\\_750um\\_GRAPPA-GRE\\_FA75\\_Band\\_T1](#)  
[rslh\\_ep3d\\_vaso\\_4e\\_axial](#)  
[rslh\\_ep3d\\_bold\\_4e\\_axial](#)  
[rslh\\_ep3d\\_vaso\\_4e\\_sagital](#)  
[rslh\\_ep3d\\_bold\\_4e\\_sagital](#)  
[t1\\_mp2rage\\_ax\\_p2\\_0.5mm\\_slab](#)  
[rslh\\_ep3d\\_vaso\\_4e\\_axial\\_p05](#)

\\CRC\protocols\studies\v5motion\localizer

TA: 0:15 PM: FIX Voxel size: 0.5×0.5×5.0 mmPAT: Off Rel. SNR: 1.00 : qfl

**Properties**

Prio recon	On
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	On
Load images to graphic segments	On
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	Single measurement

**Routine**

Slice group	1
Slices	1
Dist. factor	20 %
Position	Isocenter
Orientation	Sagittal
Phase enc. dir.	A >> P
Slice group	2
Slices	1
Dist. factor	20 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
Slice group	3
Slices	1
Dist. factor	20 %
Position	Isocenter
Orientation	Coronal
Phase enc. dir.	R >> L
AutoAlign	---
Phase oversampling	0 %
FoV read	250 mm
FoV phase	100.0 %
Slice thickness	5.0 mm
TR	8.6 ms
TE	3.69 ms
Averages	2
Concatenations	3
Filter	Elliptical filter
Coil elements	A32

**Contrast - Common**

TR	8.6 ms
TE	3.69 ms
TD	0 ms
MTC	Off
Magn. preparation	None
Flip angle	20 deg
Fat suppr.	None
Water suppr.	None
SWI	Off

**Contrast - Dynamic**

Averages	2
Averaging mode	Short term
Reconstruction	Magnitude
Measurements	1

**Contrast - Dynamic**

Multiple series	Each measurement
-----------------	------------------

**Resolution - Common**

FoV read	250 mm
FoV phase	100.0 %
Slice thickness	5.0 mm
Base resolution	256
Phase resolution	100 %
Phase partial Fourier	Off
Interpolation	On

**Resolution - iPAT**

PAT mode	None
----------	------

**Resolution - Filter Image**

Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off

**Resolution - Filter Rawdata**

Raw filter	Off
Elliptical filter	On

**Geometry - Common**

Slice group	1
Slices	1
Dist. factor	20 %
Position	Isocenter
Orientation	Sagittal
Phase enc. dir.	A >> P
Slice group	2
Slices	1
Dist. factor	20 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
Slice group	3
Slices	1
Dist. factor	20 %
Position	Isocenter
Orientation	Coronal
Phase enc. dir.	R >> L
FoV read	250 mm
FoV phase	100.0 %
Slice thickness	5.0 mm
TR	8.6 ms
Multi-slice mode	Sequential
Series	Interleaved
Concatenations	3

**Geometry - AutoAlign**

Slice group	1
Position	Isocenter
Orientation	Sagittal
Phase enc. dir.	A >> P
Slice group	2
Position	Isocenter
Orientation	Transversal

**Geometry - AutoAlign**

Phase enc. dir.	A >> P
Slice group	3
Position	Isocenter
Orientation	Coronal
Phase enc. dir.	R >> L
AutoAlign	---
Initial Position	Isocenter
L	0.0 mm
P	0.0 mm
H	0.0 mm
Initial Rotation	0.00 deg
Initial Orientation	Sagittal

**Geometry - Saturation**

Saturation mode	Standard
Fat suppr.	None
Water suppr.	None
Special sat.	None

**Geometry - Tim CT**

Tim CT mode	Off
Slices	1
Slice thickness	5.0 mm
Dist. factor	20 %
FoV read	250 mm
FoV phase	100.0 %
Segments	1

**System - Miscellaneous**

Positioning mode	FIX
Table position	F
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
Save uncombined	Off
Matrix Optimization	Off
AutoAlign	---
Coil Select Mode	Off - AutoCoilSelect

**System - Adjustments**

B0 Shim mode	Tune up
B1 Shim mode	TrueForm
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

**System - Adjust Volume**

Position	Isocenter
Orientation	Transversal
Rotation	0.00 deg
A >> P	263 mm
R >> L	350 mm
F >> H	350 mm
Reset	Off

**System - Tx/Rx**

Frequency 1H	297.162474 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000

**System - Tx/Rx**

Reset	Off
? Ref. amplitude 1H	0.000 V

**Physio - Signal1**

1st Signal/Mode	None
TR	8.6 ms
Concatenations	3
Segments	1

**Physio - Cardiac**

Magn. preparation	None
Fat suppr.	None
Dark blood	Off
FoV read	250 mm
FoV phase	100.0 %
Phase resolution	100 %

**Physio - PACE**

Resp. control	Off
Concatenations	3

**Inline - Common**

Subtract	Off
Measurements	1
StdDev	Off
Liver registration	Off
Save original images	On

**Inline - MIP**

MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On

**Inline - Soft Tissue**

Wash - In	Off
Wash - Out	Off
TTP	Off
PEI	Off
MIP - time	Off
Measurements	1

**Inline - Composing**

Distortion Corr.	Off
------------------	-----

**Sequence - Part 1**

Introduction	On
Dimension	2D
Phase stabilisation	Off
Asymmetric echo	Allowed
Contrasts	1
Flow comp.	No
Multi-slice mode	Sequential
Bandwidth	320 Hz/Px

**Sequence - Part 2**

Segments	1
Acoustic noise reduction	Active
RF pulse type	Normal
Gradient mode	Normal
Excitation	Slice-sel.
RF spoiling	On

**Sequence - Nuclei**

TX/RX Nucleus	1H
TX/RX delta frequency	0 Hz
TX Nucleus	None
TX delta frequency	0 Hz
Coil elements	A32

**Sequence - Assistant**

Mode	Off
------	-----

\\CRC\protocols\studies\v5motion\AAHead\_Scout\_32ch-head-coil

TA: 0:14 PM: REF Voxel size: 1.6×1.6×1.6 mmPAT: 3 Rel. SNR: 1.00 : fl

**Properties**

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	On
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	Single measurement

**Routine**

Slab group	1
Slabs	1
Dist. factor	20 %
Position	L0.0 A30.0 H0.0 mm
Orientation	Sagittal
Phase enc. dir.	A >> P
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	128
FoV read	260 mm
FoV phase	100.0 %
Slice thickness	1.6 mm
TR	3.25 ms
TE	1.53 ms
Averages	1
Concatenations	1
Filter	B1 filter
Coil elements	A32

**Contrast - Common**

TR	3.25 ms
TE	1.53 ms
Flip angle	16 deg

**Contrast - Dynamic**

Averages	1
Averaging mode	Short term
Reconstruction	Magnitude
Measurements	1

**Resolution - Common**

FoV read	260 mm
FoV phase	100.0 %
Slice thickness	1.6 mm
Base resolution	160
Phase resolution	100 %
Slice resolution	69 %
Phase partial Fourier	6/8
Slice partial Fourier	6/8
Trajectory	Cartesian

**Resolution - iPAT**

PAT mode	GRAPPA
Accel. factor PE	3
Ref. lines PE	24
Accel. factor 3D	1

**Resolution - iPAT**

Reference scan mode	Integrated
---------------------	------------

**Resolution - Filter Image**

Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	On
Unfiltered images	Off

**Resolution - Filter Rawdata**

Raw filter	Off
Elliptical filter	Off

**Geometry - Common**

Slab group	1
Slabs	1
Dist. factor	20 %
Position	L0.0 A30.0 H0.0 mm
Orientation	Sagittal
Phase enc. dir.	A >> P
Slice oversampling	0.0 %
Slices per slab	128
FoV read	260 mm
FoV phase	100.0 %
Slice thickness	1.6 mm
TR	3.25 ms
Multi-slice mode	Sequential
Series	Ascending
Concatenations	1

**Geometry - AutoAlign**

Slab group	1
Position	L0.0 A30.0 H0.0 mm
Orientation	Sagittal
Phase enc. dir.	A >> P
Initial Position	Isocenter
L	0.0 mm
P	0.0 mm
H	0.0 mm
Initial Rotation	0.00 deg
Initial Orientation	Transversal

**System - Miscellaneous**

Positioning mode	REF
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
Save uncombined	Off
Matrix Optimization	Off
Coil Select Mode	Default

**System - Adjustments**

B0 Shim mode	Tune up
B1 Shim mode	TrueForm
Confirm freq. adjustment	Off

**System - Adjustments**

Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

**System - Adjust Volume**

Position	Isocenter
Orientation	Transversal
Rotation	0.00 deg
A >> P	263 mm
R >> L	350 mm
F >> H	350 mm
Reset	Off

**System - Tx/Rx**

Frequency 1H	297.162474 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

**Physio - PACE**

Resp. control	Off
Concatenations	1

**Inline - Common**

Flip angle	16 deg
Measurements	1
Time to center	6.3 s

**Inline - Inline**

Subtract	Off
Measurements	1
StdDev	Off
Save original images	On

**Inline - MIP**

MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On

**Inline - Composing**

Distortion Corr.	Off
------------------	-----

**Sequence - Part 1**

Introduction	On
Dimension	3D
Asymmetric echo	Weak
Contrasts	1
Multi-slice mode	Sequential
Bandwidth	540 Hz/Px

**Sequence - Part 2**

RF pulse type	Fast
Gradient mode	Normal
Excitation	Non-sel.
RF spoiling	On

**Sequence - Nuclei**

TX/RX Nucleus	1H
TX/RX delta frequency	0 Hz

**Sequence - Nuclei**

TX Nucleus	None
TX delta frequency	0 Hz
Coil elements	A32

**Sequence - Assistant**

Mode	Off
------	-----

## \\CRC\protocols\studies\v5motion\b1map\_sag\_p2

TA: 9.2 s PM: FIX Voxel size: 4.0×4.0×4.0 mmPAT: 2 Rel. SNR: 1.00 : tfl

**Properties**

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	Off
Start measurements	Single measurement

**Routine**

Slice group	1
Slices	25
Dist. factor	100 %
Position	Isocenter
Orientation	Sagittal
Phase enc. dir.	A >> P
AutoAlign	---
Phase oversampling	0 %
FoV read	256 mm
FoV phase	100.0 %
Slice thickness	4.0 mm
TR	4000.0 ms
TE	1.72 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	A32

**Contrast - Common**

TR	4000.0 ms
TE	1.72 ms
Magn. preparation	None
Flip angle	10 deg
Fat suppr.	None
Water suppr.	None

**Contrast - Dynamic**

Averages	1
Reconstruction	Magnitude
Measurements	1
Multiple series	Each measurement

**Resolution - Common**

FoV read	256 mm
FoV phase	100.0 %
Slice thickness	4.0 mm
Base resolution	64
Phase resolution	100 %
Phase partial Fourier	Off
Interpolation	Off

**Resolution - iPAT**

PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	16
Reference scan mode	Integrated

**Resolution - Filter Image**

Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off

**Resolution - Filter Rawdata**

Raw filter	Off
Elliptical filter	Off

**Geometry - Common**

Slice group	1
Slices	25
Dist. factor	100 %
Position	Isocenter
Orientation	Sagittal
Phase enc. dir.	A >> P
FoV read	256 mm
FoV phase	100.0 %
Slice thickness	4.0 mm
TR	4000.0 ms
Multi-slice mode	Interleaved
Series	Interleaved
Concatenations	1

**Geometry - AutoAlign**

Slice group	1
Position	Isocenter
Orientation	Sagittal
Phase enc. dir.	A >> P
AutoAlign	---
Initial Position	Isocenter
L	0.0 mm
P	0.0 mm
H	0.0 mm
Initial Rotation	0.00 deg
Initial Orientation	Sagittal

**System - Miscellaneous**

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
Save uncombined	Off
Matrix Optimization	Off
AutoAlign	---
Coil Select Mode	Default

**System - Adjustments**

B0 Shim mode	Tune up
B1 Shim mode	TrueForm
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

**System - Adjust Volume**

Position	Isocenter
Orientation	Transversal
Rotation	0.00 deg
A >> P	263 mm
R >> L	350 mm
F >> H	350 mm
Reset	Off

**System - Tx/Rx**

Frequency 1H	297.162474 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

**Inline - Common**

Subtract	Off
Measurements	1
StdDev	Off
Save original images	On

**Inline - MIP**

MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On

**Inline - Composing**

Distortion Corr.	Off
------------------	-----

**Sequence - Part 1**

Introduction	Off
Dimension	2D
Asymmetric echo	Allowed
Flow comp.	No
Multi-slice mode	Interleaved
Echo spacing	3.9 ms
Bandwidth	490 Hz/Px

**Sequence - Part 2**

RF pulse type	Fast
Gradient mode	Normal
Excitation	Slice-sel.
RF spoiling	On
Turbo factor	64

**Sequence - Assistant**

Mode	Off
------	-----



\\CRC\protocols\studies\v5motion\gre\_field\_mapping\_2mm\_AAbra

TA: 1:38 PM: FIX Voxel size: 2.0×2.0×2.0 mmRel. SNR: 1.00 : fm

**Properties**

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	Off
Start measurements	Single measurement

**Routine**

Slice group	1
Slices	96
Dist. factor	0 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	R >> L
AutoAlign	Head > Brain
Phase oversampling	0 %
FoV read	256 mm
FoV phase	75.0 %
Slice thickness	2.0 mm
TR	5.2 ms
TE 1	2.26 ms
TE 2	3.28 ms
Averages	1
Concatenations	96
Filter	None
Coil elements	A32

**Contrast - Common**

TR	5.2 ms
TE 1	2.26 ms
TE 2	3.28 ms
MTC	Off
Flip angle	15 deg
Fat suppr.	None

**Contrast - Dynamic**

Averages	1
Averaging mode	Long term
Reconstruction	Magn./Phase
Measurements	1
Multiple series	Off

**Resolution - Common**

FoV read	256 mm
FoV phase	75.0 %
Slice thickness	2.0 mm
Base resolution	128
Phase resolution	100 %
Phase partial Fourier	Off
Interpolation	Off

**Resolution - Filter Image**

Image Filter	Off
Distortion Corr.	Off

**Resolution - Filter Image**

Prescan Normalize	Off
Normalize	Off
B1 filter	Off

**Resolution - Filter Rawdata**

Raw filter	Off
Elliptical filter	Off

**Geometry - Common**

Slice group	1
Slices	96
Dist. factor	0 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	R >> L
FoV read	256 mm
FoV phase	75.0 %
Slice thickness	2.0 mm
TR	5.2 ms
Multi-slice mode	Sequential
Series	Ascending
Concatenations	96

**Geometry - AutoAlign**

Slice group	1
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	R >> L
AutoAlign	Head > Brain
Initial Position	Isocenter
L	0.0 mm
P	0.0 mm
H	0.0 mm
Initial Rotation	90.00 deg
Initial Orientation	Transversal

**Geometry - Saturation**

Fat suppr.	None
Special sat.	None

**System - Miscellaneous**

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
Save uncombined	Off
Matrix Optimization	Off
AutoAlign	Head > Brain
Coil Select Mode	Default

**System - Adjustments**

B0 Shim mode	Brain
B1 Shim mode	TrueForm
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off

**System - Adjustments**

Adjustment Tolerance	Auto
----------------------	------

**System - Adjust Volume**

Position	Isocenter
Orientation	Transversal
Rotation	90.00 deg
R >> L	192 mm
A >> P	256 mm
F >> H	192 mm
Reset	Off

**System - Tx/Rx**

Frequency 1H	297.162474 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

**Sequence - Part 1**

Introduction	On
Dimension	2D
Asymmetric echo	Off
Contrasts	2
Flow comp.	No
Multi-slice mode	Sequential
Bandwidth	737 Hz/Px

**Sequence - Part 2**

RF pulse type	Fast
Gradient mode	Fast
RF spoiling	On

**Sequence - Assistant**

Mode	Off
------	-----

\\CRC\protocols\studies\v5motion\cmrr\_mbep2d\_p3\_mb2\_1.6iso\_AABrain

TA: 10:19 PM: FIX Voxel size: 1.6×1.6×1.6 mmPAT: 3 Rel. SNR: 1.00 : efpid

**Properties**

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	Off
Start measurements	Single measurement

**Routine**

Slice group	1
Slices	72
Dist. factor	0 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	Head > Brain
Phase oversampling	0 %
FoV read	216 mm
FoV phase	101.5 %
Slice thickness	1.60 mm
TR	1800 ms
TE	22.00 ms
Multi-band accel. factor	2
Filter	None
Coil elements	A32

**Contrast - Common**

TR	1800 ms
TE	22.00 ms
MTC	Off
Magn. preparation	None
Flip angle	70 deg
Fat suppr.	Fat sat.

**Contrast - Dynamic**

Averaging mode	Long term
Reconstruction	Magnitude
Measurements	330
Delay in TR	0 ms
Multiple series	Off

**Resolution - Common**

FoV read	216 mm
FoV phase	101.5 %
Slice thickness	1.60 mm
Base resolution	134
Phase resolution	100 %
Phase partial Fourier	Off
Interpolation	Off

**Resolution - iPAT**

PAT mode	GRAPPA
Accel. factor PE	3
Ref. lines PE	36
Reference scan mode	Segmented

**Resolution - Filter Image**

Distortion Corr.	Off
Prescan Normalize	Off

**Resolution - Filter Rawdata**

Raw filter	Off
Elliptical filter	Off
Hamming	Off

**Geometry - Common**

Slice group	1
Slices	72
Dist. factor	0 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
FoV read	216 mm
FoV phase	101.5 %
Slice thickness	1.60 mm
TR	1800 ms
Multi-slice mode	Interleaved
Series	Interleaved
Multi-band accel. factor	2

**Geometry - AutoAlign**

Slice group	1
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	Head > Brain
Initial Position	Isocenter
L	0.0 mm
P	0.0 mm
H	0.0 mm
Initial Rotation	-0.01 deg
Initial Orientation	Transversal

**Geometry - Saturation**

Fat suppr.	Fat sat.
Special sat.	None

**System - Miscellaneous**

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
Matrix Optimization	Off
AutoAlign	Head > Brain
Coil Select Mode	Default

**System - Adjustments**

B0 Shim mode	Standard
B1 Shim mode	TrueForm
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

**System - Adjust Volume**

Position	Isocenter
Orientation	Transversal
Rotation	89.99 deg
R >> L	216 mm
A >> P	220 mm
F >> H	116 mm
Reset	Off

**System - Tx/Rx**

Frequency 1H	297.162474 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

**Physio - Signal1**

1st Signal/Mode	None
TR	1800 ms
Multi-band accel. factor	2

**BOLD**

GLM Statistics	Off
Dynamic t-maps	Off
Ignore meas. at start	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	20
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active
Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Meas[16]	Active
Meas[17]	Active
Meas[18]	Active
Meas[19]	Active
Meas[20]	Active
Motion correction	Off
Spatial filter	Off
Measurements	330
Delay in TR	0 ms
Multiple series	Off

**Sequence - Part 1**

Introduction	Off
Contrasts	1
Flow comp.	No
Multi-slice mode	Interleaved
Free echo spacing	Off
Echo spacing	0.59 ms
Bandwidth	2072 Hz/Px

**Sequence - Part 2**

EPI factor	136
Gradient mode	Fast
RF spoiling	Off

**Sequence - Special**

Excite pulse duration	7200 us
Single-band images	Off
MB LeakBlock kernel	On
MB dual kernel	Off
MB RF phase scramble	Off
SENSE1 coil combine	Off
Invert RO/PE polarity	Off
Disable freq. update	Off
Force equal slice timing	Off
Online multi-band recon.	Online
FFT scale factor	1.00
Fat saturation FA	110.0 deg
Physio recording	DICOM
Triggering scheme	Standard

\\CRC\protocols\studies\v5motion\cmrr\_mbep2d\_p3\_mb2\_1.6iso\_reversePE\_AABrain

TA: 0:34 PM: FIX Voxel size: 1.6×1.6×1.6 mmPAT: 3 Rel. SNR: 1.00 : epfid

**Properties**

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	Off
Start measurements	Single measurement

**Routine**

Slice group	1
Slices	72
Dist. factor	0 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
Phase oversampling	0 %
FoV read	216 mm
FoV phase	101.5 %
Slice thickness	1.60 mm
TR	1800 ms
TE	22.00 ms
Multi-band accel. factor	2
Filter	None
Coil elements	A32

**Contrast - Common**

TR	1800 ms
TE	22.00 ms
MTC	Off
Magn. preparation	None
Flip angle	70 deg
Fat suppr.	Fat sat.

**Contrast - Dynamic**

Averaging mode	Long term
Reconstruction	Magnitude
Measurements	5
Delay in TR	0 ms
Multiple series	Off

**Resolution - Common**

FoV read	216 mm
FoV phase	101.5 %
Slice thickness	1.60 mm
Base resolution	134
Phase resolution	100 %
Phase partial Fourier	Off
Interpolation	Off

**Resolution - iPAT**

PAT mode	GRAPPA
Accel. factor PE	3
Ref. lines PE	36
Reference scan mode	Segmented

**Resolution - Filter Image**

Distortion Corr.	Off
Prescan Normalize	Off

**Resolution - Filter Rawdata**

Raw filter	Off
Elliptical filter	Off
Hamming	Off

**Geometry - Common**

Slice group	1
Slices	72
Dist. factor	0 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
FoV read	216 mm
FoV phase	101.5 %
Slice thickness	1.60 mm
TR	1800 ms
Multi-slice mode	Interleaved
Series	Interleaved
Multi-band accel. factor	2

**Geometry - AutoAlign**

Slice group	1
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
Initial Position	Isocenter
L	0.0 mm
P	0.0 mm
H	0.0 mm
Initial Rotation	-0.01 deg
Initial Orientation	Transversal

**Geometry - Saturation**

Fat suppr.	Fat sat.
Special sat.	None

**System - Miscellaneous**

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
Matrix Optimization	Off
AutoAlign	---
Coil Select Mode	Default

**System - Adjustments**

B0 Shim mode	Standard
B1 Shim mode	TrueForm
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

**System - Adjust Volume**

Position	Isocenter
Orientation	Transversal
Rotation	89.99 deg
R >> L	216 mm
A >> P	220 mm
F >> H	116 mm
Reset	Off

**System - Tx/Rx**

Frequency 1H	297.162474 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

**Physio - Signal1**

1st Signal/Mode	None
TR	1800 ms
Multi-band accel. factor	2

**BOLD**

GLM Statistics	Off
Dynamic t-maps	Off
Ignore meas. at start	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	20
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active
Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Meas[16]	Active
Meas[17]	Active
Meas[18]	Active
Meas[19]	Active
Meas[20]	Active
Motion correction	Off
Spatial filter	Off
Measurements	5
Delay in TR	0 ms
Multiple series	Off

**Sequence - Part 1**

Introduction	Off
Contrasts	1
Flow comp.	No
Multi-slice mode	Interleaved
Free echo spacing	Off
Echo spacing	0.59 ms
Bandwidth	2072 Hz/Px

**Sequence - Part 2**

EPI factor	136
Gradient mode	Fast
RF spoiling	Off

**Sequence - Special**

Excite pulse duration	7200 us
Single-band images	Off
MB LeakBlock kernel	On
MB dual kernel	Off
MB RF phase scramble	Off
SENSE1 coil combine	Off
Invert RO/PE polarity	On
Disable freq. update	Off
Force equal slice timing	Off
Online multi-band recon.	Online
FFT scale factor	1.00
Fat saturation FA	110.0 deg
Physio recording	Off
Triggering scheme	Standard

\\CRC\protocols\studies\v5motion\t1\_mp2rage\_sag\_p3\_0.75mm\_AAbasis

TA: 8:50 PM: REF Voxel size: 0.8×0.8×0.8 mmPAT: 3 Rel. SNR: 1.00 : tfl

**Properties**

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	Off
Start measurements	Single measurement

**Routine**

Slab group	1
Slabs	1
Dist. factor	50 %
Position	Isocenter
Orientation	Sagittal
Phase enc. dir.	A >> P
AutoAlign	Head > Basis
Phase oversampling	0 %
Slice oversampling	8.3 %
Slices per slab	192
FoV read	240 mm
FoV phase	93.8 %
Slice thickness	0.75 mm
TR	4300.0 ms
TE	2.27 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	A32

**Contrast - Common**

TR	4300.0 ms
TE	2.27 ms
Magn. preparation	Non-sel. IR
TI 1	1000 ms
TI 2	3200 ms
Flip angle 1	4.0 deg
Flip angle 2	4.0 deg
Fat suppr.	Water excit. fast
Water suppr.	None

**Contrast - Dynamic**

Averages	1
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	1
Multiple series	Each measurement

**Resolution - Common**

FoV read	240 mm
FoV phase	93.8 %
Slice thickness	0.75 mm
Base resolution	320
Phase resolution	100 %
Slice resolution	100 %
Phase partial Fourier	Off

**Resolution - Common**

Slice partial Fourier	6/8
Interpolation	Off

**Resolution - iPAT**

PAT mode	GRAPPA
Accel. factor PE	3
Ref. lines PE	36
Accel. factor 3D	1
Reference scan mode	Integrated

**Resolution - Filter Image**

Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off

**Resolution - Filter Rawdata**

Raw filter	Off
Elliptical filter	Off

**Geometry - Common**

Slab group	1
Slabs	1
Dist. factor	50 %
Position	Isocenter
Orientation	Sagittal
Phase enc. dir.	A >> P
Slice oversampling	8.3 %
Slices per slab	192
FoV read	240 mm
FoV phase	93.8 %
Slice thickness	0.75 mm
TR	4300.0 ms
Multi-slice mode	Single shot
Series	Interleaved
Concatenations	1

**Geometry - AutoAlign**

Slab group	1
Position	Isocenter
Orientation	Sagittal
Phase enc. dir.	A >> P
AutoAlign	Head > Basis
Initial Position	Isocenter
L	0.0 mm
P	0.0 mm
H	0.0 mm
Initial Rotation	0.00 deg
Initial Orientation	Sagittal

**Geometry - Navigator****System - Miscellaneous**

Positioning mode	REF
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P

**System - Miscellaneous**

Transversal	F >> H
Coil Combine Mode	Sum of Squares
Save uncombined	Off
Matrix Optimization	Off
AutoAlign	Head > Basis
Coil Select Mode	Default

**System - Adjustments**

B0 Shim mode	Brain
B1 Shim mode	TrueForm
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

**System - Adjust Volume**

Position	Isocenter
Orientation	Sagittal
Rotation	0.00 deg
A >> P	225 mm
F >> H	240 mm
R >> L	144 mm
Reset	Off

**System - Tx/Rx**

Frequency 1H	297.162474 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

**Physio - Signal1**

1st Signal/Mode	None
TR	4300.0 ms
Concatenations	1

**Physio - Cardiac**

Magn. preparation	Non-sel. IR
TI 1	1000 ms
TI 2	3200 ms
Fat suppr.	Water excit. fast
Dark blood	Off
FoV read	240 mm
FoV phase	93.8 %
Phase resolution	100 %

**Physio - PACE**

Resp. control	Off
Concatenations	1

**Inline - Common**

Subtract	Off
Measurements	1
StdDev	Off
Save original images	On

**Inline - MIP**

MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On

**Inline - Composing**

Distortion Corr.	Off
------------------	-----

**Sequence - Part 1**

Introduction	On
Dimension	3D
Elliptical scanning	Off
Reordering	Linear
Asymmetric echo	Allowed
Flow comp.	No
Multi-slice mode	Single shot
Echo spacing	7.2 ms
Bandwidth	200 Hz/Px

**Sequence - Part 2**

RF pulse type	Fast
Gradient mode	Fast*
Excitation	Non-sel.
RF spoiling	On
Incr. Gradient spoiling	Off
Turbo factor	156

**Sequence - Nuclei**

TX/RX Nucleus	1H
TX/RX delta frequency	0 Hz
TX Nucleus	None
TX delta frequency	0 Hz
Coil elements	A32

**Sequence - Assistant**

Mode	Off
------	-----



\\CRC\protocols\studies\v5motion\mp2rage\_siemens\_ipat3\_0.69mm\_AAbasis

TA: 9:05 PM: FIX Voxel size: 0.7×0.7×0.7 mmPAT: 3 Rel. SNR: 1.00 : tfl

**Properties**

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	Off
Start measurements	Single measurement

**Routine**

Slab group	1
Slabs	1
Dist. factor	50 %
Position	R1.8 A29.6 F44.8 mm
Orientation	Sagittal
Phase enc. dir.	A >> P
AutoAlign	Head > Basis
Phase oversampling	0 %
Slice oversampling	7.1 %
Slices per slab	224
FoV read	242 mm
FoV phase	100.0 %
Slice thickness	0.69 mm
TR	5000.0 ms
TE	2.58 ms
Averages	1
Concatenations	1
Filter	Distortion Corr.(3D)
Coil elements	A32

**Contrast - Common**

TR	5000.0 ms
TE	2.58 ms
Magn. preparation	Non-sel. IR
T1 1	700 ms
T1 2	2700 ms
Flip angle 1	4.0 deg
Flip angle 2	5.0 deg
Fat suppr.	Water excit. normal
Water suppr.	None

**Contrast - Dynamic**

Averages	1
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	1
Multiple series	Each measurement

**Resolution - Common**

FoV read	242 mm
FoV phase	100.0 %
Slice thickness	0.69 mm
Base resolution	352
Phase resolution	100 %
Slice resolution	100 %
Phase partial Fourier	6/8

**Resolution - Common**

Slice partial Fourier	6/8
Interpolation	Off

**Resolution - iPAT**

PAT mode	GRAPPA
Accel. factor PE	3
Ref. lines PE	32
Accel. factor 3D	1
Reference scan mode	Integrated

**Resolution - Filter Image**

Image Filter	Off
Distortion Corr.	On
Mode	3D
Unfiltered images	On
Prescan Normalize	Off
Normalize	Off
B1 filter	Off

**Resolution - Filter Rawdata**

Raw filter	Off
Elliptical filter	Off

**Geometry - Common**

Slab group	1
Slabs	1
Dist. factor	50 %
Position	R1.8 A29.6 F44.8 mm
Orientation	Sagittal
Phase enc. dir.	A >> P
Slice oversampling	7.1 %
Slices per slab	224
FoV read	242 mm
FoV phase	100.0 %
Slice thickness	0.69 mm
TR	5000.0 ms
Multi-slice mode	Single shot
Series	Interleaved
Concatenations	1

**Geometry - AutoAlign**

Slab group	1
Position	R1.8 A29.6 F44.8 mm
Orientation	Sagittal
Phase enc. dir.	A >> P
AutoAlign	Head > Basis
Initial Position	R1.8 A29.6 F44.8
R	1.8 mm
A	29.6 mm
F	44.8 mm
Initial Rotation	0.00 deg
Initial Orientation	Sagittal

**Geometry - Navigator****System - Miscellaneous**

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T

**System - Miscellaneous**

Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
Save uncombined	Off
Matrix Optimization	Off
AutoAlign	Head > Basis
Coil Select Mode	Default

**System - Adjustments**

B0 Shim mode	Brain
B1 Shim mode	TrueForm
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

**System - Adjust Volume**

Position	R1.8 A29.6 F44.8 mm
Orientation	Sagittal
Rotation	0.00 deg
A >> P	242 mm
F >> H	242 mm
R >> L	155 mm
Reset	Off

**System - Tx/Rx**

Frequency 1H	297.162474 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	3.000
Reset	Off
? Ref. amplitude 1H	0.000 V

**Physio - Signal1**

1st Signal/Mode	None
TR	5000.0 ms
Concatenations	1

**Physio - Cardiac**

Magn. preparation	Non-sel. IR
TI 1	700 ms
TI 2	2700 ms
Fat suppr.	Water excit. normal
Dark blood	Off
FoV read	242 mm
FoV phase	100.0 %
Phase resolution	100 %

**Physio - PACE**

Resp. control	Off
Concatenations	1

**Inline - Common**

Subtract	Off
Measurements	1
StdDev	Off
Save original images	On

**Inline - MIP**

MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off

**Inline - MIP**

MIP-Time	Off
Save original images	On

**Inline - Composing**

Distortion Corr.	On
Mode	3D
Unfiltered images	On

**Sequence - Part 1**

Introduction	Off
Dimension	3D
Elliptical scanning	Off
Reordering	Linear
Asymmetric echo	Off
Flow comp.	No
Multi-slice mode	Single shot
Echo spacing	6.1 ms
Bandwidth	490 Hz/Px

**Sequence - Part 2**

RF pulse type	Fast
Gradient mode	Fast*
Excitation	Non-sel.
RF spoiling	On
Incr. Gradient spoiling	Off
Turbo factor	180

**Sequence - Nuclei**

TX/RX Nucleus	1H
TX/RX delta frequency	0 Hz
TX Nucleus	None
TX delta frequency	0 Hz
Coil elements	A32

**Sequence - Assistant**

Mode	Off
------	-----

\\CRC\protocols\studies\v5motion\cmrr\_mbep2d\_p3\_mb2\_1.6iso\_Neuro3D\_Visual\_AABrain

TA: 15:25 PM: FIX Voxel size: 1.6×1.6×1.6 mmPAT: 3 Rel. SNR: 1.00 : epfid

**Properties**

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	Off
Start measurements	Single measurement

**Routine**

Slice group	1
Slices	72
Dist. factor	0 %
Position	L1.3 A5.2 H4.5 mm
Orientation	T > C-8.0
Phase enc. dir.	A >> P
AutoAlign	Head > Brain
Phase oversampling	0 %
FoV read	216 mm
FoV phase	101.5 %
Slice thickness	1.60 mm
TR	1800 ms
TE	22.00 ms
Multi-band accel. factor	2
Filter	None
Coil elements	A32

**Contrast - Common**

TR	1800 ms
TE	22.00 ms
MTC	Off
Magn. preparation	None
Flip angle	70 deg
Fat suppr.	Fat sat.

**Contrast - Dynamic**

Averaging mode	Long term
Reconstruction	Magnitude
Measurements	500
Delay in TR	0 ms
Multiple series	Off

**Resolution - Common**

FoV read	216 mm
FoV phase	101.5 %
Slice thickness	1.60 mm
Base resolution	134
Phase resolution	100 %
Phase partial Fourier	Off
Interpolation	Off

**Resolution - iPAT**

PAT mode	GRAPPA
Accel. factor PE	3
Ref. lines PE	36
Reference scan mode	Segmented

**Resolution - Filter Image**

Distortion Corr.	Off
Prescan Normalize	Off

**Resolution - Filter Rawdata**

Raw filter	Off
Elliptical filter	Off
Hamming	Off

**Geometry - Common**

Slice group	1
Slices	72
Dist. factor	0 %
Position	L1.3 A5.2 H4.5 mm
Orientation	T > C-8.0
Phase enc. dir.	A >> P
FoV read	216 mm
FoV phase	101.5 %
Slice thickness	1.60 mm
TR	1800 ms
Multi-slice mode	Interleaved
Series	Interleaved
Multi-band accel. factor	2

**Geometry - AutoAlign**

Slice group	1
Position	L1.3 A5.2 H4.5 mm
Orientation	T > C-8.0
Phase enc. dir.	A >> P
AutoAlign	Head > Brain
Initial Position	L1.3 A5.2 H4.5
L	1.3 mm
A	5.2 mm
H	4.5 mm
Initial Rotation	-0.25 deg
Initial Orientation	T > C
T > C	-8.0
> S	0.0

**Geometry - Saturation**

Fat suppr.	Fat sat.
Special sat.	None

**System - Miscellaneous**

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
Matrix Optimization	Off
AutoAlign	Head > Brain
Coil Select Mode	Default

**System - Adjustments**

B0 Shim mode	Standard
B1 Shim mode	TrueForm
Confirm freq. adjustment	Off
Assume Dominant Fat	Off

**System - Adjustments**

Assume Silicone	Off
Adjustment Tolerance	Auto

**System - Adjust Volume**

Position	L1.3 A5.2 H4.5 mm
Orientation	T > C-8.0
Rotation	89.75 deg
R >> L	216 mm
A >> P	220 mm
F >> H	116 mm
Reset	Off

**System - Tx/Rx**

Frequency 1H	297.162474 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

**Physio - Signal1**

1st Signal/Mode	None
TR	1800 ms
Multi-band accel. factor	2

**BOLD**

GLM Statistics	On
Dynamic t-maps	On
Ignore meas. at start	5
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	1.00
Paradigm size	192
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Baseline
Meas[12]	Baseline
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Meas[16]	Active
Meas[17]	Active
Meas[18]	Active
Meas[19]	Active
Meas[20]	Active
Meas[21]	Active
Meas[22]	Active
Meas[23]	Active
Meas[24]	Active
Meas[25]	Baseline
Meas[26]	Baseline
Meas[27]	Baseline
Meas[28]	Baseline
Meas[29]	Baseline
Meas[30]	Baseline

**BOLD**

Meas[31]	Baseline
Meas[32]	Baseline
Meas[33]	Baseline
Meas[34]	Baseline
Meas[35]	Baseline
Meas[36]	Baseline
Meas[37]	Active
Meas[38]	Active
Meas[39]	Active
Meas[40]	Active
Meas[41]	Active
Meas[42]	Active
Meas[43]	Active
Meas[44]	Active
Meas[45]	Active
Meas[46]	Active
Meas[47]	Active
Meas[48]	Active
Meas[49]	Baseline
Meas[50]	Baseline
Meas[51]	Baseline
Meas[52]	Baseline
Meas[53]	Baseline
Meas[54]	Baseline
Meas[55]	Baseline
Meas[56]	Baseline
Meas[57]	Baseline
Meas[58]	Baseline
Meas[59]	Baseline
Meas[60]	Baseline
Meas[61]	Active
Meas[62]	Active
Meas[63]	Active
Meas[64]	Active
Meas[65]	Active
Meas[66]	Active
Meas[67]	Active
Meas[68]	Active
Meas[69]	Active
Meas[70]	Active
Meas[71]	Active
Meas[72]	Active
Meas[73]	Baseline
Meas[74]	Baseline
Meas[75]	Baseline
Meas[76]	Baseline
Meas[77]	Baseline
Meas[78]	Baseline
Meas[79]	Baseline
Meas[80]	Baseline
Meas[81]	Baseline
Meas[82]	Baseline
Meas[83]	Baseline
Meas[84]	Baseline
Meas[85]	Active
Meas[86]	Active
Meas[87]	Active
Meas[88]	Active
Meas[89]	Active
Meas[90]	Active
Meas[91]	Active
Meas[92]	Active
Meas[93]	Active
Meas[94]	Active
Meas[95]	Active

**BOLD**

Meas[96]	Active
Meas[97]	Baseline
Meas[98]	Baseline
Meas[99]	Baseline
Meas[100]	Baseline
Meas[101]	Baseline
Meas[102]	Baseline
Meas[103]	Baseline
Meas[104]	Baseline
Meas[105]	Baseline
Meas[106]	Baseline
Meas[107]	Baseline
Meas[108]	Baseline
Meas[109]	Active
Meas[110]	Active
Meas[111]	Active
Meas[112]	Active
Meas[113]	Active
Meas[114]	Active
Meas[115]	Active
Meas[116]	Active
Meas[117]	Active
Meas[118]	Active
Meas[119]	Active
Meas[120]	Active
Meas[121]	Baseline
Meas[122]	Baseline
Meas[123]	Baseline
Meas[124]	Baseline
Meas[125]	Baseline
Meas[126]	Baseline
Meas[127]	Baseline
Meas[128]	Baseline
Meas[129]	Baseline
Meas[130]	Baseline
Meas[131]	Baseline
Meas[132]	Baseline
Meas[133]	Active
Meas[134]	Active
Meas[135]	Active
Meas[136]	Active
Meas[137]	Active
Meas[138]	Active
Meas[139]	Active
Meas[140]	Active
Meas[141]	Active
Meas[142]	Active
Meas[143]	Active
Meas[144]	Active
Meas[145]	Baseline
Meas[146]	Baseline
Meas[147]	Baseline
Meas[148]	Baseline
Meas[149]	Baseline
Meas[150]	Baseline
Meas[151]	Baseline
Meas[152]	Baseline
Meas[153]	Baseline
Meas[154]	Baseline
Meas[155]	Baseline
Meas[156]	Baseline
Meas[157]	Active
Meas[158]	Active
Meas[159]	Active
Meas[160]	Active

**BOLD**

Meas[161]	Active
Meas[162]	Active
Meas[163]	Active
Meas[164]	Active
Meas[165]	Active
Meas[166]	Active
Meas[167]	Active
Meas[168]	Active
Meas[169]	Baseline
Meas[170]	Baseline
Meas[171]	Baseline
Meas[172]	Baseline
Meas[173]	Baseline
Meas[174]	Baseline
Meas[175]	Baseline
Meas[176]	Baseline
Meas[177]	Baseline
Meas[178]	Baseline
Meas[179]	Baseline
Meas[180]	Baseline
Meas[181]	Active
Meas[182]	Active
Meas[183]	Active
Meas[184]	Active
Meas[185]	Active
Meas[186]	Active
Meas[187]	Active
Meas[188]	Active
Meas[189]	Active
Meas[190]	Active
Meas[191]	Active
Meas[192]	Active
Motion correction	On
Spatial filter	Off
Measurements	500
Delay in TR	0 ms
Multiple series	Off

**Sequence - Part 1**

Introduction	Off
Contrasts	1
Flow comp.	No
Multi-slice mode	Interleaved
Free echo spacing	Off
Echo spacing	0.59 ms
Bandwidth	2072 Hz/Px

**Sequence - Part 2**

EPI factor	136
Gradient mode	Fast
RF spoiling	Off

**Sequence - Special**

Excite pulse duration	7200 us
Single-band images	Off
MB LeakBlock kernel	On
MB dual kernel	Off
MB RF phase scramble	Off
SENSE1 coil combine	Off
Invert RO/PE polarity	Off
Disable freq. update	Off
Force equal slice timing	Off
Online multi-band recon.	Online
FFT scale factor	1.00

**Sequence - Special**

Fat saturation FA	110.0 deg
Physio recording	DICOM
Triggering scheme	Standard

\\CRC\protocols\studies\v5motion\cmrr\_mbep2d\_p3\_mb2\_1.6iso\_Neuro3D\_Audio\_AABrain

TA: 15:25 PM: FIX Voxel size: 1.6×1.6×1.6 mmPAT: 3 Rel. SNR: 1.00 : epfid

**Properties**

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	Off
Start measurements	Single measurement

**Routine**

Slice group	1
Slices	72
Dist. factor	0 %
Position	L1.3 A5.2 H4.5 mm
Orientation	T > C-8.0
Phase enc. dir.	A >> P
AutoAlign	Head > Brain
Phase oversampling	0 %
FoV read	216 mm
FoV phase	101.5 %
Slice thickness	1.60 mm
TR	1800 ms
TE	22.00 ms
Multi-band accel. factor	2
Filter	None
Coil elements	A32

**Contrast - Common**

TR	1800 ms
TE	22.00 ms
MTC	Off
Magn. preparation	None
Flip angle	70 deg
Fat suppr.	Fat sat.

**Contrast - Dynamic**

Averaging mode	Long term
Reconstruction	Magnitude
Measurements	500
Delay in TR	0 ms
Multiple series	Off

**Resolution - Common**

FoV read	216 mm
FoV phase	101.5 %
Slice thickness	1.60 mm
Base resolution	134
Phase resolution	100 %
Phase partial Fourier	Off
Interpolation	Off

**Resolution - iPAT**

PAT mode	GRAPPA
Accel. factor PE	3
Ref. lines PE	36
Reference scan mode	Segmented

**Resolution - Filter Image**

Distortion Corr.	Off
Prescan Normalize	Off

**Resolution - Filter Rawdata**

Raw filter	Off
Elliptical filter	Off
Hamming	Off

**Geometry - Common**

Slice group	1
Slices	72
Dist. factor	0 %
Position	L1.3 A5.2 H4.5 mm
Orientation	T > C-8.0
Phase enc. dir.	A >> P
FoV read	216 mm
FoV phase	101.5 %
Slice thickness	1.60 mm
TR	1800 ms
Multi-slice mode	Interleaved
Series	Interleaved
Multi-band accel. factor	2

**Geometry - AutoAlign**

Slice group	1
Position	L1.3 A5.2 H4.5 mm
Orientation	T > C-8.0
Phase enc. dir.	A >> P
AutoAlign	Head > Brain
Initial Position	L1.3 A5.2 H4.5
L	1.3 mm
A	5.2 mm
H	4.5 mm
Initial Rotation	-0.25 deg
Initial Orientation	T > C
T > C	-8.0
> S	0.0

**Geometry - Saturation**

Fat suppr.	Fat sat.
Special sat.	None

**System - Miscellaneous**

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
Matrix Optimization	Off
AutoAlign	Head > Brain
Coil Select Mode	Default

**System - Adjustments**

B0 Shim mode	Standard
B1 Shim mode	TrueForm
Confirm freq. adjustment	Off
Assume Dominant Fat	Off

**System - Adjustments**

Assume Silicone	Off
Adjustment Tolerance	Auto

**System - Adjust Volume**

Position	L1.3 A5.2 H4.5 mm
Orientation	T > C-8.0
Rotation	89.75 deg
R >> L	216 mm
A >> P	220 mm
F >> H	116 mm
Reset	Off

**System - Tx/Rx**

Frequency 1H	297.162474 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

**Physio - Signal1**

1st Signal/Mode	None
TR	1800 ms
Multi-band accel. factor	2

**BOLD**

GLM Statistics	On
Dynamic t-maps	On
Ignore meas. at start	5
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	1.00
Paradigm size	336
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Baseline
Meas[12]	Baseline
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Meas[16]	Active
Meas[17]	Active
Meas[18]	Active
Meas[19]	Active
Meas[20]	Active
Meas[21]	Active
Meas[22]	Active
Meas[23]	Active
Meas[24]	Active
Meas[25]	Baseline
Meas[26]	Baseline
Meas[27]	Baseline
Meas[28]	Baseline
Meas[29]	Baseline
Meas[30]	Baseline

**BOLD**

Meas[31]	Baseline
Meas[32]	Baseline
Meas[33]	Baseline
Meas[34]	Baseline
Meas[35]	Baseline
Meas[36]	Baseline
Meas[37]	Active
Meas[38]	Active
Meas[39]	Active
Meas[40]	Active
Meas[41]	Active
Meas[42]	Active
Meas[43]	Active
Meas[44]	Active
Meas[45]	Active
Meas[46]	Active
Meas[47]	Active
Meas[48]	Active
Meas[49]	Baseline
Meas[50]	Baseline
Meas[51]	Baseline
Meas[52]	Baseline
Meas[53]	Baseline
Meas[54]	Baseline
Meas[55]	Baseline
Meas[56]	Baseline
Meas[57]	Baseline
Meas[58]	Baseline
Meas[59]	Baseline
Meas[60]	Baseline
Meas[61]	Active
Meas[62]	Active
Meas[63]	Active
Meas[64]	Active
Meas[65]	Active
Meas[66]	Active
Meas[67]	Active
Meas[68]	Active
Meas[69]	Active
Meas[70]	Active
Meas[71]	Active
Meas[72]	Active
Meas[73]	Baseline
Meas[74]	Baseline
Meas[75]	Baseline
Meas[76]	Baseline
Meas[77]	Baseline
Meas[78]	Baseline
Meas[79]	Baseline
Meas[80]	Baseline
Meas[81]	Baseline
Meas[82]	Baseline
Meas[83]	Baseline
Meas[84]	Baseline
Meas[85]	Active
Meas[86]	Active
Meas[87]	Active
Meas[88]	Active
Meas[89]	Active
Meas[90]	Active
Meas[91]	Active
Meas[92]	Active
Meas[93]	Active
Meas[94]	Active
Meas[95]	Active



## SIEMENS MAGNETOM Terra

**BOLD**

Meas[96]	Active
Meas[97]	Baseline
Meas[98]	Baseline
Meas[99]	Baseline
Meas[100]	Baseline
Meas[101]	Baseline
Meas[102]	Baseline
Meas[103]	Baseline
Meas[104]	Baseline
Meas[105]	Baseline
Meas[106]	Baseline
Meas[107]	Baseline
Meas[108]	Baseline
Meas[109]	Active
Meas[110]	Active
Meas[111]	Active
Meas[112]	Active
Meas[113]	Active
Meas[114]	Active
Meas[115]	Active
Meas[116]	Active
Meas[117]	Active
Meas[118]	Active
Meas[119]	Active
Meas[120]	Active
Meas[121]	Baseline
Meas[122]	Baseline
Meas[123]	Baseline
Meas[124]	Baseline
Meas[125]	Baseline
Meas[126]	Baseline
Meas[127]	Baseline
Meas[128]	Baseline
Meas[129]	Baseline
Meas[130]	Baseline
Meas[131]	Baseline
Meas[132]	Baseline
Meas[133]	Active
Meas[134]	Active
Meas[135]	Active
Meas[136]	Active
Meas[137]	Active
Meas[138]	Active
Meas[139]	Active
Meas[140]	Active
Meas[141]	Active
Meas[142]	Active
Meas[143]	Active
Meas[144]	Active
Meas[145]	Baseline
Meas[146]	Baseline
Meas[147]	Baseline
Meas[148]	Baseline
Meas[149]	Baseline
Meas[150]	Baseline
Meas[151]	Baseline
Meas[152]	Baseline
Meas[153]	Baseline
Meas[154]	Baseline
Meas[155]	Baseline
Meas[156]	Baseline
Meas[157]	Active
Meas[158]	Active
Meas[159]	Active
Meas[160]	Active

**BOLD**

Meas[161]	Active
Meas[162]	Active
Meas[163]	Active
Meas[164]	Active
Meas[165]	Active
Meas[166]	Active
Meas[167]	Active
Meas[168]	Active
Meas[169]	Baseline
Meas[170]	Baseline
Meas[171]	Baseline
Meas[172]	Baseline
Meas[173]	Baseline
Meas[174]	Baseline
Meas[175]	Baseline
Meas[176]	Baseline
Meas[177]	Baseline
Meas[178]	Baseline
Meas[179]	Baseline
Meas[180]	Baseline
Meas[181]	Active
Meas[182]	Active
Meas[183]	Active
Meas[184]	Active
Meas[185]	Active
Meas[186]	Active
Meas[187]	Active
Meas[188]	Active
Meas[189]	Active
Meas[190]	Active
Meas[191]	Active
Meas[192]	Active
Meas[193]	Baseline
Meas[194]	Baseline
Meas[195]	Baseline
Meas[196]	Baseline
Meas[197]	Baseline
Meas[198]	Baseline
Meas[199]	Baseline
Meas[200]	Baseline
Meas[201]	Baseline
Meas[202]	Baseline
Meas[203]	Baseline
Meas[204]	Baseline
Meas[205]	Active
Meas[206]	Active
Meas[207]	Active
Meas[208]	Active
Meas[209]	Active
Meas[210]	Active
Meas[211]	Active
Meas[212]	Active
Meas[213]	Active
Meas[214]	Active
Meas[215]	Active
Meas[216]	Active
Meas[217]	Baseline
Meas[218]	Baseline
Meas[219]	Baseline
Meas[220]	Baseline
Meas[221]	Baseline
Meas[222]	Baseline
Meas[223]	Baseline
Meas[224]	Baseline
Meas[225]	Baseline

**BOLD**

Meas[226]	Baseline
Meas[227]	Baseline
Meas[228]	Baseline
Meas[229]	Active
Meas[230]	Active
Meas[231]	Active
Meas[232]	Active
Meas[233]	Active
Meas[234]	Active
Meas[235]	Active
Meas[236]	Active
Meas[237]	Active
Meas[238]	Active
Meas[239]	Active
Meas[240]	Active
Meas[241]	Baseline
Meas[242]	Baseline
Meas[243]	Baseline
Meas[244]	Baseline
Meas[245]	Baseline
Meas[246]	Baseline
Meas[247]	Baseline
Meas[248]	Baseline
Meas[249]	Baseline
Meas[250]	Baseline
Meas[251]	Baseline
Meas[252]	Baseline
Meas[253]	Active
Meas[254]	Active
Meas[255]	Active
Meas[256]	Active
Meas[257]	Active
Meas[258]	Active
Meas[259]	Active
Meas[260]	Active
Meas[261]	Active
Meas[262]	Active
Meas[263]	Active
Meas[264]	Active
Meas[265]	Baseline
Meas[266]	Baseline
Meas[267]	Baseline
Meas[268]	Baseline
Meas[269]	Baseline
Meas[270]	Baseline
Meas[271]	Baseline
Meas[272]	Baseline
Meas[273]	Baseline
Meas[274]	Baseline
Meas[275]	Baseline
Meas[276]	Baseline
Meas[277]	Active
Meas[278]	Active
Meas[279]	Active
Meas[280]	Active
Meas[281]	Active
Meas[282]	Active
Meas[283]	Active
Meas[284]	Active
Meas[285]	Active
Meas[286]	Active
Meas[287]	Active
Meas[288]	Active
Meas[289]	Baseline
Meas[290]	Baseline

**BOLD**

Meas[291]	Baseline
Meas[292]	Baseline
Meas[293]	Baseline
Meas[294]	Baseline
Meas[295]	Baseline
Meas[296]	Baseline
Meas[297]	Baseline
Meas[298]	Baseline
Meas[299]	Baseline
Meas[300]	Baseline
Meas[301]	Active
Meas[302]	Active
Meas[303]	Active
Meas[304]	Active
Meas[305]	Active
Meas[306]	Active
Meas[307]	Active
Meas[308]	Active
Meas[309]	Active
Meas[310]	Active
Meas[311]	Active
Meas[312]	Active
Meas[313]	Baseline
Meas[314]	Baseline
Meas[315]	Baseline
Meas[316]	Baseline
Meas[317]	Baseline
Meas[318]	Baseline
Meas[319]	Baseline
Meas[320]	Baseline
Meas[321]	Baseline
Meas[322]	Baseline
Meas[323]	Baseline
Meas[324]	Baseline
Meas[325]	Active
Meas[326]	Active
Meas[327]	Active
Meas[328]	Active
Meas[329]	Active
Meas[330]	Active
Meas[331]	Active
Meas[332]	Active
Meas[333]	Active
Meas[334]	Active
Meas[335]	Active
Meas[336]	Active
Motion correction	On
Spatial filter	Off
Measurements	500
Delay in TR	0 ms
Multiple series	Off

**Sequence - Part 1**

Introduction	Off
Contrasts	1
Flow comp.	No
Multi-slice mode	Interleaved
Free echo spacing	Off
Echo spacing	0.59 ms
Bandwidth	2072 Hz/Px

**Sequence - Part 2**

EPI factor	136
Gradient mode	Fast

**Sequence - Part 2**

RF spoiling	Off
-------------	-----

**Sequence - Special**

Excite pulse duration	7200 us
Single-band images	Off
MB LeakBlock kernel	On
MB dual kernel	Off
MB RF phase scramble	Off
SENSE1 coil combine	Off
Invert RO/PE polarity	Off
Disable freq. update	Off
Force equal slice timing	Off
Online multi-band recon.	Online
FFT scale factor	1.00
Fat saturation FA	110.0 deg
Physio recording	DICOM
Triggering scheme	Standard

\\CRC\protocols\studies\v5motion\cmrr\_mbep2d\_p3\_mb2\_1.6iso\_Neuro3D\_Visual\_v2\_AABrain

TA: 15:25 PM: FIX Voxel size: 1.6×1.6×1.6 mmPAT: 3 Rel. SNR: 1.00 : epfid

**Properties**

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	Off
Start measurements	Single measurement

**Routine**

Slice group	1
Slices	72
Dist. factor	0 %
Position	L1.3 A5.2 H4.5 mm
Orientation	T > C-8.0
Phase enc. dir.	A >> P
AutoAlign	Head > Brain
Phase oversampling	0 %
FoV read	216 mm
FoV phase	101.5 %
Slice thickness	1.60 mm
TR	1800 ms
TE	22.00 ms
Multi-band accel. factor	2
Filter	None
Coil elements	A32

**Contrast - Common**

TR	1800 ms
TE	22.00 ms
MTC	Off
Magn. preparation	None
Flip angle	70 deg
Fat suppr.	Fat sat.

**Contrast - Dynamic**

Averaging mode	Long term
Reconstruction	Magnitude
Measurements	500
Delay in TR	0 ms
Multiple series	Off

**Resolution - Common**

FoV read	216 mm
FoV phase	101.5 %
Slice thickness	1.60 mm
Base resolution	134
Phase resolution	100 %
Phase partial Fourier	Off
Interpolation	Off

**Resolution - iPAT**

PAT mode	GRAPPA
Accel. factor PE	3
Ref. lines PE	36
Reference scan mode	Segmented

**Resolution - Filter Image**

Distortion Corr.	Off
Prescan Normalize	Off

**Resolution - Filter Rawdata**

Raw filter	Off
Elliptical filter	Off
Hamming	Off

**Geometry - Common**

Slice group	1
Slices	72
Dist. factor	0 %
Position	L1.3 A5.2 H4.5 mm
Orientation	T > C-8.0
Phase enc. dir.	A >> P
FoV read	216 mm
FoV phase	101.5 %
Slice thickness	1.60 mm
TR	1800 ms
Multi-slice mode	Interleaved
Series	Interleaved
Multi-band accel. factor	2

**Geometry - AutoAlign**

Slice group	1
Position	L1.3 A5.2 H4.5 mm
Orientation	T > C-8.0
Phase enc. dir.	A >> P
AutoAlign	Head > Brain
Initial Position	L1.3 A5.2 H4.5
L	1.3 mm
A	5.2 mm
H	4.5 mm
Initial Rotation	-0.25 deg
Initial Orientation	T > C
T > C	-8.0
> S	0.0

**Geometry - Saturation**

Fat suppr.	Fat sat.
Special sat.	None

**System - Miscellaneous**

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
Matrix Optimization	Off
AutoAlign	Head > Brain
Coil Select Mode	Default

**System - Adjustments**

B0 Shim mode	Standard
B1 Shim mode	TrueForm
Confirm freq. adjustment	Off
Assume Dominant Fat	Off

**System - Adjustments**

Assume Silicone	Off
Adjustment Tolerance	Auto

**System - Adjust Volume**

Position	L1.3 A5.2 H4.5 mm
Orientation	T > C-8.0
Rotation	89.75 deg
R >> L	216 mm
A >> P	220 mm
F >> H	116 mm
Reset	Off

**System - Tx/Rx**

Frequency 1H	297.162474 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

**Physio - Signal1**

1st Signal/Mode	None
TR	1800 ms
Multi-band accel. factor	2

**BOLD**

GLM Statistics	On
Dynamic t-maps	On
Ignore meas. at start	5
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	1.00
Paradigm size	270
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Active
Meas[11]	Active
Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Meas[16]	Active
Meas[17]	Active
Meas[18]	Active
Meas[19]	Baseline
Meas[20]	Baseline
Meas[21]	Baseline
Meas[22]	Baseline
Meas[23]	Baseline
Meas[24]	Baseline
Meas[25]	Baseline
Meas[26]	Baseline
Meas[27]	Baseline
Meas[28]	Active
Meas[29]	Active
Meas[30]	Active

**BOLD**

Meas[31]	Active
Meas[32]	Active
Meas[33]	Active
Meas[34]	Active
Meas[35]	Active
Meas[36]	Active
Meas[37]	Baseline
Meas[38]	Baseline
Meas[39]	Baseline
Meas[40]	Baseline
Meas[41]	Baseline
Meas[42]	Baseline
Meas[43]	Baseline
Meas[44]	Baseline
Meas[45]	Baseline
Meas[46]	Active
Meas[47]	Active
Meas[48]	Active
Meas[49]	Active
Meas[50]	Active
Meas[51]	Active
Meas[52]	Active
Meas[53]	Active
Meas[54]	Active
Meas[55]	Baseline
Meas[56]	Baseline
Meas[57]	Baseline
Meas[58]	Baseline
Meas[59]	Baseline
Meas[60]	Baseline
Meas[61]	Baseline
Meas[62]	Baseline
Meas[63]	Baseline
Meas[64]	Active
Meas[65]	Active
Meas[66]	Active
Meas[67]	Active
Meas[68]	Active
Meas[69]	Active
Meas[70]	Active
Meas[71]	Active
Meas[72]	Active
Meas[73]	Baseline
Meas[74]	Baseline
Meas[75]	Baseline
Meas[76]	Baseline
Meas[77]	Baseline
Meas[78]	Baseline
Meas[79]	Baseline
Meas[80]	Baseline
Meas[81]	Baseline
Meas[82]	Active
Meas[83]	Active
Meas[84]	Active
Meas[85]	Active
Meas[86]	Active
Meas[87]	Active
Meas[88]	Active
Meas[89]	Active
Meas[90]	Active
Meas[91]	Baseline
Meas[92]	Baseline
Meas[93]	Baseline
Meas[94]	Baseline
Meas[95]	Baseline

## SIEMENS MAGNETOM Terra

**BOLD**

Meas[96]	Baseline
Meas[97]	Baseline
Meas[98]	Baseline
Meas[99]	Baseline
Meas[100]	Active
Meas[101]	Active
Meas[102]	Active
Meas[103]	Active
Meas[104]	Active
Meas[105]	Active
Meas[106]	Active
Meas[107]	Active
Meas[108]	Active
Meas[109]	Baseline
Meas[110]	Baseline
Meas[111]	Baseline
Meas[112]	Baseline
Meas[113]	Baseline
Meas[114]	Baseline
Meas[115]	Baseline
Meas[116]	Baseline
Meas[117]	Baseline
Meas[118]	Active
Meas[119]	Active
Meas[120]	Active
Meas[121]	Active
Meas[122]	Active
Meas[123]	Active
Meas[124]	Active
Meas[125]	Active
Meas[126]	Active
Meas[127]	Baseline
Meas[128]	Baseline
Meas[129]	Baseline
Meas[130]	Baseline
Meas[131]	Baseline
Meas[132]	Baseline
Meas[133]	Baseline
Meas[134]	Baseline
Meas[135]	Baseline
Meas[136]	Active
Meas[137]	Active
Meas[138]	Active
Meas[139]	Active
Meas[140]	Active
Meas[141]	Active
Meas[142]	Active
Meas[143]	Active
Meas[144]	Active
Meas[145]	Baseline
Meas[146]	Baseline
Meas[147]	Baseline
Meas[148]	Baseline
Meas[149]	Baseline
Meas[150]	Baseline
Meas[151]	Baseline
Meas[152]	Baseline
Meas[153]	Baseline
Meas[154]	Active
Meas[155]	Active
Meas[156]	Active
Meas[157]	Active
Meas[158]	Active
Meas[159]	Active
Meas[160]	Active

**BOLD**

Meas[161]	Active
Meas[162]	Active
Meas[163]	Baseline
Meas[164]	Baseline
Meas[165]	Baseline
Meas[166]	Baseline
Meas[167]	Baseline
Meas[168]	Baseline
Meas[169]	Baseline
Meas[170]	Baseline
Meas[171]	Baseline
Meas[172]	Active
Meas[173]	Active
Meas[174]	Active
Meas[175]	Active
Meas[176]	Active
Meas[177]	Active
Meas[178]	Active
Meas[179]	Active
Meas[180]	Active
Meas[181]	Baseline
Meas[182]	Baseline
Meas[183]	Baseline
Meas[184]	Baseline
Meas[185]	Baseline
Meas[186]	Baseline
Meas[187]	Baseline
Meas[188]	Baseline
Meas[189]	Baseline
Meas[190]	Active
Meas[191]	Active
Meas[192]	Active
Meas[193]	Active
Meas[194]	Active
Meas[195]	Active
Meas[196]	Active
Meas[197]	Active
Meas[198]	Active
Meas[199]	Baseline
Meas[200]	Baseline
Meas[201]	Baseline
Meas[202]	Baseline
Meas[203]	Baseline
Meas[204]	Baseline
Meas[205]	Baseline
Meas[206]	Baseline
Meas[207]	Baseline
Meas[208]	Active
Meas[209]	Active
Meas[210]	Active
Meas[211]	Active
Meas[212]	Active
Meas[213]	Active
Meas[214]	Active
Meas[215]	Active
Meas[216]	Active
Meas[217]	Baseline
Meas[218]	Baseline
Meas[219]	Baseline
Meas[220]	Baseline
Meas[221]	Baseline
Meas[222]	Baseline
Meas[223]	Baseline
Meas[224]	Baseline
Meas[225]	Baseline

**BOLD**

Meas[226]	Active
Meas[227]	Active
Meas[228]	Active
Meas[229]	Active
Meas[230]	Active
Meas[231]	Active
Meas[232]	Active
Meas[233]	Active
Meas[234]	Active
Meas[235]	Baseline
Meas[236]	Baseline
Meas[237]	Baseline
Meas[238]	Baseline
Meas[239]	Baseline
Meas[240]	Baseline
Meas[241]	Baseline
Meas[242]	Baseline
Meas[243]	Baseline
Meas[244]	Active
Meas[245]	Active
Meas[246]	Active
Meas[247]	Active
Meas[248]	Active
Meas[249]	Active
Meas[250]	Active
Meas[251]	Active
Meas[252]	Active
Meas[253]	Baseline
Meas[254]	Baseline
Meas[255]	Baseline
Meas[256]	Baseline
Meas[257]	Baseline
Meas[258]	Baseline
Meas[259]	Baseline
Meas[260]	Baseline
Meas[261]	Baseline
Meas[262]	Active
Meas[263]	Active
Meas[264]	Active
Meas[265]	Active
Meas[266]	Active
Meas[267]	Active
Meas[268]	Active
Meas[269]	Active
Meas[270]	Active
Motion correction	On
Spatial filter	Off
Measurements	500
Delay in TR	0 ms
Multiple series	Off

**Sequence - Special**

Excite pulse duration	7200 us
Single-band images	Off
MB LeakBlock kernel	On
MB dual kernel	Off
MB RF phase scramble	Off
SENSE1 coil combine	Off
Invert RO/PE polarity	Off
Disable freq. update	Off
Force equal slice timing	Off
Online multi-band recon.	Online
FFT scale factor	1.00
Fat saturation FA	110.0 deg
Physio recording	DICOM
Triggering scheme	Standard

**Sequence - Part 1**

Introduction	Off
Contrasts	1
Flow comp.	No
Multi-slice mode	Interleaved
Free echo spacing	Off
Echo spacing	0.59 ms
Bandwidth	2072 Hz/Px

**Sequence - Part 2**

EPI factor	136
Gradient mode	Fast
RF spoiling	Off

\\CRC\protocols\studies\v5motion\cmrr\_mbep2d\_p3\_mb2\_1.6iso\_Neuro3D\_Audio\_v2\_AABrain

TA: 15:25 PM: FIX Voxel size: 1.6×1.6×1.6 mmPAT: 3 Rel. SNR: 1.00 : epfid

**Properties**

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	Off
Start measurements	Single measurement

**Routine**

Slice group	1
Slices	72
Dist. factor	0 %
Position	L1.3 A5.2 H4.5 mm
Orientation	T > C-8.0
Phase enc. dir.	A >> P
AutoAlign	Head > Brain
Phase oversampling	0 %
FoV read	216 mm
FoV phase	101.5 %
Slice thickness	1.60 mm
TR	1800 ms
TE	22.00 ms
Multi-band accel. factor	2
Filter	None
Coil elements	A32

**Contrast - Common**

TR	1800 ms
TE	22.00 ms
MTC	Off
Magn. preparation	None
Flip angle	70 deg
Fat suppr.	Fat sat.

**Contrast - Dynamic**

Averaging mode	Long term
Reconstruction	Magnitude
Measurements	500
Delay in TR	0 ms
Multiple series	Off

**Resolution - Common**

FoV read	216 mm
FoV phase	101.5 %
Slice thickness	1.60 mm
Base resolution	134
Phase resolution	100 %
Phase partial Fourier	Off
Interpolation	Off

**Resolution - iPAT**

PAT mode	GRAPPA
Accel. factor PE	3
Ref. lines PE	36
Reference scan mode	Segmented

**Resolution - Filter Image**

Distortion Corr.	Off
Prescan Normalize	Off

**Resolution - Filter Rawdata**

Raw filter	Off
Elliptical filter	Off
Hamming	Off

**Geometry - Common**

Slice group	1
Slices	72
Dist. factor	0 %
Position	L1.3 A5.2 H4.5 mm
Orientation	T > C-8.0
Phase enc. dir.	A >> P
FoV read	216 mm
FoV phase	101.5 %
Slice thickness	1.60 mm
TR	1800 ms
Multi-slice mode	Interleaved
Series	Interleaved
Multi-band accel. factor	2

**Geometry - AutoAlign**

Slice group	1
Position	L1.3 A5.2 H4.5 mm
Orientation	T > C-8.0
Phase enc. dir.	A >> P
AutoAlign	Head > Brain
Initial Position	L1.3 A5.2 H4.5
L	1.3 mm
A	5.2 mm
H	4.5 mm
Initial Rotation	-0.25 deg
Initial Orientation	T > C
T > C	-8.0
> S	0.0

**Geometry - Saturation**

Fat suppr.	Fat sat.
Special sat.	None

**System - Miscellaneous**

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
Matrix Optimization	Off
AutoAlign	Head > Brain
Coil Select Mode	Default

**System - Adjustments**

B0 Shim mode	Standard
B1 Shim mode	TrueForm
Confirm freq. adjustment	Off
Assume Dominant Fat	Off



**System - Adjustments**

Assume Silicone	Off
Adjustment Tolerance	Auto

**System - Adjust Volume**

Position	L1.3 A5.2 H4.5 mm
Orientation	T > C-8.0
Rotation	89.75 deg
R >> L	216 mm
A >> P	220 mm
F >> H	116 mm
Reset	Off

**System - Tx/Rx**

Frequency 1H	297.162474 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

**Physio - Signal1**

1st Signal/Mode	None
TR	1800 ms
Multi-band accel. factor	2

**BOLD**

GLM Statistics	On
Dynamic t-maps	On
Ignore meas. at start	5
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	1.00
Paradigm size	400
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Active
Meas[11]	Active
Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Meas[16]	Active
Meas[17]	Active
Meas[18]	Active
Meas[19]	Baseline
Meas[20]	Baseline
Meas[21]	Baseline
Meas[22]	Baseline
Meas[23]	Baseline
Meas[24]	Baseline
Meas[25]	Baseline
Meas[26]	Baseline
Meas[27]	Baseline
Meas[28]	Active
Meas[29]	Active
Meas[30]	Active

**BOLD**

Meas[31]	Active
Meas[32]	Active
Meas[33]	Active
Meas[34]	Active
Meas[35]	Active
Meas[36]	Active
Meas[37]	Baseline
Meas[38]	Baseline
Meas[39]	Baseline
Meas[40]	Baseline
Meas[41]	Baseline
Meas[42]	Baseline
Meas[43]	Baseline
Meas[44]	Baseline
Meas[45]	Baseline
Meas[46]	Active
Meas[47]	Active
Meas[48]	Active
Meas[49]	Active
Meas[50]	Active
Meas[51]	Active
Meas[52]	Active
Meas[53]	Active
Meas[54]	Active
Meas[55]	Baseline
Meas[56]	Baseline
Meas[57]	Baseline
Meas[58]	Baseline
Meas[59]	Baseline
Meas[60]	Baseline
Meas[61]	Baseline
Meas[62]	Baseline
Meas[63]	Baseline
Meas[64]	Active
Meas[65]	Active
Meas[66]	Active
Meas[67]	Active
Meas[68]	Active
Meas[69]	Active
Meas[70]	Active
Meas[71]	Active
Meas[72]	Active
Meas[73]	Baseline
Meas[74]	Baseline
Meas[75]	Baseline
Meas[76]	Baseline
Meas[77]	Baseline
Meas[78]	Baseline
Meas[79]	Baseline
Meas[80]	Baseline
Meas[81]	Baseline
Meas[82]	Active
Meas[83]	Active
Meas[84]	Active
Meas[85]	Active
Meas[86]	Active
Meas[87]	Active
Meas[88]	Active
Meas[89]	Active
Meas[90]	Active
Meas[91]	Baseline
Meas[92]	Baseline
Meas[93]	Baseline
Meas[94]	Baseline
Meas[95]	Baseline

# SIEMENS MAGNETOM Terra

## BOLD

Meas[96]	Baseline
Meas[97]	Baseline
Meas[98]	Baseline
Meas[99]	Baseline
Meas[100]	Active
Meas[101]	Active
Meas[102]	Active
Meas[103]	Active
Meas[104]	Active
Meas[105]	Active
Meas[106]	Active
Meas[107]	Active
Meas[108]	Active
Meas[109]	Baseline
Meas[110]	Baseline
Meas[111]	Baseline
Meas[112]	Baseline
Meas[113]	Baseline
Meas[114]	Baseline
Meas[115]	Baseline
Meas[116]	Baseline
Meas[117]	Baseline
Meas[118]	Active
Meas[119]	Active
Meas[120]	Active
Meas[121]	Active
Meas[122]	Active
Meas[123]	Active
Meas[124]	Active
Meas[125]	Active
Meas[126]	Active
Meas[127]	Baseline
Meas[128]	Baseline
Meas[129]	Baseline
Meas[130]	Baseline
Meas[131]	Baseline
Meas[132]	Baseline
Meas[133]	Baseline
Meas[134]	Baseline
Meas[135]	Baseline
Meas[136]	Active
Meas[137]	Active
Meas[138]	Active
Meas[139]	Active
Meas[140]	Active
Meas[141]	Active
Meas[142]	Active
Meas[143]	Active
Meas[144]	Active
Meas[145]	Baseline
Meas[146]	Baseline
Meas[147]	Baseline
Meas[148]	Baseline
Meas[149]	Baseline
Meas[150]	Baseline
Meas[151]	Baseline
Meas[152]	Baseline
Meas[153]	Baseline
Meas[154]	Active
Meas[155]	Active
Meas[156]	Active
Meas[157]	Active
Meas[158]	Active
Meas[159]	Active
Meas[160]	Active

## BOLD

Meas[161]	Active
Meas[162]	Active
Meas[163]	Baseline
Meas[164]	Baseline
Meas[165]	Baseline
Meas[166]	Baseline
Meas[167]	Baseline
Meas[168]	Baseline
Meas[169]	Baseline
Meas[170]	Baseline
Meas[171]	Baseline
Meas[172]	Active
Meas[173]	Active
Meas[174]	Active
Meas[175]	Active
Meas[176]	Active
Meas[177]	Active
Meas[178]	Active
Meas[179]	Active
Meas[180]	Active
Meas[181]	Baseline
Meas[182]	Baseline
Meas[183]	Baseline
Meas[184]	Baseline
Meas[185]	Baseline
Meas[186]	Baseline
Meas[187]	Baseline
Meas[188]	Baseline
Meas[189]	Baseline
Meas[190]	Active
Meas[191]	Active
Meas[192]	Active
Meas[193]	Active
Meas[194]	Active
Meas[195]	Active
Meas[196]	Active
Meas[197]	Active
Meas[198]	Active
Meas[199]	Baseline
Meas[200]	Baseline
Meas[201]	Baseline
Meas[202]	Baseline
Meas[203]	Baseline
Meas[204]	Baseline
Meas[205]	Baseline
Meas[206]	Baseline
Meas[207]	Baseline
Meas[208]	Active
Meas[209]	Active
Meas[210]	Active
Meas[211]	Active
Meas[212]	Active
Meas[213]	Active
Meas[214]	Active
Meas[215]	Active
Meas[216]	Active
Meas[217]	Baseline
Meas[218]	Baseline
Meas[219]	Baseline
Meas[220]	Baseline
Meas[221]	Baseline
Meas[222]	Baseline
Meas[223]	Baseline
Meas[224]	Baseline
Meas[225]	Baseline

# SIEMENS MAGNETOM Terra

## BOLD

Meas[226]	Active
Meas[227]	Active
Meas[228]	Active
Meas[229]	Active
Meas[230]	Active
Meas[231]	Active
Meas[232]	Active
Meas[233]	Active
Meas[234]	Active
Meas[235]	Baseline
Meas[236]	Baseline
Meas[237]	Baseline
Meas[238]	Baseline
Meas[239]	Baseline
Meas[240]	Baseline
Meas[241]	Baseline
Meas[242]	Baseline
Meas[243]	Baseline
Meas[244]	Active
Meas[245]	Active
Meas[246]	Active
Meas[247]	Active
Meas[248]	Active
Meas[249]	Active
Meas[250]	Active
Meas[251]	Active
Meas[252]	Active
Meas[253]	Baseline
Meas[254]	Baseline
Meas[255]	Baseline
Meas[256]	Baseline
Meas[257]	Baseline
Meas[258]	Baseline
Meas[259]	Baseline
Meas[260]	Baseline
Meas[261]	Baseline
Meas[262]	Active
Meas[263]	Active
Meas[264]	Active
Meas[265]	Active
Meas[266]	Active
Meas[267]	Active
Meas[268]	Active
Meas[269]	Active
Meas[270]	Active
Meas[271]	Baseline
Meas[272]	Baseline
Meas[273]	Baseline
Meas[274]	Baseline
Meas[275]	Baseline
Meas[276]	Baseline
Meas[277]	Baseline
Meas[278]	Baseline
Meas[279]	Baseline
Meas[280]	Active
Meas[281]	Active
Meas[282]	Active
Meas[283]	Active
Meas[284]	Active
Meas[285]	Active
Meas[286]	Active
Meas[287]	Active
Meas[288]	Active
Meas[289]	Baseline
Meas[290]	Baseline

## BOLD

Meas[291]	Baseline
Meas[292]	Baseline
Meas[293]	Baseline
Meas[294]	Baseline
Meas[295]	Baseline
Meas[296]	Baseline
Meas[297]	Baseline
Meas[298]	Active
Meas[299]	Active
Meas[300]	Active
Meas[301]	Active
Meas[302]	Active
Meas[303]	Active
Meas[304]	Active
Meas[305]	Active
Meas[306]	Active
Meas[307]	Baseline
Meas[308]	Baseline
Meas[309]	Baseline
Meas[310]	Baseline
Meas[311]	Baseline
Meas[312]	Baseline
Meas[313]	Baseline
Meas[314]	Baseline
Meas[315]	Baseline
Meas[316]	Active
Meas[317]	Active
Meas[318]	Active
Meas[319]	Active
Meas[320]	Active
Meas[321]	Active
Meas[322]	Active
Meas[323]	Active
Meas[324]	Active
Meas[325]	Baseline
Meas[326]	Baseline
Meas[327]	Baseline
Meas[328]	Baseline
Meas[329]	Baseline
Meas[330]	Baseline
Meas[331]	Baseline
Meas[332]	Baseline
Meas[333]	Baseline
Meas[334]	Active
Meas[335]	Active
Meas[336]	Active
Meas[337]	Active
Meas[338]	Active
Meas[339]	Active
Meas[340]	Active
Meas[341]	Active
Meas[342]	Active
Meas[343]	Baseline
Meas[344]	Baseline
Meas[345]	Baseline
Meas[346]	Baseline
Meas[347]	Baseline
Meas[348]	Baseline
Meas[349]	Baseline
Meas[350]	Baseline
Meas[351]	Baseline
Meas[352]	Active
Meas[353]	Active
Meas[354]	Active
Meas[355]	Active

**BOLD**

Meas[356]	Active
Meas[357]	Active
Meas[358]	Active
Meas[359]	Active
Meas[360]	Active
Meas[361]	Active
Meas[362]	Active
Meas[363]	Active
Meas[364]	Active
Meas[365]	Active
Meas[366]	Active
Meas[367]	Active
Meas[368]	Active
Meas[369]	Active
Meas[370]	Active
Meas[371]	Active
Meas[372]	Active
Meas[373]	Active
Meas[374]	Active
Meas[375]	Active
Meas[376]	Active
Meas[377]	Active
Meas[378]	Active
Meas[379]	Baseline
Meas[380]	Baseline
Meas[381]	Baseline
Meas[382]	Baseline
Meas[383]	Baseline
Meas[384]	Baseline
Meas[385]	Baseline
Meas[386]	Baseline
Meas[387]	Baseline
Meas[388]	Active
Meas[389]	Active
Meas[390]	Active
Meas[391]	Active
Meas[392]	Active
Meas[393]	Active
Meas[394]	Active
Meas[395]	Active
Meas[396]	Active
Meas[397]	Baseline
Meas[398]	Baseline
Meas[399]	Baseline
Meas[400]	Baseline
Motion correction	On
Spatial filter	Off
Measurements	500
Delay in TR	0 ms
Multiple series	Off

**Sequence - Special**

Excite pulse duration	7200 us
Single-band images	Off
MB LeakBlock kernel	On
MB dual kernel	Off
MB RF phase scramble	Off
SENSE1 coil combine	Off
Invert RO/PE polarity	Off
Disable freq. update	Off
Force equal slice timing	Off
Online multi-band recon.	Online
FFT scale factor	1.00
Fat saturation FA	110.0 deg
Physio recording	DICOM
Triggering scheme	Standard

**Sequence - Part 1**

Introduction	Off
Contrasts	1
Flow comp.	No
Multi-slice mode	Interleaved
Free echo spacing	Off
Echo spacing	0.59 ms
Bandwidth	2072 Hz/Px

**Sequence - Part 2**

EPI factor	136
Gradient mode	Fast
RF spoiling	Off

\\CRC\protocols\studies\v5motion\cmrr\_mbep2d\_p3\_mb2\_1iso\_AAbrain

TA: 10:42 PM: REF Voxel size: 1.0×1.0×1.0 mmPAT: 3 Rel. SNR: 1.00 : epdf

**Properties**

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	Off
Start measurements	Single measurement

**Routine**

Slice group	1
Slices	60
Dist. factor	0 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	Head > Brain
Phase oversampling	0 %
FoV read	220 mm
FoV phase	98.2 %
Slice thickness	1.00 mm
TR	2220 ms
TE	30.20 ms
Multi-band accel. factor	2
Filter	None
Coil elements	A32

**Contrast - Common**

TR	2220 ms
TE	30.20 ms
MTC	Off
Magn. preparation	None
Flip angle	90 deg
Fat suppr.	Fat sat.

**Contrast - Dynamic**

Averaging mode	Long term
Reconstruction	Magnitude
Measurements	275
Delay in TR	0 ms
Multiple series	Off

**Resolution - Common**

FoV read	220 mm
FoV phase	98.2 %
Slice thickness	1.00 mm
Base resolution	220
Phase resolution	100 %
Phase partial Fourier	7/8
Interpolation	Off

**Resolution - iPAT**

PAT mode	GRAPPA
Accel. factor PE	3
Ref. lines PE	36
Reference scan mode	Segmented

**Resolution - Filter Image**

Distortion Corr.	Off
Prescan Normalize	Off

**Resolution - Filter Rawdata**

Raw filter	Off
Elliptical filter	Off
Hamming	Off

**Geometry - Common**

Slice group	1
Slices	60
Dist. factor	0 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
FoV read	220 mm
FoV phase	98.2 %
Slice thickness	1.00 mm
TR	2220 ms
Multi-slice mode	Interleaved
Series	Interleaved
Multi-band accel. factor	2

**Geometry - AutoAlign**

Slice group	1
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	Head > Brain
Initial Position	Isocenter
L	0.0 mm
P	0.0 mm
H	0.0 mm
Initial Rotation	0.00 deg
Initial Orientation	Transversal

**Geometry - Saturation**

Fat suppr.	Fat sat.
Special sat.	None

**System - Miscellaneous**

Positioning mode	REF
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
Matrix Optimization	Off
AutoAlign	Head > Brain
Coil Select Mode	Default

**System - Adjustments**

B0 Shim mode	Standard
B1 Shim mode	TrueForm
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

**System - Adjust Volume**

Position	Isocenter
Orientation	Transversal
Rotation	0.00 deg
A >> P	216 mm
R >> L	220 mm
F >> H	60 mm
Reset	Off

**System - Tx/Rx**

Frequency 1H	297.162474 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

**Physio - Signal1**

1st Signal/Mode	None
TR	2220 ms
Multi-band accel. factor	2

**BOLD**

GLM Statistics	Off
Dynamic t-maps	Off
Ignore meas. at start	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	20
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active
Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Meas[16]	Active
Meas[17]	Active
Meas[18]	Active
Meas[19]	Active
Meas[20]	Active
Motion correction	Off
Spatial filter	Off
Measurements	275
Delay in TR	0 ms
Multiple series	Off

**Sequence - Part 1**

Introduction	Off
Contrasts	1
Flow comp.	No
Multi-slice mode	Interleaved
Free echo spacing	Off
Echo spacing	0.81 ms
Bandwidth	1420 Hz/Px

**Sequence - Part 2**

EPI factor	216
Gradient mode	Fast
RF spoiling	Off

**Sequence - Special**

Excite pulse duration	6000 us
Single-band images	Off
MB LeakBlock kernel	On
MB dual kernel	Off
MB RF phase scramble	Off
SENSE1 coil combine	Off
Invert RO/PE polarity	Off
PF omits higher k-space	Off
Disable freq. update	Off
Force equal slice timing	Off
Online multi-band recon.	Online
FFT scale factor	1.00
Fat saturation FA	110.0 deg
Physio recording	Off
Triggering scheme	Standard

\\CRC\protocols\studies\v5motion\cmrr\_mbep2d\_p4\_mb2\_750um\_AAbrain

TA: 10:35 PM: REF Voxel size: 0.7×0.7×0.8 mmPAT: 4 Rel. SNR: 1.00 : epdf

**Properties**

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	Off
Start measurements	Single measurement

**Routine**

Slice group	1
Slices	60
Dist. factor	0 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	Head > Brain
Phase oversampling	0 %
FoV read	220 mm
FoV phase	89.8 %
Slice thickness	0.75 mm
TR	2181 ms
TE	25.80 ms
Multi-band accel. factor	2
Filter	None
Coil elements	A32

**Contrast - Common**

TR	2181 ms
TE	25.80 ms
MTC	Off
Magn. preparation	None
Flip angle	90 deg
Fat suppr.	Fat sat.

**Contrast - Dynamic**

Averaging mode	Long term
Reconstruction	Magnitude
Measurements	275
Delay in TR	0 ms
Multiple series	Off

**Resolution - Common**

FoV read	220 mm
FoV phase	89.8 %
Slice thickness	0.75 mm
Base resolution	294
Phase resolution	100 %
Phase partial Fourier	6/8
Interpolation	Off

**Resolution - iPAT**

PAT mode	GRAPPA
Accel. factor PE	4
Ref. lines PE	48
Reference scan mode	Segmented

**Resolution - Filter Image**

Distortion Corr.	Off
Prescan Normalize	Off

**Resolution - Filter Rawdata**

Raw filter	Off
Elliptical filter	Off
Hamming	Off

**Geometry - Common**

Slice group	1
Slices	60
Dist. factor	0 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
FoV read	220 mm
FoV phase	89.8 %
Slice thickness	0.75 mm
TR	2181 ms
Multi-slice mode	Interleaved
Series	Interleaved
Multi-band accel. factor	2

**Geometry - AutoAlign**

Slice group	1
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	Head > Brain
Initial Position	Isocenter
L	0.0 mm
P	0.0 mm
H	0.0 mm
Initial Rotation	0.00 deg
Initial Orientation	Transversal

**Geometry - Saturation**

Fat suppr.	Fat sat.
Special sat.	None

**System - Miscellaneous**

Positioning mode	REF
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
Matrix Optimization	Off
AutoAlign	Head > Brain
Coil Select Mode	Default

**System - Adjustments**

B0 Shim mode	Standard
B1 Shim mode	TrueForm
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

**System - Adjust Volume**

Position	Isocenter
Orientation	Transversal
Rotation	0.00 deg
A >> P	198 mm
R >> L	220 mm
F >> H	45 mm
Reset	Off

**System - Tx/Rx**

Frequency 1H	297.162474 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

**Physio - Signal1**

1st Signal/Mode	None
TR	2181 ms
Multi-band accel. factor	2

**BOLD**

GLM Statistics	Off
Dynamic t-maps	Off
Ignore meas. at start	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	20
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active
Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Meas[16]	Active
Meas[17]	Active
Meas[18]	Active
Meas[19]	Active
Meas[20]	Active
Motion correction	Off
Spatial filter	Off
Measurements	275
Delay in TR	0 ms
Multiple series	Off

**Sequence - Part 1**

Introduction	Off
Contrasts	1
Flow comp.	No
Multi-slice mode	Interleaved
Free echo spacing	Off
Echo spacing	1.01 ms
Bandwidth	1214 Hz/Px

**Sequence - Part 2**

EPI factor	264
Gradient mode	Fast
RF spoiling	Off

**Sequence - Special**

Excite pulse duration	6000 us
Single-band images	Off
MB LeakBlock kernel	On
MB dual kernel	Off
MB RF phase scramble	Off
SENSE1 coil combine	On
Invert RO/PE polarity	Off
PF omits higher k-space	Off
Disable freq. update	Off
Force equal slice timing	Off
Online multi-band recon.	Online
FFT scale factor	1.00
Fat saturation FA	110.0 deg
Physio recording	Off
Triggering scheme	Standard



\\CRC\protocols\studies\v5motion\cmrr\_mbep2d\_p3\_mb2\_1.6iso\_AABrain

TA: 10:19 PM: FIX Voxel size: 1.6×1.6×1.6 mmPAT: 3 Rel. SNR: 1.00 : epfid

**Properties**

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	Off
Start measurements	Single measurement

**Routine**

Slice group	1
Slices	72
Dist. factor	0 %
Position	L0.5 P2.1 F1.0 mm
Orientation	T > C-32.2 > S-0.6
Phase enc. dir.	A >> P
AutoAlign	Head > Brain
Phase oversampling	0 %
FoV read	216 mm
FoV phase	101.5 %
Slice thickness	1.60 mm
TR	1800 ms
TE	22.00 ms
Multi-band accel. factor	2
Filter	None
Coil elements	A32

**Contrast - Common**

TR	1800 ms
TE	22.00 ms
MTC	Off
Magn. preparation	None
Flip angle	70 deg
Fat suppr.	Fat sat.

**Contrast - Dynamic**

Averaging mode	Long term
Reconstruction	Magnitude
Measurements	330
Delay in TR	0 ms
Multiple series	Off

**Resolution - Common**

FoV read	216 mm
FoV phase	101.5 %
Slice thickness	1.60 mm
Base resolution	134
Phase resolution	100 %
Phase partial Fourier	Off
Interpolation	Off

**Resolution - iPAT**

PAT mode	GRAPPA
Accel. factor PE	3
Ref. lines PE	36
Reference scan mode	Segmented

**Resolution - Filter Image**

Distortion Corr.	Off
Prescan Normalize	Off

**Resolution - Filter Rawdata**

Raw filter	Off
Elliptical filter	Off
Hamming	Off

**Geometry - Common**

Slice group	1
Slices	72
Dist. factor	0 %
Position	L0.5 P2.1 F1.0 mm
Orientation	T > C-32.2 > S-0.6
Phase enc. dir.	A >> P
FoV read	216 mm
FoV phase	101.5 %
Slice thickness	1.60 mm
TR	1800 ms
Multi-slice mode	Interleaved
Series	Interleaved
Multi-band accel. factor	2

**Geometry - AutoAlign**

Slice group	1
Position	L0.5 P2.1 F1.0 mm
Orientation	T > C-32.2 > S-0.6
Phase enc. dir.	A >> P
AutoAlign	Head > Brain
Initial Position	L0.0 P2.2 H7.1
R	0.0 mm
P	2.2 mm
H	7.1 mm
Initial Rotation	0.31 deg
Initial Orientation	T > C
T > C	-21.4
> S	0.0

**Geometry - Saturation**

Fat suppr.	Fat sat.
Special sat.	None

**System - Miscellaneous**

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
Matrix Optimization	Off
AutoAlign	Head > Brain
Coil Select Mode	Default

**System - Adjustments**

B0 Shim mode	Standard
B1 Shim mode	TrueForm
Confirm freq. adjustment	Off
Assume Dominant Fat	Off

**System - Adjustments**

Assume Silicone	Off
Adjustment Tolerance	Auto

**System - Adjust Volume**

Position	L0.5 P2.1 F1.0 mm
Orientation	T > C-32.2 > S-0.6
Rotation	88.70 deg
R >> L	216 mm
A >> P	220 mm
F >> H	116 mm
Reset	Off

**System - Tx/Rx**

Frequency 1H	297.162474 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

**Physio - Signal1**

1st Signal/Mode	None
TR	1800 ms
Multi-band accel. factor	2

**BOLD**

GLM Statistics	Off
Dynamic t-maps	Off
Ignore meas. at start	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	20
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active
Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Meas[16]	Active
Meas[17]	Active
Meas[18]	Active
Meas[19]	Active
Meas[20]	Active
Motion correction	Off
Spatial filter	Off
Measurements	330
Delay in TR	0 ms
Multiple series	Off

**Sequence - Part 1**

Introduction	Off
Contrasts	1
Flow comp.	No

**Sequence - Part 1**

Multi-slice mode	Interleaved
Free echo spacing	Off
Echo spacing	0.59 ms
Bandwidth	2072 Hz/Px

**Sequence - Part 2**

EPI factor	136
Gradient mode	Fast
RF spoiling	Off

**Sequence - Special**

Excite pulse duration	7200 us
Single-band images	Off
MB LeakBlock kernel	On
MB dual kernel	Off
MB RF phase scramble	Off
SENSE1 coil combine	Off
Invert RO/PE polarity	Off
Disable freq. update	Off
Force equal slice timing	Off
Online multi-band recon.	Online
FFT scale factor	1.00
Fat saturation FA	110.0 deg
Physio recording	DICOM
Triggering scheme	Standard

\\CRC\protocols\studies\lv5motion\cmrr\_mbep2d\_p3\_mb2\_1.6iso\_reversePE\_AABrain

TA: 0:34 PM: FIX Voxel size: 1.6×1.6×1.6 mmPAT: 3 Rel. SNR: 1.00 : epfid

**Properties**

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	Off
Start measurements	Single measurement

**Routine**

Slice group	1
Slices	72
Dist. factor	0 %
Position	L0.5 P2.1 F1.0 mm
Orientation	T > C-32.2 > S-0.6
Phase enc. dir.	A >> P
AutoAlign	Head > Brain
Phase oversampling	0 %
FoV read	216 mm
FoV phase	101.5 %
Slice thickness	1.60 mm
TR	1800 ms
TE	22.00 ms
Multi-band accel. factor	2
Filter	None
Coil elements	A32

**Contrast - Common**

TR	1800 ms
TE	22.00 ms
MTC	Off
Magn. preparation	None
Flip angle	70 deg
Fat suppr.	Fat sat.

**Contrast - Dynamic**

Averaging mode	Long term
Reconstruction	Magnitude
Measurements	5
Delay in TR	0 ms
Multiple series	Off

**Resolution - Common**

FoV read	216 mm
FoV phase	101.5 %
Slice thickness	1.60 mm
Base resolution	134
Phase resolution	100 %
Phase partial Fourier	Off
Interpolation	Off

**Resolution - iPAT**

PAT mode	GRAPPA
Accel. factor PE	3
Ref. lines PE	36
Reference scan mode	Segmented

**Resolution - Filter Image**

Distortion Corr.	Off
Prescan Normalize	Off

**Resolution - Filter Rawdata**

Raw filter	Off
Elliptical filter	Off
Hamming	Off

**Geometry - Common**

Slice group	1
Slices	72
Dist. factor	0 %
Position	L0.5 P2.1 F1.0 mm
Orientation	T > C-32.2 > S-0.6
Phase enc. dir.	A >> P
FoV read	216 mm
FoV phase	101.5 %
Slice thickness	1.60 mm
TR	1800 ms
Multi-slice mode	Interleaved
Series	Interleaved
Multi-band accel. factor	2

**Geometry - AutoAlign**

Slice group	1
Position	L0.5 P2.1 F1.0 mm
Orientation	T > C-32.2 > S-0.6
Phase enc. dir.	A >> P
AutoAlign	Head > Brain
Initial Position	L0.0 P2.2 H7.1
R	0.0 mm
P	2.2 mm
H	7.1 mm
Initial Rotation	0.31 deg
Initial Orientation	T > C
T > C	-21.4
> S	0.0

**Geometry - Saturation**

Fat suppr.	Fat sat.
Special sat.	None

**System - Miscellaneous**

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
Matrix Optimization	Off
AutoAlign	Head > Brain
Coil Select Mode	Default

**System - Adjustments**

B0 Shim mode	Standard
B1 Shim mode	TrueForm
Confirm freq. adjustment	Off
Assume Dominant Fat	Off

**System - Adjustments**

Assume Silicone	Off
Adjustment Tolerance	Auto

**System - Adjust Volume**

Position	L0.5 P2.1 F1.0 mm
Orientation	T > C-32.2 > S-0.6
Rotation	88.70 deg
R >> L	216 mm
A >> P	220 mm
F >> H	116 mm
Reset	Off

**System - Tx/Rx**

Frequency 1H	297.162474 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

**Physio - Signal1**

1st Signal/Mode	None
TR	1800 ms
Multi-band accel. factor	2

**BOLD**

GLM Statistics	Off
Dynamic t-maps	Off
Ignore meas. at start	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	20
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active
Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Meas[16]	Active
Meas[17]	Active
Meas[18]	Active
Meas[19]	Active
Meas[20]	Active
Motion correction	Off
Spatial filter	Off
Measurements	5
Delay in TR	0 ms
Multiple series	Off

**Sequence - Part 1**

Introduction	Off
Contrasts	1
Flow comp.	No

**Sequence - Part 1**

Multi-slice mode	Interleaved
Free echo spacing	Off
Echo spacing	0.59 ms
Bandwidth	2072 Hz/Px

**Sequence - Part 2**

EPI factor	136
Gradient mode	Fast
RF spoiling	Off

**Sequence - Special**

Excite pulse duration	7200 us
Single-band images	Off
MB LeakBlock kernel	On
MB dual kernel	Off
MB RF phase scramble	Off
SENSE1 coil combine	Off
Invert RO/PE polarity	On
Disable freq. update	Off
Force equal slice timing	Off
Online multi-band recon.	Online
FFT scale factor	1.00
Fat saturation FA	110.0 deg
Physio recording	Off
Triggering scheme	Standard

\\CRC\protocols\studies\v5motion\gre\_field\_mapping\_2mm\_eufind

TA: 1:38 PM: FIX Voxel size: 2.0×2.0×2.0 mmRel. SNR: 1.00 : fm

**Properties**

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	Off
Start measurements	Single measurement

**Routine**

Slice group	1
Slices	96
Dist. factor	0 %
Position	L0.5 P2.1 F1.0 mm
Orientation	T > C-32.2 > S-0.6
Phase enc. dir.	R >> L
AutoAlign	Head > Brain
Phase oversampling	0 %
FoV read	256 mm
FoV phase	75.0 %
Slice thickness	2.0 mm
TR	5.2 ms
TE 1	2.26 ms
TE 2	3.28 ms
Averages	1
Concatenations	96
Filter	None
Coil elements	A32

**Contrast - Common**

TR	5.2 ms
TE 1	2.26 ms
TE 2	3.28 ms
MTC	Off
Flip angle	15 deg
Fat suppr.	None

**Contrast - Dynamic**

Averages	1
Averaging mode	Long term
Reconstruction	Magn./Phase
Measurements	1
Multiple series	Off

**Resolution - Common**

FoV read	256 mm
FoV phase	75.0 %
Slice thickness	2.0 mm
Base resolution	128
Phase resolution	100 %
Phase partial Fourier	Off
Interpolation	Off

**Resolution - Filter Image**

Image Filter	Off
Distortion Corr.	Off

**Resolution - Filter Image**

Prescan Normalize	Off
Normalize	Off
B1 filter	Off

**Resolution - Filter Rawdata**

Raw filter	Off
Elliptical filter	Off

**Geometry - Common**

Slice group	1
Slices	96
Dist. factor	0 %
Position	L0.5 P2.1 F1.0 mm
Orientation	T > C-32.2 > S-0.6
Phase enc. dir.	R >> L
FoV read	256 mm
FoV phase	75.0 %
Slice thickness	2.0 mm
TR	5.2 ms
Multi-slice mode	Sequential
Series	Ascending
Concatenations	96

**Geometry - AutoAlign**

Slice group	1
Position	L0.5 P2.1 F1.0 mm
Orientation	T > C-32.2 > S-0.6
Phase enc. dir.	R >> L
AutoAlign	Head > Brain
Initial Position	L0.0 P2.2 H7.1
R	0.0 mm
P	2.2 mm
H	7.1 mm
Initial Rotation	90.31 deg
Initial Orientation	T > C
T > C	-21.4
> S	0.0

**Geometry - Saturation**

Fat suppr.	None
Special sat.	None

**System - Miscellaneous**

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
Save uncombined	Off
Matrix Optimization	Off
AutoAlign	Head > Brain
Coil Select Mode	Default

**System - Adjustments**

B0 Shim mode	Standard
B1 Shim mode	TrueForm
Confirm freq. adjustment	Off

**System - Adjustments**

Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

**System - Adjust Volume**

! Position	L0.5 P2.1 F1.0 mm
! Orientation	T > C-32.2 > S-0.6
! Rotation	88.70 deg
! R >> L	216 mm
! A >> P	220 mm
! F >> H	116 mm
Reset	Off

**System - Tx/Rx**

Frequency 1H	297.162474 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

**Sequence - Part 1**

Introduction	On
Dimension	2D
Asymmetric echo	Off
Contrasts	2
Flow comp.	No
Multi-slice mode	Sequential
Bandwidth	737 Hz/Px

**Sequence - Part 2**

RF pulse type	Fast
Gradient mode	Fast
RF spoiling	On

**Sequence - Assistant**

Mode	Off
------	-----

\\CRC\protocols\studies\v5motion\cmrr\_mbep2d\_p4\_mb2\_750um\_AAbrain

TA: 10:35 PM: REF Voxel size: 0.7×0.7×0.8 mmPAT: 4 Rel. SNR: 1.00 : epdf

**Properties**

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	Off
Start measurements	Single measurement

**Routine**

Slice group	1
Slices	60
Dist. factor	0 %
Position	L0.4 A5.3 F10.7 mm
Orientation	T > C-27.9 > S-0.7
Phase enc. dir.	A >> P
AutoAlign	Head > Brain
Phase oversampling	0 %
FoV read	220 mm
FoV phase	89.8 %
Slice thickness	0.75 mm
TR	2181 ms
TE	25.80 ms
Multi-band accel. factor	2
Filter	None
Coil elements	A32

**Contrast - Common**

TR	2181 ms
TE	25.80 ms
MTC	Off
Magn. preparation	None
Flip angle	90 deg
Fat suppr.	Fat sat.

**Contrast - Dynamic**

Averaging mode	Long term
Reconstruction	Magnitude
Measurements	275
Delay in TR	0 ms
Multiple series	Off

**Resolution - Common**

FoV read	220 mm
FoV phase	89.8 %
Slice thickness	0.75 mm
Base resolution	294
Phase resolution	100 %
Phase partial Fourier	6/8
Interpolation	Off

**Resolution - iPAT**

PAT mode	GRAPPA
Accel. factor PE	4
Ref. lines PE	48
Reference scan mode	Segmented

**Resolution - Filter Image**

Distortion Corr.	Off
Prescan Normalize	Off

**Resolution - Filter Rawdata**

Raw filter	Off
Elliptical filter	Off
Hamming	Off

**Geometry - Common**

Slice group	1
Slices	60
Dist. factor	0 %
Position	L0.4 A5.3 F10.7 mm
Orientation	T > C-27.9 > S-0.7
Phase enc. dir.	A >> P
FoV read	220 mm
FoV phase	89.8 %
Slice thickness	0.75 mm
TR	2181 ms
Multi-slice mode	Interleaved
Series	Interleaved
Multi-band accel. factor	2

**Geometry - AutoAlign**

Slice group	1
Position	L0.4 A5.3 F10.7 mm
Orientation	T > C-27.9 > S-0.7
Phase enc. dir.	A >> P
AutoAlign	Head > Brain
Initial Position	L0.0 A3.3 F3.8
R	0.0 mm
A	3.3 mm
F	3.8 mm
Initial Rotation	0.28 deg
Initial Orientation	T > C
T > C	-17.1
> S	0.0

**Geometry - Saturation**

Fat suppr.	Fat sat.
Special sat.	None

**System - Miscellaneous**

Positioning mode	REF
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
Matrix Optimization	Off
AutoAlign	Head > Brain
Coil Select Mode	Default

**System - Adjustments**

B0 Shim mode	Standard
B1 Shim mode	TrueForm
Confirm freq. adjustment	Off
Assume Dominant Fat	Off

**System - Adjustments**

Assume Silicone	Off
Adjustment Tolerance	Auto

**System - Adjust Volume**

Position	L0.4 A5.3 F10.7 mm
Orientation	T > C-27.9 > S-0.7
Rotation	-1.29 deg
A >> P	198 mm
R >> L	220 mm
F >> H	45 mm
Reset	Off

**System - Tx/Rx**

Frequency 1H	297.162474 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

**Physio - Signal1**

1st Signal/Mode	None
TR	2181 ms
Multi-band accel. factor	2

**BOLD**

GLM Statistics	Off
Dynamic t-maps	Off
Ignore meas. at start	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	20
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active
Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Meas[16]	Active
Meas[17]	Active
Meas[18]	Active
Meas[19]	Active
Meas[20]	Active
Motion correction	Off
Spatial filter	Off
Measurements	275
Delay in TR	0 ms
Multiple series	Off

**Sequence - Part 1**

Introduction	Off
Contrasts	1
Flow comp.	No

**Sequence - Part 1**

Multi-slice mode	Interleaved
Free echo spacing	Off
Echo spacing	1.01 ms
Bandwidth	1214 Hz/Px

**Sequence - Part 2**

EPI factor	264
Gradient mode	Fast
RF spoiling	Off

**Sequence - Special**

Excite pulse duration	6000 us
Single-band images	Off
MB LeakBlock kernel	On
MB dual kernel	Off
MB RF phase scramble	Off
SENSE1 coil combine	On
Invert RO/PE polarity	Off
PF omits higher k-space	Off
Disable freq. update	Off
Force equal slice timing	Off
Online multi-band recon.	Online
FFT scale factor	1.00
Fat saturation FA	110.0 deg
Physio recording	DICOM
Triggering scheme	Standard



\\CRC\protocols\studies\v5motion\cmrr\_mbep2d\_p4\_mb2\_750um\_reversePE\_AAbra

TA: 0:46 PM: FIX Voxel size: 0.7×0.7×0.8 mmPAT: 4 Rel. SNR: 1.00 : epfid

**Properties**

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	Off
Start measurements	Single measurement

**Routine**

Slice group	1
Slices	60
Dist. factor	0 %
Position	L0.4 A5.3 F10.7 mm
Orientation	T > C-27.9 > S-0.7
Phase enc. dir.	A >> P
AutoAlign	Head > Brain
Phase oversampling	0 %
FoV read	220 mm
FoV phase	89.8 %
Slice thickness	0.75 mm
TR	2181 ms
TE	25.80 ms
Multi-band accel. factor	2
Filter	None
Coil elements	A32

**Contrast - Common**

TR	2181 ms
TE	25.80 ms
MTC	Off
Magn. preparation	None
Flip angle	90 deg
Fat suppr.	Fat sat.

**Contrast - Dynamic**

Averaging mode	Long term
Reconstruction	Magnitude
Measurements	5
Delay in TR	0 ms
Multiple series	Off

**Resolution - Common**

FoV read	220 mm
FoV phase	89.8 %
Slice thickness	0.75 mm
Base resolution	294
Phase resolution	100 %
Phase partial Fourier	6/8
Interpolation	Off

**Resolution - iPAT**

PAT mode	GRAPPA
Accel. factor PE	4
Ref. lines PE	48
Reference scan mode	Segmented

**Resolution - Filter Image**

Distortion Corr.	Off
Prescan Normalize	Off

**Resolution - Filter Rawdata**

Raw filter	Off
Elliptical filter	Off
Hamming	Off

**Geometry - Common**

Slice group	1
Slices	60
Dist. factor	0 %
Position	L0.4 A5.3 F10.7 mm
Orientation	T > C-27.9 > S-0.7
Phase enc. dir.	A >> P
FoV read	220 mm
FoV phase	89.8 %
Slice thickness	0.75 mm
TR	2181 ms
Multi-slice mode	Interleaved
Series	Interleaved
Multi-band accel. factor	2

**Geometry - AutoAlign**

Slice group	1
Position	L0.4 A5.3 F10.7 mm
Orientation	T > C-27.9 > S-0.7
Phase enc. dir.	A >> P
AutoAlign	Head > Brain
Initial Position	L0.0 A3.3 F3.8
R	0.0 mm
A	3.3 mm
F	3.8 mm
Initial Rotation	0.28 deg
Initial Orientation	T > C
T > C	-17.1
> S	0.0

**Geometry - Saturation**

Fat suppr.	Fat sat.
Special sat.	None

**System - Miscellaneous**

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
Matrix Optimization	Off
AutoAlign	Head > Brain
Coil Select Mode	Default

**System - Adjustments**

B0 Shim mode	Standard
B1 Shim mode	TrueForm
Confirm freq. adjustment	Off
Assume Dominant Fat	Off

**System - Adjustments**

Assume Silicone	Off
Adjustment Tolerance	Auto

**System - Adjust Volume**

Position	L0.4 A5.3 F10.7 mm
Orientation	T > C-27.9 > S-0.7
Rotation	-1.29 deg
A >> P	198 mm
R >> L	220 mm
F >> H	45 mm
Reset	Off

**System - Tx/Rx**

Frequency 1H	297.162474 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

**Physio - Signal1**

1st Signal/Mode	None
TR	2181 ms
Multi-band accel. factor	2

**BOLD**

GLM Statistics	Off
Dynamic t-maps	Off
Ignore meas. at start	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	20
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active
Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Meas[16]	Active
Meas[17]	Active
Meas[18]	Active
Meas[19]	Active
Meas[20]	Active
Motion correction	Off
Spatial filter	Off
Measurements	5
Delay in TR	0 ms
Multiple series	Off

**Sequence - Part 1**

Introduction	Off
Contrasts	1
Flow comp.	No

**Sequence - Part 1**

Multi-slice mode	Interleaved
Free echo spacing	Off
Echo spacing	1.01 ms
Bandwidth	1214 Hz/Px

**Sequence - Part 2**

EPI factor	264
Gradient mode	Fast
RF spoiling	Off

**Sequence - Special**

Excite pulse duration	6000 us
Single-band images	Off
MB LeakBlock kernel	On
MB dual kernel	Off
MB RF phase scramble	Off
SENSE1 coil combine	On
Invert RO/PE polarity	On
PF omits higher k-space	Off
Disable freq. update	Off
Force equal slice timing	Off
Online multi-band recon.	Online
FFT scale factor	1.00
Fat saturation FA	110.0 deg
Physio recording	Off
Triggering scheme	Standard

\\CRC\protocols\studies\v5motion\cmrr\_mbep2d\_p4\_mb2\_750um\_GRAPPA-GRE\_FA75

TA: 1:07 PM: FIX Voxel size: 0.8×0.8×0.8 mmPAT: 4 Rel. SNR: 1.00 : epfid

**Properties**

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	Off
Start measurements	Single measurement

**Routine**

Slice group	1
Slices	48
Dist. factor	0 %
Position	L0.0 A12.1 F4.8 mm
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
Phase oversampling	0 %
FoV read	186 mm
FoV phase	94.4 %
Slice thickness	0.75 mm
TR	3000 ms
TE	25.00 ms
Multi-band accel. factor	2
Filter	Raw filter
Coil elements	A32

**Contrast - Common**

TR	3000 ms
TE	25.00 ms
MTC	Off
Magn. preparation	None
Flip angle	75 deg
Fat suppr.	Fat sat.

**Contrast - Dynamic**

Averaging mode	Long term
Reconstruction	Magnitude
Measurements	10
Delay in TR	0 ms
Multiple series	Off

**Resolution - Common**

FoV read	186 mm
FoV phase	94.4 %
Slice thickness	0.75 mm
Base resolution	248
Phase resolution	100 %
Phase partial Fourier	6/8
Interpolation	Off

**Resolution - iPAT**

PAT mode	GRAPPA
Accel. factor PE	4
Ref. lines PE	48
Reference scan mode	GRE

**Resolution - Filter Image**

Distortion Corr.	Off
Prescan Normalize	Off

**Resolution - Filter Rawdata**

Raw filter	On
Elliptical filter	Off
Hamming	Off

**Geometry - Common**

Slice group	1
Slices	48
Dist. factor	0 %
Position	L0.0 A12.1 F4.8 mm
Orientation	Transversal
Phase enc. dir.	A >> P
FoV read	186 mm
FoV phase	94.4 %
Slice thickness	0.75 mm
TR	3000 ms
Multi-slice mode	Interleaved
Series	Descending
Multi-band accel. factor	2

**Geometry - AutoAlign**

Slice group	1
Position	L0.0 A12.1 F4.8 mm
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
Initial Position	L0.0 A12.1 F4.8
L	0.0 mm
A	12.1 mm
F	4.8 mm
Initial Rotation	0.00 deg
Initial Orientation	Transversal

**Geometry - Saturation**

Fat suppr.	Fat sat.
Special sat.	None

**System - Miscellaneous**

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
Matrix Optimization	Off
AutoAlign	---
Coil Select Mode	Default

**System - Adjustments**

B0 Shim mode	Standard
B1 Shim mode	TrueForm
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

**System - Adjust Volume**

Position	L0.0 A12.1 F4.8 mm
Orientation	Transversal
Rotation	0.00 deg
A >> P	176 mm
R >> L	186 mm
F >> H	36 mm
Reset	Off

**System - Tx/Rx**

Frequency 1H	297.162474 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

**Physio - Signal1**

1st Signal/Mode	None
TR	3000 ms
Multi-band accel. factor	2

**BOLD**

GLM Statistics	Off
Dynamic t-maps	Off
Ignore meas. at start	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	20
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active
Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Meas[16]	Active
Meas[17]	Active
Meas[18]	Active
Meas[19]	Active
Meas[20]	Active
Motion correction	Off
Spatial filter	Off
Measurements	10
Delay in TR	0 ms
Multiple series	Off

**Sequence - Part 1**

Introduction	On
Contrasts	1
Flow comp.	No
Multi-slice mode	Interleaved
Free echo spacing	Off
Echo spacing	1.04 ms
Bandwidth	1120 Hz/Px

**Sequence - Part 2**

EPI factor	234
Gradient mode	Fast
RF spoiling	Off

**Sequence - Special**

Excite pulse duration	6000 us
Single-band images	Off
MB LeakBlock kernel	On
MB dual kernel	Off
MB RF phase scramble	Off
SENSE1 coil combine	On
Invert RO/PE polarity	Off
PF omits higher k-space	Off
Disable freq. update	Off
Force equal slice timing	Off
Online multi-band recon.	Online
FFT scale factor	1.00
Fat saturation FA	110.0 deg
GRE iPAT ref. FA	12.0 deg
Physio recording	Off
Triggering scheme	Standard

\\CRC\protocols\studies\v5motion\cmrr\_mbep2d\_p4\_mb2\_750um\_AAbrain

TA: 10:35 PM: REF Voxel size: 0.7×0.7×0.8 mmPAT: 4 Rel. SNR: 1.00 : epdf

**Properties**

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	Off
Start measurements	Single measurement

**Routine**

Slice group	1
Slices	60
Dist. factor	0 %
Position	L1.7 P11.0 H4.6 mm
Orientation	T > C-38.4 > S-0.7
Phase enc. dir.	A >> P
AutoAlign	Head > Brain
Phase oversampling	0 %
FoV read	220 mm
FoV phase	89.8 %
Slice thickness	0.75 mm
TR	2181 ms
TE	25.80 ms
Multi-band accel. factor	2
Filter	None
Coil elements	A32

**Contrast - Common**

TR	2181 ms
TE	25.80 ms
MTC	Off
Magn. preparation	None
Flip angle	90 deg
Fat suppr.	Fat sat.

**Contrast - Dynamic**

Averaging mode	Long term
Reconstruction	Magnitude
Measurements	275
Delay in TR	0 ms
Multiple series	Off

**Resolution - Common**

FoV read	220 mm
FoV phase	89.8 %
Slice thickness	0.75 mm
Base resolution	294
Phase resolution	100 %
Phase partial Fourier	6/8
Interpolation	Off

**Resolution - iPAT**

PAT mode	GRAPPA
Accel. factor PE	4
Ref. lines PE	48
Reference scan mode	Segmented

**Resolution - Filter Image**

Distortion Corr.	Off
Prescan Normalize	Off

**Resolution - Filter Rawdata**

Raw filter	Off
Elliptical filter	Off
Hamming	Off

**Geometry - Common**

Slice group	1
Slices	60
Dist. factor	0 %
Position	L1.7 P11.0 H4.6 mm
Orientation	T > C-38.4 > S-0.7
Phase enc. dir.	A >> P
FoV read	220 mm
FoV phase	89.8 %
Slice thickness	0.75 mm
TR	2181 ms
Multi-slice mode	Interleaved
Series	Interleaved
Multi-band accel. factor	2

**Geometry - AutoAlign**

Slice group	1
Position	L1.7 P11.0 H4.6 mm
Orientation	T > C-38.4 > S-0.7
Phase enc. dir.	A >> P
AutoAlign	Head > Brain
Initial Position	L0.0 P7.8 H14.2
L	0.0 mm
P	7.8 mm
H	14.2 mm
Initial Rotation	0.35 deg
Initial Orientation	T > C
T > C	-24.0
> S	0.0

**Geometry - Saturation**

Fat suppr.	Fat sat.
Special sat.	None

**System - Miscellaneous**

Positioning mode	REF
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
Matrix Optimization	Off
AutoAlign	Head > Brain
Coil Select Mode	Default

**System - Adjustments**

B0 Shim mode	Standard
B1 Shim mode	TrueForm
Confirm freq. adjustment	Off
Assume Dominant Fat	Off

**System - Adjustments**

Assume Silicone	Off
Adjustment Tolerance	Auto

**System - Adjust Volume**

Position	L1.7 P11.0 H4.6 mm
Orientation	T > C-38.4 > S-0.7
Rotation	2.74 deg
A >> P	198 mm
R >> L	220 mm
F >> H	45 mm
Reset	Off

**System - Tx/Rx**

Frequency 1H	297.162474 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

**Physio - Signal1**

1st Signal/Mode	None
TR	2181 ms
Multi-band accel. factor	2

**BOLD**

GLM Statistics	Off
Dynamic t-maps	Off
Ignore meas. at start	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	20
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active
Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Meas[16]	Active
Meas[17]	Active
Meas[18]	Active
Meas[19]	Active
Meas[20]	Active
Motion correction	Off
Spatial filter	Off
Measurements	275
Delay in TR	0 ms
Multiple series	Off

**Sequence - Part 1**

Introduction	Off
Contrasts	1
Flow comp.	No

**Sequence - Part 1**

Multi-slice mode	Interleaved
Free echo spacing	Off
Echo spacing	1.01 ms
Bandwidth	1214 Hz/Px

**Sequence - Part 2**

EPI factor	264
Gradient mode	Fast
RF spoiling	Off

**Sequence - Special**

Excite pulse duration	6000 us
Single-band images	Off
MB LeakBlock kernel	On
MB dual kernel	Off
MB RF phase scramble	Off
SENSE1 coil combine	On
Invert RO/PE polarity	Off
PF omits higher k-space	Off
Disable freq. update	Off
Force equal slice timing	Off
Online multi-band recon.	Online
FFT scale factor	1.00
Fat saturation FA	110.0 deg
Physio recording	DICOM
Triggering scheme	Standard

\\CRC\protocols\studies\v5motion\cmrr\_mbep2d\_p4\_mb2\_750um\_GRAPPA\_GRE\_FA90\_B1148

TA: 15:37 PM: FIX Voxel size: 0.8×0.8×0.8 mmPAT: 4 Rel. SNR: 1.00 : epfid

**Properties**

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	Off
Start measurements	Single measurement

**Routine**

Slice group	1
Slices	48
Dist. factor	0 %
Position	L1.7 P11.0 H4.6 mm
Orientation	T > C-38.4 > S-0.7
Phase enc. dir.	A >> P
AutoAlign	Head > Brain
Phase oversampling	0 %
FoV read	192 mm
FoV phase	93.8 %
Slice thickness	0.75 mm
TR	3000 ms
TE	25.80 ms
Multi-band accel. factor	2
Filter	Raw filter
Coil elements	A32

**Contrast - Common**

TR	3000 ms
TE	25.80 ms
MTC	Off
Magn. preparation	None
Flip angle	90 deg
Fat suppr.	Fat sat.

**Contrast - Dynamic**

Averaging mode	Long term
Reconstruction	Magnitude
Measurements	300
Delay in TR	0 ms
Multiple series	Off

**Resolution - Common**

FoV read	192 mm
FoV phase	93.8 %
Slice thickness	0.75 mm
Base resolution	256
Phase resolution	100 %
Phase partial Fourier	6/8
Interpolation	Off

**Resolution - iPAT**

PAT mode	GRAPPA
Accel. factor PE	4
Ref. lines PE	48
Reference scan mode	GRE

**Resolution - Filter Image**

Distortion Corr.	Off
Prescan Normalize	Off

**Resolution - Filter Rawdata**

Raw filter	On
Elliptical filter	Off
Hamming	Off

**Geometry - Common**

Slice group	1
Slices	48
Dist. factor	0 %
Position	L1.7 P11.0 H4.6 mm
Orientation	T > C-38.4 > S-0.7
Phase enc. dir.	A >> P
FoV read	192 mm
FoV phase	93.8 %
Slice thickness	0.75 mm
TR	3000 ms
Multi-slice mode	Interleaved
Series	Descending
Multi-band accel. factor	2

**Geometry - AutoAlign**

Slice group	1
Position	L1.7 P11.0 H4.6 mm
Orientation	T > C-38.4 > S-0.7
Phase enc. dir.	A >> P
AutoAlign	Head > Brain
Initial Position	L0.0 P7.8 H14.2
L	0.0 mm
P	7.8 mm
H	14.2 mm
Initial Rotation	0.35 deg
Initial Orientation	T > C
T > C	-24.0
> S	0.0

**Geometry - Saturation**

Fat suppr.	Fat sat.
Special sat.	None

**System - Miscellaneous**

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
Matrix Optimization	Off
AutoAlign	Head > Brain
Coil Select Mode	Default

**System - Adjustments**

B0 Shim mode	Standard
B1 Shim mode	TrueForm
Confirm freq. adjustment	Off
Assume Dominant Fat	Off

**System - Adjustments**

Assume Silicone	Off
Adjustment Tolerance	Auto

**System - Adjust Volume**

Position	L1.7 P11.0 H4.6 mm
Orientation	T > C-38.4 > S-0.7
Rotation	2.74 deg
A >> P	180 mm
R >> L	192 mm
F >> H	36 mm
Reset	Off

**System - Tx/Rx**

Frequency 1H	297.162474 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

**Physio - Signal1**

1st Signal/Mode	None
TR	3000 ms
Multi-band accel. factor	2

**BOLD**

GLM Statistics	Off
Dynamic t-maps	Off
Ignore meas. at start	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	20
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active
Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Meas[16]	Active
Meas[17]	Active
Meas[18]	Active
Meas[19]	Active
Meas[20]	Active
Motion correction	Off
Spatial filter	Off
Measurements	300
Delay in TR	0 ms
Multiple series	Off

**Sequence - Part 1**

Introduction	On
Contrasts	1
Flow comp.	No

**Sequence - Part 1**

Multi-slice mode	Interleaved
Free echo spacing	Off
Echo spacing	1.02 ms
Bandwidth	1148 Hz/Px

**Sequence - Part 2**

EPI factor	240
Gradient mode	Fast
RF spoiling	Off

**Sequence - Special**

Excite pulse duration	6000 us
Single-band images	Off
MB LeakBlock kernel	On
MB dual kernel	Off
MB RF phase scramble	Off
SENSE1 coil combine	On
Invert RO/PE polarity	Off
PF omits higher k-space	Off
Disable freq. update	Off
Force equal slice timing	Off
Online multi-band recon.	Online
FFT scale factor	1.00
Fat saturation FA	110.0 deg
GRE iPAT ref. FA	12.0 deg
Physio recording	DICOM
Triggering scheme	Standard



\\CRC\protocols\studies\v5motion\cmrr\_mbep2d\_p4\_mb2\_750um\_GRAPPA-GRE\_FA75\_Band

TA: 15:37 PM: FIX Voxel size: 0.8×0.8×0.8 mmPAT: 4 Rel. SNR: 1.00 : epfid

**Properties**

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	Off
Start measurements	Single measurement

**Routine**

Slice group	1
Slices	48
Dist. factor	0 %
Position	L1.7 P11.0 H4.6 mm
Orientation	T > C-38.4 > S-0.7
Phase enc. dir.	A >> P
AutoAlign	Head > Brain
Phase oversampling	0 %
FoV read	186 mm
FoV phase	94.4 %
Slice thickness	0.75 mm
TR	3000 ms
TE	25.00 ms
Multi-band accel. factor	2
Filter	Raw filter
Coil elements	A32

**Contrast - Common**

TR	3000 ms
TE	25.00 ms
MTC	Off
Magn. preparation	None
Flip angle	75 deg
Fat suppr.	Fat sat.

**Contrast - Dynamic**

Averaging mode	Long term
Reconstruction	Magnitude
Measurements	300
Delay in TR	0 ms
Multiple series	Off

**Resolution - Common**

FoV read	186 mm
FoV phase	94.4 %
Slice thickness	0.75 mm
Base resolution	248
Phase resolution	100 %
Phase partial Fourier	6/8
Interpolation	Off

**Resolution - iPAT**

PAT mode	GRAPPA
Accel. factor PE	4
Ref. lines PE	48
Reference scan mode	GRE

**Resolution - Filter Image**

Distortion Corr.	Off
Prescan Normalize	Off

**Resolution - Filter Rawdata**

Raw filter	On
Elliptical filter	Off
Hamming	Off

**Geometry - Common**

Slice group	1
Slices	48
Dist. factor	0 %
Position	L1.7 P11.0 H4.6 mm
Orientation	T > C-38.4 > S-0.7
Phase enc. dir.	A >> P
FoV read	186 mm
FoV phase	94.4 %
Slice thickness	0.75 mm
TR	3000 ms
Multi-slice mode	Interleaved
Series	Descending
Multi-band accel. factor	2

**Geometry - AutoAlign**

Slice group	1
Position	L1.7 P11.0 H4.6 mm
Orientation	T > C-38.4 > S-0.7
Phase enc. dir.	A >> P
AutoAlign	Head > Brain
Initial Position	L0.0 P7.8 H14.2
L	0.0 mm
P	7.8 mm
H	14.2 mm
Initial Rotation	0.35 deg
Initial Orientation	T > C
T > C	-24.0
> S	0.0

**Geometry - Saturation**

Fat suppr.	Fat sat.
Special sat.	None

**System - Miscellaneous**

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
Matrix Optimization	Off
AutoAlign	Head > Brain
Coil Select Mode	Default

**System - Adjustments**

B0 Shim mode	Standard
B1 Shim mode	TrueForm
Confirm freq. adjustment	Off
Assume Dominant Fat	Off

**System - Adjustments**

Assume Silicone	Off
Adjustment Tolerance	Auto

**System - Adjust Volume**

Position	L1.7 P11.0 H4.6 mm
Orientation	T > C-38.4 > S-0.7
Rotation	2.74 deg
A >> P	176 mm
R >> L	186 mm
F >> H	36 mm
Reset	Off

**System - Tx/Rx**

Frequency 1H	297.162474 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

**Physio - Signal1**

1st Signal/Mode	None
TR	3000 ms
Multi-band accel. factor	2

**BOLD**

GLM Statistics	Off
Dynamic t-maps	Off
Ignore meas. at start	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	20
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active
Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Meas[16]	Active
Meas[17]	Active
Meas[18]	Active
Meas[19]	Active
Meas[20]	Active
Motion correction	Off
Spatial filter	Off
Measurements	300
Delay in TR	0 ms
Multiple series	Off

**Sequence - Part 1**

Introduction	On
Contrasts	1
Flow comp.	No

**Sequence - Part 1**

Multi-slice mode	Interleaved
Free echo spacing	Off
Echo spacing	1.03 ms
Bandwidth	1260 Hz/Px

**Sequence - Part 2**

EPI factor	234
Gradient mode	Fast
RF spoiling	Off

**Sequence - Special**

Excite pulse duration	6000 us
Single-band images	Off
MB LeakBlock kernel	On
MB dual kernel	Off
MB RF phase scramble	Off
SENSE1 coil combine	On
Invert RO/PE polarity	Off
PF omits higher k-space	Off
Disable freq. update	Off
Force equal slice timing	Off
Online multi-band recon.	Online
FFT scale factor	1.00
Fat saturation FA	110.0 deg
GRE iPAT ref. FA	12.0 deg
Physio recording	DICOM
Triggering scheme	Standard

\\CRC\protocols\studies\v5motion\cmrr_mbep2d_p4_mb2_750um_GRAPPA-GRE_FA75_Band_TR21
81
TA: 18:45 PM: FIX Voxel size: 0.8×0.8×0.8 mmPAT: 4 Rel. SNR: 1.00 : epfid

**Properties**

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	Off
Start measurements	Single measurement

**Routine**

Slice group	1
Slices	48
Dist. factor	0 %
Position	L1.7 P11.0 H4.6 mm
Orientation	T > C-38.4 > S-0.7
Phase enc. dir.	A >> P
AutoAlign	Head > Brain
Phase oversampling	0 %
FoV read	186 mm
FoV phase	94.4 %
Slice thickness	0.75 mm
TR	2181 ms
TE	25.00 ms
Multi-band accel. factor	2
Filter	Raw filter
Coil elements	A32

**Contrast - Common**

TR	2181 ms
TE	25.00 ms
MTC	Off
Magn. preparation	None
Flip angle	75 deg
Fat suppr.	Fat sat.

**Contrast - Dynamic**

Averaging mode	Long term
Reconstruction	Magnitude
Measurements	500
Delay in TR	0 ms
Multiple series	Off

**Resolution - Common**

FoV read	186 mm
FoV phase	94.4 %
Slice thickness	0.75 mm
Base resolution	248
Phase resolution	100 %
Phase partial Fourier	6/8
Interpolation	Off

**Resolution - iPAT**

PAT mode	GRAPPA
Accel. factor PE	4
Ref. lines PE	48

**Resolution - iPAT**

Reference scan mode	GRE
---------------------	-----

**Resolution - Filter Image**

Distortion Corr.	Off
Prescan Normalize	Off

**Resolution - Filter Rawdata**

Raw filter	On
Elliptical filter	Off
Hamming	Off

**Geometry - Common**

Slice group	1
Slices	48
Dist. factor	0 %
Position	L1.7 P11.0 H4.6 mm
Orientation	T > C-38.4 > S-0.7
Phase enc. dir.	A >> P
FoV read	186 mm
FoV phase	94.4 %
Slice thickness	0.75 mm
TR	2181 ms
Multi-slice mode	Interleaved
Series	Descending
Multi-band accel. factor	2

**Geometry - AutoAlign**

Slice group	1
Position	L1.7 P11.0 H4.6 mm
Orientation	T > C-38.4 > S-0.7
Phase enc. dir.	A >> P
AutoAlign	Head > Brain
Initial Position	L0.0 P7.8 H14.2
L	0.0 mm
P	7.8 mm
H	14.2 mm
Initial Rotation	0.35 deg
Initial Orientation	T > C
T > C	-24.0
> S	0.0

**Geometry - Saturation**

Fat suppr.	Fat sat.
Special sat.	None

**System - Miscellaneous**

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
Matrix Optimization	Off
AutoAlign	Head > Brain
Coil Select Mode	Default

**System - Adjustments**

B0 Shim mode	Standard
B1 Shim mode	TrueForm
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

**System - Adjust Volume**

Position	L1.7 P11.0 H4.6 mm
Orientation	T > C-38.4 > S-0.7
Rotation	2.74 deg
A >> P	176 mm
R >> L	186 mm
F >> H	36 mm
Reset	Off

**System - Tx/Rx**

Frequency 1H	297.162474 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

**Physio - Signal1**

1st Signal/Mode	None
TR	2181 ms
Multi-band accel. factor	2

**BOLD**

GLM Statistics	Off
Dynamic t-maps	Off
Ignore meas. at start	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	20
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active
Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Meas[16]	Active
Meas[17]	Active
Meas[18]	Active
Meas[19]	Active
Meas[20]	Active
Motion correction	Off
Spatial filter	Off
Measurements	500
Delay in TR	0 ms
Multiple series	Off

**Sequence - Part 1**

Introduction	On
Contrasts	1
Flow comp.	No
Multi-slice mode	Interleaved
Free echo spacing	Off
Echo spacing	1.03 ms
Bandwidth	1260 Hz/Px

**Sequence - Part 2**

EPI factor	234
Gradient mode	Fast
RF spoiling	Off

**Sequence - Special**

Excite pulse duration	6000 us
Single-band images	Off
MB LeakBlock kernel	On
MB dual kernel	Off
MB RF phase scramble	Off
SENSE1 coil combine	On
Invert RO/PE polarity	Off
PF omits higher k-space	Off
Disable freq. update	Off
Force equal slice timing	Off
Online multi-band recon.	Online
FFT scale factor	1.00
Fat saturation FA	110.0 deg
GRE iPAT ref. FA	12.0 deg
Physio recording	DICOM
Triggering scheme	Standard

\\CRC\protocols\studies\v5motion\cmrr\_mbep2d\_p4\_mb2\_750um\_GRAPPA-GRE\_FA75\_Band\_TR21  
81\_rengo

TA: 18:45 PM: FIX Voxel size: 0.8×0.8×0.8 mmPAT: 4 Rel. SNR: 1.00 : epfid

### Properties

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	Off
Start measurements	Single measurement

### Routine

Slice group	1
Slices	48
Dist. factor	0 %
Position	L0.0 P7.8 H14.2 mm
Orientation	T > C-24.0
Phase enc. dir.	A >> P
AutoAlign	Head > Brain
Phase oversampling	0 %
FoV read	186 mm
FoV phase	94.4 %
Slice thickness	0.75 mm
TR	2181 ms
TE	25.00 ms
Multi-band accel. factor	2
Filter	Raw filter
Coil elements	A32

### Contrast - Common

TR	2181 ms
TE	25.00 ms
MTC	Off
Magn. preparation	None
Flip angle	75 deg
Fat suppr.	Fat sat.

### Contrast - Dynamic

Averaging mode	Long term
Reconstruction	Magnitude
Measurements	500
Delay in TR	0 ms
Multiple series	Off

### Resolution - Common

FoV read	186 mm
FoV phase	94.4 %
Slice thickness	0.75 mm
Base resolution	248
Phase resolution	100 %
Phase partial Fourier	6/8
Interpolation	Off

### Resolution - iPAT

PAT mode	GRAPPA
Accel. factor PE	4
Ref. lines PE	48

### Resolution - iPAT

Reference scan mode	GRE
---------------------	-----

### Resolution - Filter Image

Distortion Corr.	Off
Prescan Normalize	Off

### Resolution - Filter Rawdata

Raw filter	On
Elliptical filter	Off
Hamming	Off

### Geometry - Common

Slice group	1
Slices	48
Dist. factor	0 %
Position	L0.0 P7.8 H14.2 mm
Orientation	T > C-24.0
Phase enc. dir.	A >> P
FoV read	186 mm
FoV phase	94.4 %
Slice thickness	0.75 mm
TR	2181 ms
Multi-slice mode	Interleaved
Series	Descending
Multi-band accel. factor	2

### Geometry - AutoAlign

Slice group	1
Position	L0.0 P7.8 H14.2 mm
Orientation	T > C-24.0
Phase enc. dir.	A >> P
AutoAlign	Head > Brain
Initial Position	L0.0 P7.8 H14.2
L	0.0 mm
P	7.8 mm
H	14.2 mm
Initial Rotation	0.35 deg
Initial Orientation	T > C
T > C	-24.0
> S	0.0

### Geometry - Saturation

Fat suppr.	Fat sat.
Special sat.	None

### System - Miscellaneous

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
Matrix Optimization	Off
AutoAlign	Head > Brain
Coil Select Mode	Default

**System - Adjustments**

B0 Shim mode	Standard
B1 Shim mode	TrueForm
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

**System - Adjust Volume**

Position	L0.0 P7.8 H14.2 mm
Orientation	T > C-24.0
Rotation	0.35 deg
A >> P	176 mm
R >> L	186 mm
F >> H	36 mm
Reset	Off

**System - Tx/Rx**

Frequency 1H	297.162474 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

**Physio - Signal1**

1st Signal/Mode	None
TR	2181 ms
Multi-band accel. factor	2

**BOLD**

GLM Statistics	Off
Dynamic t-maps	Off
Ignore meas. at start	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	20
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active
Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Meas[16]	Active
Meas[17]	Active
Meas[18]	Active
Meas[19]	Active
Meas[20]	Active
Motion correction	Off
Spatial filter	Off
Measurements	500
Delay in TR	0 ms
Multiple series	Off

**Sequence - Part 1**

Introduction	On
Contrasts	1
Flow comp.	No
Multi-slice mode	Interleaved
Free echo spacing	Off
Echo spacing	1.03 ms
Bandwidth	1260 Hz/Px

**Sequence - Part 2**

EPI factor	234
Gradient mode	Fast
RF spoiling	Off

**Sequence - Special**

Excite pulse duration	6000 us
Single-band images	Off
MB LeakBlock kernel	On
MB dual kernel	Off
MB RF phase scramble	Off
SENSE1 coil combine	On
Invert RO/PE polarity	Off
PF omits higher k-space	Off
Disable freq. update	Off
Force equal slice timing	Off
Online multi-band recon.	Online
FFT scale factor	1.00
Fat saturation FA	110.0 deg
GRE iPAT ref. FA	12.0 deg
Physio recording	DICOM
Triggering scheme	Standard

\\CRC\protocols\studies\v5motion\cmrr\_mbep2d\_p4\_mb2\_750um\_GRAPPA-GRE\_FA75\_Band\_TR21  
81\_reference

TA: 18:45 PM: FIX Voxel size: 0.8×0.8×0.8 mmPAT: 4 Rel. SNR: 1.00 : epfid

### Properties

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	Off
Start measurements	Single measurement

### Routine

Slice group	1
Slices	48
Dist. factor	0 %
Position	L1.7 P11.0 H4.6 mm
Orientation	T > C-38.4 > S-0.7
Phase enc. dir.	A >> P
AutoAlign	Head > Brain
Phase oversampling	0 %
FoV read	186 mm
FoV phase	94.4 %
Slice thickness	0.75 mm
TR	2181 ms
TE	25.00 ms
Multi-band accel. factor	2
Filter	Raw filter
Coil elements	A32

### Contrast - Common

TR	2181 ms
TE	25.00 ms
MTC	Off
Magn. preparation	None
Flip angle	75 deg
Fat suppr.	Fat sat.

### Contrast - Dynamic

Averaging mode	Long term
Reconstruction	Magnitude
Measurements	500
Delay in TR	0 ms
Multiple series	Off

### Resolution - Common

FoV read	186 mm
FoV phase	94.4 %
Slice thickness	0.75 mm
Base resolution	248
Phase resolution	100 %
Phase partial Fourier	6/8
Interpolation	Off

### Resolution - iPAT

PAT mode	GRAPPA
Accel. factor PE	4
Ref. lines PE	48

### Resolution - iPAT

Reference scan mode	GRE
---------------------	-----

### Resolution - Filter Image

Distortion Corr.	Off
Prescan Normalize	Off

### Resolution - Filter Rawdata

Raw filter	On
Elliptical filter	Off
Hamming	Off

### Geometry - Common

Slice group	1
Slices	48
Dist. factor	0 %
Position	L1.7 P11.0 H4.6 mm
Orientation	T > C-38.4 > S-0.7
Phase enc. dir.	A >> P
FoV read	186 mm
FoV phase	94.4 %
Slice thickness	0.75 mm
TR	2181 ms
Multi-slice mode	Interleaved
Series	Descending
Multi-band accel. factor	2

### Geometry - AutoAlign

Slice group	1
Position	L1.7 P11.0 H4.6 mm
Orientation	T > C-38.4 > S-0.7
Phase enc. dir.	A >> P
AutoAlign	Head > Brain
Initial Position	L0.0 P7.8 H14.2
L	0.0 mm
P	7.8 mm
H	14.2 mm
Initial Rotation	0.35 deg
Initial Orientation	T > C
T > C	-24.0
> S	0.0

### Geometry - Saturation

Fat suppr.	Fat sat.
Special sat.	None

### System - Miscellaneous

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
Matrix Optimization	Off
AutoAlign	Head > Brain
Coil Select Mode	Default

**System - Adjustments**

B0 Shim mode	Standard
B1 Shim mode	TrueForm
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

**System - Adjust Volume**

Position	L1.7 P11.0 H4.6 mm
Orientation	T > C-38.4 > S-0.7
Rotation	2.74 deg
A >> P	176 mm
R >> L	186 mm
F >> H	36 mm
Reset	Off

**System - Tx/Rx**

Frequency 1H	297.162474 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

**Physio - Signal1**

1st Signal/Mode	None
TR	2181 ms
Multi-band accel. factor	2

**BOLD**

GLM Statistics	Off
Dynamic t-maps	Off
Ignore meas. at start	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	20
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active
Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Meas[16]	Active
Meas[17]	Active
Meas[18]	Active
Meas[19]	Active
Meas[20]	Active
Motion correction	Off
Spatial filter	Off
Measurements	500
Delay in TR	0 ms
Multiple series	Off

**Sequence - Part 1**

Introduction	On
Contrasts	1
Flow comp.	No
Multi-slice mode	Interleaved
Free echo spacing	Off
Echo spacing	1.03 ms
Bandwidth	1260 Hz/Px

**Sequence - Part 2**

EPI factor	234
Gradient mode	Fast
RF spoiling	Off

**Sequence - Special**

Excite pulse duration	6000 us
Single-band images	Off
MB LeakBlock kernel	On
MB dual kernel	Off
MB RF phase scramble	Off
SENSE1 coil combine	On
Invert RO/PE polarity	Off
PF omits higher k-space	Off
Disable freq. update	Off
Force equal slice timing	Off
Online multi-band recon.	Online
FFT scale factor	1.00
Fat saturation FA	110.0 deg
GRE iPAT ref. FA	12.0 deg
Physio recording	DICOM
Triggering scheme	Standard



\\CRC\protocols\studies\v5motion\rslh\_ep3d\_vaso\_4e\_axial

TA: 7:46 PM: REF Voxel size: 0.8×0.8×0.8 mmPAT: 3 Rel. SNR: 1.00 : 684e6dda

**Properties**

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	Off
Start measurements	Single measurement

**Resolution - Common**

Slice partial Fourier	Off
Interpolation	Off

**Resolution - iPAT**

PAT mode	CAIPIRINHA
Acc. factor PE	3
Ref. lines PE	75
Acc. factor 3D	1
Ref. lines 3D	22
CAIPI 3D Shift	0
Reference Scan Mode	GRE/separate
CAIPIRINHA mode	Free

**Routine**

Slab group	1
Slabs	1
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
Slab Scale	-10 %
Slices per slab	24
FoV read	186 mm
FoV phase	100.0 %
Slice thickness	0.75 mm
TR 1	71.3 ms
TR 2	4529 ms
TE 1	24.70 ms
Averages	1
Multi-echo Shots	1
Filter	Distortion Corr.(3D)
Coil elements	A32

**Resolution - Filter Image**

Image Filter	Off
Distortion Corr.	On
Mode	3D
Unfiltered images	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off

**Resolution - Filter Rawdata**

Raw filter	Off
Elliptical filter	Off

**Geometry - Common**

Slab group	1
Slabs	1
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
Slab Scale	-10 %
Slices per slab	24
FoV read	186 mm
FoV phase	100.0 %
Slice thickness	0.75 mm
TR 1	71.3 ms
TR 2	4529 ms
Multi-slice mode	Interleaved
Series	Ascending
Multi-echo Shots	1

**Contrast - Common**

TR 1	71.3 ms
TR 2	4529 ms
TE 1	24.70 ms
Multi-echo spacing	66.03 ms
Magn. preparation	Non-sel. HSN IR
TI 1	1515.6 ms
TI 2	3226.8 ms
Flip angle	60 deg
Fat suppr.	Fat sat.
Magn. Prep. Shots	1

**Contrast - Dynamic**

Averages	1
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	100
Pause after meas.	0.0 s

**Resolution - Common**

FoV read	186 mm
FoV phase	100.0 %
Slice thickness	0.75 mm
Base resolution	248
Phase resolution	100 %
Slice resolution	100 %
Phase partial Fourier	6/8

**Geometry - AutoAlign**

Slab group	1
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
Initial Position	Isocenter
L	0.0 mm
P	0.0 mm
H	0.0 mm
Initial Rotation	0.00 deg
Initial Orientation	Transversal

**Geometry - Saturation**

Saturation mode	Standard
Fat suppr.	Fat sat.

**System - Miscellaneous**

Positioning mode	REF
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
Save uncombined	Off
Matrix Optimization	Off
AutoAlign	---
Coil Select Mode	Default

**System - Adjustments**

B0 Shim mode	Standard
B1 Shim mode	TrueForm
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

**System - Adjust Volume**

Position	Isocenter
Orientation	Transversal
Rotation	0.00 deg
A >> P	186 mm
R >> L	186 mm
F >> H	18 mm
Reset	Off

**System - Tx/Rx**

Frequency 1H	297.162474 MHz
Correction factor	1
Gain	Low
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

**Sequence - Part 1**

Introduction	On
Dimension	3D
Reordering	Linear
Asymmetric echo	Off
Contrasts	1
Multi-slice mode	Interleaved
Echo spacing	1.05 ms
Bandwidth	1062 Hz/Px

**Sequence - Part 2**

EPI factor	62
Segmentation	1
RF pulse type	Normal
Gradient mode	Normal
Excitation	Slab-sel.
RF spoiling	On
Turbo factor	24

**Sequence - Special**

PATRef FA	3 deg
RF duration	3500 us
RF BWT product	15
Ernst T1	1200 ms
PATRef prep. shots	10
Volume dummy shots	0

**Sequence - Special**

Dummy Measurements	0
Invert PE	Off
Min. TE if PF	On
Echo Time Shift	On
Ramp Sampling	On
NORDIC	Off
Water Exc.	-none-
External PC	per Series
Saturation RF	per Shot
EPI rise time factor	1.10
Mosaic DICOMs	On
Modify IcePAT	On
HSN RF power scale	3.00
Inversion Delay	650 ms
Relaxation Delay	0 ms
Var. FA /MAGEC	4

**Sequence - Assistant**

Mode	Off
------	-----

\\CRC\protocols\studies\v5motion\rslh\_ep3d\_bold\_4e\_axial

TA: 3:33 PM: REF Voxel size: 0.8×0.8×0.8 mmPAT: 3 Rel. SNR: 1.00 : 684e6dda

**Properties**

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	Off
Start measurements	Single measurement

**Routine**

Slab group	1
Slabs	1
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
Slab Scale	-10 %
Slices per slab	24
FoV read	186 mm
FoV phase	100.0 %
Slice thickness	0.75 mm
TR 1	74.4 ms
TR 2	2004 ms
TE 1	25.70 ms
Averages	1
Multi-echo Shots	1
Filter	Distortion Corr.(3D)
Coil elements	A32

**Contrast - Common**

TR 1	74.4 ms
TR 2	2004 ms
TE 1	25.70 ms
Multi-echo spacing	69.13 ms
Magn. preparation	None
Flip angle	17 deg
Fat suppr.	Fat sat.
Magn. Prep. Shots	1

**Contrast - Dynamic**

Averages	1
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	100
Pause after meas.	0.0 s

**Resolution - Common**

FoV read	186 mm
FoV phase	100.0 %
Slice thickness	0.75 mm
Base resolution	248
Phase resolution	100 %
Slice resolution	100 %
Phase partial Fourier	6/8
Slice partial Fourier	Off
Interpolation	Off

**Resolution - iPAT**

PAT mode	CAIPIRINHA
Acc. factor PE	3
Ref. lines PE	75
Acc. factor 3D	1
Ref. lines 3D	22
CAIPI 3D Shift	0
Reference Scan Mode	GRE/separate
CAIPIRINHA mode	Free

**Resolution - Filter Image**

Image Filter	Off
Distortion Corr.	On
Mode	3D
Unfiltered images	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off

**Resolution - Filter Rawdata**

Raw filter	Off
Elliptical filter	Off

**Geometry - Common**

Slab group	1
Slabs	1
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
Slab Scale	-10 %
Slices per slab	24
FoV read	186 mm
FoV phase	100.0 %
Slice thickness	0.75 mm
TR 1	74.4 ms
TR 2	2004 ms
Multi-slice mode	Interleaved
Series	Ascending
Multi-echo Shots	1

**Geometry - AutoAlign**

Slab group	1
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
Initial Position	Isocenter
L	0.0 mm
P	0.0 mm
H	0.0 mm
Initial Rotation	0.00 deg
Initial Orientation	Transversal

**Geometry - Saturation**

Saturation mode	Standard
Fat suppr.	Fat sat.

**System - Miscellaneous**

Positioning mode	REF
Table position	H
Table position	0 mm

**System - Miscellaneous**

MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
Save uncombined	Off
Matrix Optimization	Off
AutoAlign	---
Coil Select Mode	Default

**System - Adjustments**

B0 Shim mode	Standard
B1 Shim mode	TrueForm
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

**System - Adjust Volume**

Position	Isocenter
Orientation	Transversal
Rotation	0.00 deg
A >> P	186 mm
R >> L	186 mm
F >> H	18 mm
Reset	Off

**System - Tx/Rx**

Frequency 1H	297.162474 MHz
Correction factor	1
Gain	Low
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

**Sequence - Part 1**

Introduction	On
Dimension	3D
Reordering	Linear
Asymmetric echo	Off
Contrasts	1
Multi-slice mode	Interleaved
Echo spacing	1.1 ms
Bandwidth	1008 Hz/Px

**Sequence - Part 2**

EPI factor	62
Segmentation	1
RF pulse type	Normal
Gradient mode	Normal
Excitation	Slab-sel.
RF spoiling	On
Turbo factor	24

**Sequence - Special**

PATRef FA	3 deg
RF duration	3500 us
RF BWT product	15
Ernst T1	1200 ms
PATRef prep. shots	10
Volume dummy shots	0
Dummy Measurements	0
Invert PE	Off
Min. TE if PF	On

**Sequence - Special**

Echo Time Shift	On
Ramp Sampling	On
NORDIC	Off
Water Exc.	-none-
External PC	per Series
Saturation RF	per Shot
EPI rise time factor	1.10
Mosaic DICOMs	On
Modify IcePAT	On
Var. FA /MAGEC	0

**Sequence - Assistant**

Mode	Off
------	-----

\\CRC\protocols\studies\v5motion\rslh\_ep3d\_vaso\_4e\_sagital

TA: 7:46 PM: REF Voxel size: 0.8×0.8×0.8 mmPAT: 3 Rel. SNR: 1.00 : 684e6dda

**Properties**

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	Off
Start measurements	Single measurement

**Resolution - Common**

Slice partial Fourier	Off
Interpolation	Off

**Resolution - iPAT**

PAT mode	CAIPIRINHA
Acc. factor PE	3
Ref. lines PE	75
Acc. factor 3D	1
Ref. lines 3D	22
CAIPI 3D Shift	0
Reference Scan Mode	GRE/separate
CAIPIRINHA mode	Free

**Routine**

Slab group	1
Slabs	1
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
Slab Scale	-10 %
Slices per slab	24
FoV read	186 mm
FoV phase	100.0 %
Slice thickness	0.75 mm
TR 1	71.3 ms
TR 2	4529 ms
TE 1	24.70 ms
Averages	1
Multi-echo Shots	1
Filter	Distortion Corr.(3D)
Coil elements	A32

**Resolution - Filter Image**

Image Filter	Off
Distortion Corr.	On
Mode	3D
Unfiltered images	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off

**Resolution - Filter Rawdata**

Raw filter	Off
Elliptical filter	Off

**Geometry - Common**

Slab group	1
Slabs	1
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
Slab Scale	-10 %
Slices per slab	24
FoV read	186 mm
FoV phase	100.0 %
Slice thickness	0.75 mm
TR 1	71.3 ms
TR 2	4529 ms
Multi-slice mode	Interleaved
Series	Ascending
Multi-echo Shots	1

**Contrast - Common**

TR 1	71.3 ms
TR 2	4529 ms
TE 1	24.70 ms
Multi-echo spacing	66.03 ms
Magn. preparation	Non-sel. HSN IR
TI 1	1515.6 ms
TI 2	3226.8 ms
Flip angle	60 deg
Fat suppr.	Fat sat.
Magn. Prep. Shots	1

**Contrast - Dynamic**

Averages	1
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	100
Pause after meas.	0.0 s

**Resolution - Common**

FoV read	186 mm
FoV phase	100.0 %
Slice thickness	0.75 mm
Base resolution	248
Phase resolution	100 %
Slice resolution	100 %
Phase partial Fourier	6/8

**Geometry - AutoAlign**

Slab group	1
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
Initial Position	Isocenter
L	0.0 mm
P	0.0 mm
H	0.0 mm
Initial Rotation	0.00 deg
Initial Orientation	Transversal

**Geometry - Saturation**

Saturation mode	Standard
Fat suppr.	Fat sat.

**System - Miscellaneous**

Positioning mode	REF
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
Save uncombined	Off
Matrix Optimization	Off
AutoAlign	---
Coil Select Mode	Default

**System - Adjustments**

B0 Shim mode	Standard
B1 Shim mode	TrueForm
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

**System - Adjust Volume**

Position	Isocenter
Orientation	Transversal
Rotation	0.00 deg
A >> P	186 mm
R >> L	186 mm
F >> H	18 mm
Reset	Off

**System - Tx/Rx**

Frequency 1H	297.162474 MHz
Correction factor	1
Gain	Low
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

**Sequence - Part 1**

Introduction	On
Dimension	3D
Reordering	Linear
Asymmetric echo	Off
Contrasts	1
Multi-slice mode	Interleaved
Echo spacing	1.05 ms
Bandwidth	1062 Hz/Px

**Sequence - Part 2**

EPI factor	62
Segmentation	1
RF pulse type	Normal
Gradient mode	Normal
Excitation	Slab-sel.
RF spoiling	On
Turbo factor	24

**Sequence - Special**

PATRef FA	3 deg
RF duration	3500 us
RF BWT product	15
Ernst T1	1200 ms
PATRef prep. shots	10
Volume dummy shots	0

**Sequence - Special**

Dummy Measurements	0
Invert PE	Off
Min. TE if PF	On
Echo Time Shift	On
Ramp Sampling	On
NORDIC	Off
Water Exc.	-none-
External PC	per Series
Saturation RF	per Shot
EPI rise time factor	1.10
Mosaic DICOMs	On
Modify IcePAT	On
HSN RF power scale	3.00
Inversion Delay	650 ms
Relaxation Delay	0 ms
Var. FA /MAGEC	4

**Sequence - Assistant**

Mode	Off
------	-----

\\CRC\protocols\studies\v5motion\rslh\_ep3d\_bold\_4e\_sagital

TA: 2:43 PM: REF Voxel size: 0.7×0.7×0.8 mmPAT: 3 Rel. SNR: 1.00 : 684e6dda

**Properties**

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	Off
Start measurements	Single measurement

**Routine**

Slab group	1
Slabs	1
Position	R46.6 A11.5 H0.0 mm
Orientation	Sagittal
Phase enc. dir.	A >> P
AutoAlign	---
Slab Scale	-10 %
Slices per slab	24
FoV read	130 mm
FoV phase	100.0 %
Slice thickness	0.75 mm
TR 1	53.6 ms
TR 2	1504 ms
TE 1	19.20 ms
Averages	1
Multi-echo Shots	1
Filter	Distortion Corr.(3D)
Coil elements	A32

**Contrast - Common**

TR 1	53.6 ms
TR 2	1504 ms
TE 1	19.20 ms
Multi-echo spacing	48.51 ms
Magn. preparation	None
Flip angle	17 deg
Fat suppr.	Fat sat.
Magn. Prep. Shots	1

**Contrast - Dynamic**

Averages	1
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	100
Pause after meas.	0.0 s

**Resolution - Common**

FoV read	130 mm
FoV phase	100.0 %
Slice thickness	0.75 mm
Base resolution	174
Phase resolution	100 %
Slice resolution	100 %
Phase partial Fourier	6/8
Slice partial Fourier	Off
Interpolation	Off

**Resolution - iPAT**

PAT mode	CAIPIRINHA
Acc. factor PE	3
Ref. lines PE	75
Acc. factor 3D	1
Ref. lines 3D	22
CAIPI 3D Shift	0
Reference Scan Mode	GRE/separate
CAIPIRINHA mode	Free

**Resolution - Filter Image**

Image Filter	Off
Distortion Corr.	On
Mode	3D
Unfiltered images	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off

**Resolution - Filter Rawdata**

Raw filter	Off
Elliptical filter	Off

**Geometry - Common**

Slab group	1
Slabs	1
Position	R46.6 A11.5 H0.0 mm
Orientation	Sagittal
Phase enc. dir.	A >> P
Slab Scale	-10 %
Slices per slab	24
FoV read	130 mm
FoV phase	100.0 %
Slice thickness	0.75 mm
TR 1	53.6 ms
TR 2	1504 ms
Multi-slice mode	Interleaved
Series	Ascending
Multi-echo Shots	1

**Geometry - AutoAlign**

Slab group	1
Position	R46.6 A11.5 H0.0 mm
Orientation	Sagittal
Phase enc. dir.	A >> P
AutoAlign	---
Initial Position	R46.6 A11.5 H0.0
R	46.6 mm
A	11.5 mm
H	0.0 mm
Initial Rotation	0.00 deg
Initial Orientation	Sagittal

**Geometry - Saturation**

Saturation mode	Standard
Fat suppr.	Fat sat.

**System - Miscellaneous**

Positioning mode	REF
Table position	H
Table position	0 mm

**System - Miscellaneous**

MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
Save uncombined	Off
Matrix Optimization	Off
AutoAlign	---
Coil Select Mode	Default

**System - Adjustments**

B0 Shim mode	Standard
B1 Shim mode	TrueForm
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

**System - Adjust Volume**

Position	R46.6 A11.5 H0.0 mm
Orientation	Sagittal
Rotation	0.00 deg
A >> P	130 mm
F >> H	130 mm
R >> L	18 mm
Reset	Off

**System - Tx/Rx**

Frequency 1H	297.162474 MHz
Correction factor	1
Gain	Low
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

**Sequence - Part 1**

Introduction	On
Dimension	3D
Reordering	Linear
Asymmetric echo	Off
Contrasts	1
Multi-slice mode	Interleaved
Echo spacing	1.1 ms
Bandwidth	1026 Hz/Px

**Sequence - Part 2**

EPI factor	43
Segmentation	1
RF pulse type	Normal
Gradient mode	Normal
Excitation	Slab-sel.
RF spoiling	On
Turbo factor	24

**Sequence - Special**

PATRef FA	3 deg
RF duration	3500 us
RF BWT product	15
Ernst T1	1200 ms
PATRef prep. shots	10
Volume dummy shots	0
Dummy Measurements	0
Invert PE	Off
Min. TE if PF	On

**Sequence - Special**

Echo Time Shift	On
Ramp Sampling	On
NORDIC	Off
Water Exc.	-none-
External PC	per Series
Saturation RF	per Shot
EPI rise time factor	1.10
Mosaic DICOMs	On
Modify IcePAT	On
Var. FA /MAGEC	0

**Sequence - Assistant**

Mode	Off
------	-----



\\CRC\protocols\studies\v5motion\t1\_mp2rage\_ax\_p2\_0.5mm\_slab

TA: 5:26 PM: FIX Voxel size: 0.5×0.5×0.5 mmPAT: 2 Rel. SNR: 1.00 : tfl

**Properties**

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	Off
Start measurements	Single measurement

**Routine**

Slab group	1
Slabs	1
Dist. factor	50 %
Position	L0.0 A10.3 H0.0 mm
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	72
FoV read	160 mm
FoV phase	100.0 %
Slice thickness	0.50 mm
TR	6000.0 ms
TE	2.45 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	A32

**Contrast - Common**

TR	6000.0 ms
TE	2.45 ms
Magn. preparation	Non-sel. IR
TI 1	900 ms
TI 2	2900 ms
Flip angle 1	6.0 deg
Flip angle 2	7.0 deg
Fat suppr.	None
Water suppr.	None

**Contrast - Dynamic**

Averages	1
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	1
Multiple series	Each measurement

**Resolution - Common**

FoV read	160 mm
FoV phase	100.0 %
Slice thickness	0.50 mm
Base resolution	320
Phase resolution	100 %
Slice resolution	100 %
Phase partial Fourier	Off

**Resolution - Common**

Slice partial Fourier	6/8
Interpolation	Off

**Resolution - iPAT**

PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	36
Accel. factor 3D	1
Reference scan mode	Integrated

**Resolution - Filter Image**

Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off

**Resolution - Filter Rawdata**

Raw filter	Off
Elliptical filter	Off

**Geometry - Common**

Slab group	1
Slabs	1
Dist. factor	50 %
Position	L0.0 A10.3 H0.0 mm
Orientation	Transversal
Phase enc. dir.	A >> P
Slice oversampling	0.0 %
Slices per slab	72
FoV read	160 mm
FoV phase	100.0 %
Slice thickness	0.50 mm
TR	6000.0 ms
Multi-slice mode	Single shot
Series	Interleaved
Concatenations	1

**Geometry - AutoAlign**

Slab group	1
Position	L0.0 A10.3 H0.0 mm
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
Initial Position	L0.0 A10.3 H0.0
L	0.0 mm
A	10.3 mm
H	0.0 mm
Initial Rotation	0.00 deg
Initial Orientation	Transversal

**Geometry - Navigator****System - Miscellaneous**

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P

**System - Miscellaneous**

Transversal	F >> H
Coil Combine Mode	Sum of Squares
Save uncombined	Off
Matrix Optimization	Off
AutoAlign	---
Coil Select Mode	Default

**System - Adjustments**

B0 Shim mode	Brain
B1 Shim mode	TrueForm
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

**System - Adjust Volume**

Position	L0.0 A10.3 H0.0 mm
Orientation	Transversal
Rotation	0.00 deg
A >> P	160 mm
R >> L	160 mm
F >> H	36 mm
Reset	Off

**System - Tx/Rx**

Frequency 1H	297.162474 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

**Physio - Signal1**

1st Signal/Mode	None
TR	6000.0 ms
Concatenations	1

**Physio - Cardiac**

Magn. preparation	Non-sel. IR
TI 1	900 ms
TI 2	2900 ms
Fat suppr.	None
Dark blood	Off
FoV read	160 mm
FoV phase	100.0 %
Phase resolution	100 %

**Physio - PACE**

Resp. control	Off
Concatenations	1

**Inline - Common**

Subtract	Off
Measurements	1
StdDev	Off
Save original images	On

**Inline - MIP**

MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On

**Inline - Composing**

Distortion Corr.	Off
------------------	-----

**Sequence - Part 1**

Introduction	On
Dimension	3D
Elliptical scanning	Off
Reordering	Linear rot.
Asymmetric echo	Allowed
Flow comp.	No
Multi-slice mode	Single shot
Echo spacing	7.4 ms
Bandwidth	200 Hz/Px

**Sequence - Part 2**

RF pulse type	Normal
Gradient mode	Fast*
Excitation	Non-sel.
RF spoiling	On
Incr. Gradient spoiling	Off
Turbo factor	319

**Sequence - Nuclei**

TX/RX Nucleus	1H
TX/RX delta frequency	0 Hz
TX Nucleus	None
TX delta frequency	0 Hz
Coil elements	A32

**Sequence - Assistant**

Mode	Off
------	-----

\\CRC\protocols\studies\v5motion\rslh\_ep3d\_vaso\_4e\_axial\_p05

TA: 14:14 PM: REF Voxel size: 0.5×0.5×0.5 mmPAT: 3 Rel. SNR: 1.00 : 684e6dda

**Properties**

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	Off
Start measurements	Single measurement

**Routine**

Slab group	1
Slabs	1
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
Slab Scale	-10 %
Slices per slab	24
FoV read	186 mm
FoV phase	100.0 %
Slice thickness	0.50 mm
TR 1	78.1 ms
TR 2	8391 ms
TE 1	27.60 ms
Averages	1
Multi-echo Shots	1
Filter	Distortion Corr.(3D)
Coil elements	A32

**Contrast - Common**

TR 1	78.1 ms
TR 2	8391 ms
TE 1	27.60 ms
Multi-echo spacing	71.74 ms
Magn. preparation	Non-sel. HSN IR
TI 1	1884.4 ms
TI 2	5633.2 ms
Flip angle	40 deg
Fat suppr.	Fat sat.
Magn. Prep. Shots	1

**Contrast - Dynamic**

Averages	1
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	100
Pause after meas.	0.0 s

**Resolution - Common**

FoV read	186 mm
FoV phase	100.0 %
Slice thickness	0.50 mm
Base resolution	374
Phase resolution	100 %
Slice resolution	100 %
Phase partial Fourier	6/8

**Resolution - Common**

Slice partial Fourier	Off
Interpolation	Off

**Resolution - iPAT**

PAT mode	CAIPIRINHA
Acc. factor PE	3
Ref. lines PE	75
Acc. factor 3D	1
Ref. lines 3D	22
CAIPI 3D Shift	0
Reference Scan Mode	GRE/separate
CAIPIRINHA mode	Free

**Resolution - Filter Image**

Image Filter	Off
Distortion Corr.	On
Mode	3D
Unfiltered images	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off

**Resolution - Filter Rawdata**

Raw filter	Off
Elliptical filter	Off

**Geometry - Common**

Slab group	1
Slabs	1
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
Slab Scale	-10 %
Slices per slab	24
FoV read	186 mm
FoV phase	100.0 %
Slice thickness	0.50 mm
TR 1	78.1 ms
TR 2	8391 ms
Multi-slice mode	Interleaved
Series	Ascending
Multi-echo Shots	1

**Geometry - AutoAlign**

Slab group	1
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
Initial Position	Isocenter
L	0.0 mm
P	0.0 mm
H	0.0 mm
Initial Rotation	0.00 deg
Initial Orientation	Transversal

**Geometry - Saturation**

Saturation mode	Standard
Fat suppr.	Fat sat.

**System - Miscellaneous**

Positioning mode	REF
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
Save uncombined	Off
Matrix Optimization	Off
AutoAlign	---
Coil Select Mode	Default

**System - Adjustments**

B0 Shim mode	Standard
B1 Shim mode	TrueForm
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

**System - Adjust Volume**

Position	Isocenter
Orientation	Transversal
Rotation	0.00 deg
A >> P	186 mm
R >> L	186 mm
F >> H	12 mm
Reset	Off

**System - Tx/Rx**

Frequency 1H	297.162474 MHz
Correction factor	1
Gain	Low
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

**Sequence - Part 1**

Introduction	On
Dimension	3D
Reordering	Linear
Asymmetric echo	Off
Contrasts	1
Multi-slice mode	Interleaved
Echo spacing	1.49 ms
Bandwidth	742 Hz/Px

**Sequence - Part 2**

EPI factor	47
Segmentation	2
RF pulse type	Normal
Gradient mode	Normal
Excitation	Slab-sel.
RF spoiling	On
Turbo factor	48

**Sequence - Special**

PATRef FA	3 deg
RF duration	3500 us
RF BWT product	15
Ernst T1	1200 ms
PATRef prep. shots	10
Volume dummy shots	0

**Sequence - Special**

Dummy Measurements	0
Invert PE	Off
Min. TE if PF	On
Echo Time Shift	On
Ramp Sampling	On
NORDIC	Off
Water Exc.	-none-
External PC	per Series
Saturation RF	per Shot
EPI rise time factor	1.10
Mosaic DICOMs	On
Modify IcePAT	On
HSN RF power scale	3.00
Inversion Delay	0 ms
Relaxation Delay	0 ms
Var. FA /MAGEC	0

**Sequence - Assistant**

Mode	Off
------	-----