

# SIEMENS MAGNETOM Investigational\_Device\_7T syngo MR B17

\\USER\Renzo\FAIR\neg\_BOLD\MAFI

TA: 2:17

Voxel size: 6.0x6.0x6.0 mm

Rel. SNR: 1.00

USER: Renzo\MAFI

## Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

## Routine

Slab group 1	
Slabs	1
Dist. factor	0 %
Position	L0.0 A28.5 F1.2
Orientation	Sagittal
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	28
FoV read	192 mm
FoV phase	100.0 %
Slice thickness	6.00 mm
TR	150 ms
TE 1	1.080 ms
TE 2	1.080 ms
TE 3	2.09 ms
TE 4	3.100 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	A24

## Contrast

MTC	Off
Flip angle	90 deg
Fat suppr.	None
Water suppr.	None
Averaging mode	Short term
Reconstruction	Magn./Phase
Measurements	1
Multiple series	Off

## Resolution

Base resolution	32
Phase resolution	100 %
Slice resolution	100 %
Phase partial Fourier	Off
Interpolation	Off
Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off
Raw filter	Off
Elliptical filter	Off

## Geometry

Multi-slice mode	Sequential
Series	Ascending
Special sat.	None
Table position	H
Table position	0 mm
Inline Composing	Off

## System

V24	Off
A24	On
Positioning mode	REF
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Adaptive Combine
AutoAlign	---
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
! Ref. amplitude 1H	270.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	L0.0 A28.5 F1.2
Orientation	Sagittal
Rotation	0.00 deg
F >> H	192 mm
A >> P	192 mm
R >> L	168 mm

## Physio

1st Signal/Mode	None
-----------------	------

## Inline

Subtract	Off
Std-Dev-Sag	Off
Std-Dev-Cor	Off
Std-Dev-Tra	Off
Std-Dev-Time	Off
MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On

## Sequence

Introduction	Off
Dimension	3D
Contrasts	4
Bandwidth	1560 Hz/Px
Gradient mode	Fast
RF spoiling	On
Online ICE	Off
RF pulse type	square
Pulse duration	500 us
Spoil me!	On
TR2/TR1	5

# SIEMENS MAGNETOM Investigational\_Device\_7T syngo MR B17

N dummy TRs	20
Sample T1	1800 ms
Diffusion damping d= bD	0.6000
Diffusion coefficient D	2.2000 $\mu\text{m}^2/\text{ms}$
RF spoil phase increment	129.3 deg
Number of pulse shapes	1
<hr/>	
TX/RX Nucleus	1H
TX/RX delta frequency	0 Hz
TX Nucleus	None
TX delta frequency	0 Hz