

# SIEMENS MAGNETOM Investigational\_Device\_7T syngo MR B17

\\USER\RenzoHuber\Lonike\_20201118\executed\_VASO20201118\Checklist\_ok

TA: 1:11 PAT: Off Voxel size: 1.0x1.0x5.0 mm Rel. SNR: 1.00 SIEMENS: tfl

## Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	Off
Wait for user to start	Off
Start measurements	single

## Routine

Slice group 1	
Slices	11
Dist. factor	150 %
Position	R1.8 A19.9 F14.3
Orientation	Sagittal
Phase enc. dir.	A >> P
Rotation	0.00 deg
Slice group 2	
Slices	5
Dist. factor	100 %
Position	R2.3 A15.9 H1.4
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	0.00 deg
Slice group 3	
Slices	7
Dist. factor	200 %
Position	R1.5 A21.4 F6.6
Orientation	Coronal
Phase enc. dir.	R >> L
Rotation	0.00 deg
Phase oversampling	0 %
FoV read	200 mm
FoV phase	100.0 %
Slice thickness	5.0 mm
TR	3000 ms
TE	2.24 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	A32

## Contrast

Magn. preparation	Slice-sel. IR
TI	1100 ms
Flip angle	6 deg
Fat suppr.	None
Water suppr.	None
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	1
Multiple series	Each measurement

## Resolution

Base resolution	192
Phase resolution	100 %
Phase partial Fourier	Off
Interpolation	Off

PAT mode None

Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off
Raw filter	Off
Elliptical filter	Off

## Geometry

Multi-slice mode	Single shot
Series	Interleaved

Table position	H
Table position	0 mm
Inline Composing	Off

## System

V32	Off
A32	On
Positioning mode	REF
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Adaptive Combine
AutoAlign	---
Auto Coil Select	Default
Shim mode	Tune up
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
! Ref. amplitude 1H	220.000 V
Adjustment Tolerance	Auto
Adjust volume	
! Position	L0.0 A23.4 F1.3
! Orientation	Transversal
! Rotation	0.00 deg
! R >> L	350 mm
! A >> P	213 mm
! F >> H	189 mm

## Physio

1st Signal/Mode	None
Dark blood	Off
Resp. control	Off

## Inline

Subtract	Off
Std-Dev-Sag	Off
Std-Dev-Cor	Off
Std-Dev-Tra	Off
Std-Dev-Time	Off
MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On

## Sequence

Introduction	On
--------------	----

## SIEMENS MAGNETOM Investigational\_Device\_7T syngo MR B17

Dimension	2D
Asymmetric echo	Allowed
Bandwidth	240 Hz/Px
Flow comp.	No
Echo spacing	5.5 ms
<hr/>	
RF pulse type	Normal
Gradient mode	Fast
Excitation	Slice-sel.
RF spoiling	On

# SIEMENS MAGNETOM Investigational\_Device\_7T syngo MR B17

\\USER\RenzoHuber\Lonike\_20201118\executed\_VASO20201118\cubric\_mp2rage\_fatnav\_ICE900\_0p7\_pPF

TA: 11:58 PAT: 3 Voxel size: 0.7x0.7x0.7 mm Rel. SNR: 1.00 USER: cubric\_mp2rage\_fatnav\_ICE900

## Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

## Routine

Slab group 1	
Slabs	1
Dist. factor	50 %
Position	L0.0 A27.8 F7.0
Orientation	Sagittal
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	256
FoV read	210 mm
FoV phase	100.0 %
Slice thickness	0.70 mm
TR	6000 ms
TE	2.39 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	A32

## Contrast

Magn. preparation	Non-sel. IR
TI 1	800 ms
TI 2	2750 ms
Flip angle 1	4 deg
Flip angle 2	5 deg
Fat suppr.	None
Water suppr.	None
2nd Inversion Contrast	On
Averaging mode	Long term
Reconstruction	Magn./Phase
Measurements	1
Multiple series	Each measurement

## Resolution

Base resolution	300
Phase resolution	100 %
Slice resolution	100 %
Phase partial Fourier	Off
Slice partial Fourier	7/8
Interpolation	Off
PAT mode	GRAPPA
Accel. factor PE	3
Ref. lines PE	30
Accel. factor 3D	1
Reference scan mode	Integrated

Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off
Raw filter	Off
Elliptical filter	Off

## Geometry

Multi-slice mode	Single shot
Series	Interleaved
Table position	H
Table position	0 mm
Inline Composing	Off

## System

V32	Off
A32	On
Positioning mode	REF
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Adaptive Combine
AutoAlign	---
Auto Coil Select	Default
Shim mode	Tune up
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
? Ref. amplitude 1H	0.000 V
Adjustment Tolerance	Auto
Adjust volume	
! Position	L0.0 A22.7 H15.9
! Orientation	Transversal
! Rotation	0.00 deg
! R >> L	174 mm
! A >> P	180 mm
! F >> H	49 mm

## Physio

1st Signal/Mode	None
Dark blood	Off
Resp. control	Off

## Inline

Subtract	Off
Std-Dev-Sag	Off
Std-Dev-Cor	Off
Std-Dev-Tra	Off
Std-Dev-Time	Off
MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On

## Sequence

Introduction	On
Dimension	3D
Elliptical scanning	Off

# SIEMENS MAGNETOM Investigational\_Device\_7T syngo MR B17

Asymmetric echo	Allowed
Contrasts	1
Bandwidth	180 Hz/Px
Flow comp.	No
Echo spacing	7.1 ms
<hr/>	
RF pulse type	Fast
Gradient mode	Fast
Excitation	Non-sel.
RF spoiling	On
<hr/>	
FFT Scale Factor	100 %
LIN/PAR Swap	Off
Ext. INV Pulse	On
Flip Angle	2000
Uniform Image	On
T1 Map	On
Denoise Weighting	150
Acquire FatNavs	On
FatNav resolution	2 mm
FatNav flip angle	3.0 degrees
Train-FatNav delay	20 ms

# SIEMENS MAGNETOM Investigational\_Device\_7T syngo MR B17

\\USER\RenzoHuber\Lonike\_20201118\executed\_VASO20201118\VASO\_151\_0.8mm\_sag\_FA4

TA: 3:18 PAT: 3 Voxel size: 0.8x0.8x0.8 mm Rel. SNR: 1.00 USER: VASO\_151

## Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

## Routine

Slab group 1	
Slabs	1
Dist. factor	50 %
Position	R47.2 A27.2 H8.8
Orientation	S > C3.4
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
Slice oversampling	9.1 %
Slices per slab	22
FoV read	130.0 mm
FoV phase	133.3 %
Slice thickness	0.80 mm
TR	2475.80 ms
TE	24 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	A32

## Contrast

Perfusion mode	SS-SI VASO
TI2	700 ms
TI1	50 ms
TI1s	50 ms
Flip angle	4 deg
Fat suppr.	Fat sat.
Fat sat. mode	Strong
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	80
Delay in TR	0 ms
Multiple series	Off
Perfusion mode	PICORE Q2T
Inversion time 1	50 ms
Saturation stop time	50 ms
Inversion time 2	700.0 ms
Flow limit	100 cm/s

## Resolution

Base resolution	162
Phase resolution	100 %
Slice resolution	100 %
Phase partial Fourier	6/8
Slice partial Fourier	Off
Interpolation	Off

PAT mode	GRAPPA
Accel. factor PE	3
Ref. lines PE	45
Accel. factor 3D	1
Ref. lines 3D	24
Reference scan mode	Separate

Prescan Normalize	Off
Raw filter	Off
Elliptical filter	Off
Hamming	Off

## Geometry

Multi-slice mode	Interleaved
Series	Ascending
Special sat.	Parallel R
Gap	25.0 mm
Thickness	100 mm
Table position	H
Table position	0 mm
Inline Composing	Off

## System

V32	Off
A32	On
Positioning mode	REF
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Sum of Squares
AutoAlign	Head > Brain
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
! Ref. amplitude 1H	225.000 V
Adjustment Tolerance	Auto
Adjust volume	
! Position	R45.1 A26.5 H3.3
! Orientation	Sagittal
! Rotation	90.00 deg
! A >> P	174 mm
! F >> H	152 mm
! R >> L	39 mm

## Physio

1st Signal/Mode	None
-----------------	------

## BOLD

Motion correction	Off
Spatial filter	Off

## Sequence

Introduction	On
Dimension	3D
Reordering	Linear
Contrasts	1
Bandwidth	1144 Hz/Px
Free echo spacing	Off
Echo spacing	0.98 ms
EPI factor	216

# SIEMENS MAGNETOM Investigational\_Device\_7T syngo MR B17

RF pulse type	Normal
Gradient mode	Fast
Excitation	Slab-sel.
RF spoiling	On
<hr/>	
Ampl	180
MAGEC FA	12 in deg
ph.skip 4 Robert (the one)	1
MAGEC SS-SI?	On
Maxwell Correction	Off
log physio files	On
FFT scale	2.00
dummy prepscan time	3 s
z shim	0.00 mT/m*ms
RF duration	2200 us
RF BWTP	25.0
Renzo: Delta TI	71 ms
EFFECTIVE TR	59419 ms
PatPartitions	24
EPI phase correction	local
PAT refscan mode	Flash
FlashRef BaseRes	162
FlashRef BW	100 Hz/px
FlashRef TE	7000 us
FlashRef FA	5 deg
use CAIPI	Off

## SIEMENS MAGNETOM Investigational\_Device\_7T syngo MR B17

\\USER\RenzoHuber\Lonike\_20201118\executed\_VASO20201118\VASO\_151\_0.8mm\_sag\_FA26

TA: 3:18

PAT: 3

Voxel size: 0.8x0.8x0.8 mm

Rel. SNR: 1.00

USER: VASO\_151

## Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

## Routine

Slab group 1	
Slabs	1
Dist. factor	50 %
Position	R47.2 A27.2 H8.8
Orientation	S > C3.4
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
Slice oversampling	9.1 %
Slices per slab	22
FoV read	130.0 mm
FoV phase	133.3 %
Slice thickness	0.80 mm
TR	2475.80 ms
TE	24 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	A32

## Contrast

Perfusion mode	SS-SI VASO
TI2	700 ms
TI1	50 ms
TI1s	50 ms
Flip angle	26 deg
Fat suppr.	Fat sat.
Fat sat. mode	Strong

Averaging mode	Long term
Reconstruction	Magnitude
Measurements	80
Delay in TR	0 ms
Multiple series	Off

Perfusion mode	PICORE Q2T
Inversion time 1	50 ms
Saturation stop time	50 ms
Inversion time 2	700.0 ms
Flow limit	100 cm/s

## Resolution

Base resolution	162
Phase resolution	100 %
Slice resolution	100 %
Phase partial Fourier	6/8
Slice partial Fourier	Off
Interpolation	Off

PAT mode	GRAPPA
Accel. factor PE	3
Ref. lines PE	45
Accel. factor 3D	1
Ref. lines 3D	24
Reference scan mode	Separate

Prescan Normalize	Off
Raw filter	Off
Elliptical filter	Off
Hamming	Off

## Geometry

Multi-slice mode	Interleaved
Series	Ascending
Special sat.	Parallel R
Gap	25.0 mm
Thickness	100 mm
Table position	H
Table position	0 mm
Inline Composing	Off

## System

V32	Off
A32	On
Positioning mode	REF
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Sum of Squares
AutoAlign	Head > Brain
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
! Ref. amplitude 1H	220.000 V
Adjustment Tolerance	Auto
Adjust volume	
! Position	R45.1 A26.5 H3.3
! Orientation	Sagittal
! Rotation	90.00 deg
! A >> P	174 mm
! F >> H	152 mm
! R >> L	39 mm

## Physio

1st Signal/Mode	None
-----------------	------

## BOLD

Motion correction	Off
Spatial filter	Off

## Sequence

Introduction	On
Dimension	3D
Reordering	Linear
Contrasts	1
Bandwidth	1144 Hz/Px
Free echo spacing	Off
Echo spacing	0.98 ms
EPI factor	216

# SIEMENS MAGNETOM Investigational\_Device\_7T syngo MR B17

RF pulse type	Normal
Gradient mode	Fast
Excitation	Slab-sel.
RF spoiling	On
<hr/>	
Ampl	150
MAGEC FA	12 in deg
ph.skip 4 Robert (the one)	1
MAGEC SS-SI?	On
Maxwell Correction	Off
log physio files	On
FFT scale	2.00
dummy prepscan time	3 s
z shim	0.00 mT/m*ms
RF duration	2200 us
RF BWTP	25.0
Renzo: Delta TI	71 ms
EFFECTIVE TR	59419 ms
PatPartitions	24
EPI phase correction	local
PAT refscan mode	Flash
FlashRef BaseRes	162
FlashRef BW	100 Hz/px
FlashRef TE	7000 us
FlashRef FA	5 deg
use CAIPI	Off



# SIEMENS MAGNETOM Investigational\_Device\_7T syngo MR B17

\\USER\RenzoHuber\Lonike\_20201118\executed\_VASO20201118\VASO\_151\_0.8mm\_sag\_FA26\_no\_MAGEI

TA: 3:18 PAT: 3 Voxel size: 0.8x0.8x0.8 mm Rel. SNR: 1.00 USER: VASO\_151

## Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

## Routine

Slab group 1	
Slabs	1
Dist. factor	50 %
Position	R47.2 A27.2 H8.8
Orientation	S > C3.4
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
Slice oversampling	9.1 %
Slices per slab	22
FoV read	130.0 mm
FoV phase	133.3 %
Slice thickness	0.80 mm
TR	2475.80 ms
TE	24 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	A32

## Contrast

Perfusion mode	SS-SI VASO
TI2	700 ms
TI1	50 ms
TI1s	50 ms
Flip angle	26 deg
Fat suppr.	Fat sat.
Fat sat. mode	Strong
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	80
Delay in TR	0 ms
Multiple series	Off
Perfusion mode	PICORE Q2T
Inversion time 1	50 ms
Saturation stop time	50 ms
Inversion time 2	700.0 ms
Flow limit	100 cm/s

## Resolution

Base resolution	162
Phase resolution	100 %
Slice resolution	100 %
Phase partial Fourier	6/8
Slice partial Fourier	Off
Interpolation	Off

PAT mode	GRAPPA
Accel. factor PE	3
Ref. lines PE	45
Accel. factor 3D	1
Ref. lines 3D	24
Reference scan mode	Separate

Prescan Normalize	Off
Raw filter	Off
Elliptical filter	Off
Hamming	Off

## Geometry

Multi-slice mode	Interleaved
Series	Ascending
Special sat.	Parallel R
Gap	25.0 mm
Thickness	100 mm
Table position	H
Table position	0 mm
Inline Composing	Off

## System

V32	Off
A32	On
Positioning mode	REF
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Sum of Squares
AutoAlign	Head > Brain
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
! Ref. amplitude 1H	220.000 V
Adjustment Tolerance	Auto
Adjust volume	
! Position	R45.1 A26.5 H3.3
! Orientation	Sagittal
! Rotation	90.00 deg
! A >> P	174 mm
! F >> H	152 mm
! R >> L	39 mm

## Physio

1st Signal/Mode	None
-----------------	------

## BOLD

Motion correction	Off
Spatial filter	Off

## Sequence

Introduction	On
Dimension	3D
Reordering	Linear
Contrasts	1
Bandwidth	1144 Hz/Px
Free echo spacing	Off
Echo spacing	0.98 ms
EPI factor	216

# SIEMENS MAGNETOM Investigational\_Device\_7T syngo MR B17

RF pulse type	Normal
Gradient mode	Fast
Excitation	Slab-sel.
RF spoiling	On
<hr/>	
Ampl	150
MAGEC FA	12 in deg
ph.skip 4 Robert (the one)	1
MAGEC SS-SI?	Off
Maxwell Correction	Off
log physio files	On
FFT scale	2.00
dummy prepscan time	3 s
z shim	0.00 mT/m*ms
RF duration	2200 us
RF BWTP	25.0
Renzo: Delta TI	71 ms
EFFECTIVE TR	59419 ms
PatPartitions	24
EPI phase correction	local
PAT refscan mode	Flash
FlashRef BaseRes	162
FlashRef BW	100 Hz/px
FlashRef TE	7000 us
FlashRef FA	5 deg
use CAIPI	Off

# SIEMENS MAGNETOM Investigational\_Device\_7T syngo MR B17

\\USER\RenzoHuber\Lonike\_20201118\executed\_VASO20201118\VAPER\_FAalter\_FOV133\_chai\_sag

TA: 3:10 PAT: 3 Voxel size: 0.8x0.8x0.9 mm Rel. SNR: 1.00 USER: VAPER\_FAalter\_FOV133

## Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

## Routine

Slab group 1	
Slabs	1
Dist. factor	50 %
Position	R47.2 A27.2 H8.8
Orientation	S > C3.4
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
Slice oversampling	9.1 %
Slices per slab	22
FoV read	130.0 mm
FoV phase	133.3 %
Slice thickness	0.89 mm
TR	1795.3 ms
TE	24 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	A32

## Contrast

Perfusion mode	Picore Q2TIPS
TI2	1100 ms
TI1	50 ms
TI1s	50 ms
Flip angle	27.0 deg
Fat suppr.	None
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	106
Delay in TR	0 ms
Multiple series	Off
Perfusion mode	PICORE Q2T
Inversion time 1	50 ms
Saturation stop time	50 ms
Inversion time 2	1100.0 ms
Flow limit	100 cm/s

## Resolution

Base resolution	162
Phase resolution	100 %
Slice resolution	100 %
Phase partial Fourier	6/8
Slice partial Fourier	Off
Interpolation	Off
PAT mode	GRAPPA

Accel. factor PE	3
Ref. lines PE	45
Accel. factor 3D	1
Ref. lines 3D	24
Reference scan mode	Separate
Prescan Normalize	Off
Raw filter	Off
Elliptical filter	Off
Hamming	Off

## Geometry

Multi-slice mode	Interleaved
Series	Ascending
Special sat.	Parallel R
Gap	25.0 mm
Thickness	100 mm
Table position	H
Table position	0 mm
Inline Composing	Off

## System

V32	Off
A32	On
Positioning mode	FIX
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Sum of Squares
AutoAlign	Head > Brain
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
! Ref. amplitude 1H	250.000 V
Adjustment Tolerance	Auto
Adjust volume	
! Position	R45.1 A26.5 H3.3
! Orientation	Sagittal
! Rotation	90.00 deg
! A >> P	174 mm
! F >> H	152 mm
! R >> L	39 mm

## Physio

1st Signal/Mode	None
-----------------	------

## BOLD

Motion correction	Off
Spatial filter	Off

## Sequence

Introduction	On
Dimension	3D
Reordering	Linear
Contrasts	1
Bandwidth	1144 Hz/Px
Free echo spacing	Off
Echo spacing	0.98 ms
EPI factor	216
RF pulse type	Normal

Gradient mode	Fast
Excitation	Slab-sel.
RF spoiling	On
<hr/>	
Read Diff Amp	0.0 mT/m
Phase Diff Amp	0.0 mT/m
Slice Diff Amp	0.0 mT/m
Dante puls # in 1st par	38
Dante puls # in 2nd par	38
MT puls # each DANTE	0
Pulses FA in DANTE	10.5 degree
TAU in DANTE	200 us
diff TAU in MT	0 us
DANTE-RF dur	150 us
FA diff in DANTE	-3.0 degree
use Ernst angle	Off
Maxwell Correction	Off
log physio files	Off
FFT scale	2.00
dummy prepscan time	3 s
z shim	0.00 mT/m*ms
RF duration	2200 us
RF BWTP	25.0
EFFECTIVE TR	72 ms
PatPartitions	24
EPI phase correction	local
PAT refscan mode	Flash
FlashRef BaseRes	162
FlashRef BW	100 Hz/px
FlashRef TE	7000 us
FlashRef FA	5 deg
use CAIPI	Off

# SIEMENS MAGNETOM Investigational\_Device\_7T syngo MR B17

\\USER\RenzoHuber\Lonike\_20201118\executed\_VASO20201118\VASO\_151\_0.8mm\_sag\_FA26\_no\_MAGEI

TA: 8:00 PAT: 3 Voxel size: 0.8x0.8x0.8 mm Rel. SNR: 1.00 USER: VASO\_151

## Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

## Routine

Slab group 1	
Slabs	1
Dist. factor	50 %
Position	R47.2 A27.2 H8.8
Orientation	S > C3.4
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
Slice oversampling	9.1 %
Slices per slab	22
FoV read	130.0 mm
FoV phase	133.3 %
Slice thickness	0.80 mm
TR	2475.80 ms
TE	24 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	A32

## Contrast

Perfusion mode	SS-SI VASO
TI2	700 ms
TI1	50 ms
TI1s	50 ms
Flip angle	26 deg
Fat suppr.	Fat sat.
Fat sat. mode	Strong
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	194
Delay in TR	0 ms
Multiple series	Off
Perfusion mode	PICORE Q2T
Inversion time 1	50 ms
Saturation stop time	50 ms
Inversion time 2	700.0 ms
Flow limit	100 cm/s

## Resolution

Base resolution	162
Phase resolution	100 %
Slice resolution	100 %
Phase partial Fourier	6/8
Slice partial Fourier	Off
Interpolation	Off

PAT mode	GRAPPA
Accel. factor PE	3
Ref. lines PE	45
Accel. factor 3D	1
Ref. lines 3D	24
Reference scan mode	Separate

Prescan Normalize	Off
Raw filter	Off
Elliptical filter	Off
Hamming	Off

## Geometry

Multi-slice mode	Interleaved
Series	Ascending
Special sat.	Parallel R
Gap	25.0 mm
Thickness	100 mm
Table position	H
Table position	0 mm
Inline Composing	Off

## System

V32	Off
A32	On
Positioning mode	REF
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Sum of Squares
AutoAlign	Head > Brain
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
! Ref. amplitude 1H	220.000 V
Adjustment Tolerance	Auto
Adjust volume	
! Position	R45.1 A26.5 H3.3
! Orientation	Sagittal
! Rotation	90.00 deg
! A >> P	174 mm
! F >> H	152 mm
! R >> L	39 mm

## Physio

1st Signal/Mode	None
-----------------	------

## BOLD

Motion correction	Off
Spatial filter	Off

## Sequence

Introduction	On
Dimension	3D
Reordering	Linear
Contrasts	1
Bandwidth	1144 Hz/Px
Free echo spacing	Off
Echo spacing	0.98 ms
EPI factor	216

# SIEMENS MAGNETOM Investigational\_Device\_7T syngo MR B17

RF pulse type	Normal
Gradient mode	Fast
Excitation	Slab-sel.
RF spoiling	On
<hr/>	
Ampl	150
MAGEC FA	12 in deg
ph.skip 4 Robert (the one)	1
MAGEC SS-SI?	Off
Maxwell Correction	Off
log physio files	On
FFT scale	2.00
dummy prepscan time	3 s
z shim	0.00 mT/m*ms
RF duration	2200 us
RF BWTP	25.0
Renzo: Delta TI	71 ms
EFFECTIVE TR	59419 ms
PatPartitions	24
EPI phase correction	local
PAT refscan mode	Flash
FlashRef BaseRes	162
FlashRef BW	100 Hz/px
FlashRef TE	7000 us
FlashRef FA	5 deg
use CAIPI	Off

# SIEMENS MAGNETOM Investigational\_Device\_7T syngo MR B17

\\USER\RenzoHuber\Lonike\_20201118\executed\_VASO20201118\VASO\_151\_0.8mm\_axial\_FA4

TA: 4:22 PAT: 3 Voxel size: 0.8x0.8x0.8 mm Rel. SNR: 1.00 USER: VASO\_151

## Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

## Routine

Slab group 1	
Slabs	1
Dist. factor	50 %
Position	R2.3 A17.4 H9.0
Orientation	T > C15.8
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
Slice oversampling	9.1 %
Slices per slab	22
FoV read	130.0 mm
FoV phase	133.3 %
Slice thickness	0.80 mm
TR	2622.40 ms
TE	24 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	A32

## Contrast

Perfusion mode	SS-SI VASO
TI2	700 ms
TI1	50 ms
TI1s	50 ms
Flip angle	4 deg
Fat suppr.	Fat sat.
Fat sat. mode	Strong
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	100
Delay in TR	0 ms
Multiple series	Off
Perfusion mode	PICORE Q2T
Inversion time 1	50 ms
Saturation stop time	50 ms
Inversion time 2	700.0 ms
Flow limit	100 cm/s

## Resolution

Base resolution	162
Phase resolution	100 %
Slice resolution	100 %
Phase partial Fourier	6/8
Slice partial Fourier	Off
Interpolation	Off

PAT mode	GRAPPA
Accel. factor PE	3
Ref. lines PE	45
Accel. factor 3D	1
Ref. lines 3D	24
Reference scan mode	Separate

Prescan Normalize	Off
Raw filter	Off
Elliptical filter	Off
Hamming	Off

## Geometry

Multi-slice mode	Interleaved
Series	Ascending
Special sat.	Parallel F
Gap	25.0 mm
Thickness	100 mm
Table position	H
Table position	0 mm
Inline Composing	Off

## System

V32	Off
A32	On
Positioning mode	REF
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
! Ref. amplitude 1H	225.000 V
Adjustment Tolerance	Auto
Adjust volume	
! Position	L0.0 A22.7 H11.6
! Orientation	T > C14.2
! Rotation	0.00 deg
! R >> L	174 mm
! A >> P	180 mm
! F >> H	49 mm

## Physio

1st Signal/Mode	None
-----------------	------

## BOLD

Motion correction	Off
Spatial filter	Off

## Sequence

Introduction	On
Dimension	3D
Reordering	Linear
Contrasts	1
Bandwidth	1144 Hz/Px
Free echo spacing	Off
Echo spacing	0.98 ms
EPI factor	216

# SIEMENS MAGNETOM Investigational\_Device\_7T syngo MR B17

RF pulse type	Normal
Gradient mode	Fast
Excitation	Slab-sel.
RF spoiling	On
<hr/>	
Ampl	180
MAGEC FA	12 in deg
ph.skip 4 Robert (the one)	1
MAGEC SS-SI?	On
Maxwell Correction	Off
log physio files	On
FFT scale	2.00
dummy prepscan time	3 s
z shim	0.00 mT/m*ms
RF duration	2000 us
RF BWTP	25.0
Renzo: Delta TI	71 ms
EFFECTIVE TR	62937 ms
PatPartitions	24
EPI phase correction	local
PAT refscan mode	Flash
FlashRef BaseRes	162
FlashRef BW	100 Hz/px
FlashRef TE	7000 us
FlashRef FA	5 deg
use CAIPI	Off



# SIEMENS MAGNETOM Investigational\_Device\_7T syngo MR B17

\\USER\RenzoHuber\Lonike\_20201118\executed\_VASO20201118\VASO\_151\_0.8mm\_axial\_FA26

TA: 2:11 PAT: 3 Voxel size: 0.8x0.8x0.8 mm Rel. SNR: 1.00 USER: VASO\_151

## Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

## Routine

Slab group 1	
Slabs	1
Dist. factor	50 %
Position	R2.3 A17.4 H9.0
Orientation	T > C15.8
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
Slice oversampling	9.1 %
Slices per slab	22
FoV read	130.0 mm
FoV phase	133.3 %
Slice thickness	0.80 mm
TR	2622.40 ms
TE	24 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	A32

## Contrast

Perfusion mode	SS-SI VASO
TI2	700 ms
TI1	50 ms
TI1s	50 ms
Flip angle	26 deg
Fat suppr.	Fat sat.
Fat sat. mode	Strong
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	50
Delay in TR	0 ms
Multiple series	Off
Perfusion mode	PICORE Q2T
Inversion time 1	50 ms
Saturation stop time	50 ms
Inversion time 2	700.0 ms
Flow limit	100 cm/s

## Resolution

Base resolution	162
Phase resolution	100 %
Slice resolution	100 %
Phase partial Fourier	6/8
Slice partial Fourier	Off
Interpolation	Off

PAT mode	GRAPPA
Accel. factor PE	3
Ref. lines PE	45
Accel. factor 3D	1
Ref. lines 3D	24
Reference scan mode	Separate

Prescan Normalize	Off
Raw filter	Off
Elliptical filter	Off
Hamming	Off

## Geometry

Multi-slice mode	Interleaved
Series	Ascending
Special sat.	Parallel F
Gap	25.0 mm
Thickness	100 mm
Table position	H
Table position	0 mm
Inline Composing	Off

## System

V32	Off
A32	On
Positioning mode	REF
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
! Ref. amplitude 1H	225.000 V
Adjustment Tolerance	Auto
Adjust volume	
! Position	L0.0 A22.7 H11.6
! Orientation	T > C14.2
! Rotation	0.00 deg
! R >> L	174 mm
! A >> P	180 mm
! F >> H	49 mm

## Physio

1st Signal/Mode	None
-----------------	------

## BOLD

Motion correction	Off
Spatial filter	Off

## Sequence

Introduction	On
Dimension	3D
Reordering	Linear
Contrasts	1
Bandwidth	1144 Hz/Px
Free echo spacing	Off
Echo spacing	0.98 ms
EPI factor	216

# SIEMENS MAGNETOM Investigational\_Device\_7T syngo MR B17

RF pulse type	Normal
Gradient mode	Fast
Excitation	Slab-sel.
RF spoiling	On
<hr/>	
Ampl	180
MAGEC FA	12 in deg
ph.skip 4 Robert (the one)	1
MAGEC SS-SI?	On
Maxwell Correction	Off
log physio files	On
FFT scale	2.00
dummy prepsan time	3 s
z shim	0.00 mT/m*ms
RF duration	2000 us
RF BWTP	25.0
Renzo: Delta TI	71 ms
EFFECTIVE TR	62937 ms
PatPartitions	24
EPI phase correction	local
PAT refscan mode	Flash
FlashRef BaseRes	162
FlashRef BW	100 Hz/px
FlashRef TE	7000 us
FlashRef FA	5 deg
use CAIPI	Off