

\\USER\ \_ HighRes\_MP2RAGE\t1\_mp2rage\_sag\_p3\_0.5mm\_highres

TA: 16:14 PM: FIX Voxel size: 0.5×0.5×0.5 mmPAT: 2 Rel. SNR: 1.00 : tfl

**Properties**

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	On
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	Off
Start measurements	Single measurement

**Routine**

Slab group	1
Slabs	1
Dist. factor	50 %
Position	Isocenter
Orientation	T > S34.7 > C0.1
Phase enc. dir.	A >> P
AutoAlign	---
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	72
FoV read	161 mm
FoV phase	100.0 %
Slice thickness	0.50 mm
TR	6000.0 ms
TE	4.13 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	A32

**Contrast - Common**

TR	6000.0 ms
TE	4.13 ms
Magn. preparation	Non-sel. IR
T1 1	900 ms
T1 2	2900 ms
Flip angle 1	6.0 deg
Flip angle 2	7.0 deg
Fat suppr.	None
Water suppr.	None

**Contrast - Dynamic**

Averages	1
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	3
Pause after meas. 1	0.0 s
Pause after meas. 2	0.0 s
Multiple series	Each measurement

**Resolution - Common**

FoV read	161 mm
FoV phase	100.0 %
Slice thickness	0.50 mm
Base resolution	320
Phase resolution	100 %

**Resolution - Common**

Slice resolution	100 %
Phase partial Fourier	Off
Slice partial Fourier	6/8
Interpolation	Off

**Resolution - iPAT**

PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	24
Accel. factor 3D	1
Reference scan mode	Integrated

**Resolution - Filter Image**

Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off

**Resolution - Filter Rawdata**

Raw filter	Off
Elliptical filter	Off

**Geometry - Common**

Slab group	1
Slabs	1
Dist. factor	50 %
Position	Isocenter
Orientation	T > S34.7 > C0.1
Phase enc. dir.	A >> P
Slice oversampling	0.0 %
Slices per slab	72
FoV read	161 mm
FoV phase	100.0 %
Slice thickness	0.50 mm
TR	6000.0 ms
Multi-slice mode	Single shot
Series	Interleaved
Concatenations	1

**Geometry - AutoAlign**

Slab group	1
Position	Isocenter
Orientation	T > S34.7 > C0.1
Phase enc. dir.	A >> P
AutoAlign	---
Initial Position	Isocenter
L	0.0 mm
P	0.0 mm
H	0.0 mm
Initial Rotation	0.00 deg
Initial Orientation	T > S
T > S	34.7
> C	0.1

**Geometry - Navigator****System - Miscellaneous**

Positioning mode	FIX
Table position	H

**System - Miscellaneous**

Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Adaptive Combine
Save uncombined	Off
Matrix Optimization	Off
AutoAlign	---
Coil Select Mode	Default

**System - Adjustments**

B0 Shim mode	Standard
B1 Shim mode	TrueForm
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

**System - Adjust Volume**

Position	Isocenter
Orientation	T > S34.7 > C0.1
Rotation	0.00 deg
A >> P	161 mm
R >> L	161 mm
F >> H	36 mm
Reset	Off

**System - Tx/Rx**

Frequency 1H	297.203686 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

**Physio - Signal1**

1st Signal/Mode	None
TR	6000.0 ms
Concatenations	1

**Physio - Cardiac**

Magn. preparation	Non-sel. IR
TI 1	900 ms
TI 2	2900 ms
Fat suppr.	None
Dark blood	Off
FoV read	161 mm
FoV phase	100.0 %
Phase resolution	100 %

**Physio - PACE**

Resp. control	Off
Concatenations	1

**Inline - Common**

Subtract	Off
Measurements	3
StdDev	Off
Save original images	On

**Inline - MIP**

MIP-Sag	Off
---------	-----

**Inline - MIP**

MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On

**Inline - Composing**

Distortion Corr.	Off
------------------	-----

**Sequence - Part 1**

Introduction	On
Dimension	3D
Elliptical scanning	Off
Reordering	Linear rot.
Asymmetric echo	Off
Flow comp.	No
Multi-slice mode	Single shot
Echo spacing	8.5 ms
Bandwidth	190 Hz/Px

**Sequence - Part 2**

RF pulse type	Fast
Gradient mode	Fast
Excitation	Slab-sel.
RF spoiling	On
Incr. Gradient spoiling	Off
Turbo factor	319

**Sequence - Nuclei**

TX/RX Nucleus	1H
TX/RX delta frequency	0 Hz
TX Nucleus	None
TX delta frequency	0 Hz
Coil elements	A32

**Sequence - Assistant**

Mode	TR
Max. TR	4000.0 ms