

Table of contents

\\USER

<table border="1"><tr><td>Test</td></tr></table>	Test
Test	

<table border="1"><tr><td><table border="1"><tr><td>Joseph</td></tr></table></td></tr></table>	<table border="1"><tr><td>Joseph</td></tr></table>	Joseph
<table border="1"><tr><td>Joseph</td></tr></table>	Joseph	
Joseph		

<table border="1"><tr><td><table border="1"><tr><td>20210909_PALI_spike_test</td></tr></table></td></tr></table>	<table border="1"><tr><td>20210909_PALI_spike_test</td></tr></table>	20210909_PALI_spike_test
<table border="1"><tr><td>20210909_PALI_spike_test</td></tr></table>	20210909_PALI_spike_test	
20210909_PALI_spike_test		

<table border="1"><tr><td><table border="1"><tr><td>mp2rage_tfl_renzo</td></tr></table></td></tr></table>	<table border="1"><tr><td>mp2rage_tfl_renzo</td></tr></table>	mp2rage_tfl_renzo
<table border="1"><tr><td>mp2rage_tfl_renzo</td></tr></table>	mp2rage_tfl_renzo	
mp2rage_tfl_renzo		

\\USER\Test\Joseph\20210909_PALI_spike_test\mp2rage_tfl_rengo

TA: 27:47 PM: FIX Voxel size: 0.3×0.3×0.3 mmPAT: 2 Rel. SNR: 1.00 : rengo

Properties

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	Off
Start measurements	Single measurement

Routine

Slab group	1
Slabs	1
Dist. factor	50 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	Head > Basis
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	72
FoV read	180 mm
FoV phase	98.5 %
Slice thickness	0.35 mm
TR	5140.0 ms
TE	3.97 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	R96

Contrast - Common

TR	5140.0 ms
TE	3.97 ms
Magn. preparation	Non-sel. IR
TI 1	1290 ms
TI 2	3860 ms
Flip angle 1	9.0 deg
Flip angle 2	7.0 deg
Fat suppr.	None
Water suppr.	None

Contrast - Dynamic

Averages	1
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	6
Pause after meas. 1	0.0 s
Pause after meas. 2	0.0 s
Pause after meas. 3	0.0 s
Pause after meas. 4	0.0 s
Pause after meas. 5	0.0 s
Multiple series	Each measurement

Resolution - Common

FoV read	180 mm
FoV phase	98.5 %

Resolution - Common

Slice thickness	0.35 mm
Base resolution	520
Phase resolution	100 %
Slice resolution	100 %
Phase partial Fourier	Off
Slice partial Fourier	6/8
Interpolation	Off

Resolution - iPAT

PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	24
Accel. factor 3D	1
Reference scan mode	Integrated

Resolution - Filter Image

Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off

Resolution - Filter Rawdata

Raw filter	Off
Elliptical filter	Off

Geometry - Common

Slab group	1
Slabs	1
Dist. factor	50 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
Slice oversampling	0.0 %
Slices per slab	72
FoV read	180 mm
FoV phase	98.5 %
Slice thickness	0.35 mm
TR	5140.0 ms
Multi-slice mode	Single shot
Series	Interleaved
Concatenations	1

Geometry - AutoAlign

Slab group	1
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	Head > Basis
Initial Position	Isocenter
L	0.0 mm
P	0.0 mm
H	0.0 mm
Initial Rotation	0.00 deg
Initial Orientation	Transversal

Geometry - Navigator**Geometry - Tim Planning Suite**

Set-n-Go Protocol	Off
-------------------	-----

Geometry - Tim Planning Suite

Table position	H
Table position	0 mm
Inline Composing	Off

System - Miscellaneous

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
Save uncombined	Off
Matrix Optimization	Off
AutoAlign	Head > Basis
Coil Select Mode	Default

System - Adjustments

B0 Shim mode	Standard
B1 Shim mode	TrueForm
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

System - Adjust Volume

Position	Isocenter
Orientation	Transversal
Rotation	0.00 deg
A >> P	178 mm
R >> L	180 mm
F >> H	26 mm
Reset	Off

System - pTx Volumes

B1 Shim mode	TrueForm
Excitation	Slab-sel.

System - Tx/Rx

Frequency 1H	297.210739 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

Physio - Signal1

1st Signal/Mode	None
TR	5140.0 ms
Concatenations	1

Physio - Cardiac

Magn. preparation	Non-sel. IR
TI 1	1290 ms
TI 2	3860 ms
Fat suppr.	None
Dark blood	Off
FoV read	180 mm
FoV phase	98.5 %
Phase resolution	100 %

Physio - PACE

Resp. control	Off
Concatenations	1

Inline - Common

Subtract	Off
Measurements	6
StdDev	Off
Save original images	On

Inline - MIP

MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On

Inline - Composing

Inline Composing	Off
Distortion Corr.	Off

Inline - MapIt

Save original images	On
MapIt	None
Flip angle 1	9.0 deg
Flip angle 2	7.0 deg
Measurements	6
TR	5140.0 ms
TE	3.97 ms

Sequence - Part 1

Introduction	On
Dimension	3D
Elliptical scanning	Off
Reordering	Linear rot.
Asymmetric echo	Allowed
Flow comp.	No
Multi-slice mode	Single shot
Echo spacing	9.5 ms
Bandwidth	250 Hz/Px

Sequence - Part 2

RF pulse type	Normal
Gradient mode	Fast
Excitation	Slab-sel.
RF spoiling	On
Incr. Gradient spoiling	Off
Turbo factor	511

Sequence - Assistant

Mode	Off
------	-----