

**Table of contents**

\USER

FMRIF

[XT-ID:93-M-0170]]Renzo

20230327\_ptx\_used

localizer\_irtfl\_ptx  
rslh\_ep3d\_vaso\_nih5k\_saggital\_normal  
rslh\_ep3d\_vaso\_nih5k\_saggital\_normal  
rslh\_ep3d\_vaso\_nih5k\_saggital\_normal  
rslh\_ep3d\_vaso\_nih5k\_saggital\_normal\_gfactor  
cstfl\_wip925b\_protocolfromTobi

\\USER\FMRIF\XT-ID:93-M-0170\Renzo\20230327\_ptx\_used\localizer\_irtfl\_ptx

TA: 1:00 PM: REF Voxel size: 1.0×1.0×2.0 mmPAT: Off Rel. SNR: 1.00 : tfl

**Properties**

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	On
Auto close inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	Single measurement

**Contrast - Dynamic**

Measurements	1
Multiple series	Each measurement

**Resolution - Common**

FoV read	200 mm
FoV phase	100.0 %
Slice thickness	2.0 mm
Base resolution	192
Phase resolution	100 %
Phase partial Fourier	Off
Interpolation	Off

**Resolution - iPAT**

PAT mode	None
----------	------

**Resolution - Filter Image**

Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off

**Resolution - Filter Rawdata**

Raw filter	Off
Elliptical filter	Off

**Routine**

Slice group	1
Slices	6
Dist. factor	750 %
Position	L1.9 A5.0 F8.0 mm
Orientation	Sagittal
Phase enc. dir.	A >> P
Slice group	2
Slices	4
Dist. factor	600 %
Position	L1.5 A1.1 H6.0 mm
Orientation	Transversal
Phase enc. dir.	A >> P
Slice group	3
Slices	4
Dist. factor	800 %
Position	L1.5 P11.1 F7.2 mm
Orientation	Coronal
Phase enc. dir.	R >> L
AutoAlign	---
Phase oversampling	0 %
FoV read	200 mm
FoV phase	100.0 %
Slice thickness	2.0 mm
TR	4200.0 ms
TE	3.46 ms
Averages	1
Concatenations	14
Filter	None
Coil elements	AC

**Geometry - Common**

Slice group	1
Slices	6
Dist. factor	750 %
Position	L1.9 A5.0 F8.0 mm
Orientation	Sagittal
Phase enc. dir.	A >> P
Slice group	2
Slices	4
Dist. factor	600 %
Position	L1.5 A1.1 H6.0 mm
Orientation	Transversal
Phase enc. dir.	A >> P
Slice group	3
Slices	4
Dist. factor	800 %
Position	L1.5 P11.1 F7.2 mm
Orientation	Coronal
Phase enc. dir.	R >> L
FoV read	200 mm
FoV phase	100.0 %
Slice thickness	2.0 mm
TR	4200.0 ms
Multi-slice mode	Sequential
Series	Ascending
Concatenations	14

**Contrast - Common**

TR	4200.0 ms
TE	3.46 ms
TD	0 ms
Magn. preparation	Slice-sel. IR
T1 1	840 ms
T1 2	2540 ms
Flip angle 1	5.0 deg
Flip angle 2	8.0 deg
Fat suppr.	None
Water suppr.	None

**Geometry - AutoAlign**

Slice group	1
Position	L1.9 A5.0 F8.0 mm
Orientation	Sagittal
Phase enc. dir.	A >> P
Slice group	2
Position	L1.5 A1.1 H6.0 mm

**Contrast - Dynamic**

Averages	1
Averaging mode	Long term
Reconstruction	Magnitude

**Geometry - AutoAlign**

Orientation	Transversal
Phase enc. dir.	A >> P
Slice group	3
Position	L1.5 P11.1 F7.2 mm
Orientation	Coronal
Phase enc. dir.	R >> L
AutoAlign	---
Initial Position	L1.9 A5.0 F8.0
L	1.9 mm
A	5.0 mm
F	8.0 mm
Initial Rotation	0.00 deg
Initial Orientation	Sagittal

**Geometry - Navigator****Geometry - Tim Planning Suite**

Set-n-Go Protocol	Off
Table position	H
Table position	0 mm
Inline Composing	Off

**System - Miscellaneous**

Positioning mode	REF
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Adaptive Combine
Save uncombined	Off
Matrix Optimization	Off
AutoAlign	---
Coil Select Mode	Default

**System - Adjustments**

B0 Shim mode	Tune up
B1 Shim mode	TrueForm
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

**System - Adjust Volume**

Position	Isocenter
Orientation	Transversal
Rotation	0.00 deg
A >> P	263 mm
R >> L	350 mm
F >> H	350 mm
Reset	Off

**System - pTx Volumes**

B1 Shim mode	TrueForm
Excitation	Slice-sel.

**System - Tx/Rx**

Frequency 1H	297.145077 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off

**System - Tx/Rx**

? Ref. amplitude 1H	0.000 V
---------------------	---------

**Physio - Signal1**

1st Signal/Mode	None
TR	4200.0 ms
Concatenations	14

**Physio - Cardiac**

Magn. preparation	Slice-sel. IR
TI 1	840 ms
TI 2	2540 ms
Fat suppr.	None
Dark blood	Off
FoV read	200 mm
FoV phase	100.0 %
Phase resolution	100 %

**Physio - PACE**

Resp. control	Off
Concatenations	14

**Inline - Common**

Subtract	Off
Measurements	1
StdDev	Off
Save original images	On

**Inline - MIP**

MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On

**Inline - Composing**

Inline Composing	Off
Distortion Corr.	Off

**Inline - MapIt**

Save original images	On
MapIt	None
Flip angle 1	5.0 deg
Flip angle 2	8.0 deg
Measurements	1
TR	4200.0 ms
TE	3.46 ms

**Sequence - Part 1**

Introduction	On
Dimension	2D
Asymmetric echo	Off
Flow comp.	No
Multi-slice mode	Sequential
Echo spacing	6.7 ms
Bandwidth	240 Hz/Px

**Sequence - Part 2**

RF pulse type	Normal
Gradient mode	Fast
Excitation	Slice-sel.
RF spoiling	On
Incr. Gradient spoiling	Off
Turbo factor	192

**Sequence - Assistant**

Mode	Off
------	-----

\\USER\FMRIF\XT-ID:93-M-0170\Renzo\20230327_ptx_used\rslh_ep3d_vaso_nih5k_saggital_norma I
TA: 1:12 PM: REF Voxel size: 0.8×0.8×0.8 mmPAT: 3 Rel. SNR: 1.00 : 5k

**Properties**

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	Off
Start measurements	Single measurement

**Routine**

Slab group	1
Slabs	1
Position	L26.0 A0.1 F16.4 mm
Orientation	Sagittal
Phase enc. dir.	P >> A
AutoAlign	---
Slab Scale	-10 %
Slices per slab	36
FoV read	186 mm
FoV phase	94.5 %
Slice thickness	0.84 mm
TR 1	68.9 ms
TR 2	6136 ms
TE 1	22.90 ms
Averages	1
Filter	None
Coil elements	AC

**Contrast - Common**

TR 1	68.9 ms
TR 2	6136 ms
TE 1	22.90 ms
Multi-echo spacing	61.2 ms
Magn. preparation	Non-sel. HSN IR
TI 1	1750.2 ms
TI 2	4230.6 ms
Flip angle	57 deg
Fat suppr.	Fat sat.
Magn. Prep. Shots	1

**Contrast - Dynamic**

Averages	1
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	10
Pause after meas. 1	0.0 s
Pause after meas. 2	0.0 s
Pause after meas. 3	0.0 s
Pause after meas. 4	0.0 s
Pause after meas. 5	0.0 s
Pause after meas. 6	0.0 s
Pause after meas. 7	0.0 s
Pause after meas. 8	0.0 s
Pause after meas. 9	0.0 s

**Resolution - Common**

FoV read	186 mm
FoV phase	94.5 %
Slice thickness	0.84 mm
Base resolution	220
Phase resolution	100 %
Slice resolution	100 %
Phase partial Fourier	6/8
Slice partial Fourier	Off
Interpolation	Off

**Resolution - iPAT**

PAT mode	CAIPIRINHA
Acc. factor PE	1
Ref. lines PE	48
Acc. factor 3D	3
Ref. lines 3D	36
CAIPI 3D Shift	1
Reference Scan Mode	EPI/separate
CAIPI Mode (tooltip)	Skipped-CAIPI
Total PAT factor	3

**Resolution - Filter Image**

Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off

**Resolution - Filter Rawdata**

Raw filter	Off
Elliptical filter	Off

**Geometry - Common**

Slab group	1
Slabs	1
Position	L26.0 A0.1 F16.4 mm
Orientation	Sagittal
Phase enc. dir.	P >> A
Slab Scale	-10 %
Slices per slab	36
FoV read	186 mm
FoV phase	94.5 %
Slice thickness	0.84 mm
TR 1	68.9 ms
TR 2	6136 ms

**Geometry - AutoAlign**

Slab group	1
Position	L26.0 A0.1 F16.4 mm
Orientation	Sagittal
Phase enc. dir.	P >> A
AutoAlign	---
Initial Position	L26.0 A0.1 F16.4
L	26.0 mm
A	0.1 mm
F	16.4 mm
Initial Rotation	-180.00 deg
Initial Orientation	Sagittal

**Geometry - Saturation**

Saturation mode	Standard
Fat suppr.	Fat sat.

**Geometry - Tim Planning Suite**

Set-n-Go Protocol	Off
Table position	H
Table position	0 mm
Inline Composing	Off

**System - Miscellaneous**

Positioning mode	REF
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
Save uncombined	Off
Matrix Optimization	Off
AutoAlign	---
Coil Select Mode	Default

**System - Adjustments**

B0 Shim mode	Brain
B1 Shim mode	TrueForm
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

**System - Adjust Volume**

! Position	L29.1 A2.4 H7.6 mm
! Orientation	Sagittal
! Rotation	0.00 deg
! A >> P	193 mm
! F >> H	124 mm
! R >> L	45 mm
Reset	Off

**System - pTx Volumes**

B1 Shim mode	TrueForm
Excitation	Slab-sel.

**System - Tx/Rx**

Frequency 1H	297.145077 MHz
Correction factor	1
Gain	Low
Img. Scale Cor.	1.000
Reset	Off
! Ref. amplitude 1H	220.000 V

**Sequence - Part 1**

Introduction	On
Dimension	3D
Reordering	Linear
Contrasts	1
Echo spacing	1.16 ms
Bandwidth	946 Hz/Px

**Sequence - Part 2**

EPI factor	52
Segmentation	3

**Sequence - Part 2**

RF pulse type	Normal
Gradient mode	Fast
Excitation	Slab-sel.
RF spoiling	On
Turbo factor	36

**Sequence - Special**

PATRef FA	3 deg
RF duration	2500 us
RF BWT product	8
Ernst T1	1200 ms
PATRef prep. shots	10
Volume dummy shots	0
Dummy Measurements	0
ETL per RTEB	1
CHECK FLIP ANGLE!	On
Invert PE	On
Min. TE if PF	On
Echo Time Shift	On
Ramp Sampling	On
NORDIC	Off
SVDPC	Off
Sym VASO	Off
Dual-pol. EPI	Off
Invert RO	On
Invert 3D	Off
Disable PF reco	Off
Disable PF reco	Off
Save sampling	Off
PE VComp	Off
Water Exc.	-none-
External PC	per Series
Saturation RF	per Shot
FIDNavs	-none-
EPI rise time factor	1.10
Mosaic DICOMs	On
Modify Ice Config	On
G-factor map	Off
GRAPPA Regularization	10 10 <sup>-6</sup>
HSN RF power scale	2.50
Inversion Delay	500 ms
Relaxation Delay	0 ms
Var. FA /MAGEC	4

**Sequence - Assistant**

Mode	Off
------	-----

\\USER\FMRIF\[XT-ID:93-M-0170]\Renzo\20230327_ptx_used\rslh_ep3d_vaso_nih5k_saggital_norma I
TA: 10:24 PM: REF Voxel size: 0.8×0.8×0.8 mmPAT: 3 Rel. SNR: 1.00 : 5k

**Properties**

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	Off
Start measurements	Single measurement

**Routine**

Slab group	1
Slabs	1
Position	L26.0 A0.1 F16.4 mm
Orientation	Sagittal
Phase enc. dir.	P >> A
AutoAlign	---
Slab Scale	-10 %
Slices per slab	36
FoV read	186 mm
FoV phase	94.5 %
Slice thickness	0.84 mm
TR 1	68.9 ms
TR 2	6136 ms
TE 1	22.90 ms
Averages	1
Filter	None
Coil elements	AC

**Contrast - Common**

TR 1	68.9 ms
TR 2	6136 ms
TE 1	22.90 ms
Multi-echo spacing	61.2 ms
Magn. preparation	Non-sel. HSN IR
TI 1	1750.2 ms
TI 2	4230.6 ms
Flip angle	57 deg
Fat suppr.	Fat sat.
Magn. Prep. Shots	1

**Contrast - Dynamic**

Averages	1
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	100
Pause after meas.	0.0 s

**Resolution - Common**

FoV read	186 mm
FoV phase	94.5 %
Slice thickness	0.84 mm
Base resolution	220
Phase resolution	100 %
Slice resolution	100 %
Phase partial Fourier	6/8

**Resolution - Common**

Slice partial Fourier	Off
Interpolation	Off

**Resolution - iPAT**

PAT mode	CAIPIRINHA
Acc. factor PE	1
Ref. lines PE	48
Acc. factor 3D	3
Ref. lines 3D	36
CAIPI 3D Shift	1
Reference Scan Mode	EPI/separate
CAIPI Mode (tooltip)	Skipped-CAIPI
Total PAT factor	3

**Resolution - Filter Image**

Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off

**Resolution - Filter Rawdata**

Raw filter	Off
Elliptical filter	Off

**Geometry - Common**

Slab group	1
Slabs	1
Position	L26.0 A0.1 F16.4 mm
Orientation	Sagittal
Phase enc. dir.	P >> A
Slab Scale	-10 %
Slices per slab	36
FoV read	186 mm
FoV phase	94.5 %
Slice thickness	0.84 mm
TR 1	68.9 ms
TR 2	6136 ms

**Geometry - AutoAlign**

Slab group	1
Position	L26.0 A0.1 F16.4 mm
Orientation	Sagittal
Phase enc. dir.	P >> A
AutoAlign	---
Initial Position	L26.0 A0.1 F16.4
L	26.0 mm
A	0.1 mm
F	16.4 mm
Initial Rotation	-180.00 deg
Initial Orientation	Sagittal

**Geometry - Saturation**

Saturation mode	Standard
Fat suppr.	Fat sat.

**Geometry - Tim Planning Suite**

Set-n-Go Protocol	Off
-------------------	-----

**Geometry - Tim Planning Suite**

Table position	H
Table position	0 mm
Inline Composing	Off

**System - Miscellaneous**

Positioning mode	REF
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
Save uncombined	Off
Matrix Optimization	Off
AutoAlign	---
Coil Select Mode	Default

**System - Adjustments**

B0 Shim mode	Brain
B1 Shim mode	TrueForm
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

**System - Adjust Volume**

! Position	L29.1 A2.4 H7.6 mm
! Orientation	Sagittal
! Rotation	0.00 deg
! A >> P	193 mm
! F >> H	124 mm
! R >> L	45 mm
Reset	Off

**System - pTx Volumes**

B1 Shim mode	TrueForm
Excitation	Slab-sel.

**System - Tx/Rx**

Frequency 1H	297.145077 MHz
Correction factor	1
Gain	Low
Img. Scale Cor.	1.000
Reset	Off
! Ref. amplitude 1H	220.000 V

**Sequence - Part 1**

Introduction	On
Dimension	3D
Reordering	Linear
Contrasts	1
Echo spacing	1.16 ms
Bandwidth	946 Hz/Px

**Sequence - Part 2**

EPI factor	52
Segmentation	3
RF pulse type	Normal
Gradient mode	Fast
Excitation	Slab-sel.
RF spoiling	On
Turbo factor	36

**Sequence - Special**

PATRef FA	3 deg
RF duration	2500 us
RF BWT product	8
Ernst T1	1200 ms
PATRef prep. shots	10
Volume dummy shots	0
Dummy Measurements	0
ETL per RTEB	1
CHECK FLIP ANGLE!	On
Invert PE	On
Min. TE if PF	On
Echo Time Shift	On
Ramp Sampling	On
NORDIC	Off
SVDPC	Off
Sym VASO	Off
Dual-pol. EPI	Off
Invert RO	On
Invert 3D	Off
Disable PF reco	Off
Disable PF reco	Off
Save sampling	Off
PE VComp	Off
Water Exc.	-none-
External PC	per Series
Saturation RF	per Shot
FIDNavs	-none-
EPI rise time factor	1.10
Mosaic DICOMs	On
Modify Ice Config	On
G-factor map	Off
GRAPPA Regularization	10 10 <sup>-6</sup>
HSN RF power scale	2.50
Inversion Delay	500 ms
Relaxation Delay	0 ms
Var. FA /MAGEC	4

**Sequence - Assistant**

Mode	Off
------	-----



\\USER\FMRIF\[XT-ID:93-M-0170]\Renzo\20230327_ptx_used\rslh_ep3d_vaso_nih5k_saggital_norma I
TA: 10:24 PM: REF Voxel size: 0.8×0.8×0.8 mmPAT: 3 Rel. SNR: 1.00 : 5k

**Properties**

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	Off
Start measurements	Single measurement

**Routine**

Slab group	1
Slabs	1
Position	L26.0 A0.1 F16.4 mm
Orientation	Sagittal
Phase enc. dir.	P >> A
AutoAlign	---
Slab Scale	-10 %
Slices per slab	36
FoV read	186 mm
FoV phase	94.5 %
Slice thickness	0.84 mm
TR 1	68.9 ms
TR 2	6136 ms
TE 1	22.90 ms
Averages	1
Filter	None
Coil elements	AC

**Contrast - Common**

TR 1	68.9 ms
TR 2	6136 ms
TE 1	22.90 ms
Multi-echo spacing	61.2 ms
Magn. preparation	Non-sel. HSN IR
TI 1	1750.2 ms
TI 2	4230.6 ms
Flip angle	57 deg
Fat suppr.	Fat sat.
Magn. Prep. Shots	1

**Contrast - Dynamic**

Averages	1
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	100
Pause after meas.	0.0 s

**Resolution - Common**

FoV read	186 mm
FoV phase	94.5 %
Slice thickness	0.84 mm
Base resolution	220
Phase resolution	100 %
Slice resolution	100 %
Phase partial Fourier	6/8

**Resolution - Common**

Slice partial Fourier	Off
Interpolation	Off

**Resolution - iPAT**

PAT mode	CAIPIRINHA
Acc. factor PE	1
Ref. lines PE	48
Acc. factor 3D	3
Ref. lines 3D	36
CAIPI 3D Shift	1
Reference Scan Mode	EPI/separate
CAIPI Mode (tooltip)	Skipped-CAIPI
Total PAT factor	3

**Resolution - Filter Image**

Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off

**Resolution - Filter Rawdata**

Raw filter	Off
Elliptical filter	Off

**Geometry - Common**

Slab group	1
Slabs	1
Position	L26.0 A0.1 F16.4 mm
Orientation	Sagittal
Phase enc. dir.	P >> A
Slab Scale	-10 %
Slices per slab	36
FoV read	186 mm
FoV phase	94.5 %
Slice thickness	0.84 mm
TR 1	68.9 ms
TR 2	6136 ms

**Geometry - AutoAlign**

Slab group	1
Position	L26.0 A0.1 F16.4 mm
Orientation	Sagittal
Phase enc. dir.	P >> A
AutoAlign	---
Initial Position	L26.0 A0.1 F16.4
L	26.0 mm
A	0.1 mm
F	16.4 mm
Initial Rotation	-180.00 deg
Initial Orientation	Sagittal

**Geometry - Saturation**

Saturation mode	Standard
Fat suppr.	Fat sat.

**Geometry - Tim Planning Suite**

Set-n-Go Protocol	Off
-------------------	-----

**Geometry - Tim Planning Suite**

Table position	H
Table position	0 mm
Inline Composing	Off

**System - Miscellaneous**

Positioning mode	REF
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
Save uncombined	Off
Matrix Optimization	Off
AutoAlign	---
Coil Select Mode	Default

**System - Adjustments**

B0 Shim mode	Brain
B1 Shim mode	TrueForm
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

**System - Adjust Volume**

! Position	L29.1 A2.4 H7.6 mm
! Orientation	Sagittal
! Rotation	0.00 deg
! A >> P	193 mm
! F >> H	124 mm
! R >> L	45 mm
Reset	Off

**System - pTx Volumes**

B1 Shim mode	TrueForm
Excitation	Slab-sel.

**System - Tx/Rx**

Frequency 1H	297.145077 MHz
Correction factor	1
Gain	Low
Img. Scale Cor.	1.000
Reset	Off
! Ref. amplitude 1H	220.000 V

**Sequence - Part 1**

Introduction	On
Dimension	3D
Reordering	Linear
Contrasts	1
Echo spacing	1.16 ms
Bandwidth	946 Hz/Px

**Sequence - Part 2**

EPI factor	52
Segmentation	3
RF pulse type	Normal
Gradient mode	Fast
Excitation	Slab-sel.
RF spoiling	On
Turbo factor	36

**Sequence - Special**

PATRef FA	3 deg
RF duration	2500 us
RF BWT product	8
Ernst T1	1200 ms
PATRef prep. shots	10
Volume dummy shots	0
Dummy Measurements	0
ETL per RTEB	1
CHECK FLIP ANGLE!	On
Invert PE	On
Min. TE if PF	On
Echo Time Shift	On
Ramp Sampling	On
NORDIC	Off
SVDPC	Off
Sym VASO	Off
Dual-pol. EPI	Off
Invert RO	On
Invert 3D	Off
Disable PF reco	Off
Disable PF reco	Off
Save sampling	Off
PE VComp	Off
Water Exc.	-none-
External PC	per Series
Saturation RF	per Shot
FIDNavs	-none-
EPI rise time factor	1.10
Mosaic DICOMs	On
Modify Ice Config	On
G-factor map	Off
GRAPPA Regularization	10 10 <sup>-6</sup>
HSN RF power scale	2.50
Inversion Delay	500 ms
Relaxation Delay	0 ms
Var. FA /MAGEC	4

**Sequence - Assistant**

Mode	Off
------	-----

\\USER\FMRIF\XT-ID:93-M-0170\Renzo\20230327_ptx_used\rslh_ep3d_vaso_nih5k_saggital_norma l_gfactor
TA: 0:16 PM: REF Voxel size: 0.8×0.8×0.8 mmPAT: 3 Rel. SNR: 1.00 : 5k

**Properties**

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	Off
Start measurements	Single measurement

**Routine**

Slab group	1
Slabs	1
Position	L26.0 A0.1 F16.4 mm
Orientation	Sagittal
Phase enc. dir.	P >> A
AutoAlign	---
Slab Scale	-10 %
Slices per slab	36
FoV read	186 mm
FoV phase	94.5 %
Slice thickness	0.84 mm
TR 1	68.9 ms
TR 2	6136 ms
TE 1	22.90 ms
Averages	1
Filter	None
Coil elements	AC

**Contrast - Common**

TR 1	68.9 ms
TR 2	6136 ms
TE 1	22.90 ms
Multi-echo spacing	61.2 ms
Magn. preparation	Non-sel. HSN IR
TI 1	1750.2 ms
TI 2	4230.6 ms
Flip angle	57 deg
Fat suppr.	Fat sat.
Magn. Prep. Shots	1

**Contrast - Dynamic**

Averages	1
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	1

**Resolution - Common**

FoV read	186 mm
FoV phase	94.5 %
Slice thickness	0.84 mm
Base resolution	220
Phase resolution	100 %
Slice resolution	100 %
Phase partial Fourier	6/8
Slice partial Fourier	Off

**Resolution - Common**

Interpolation	Off
---------------	-----

**Resolution - iPAT**

PAT mode	CAIPIRINHA
Acc. factor PE	1
Ref. lines PE	48
Acc. factor 3D	3
Ref. lines 3D	36
CAIPI 3D Shift	1
Reference Scan Mode	EPI/separate
CAIPI Mode (tooltip)	Skipped-CAIPI
Total PAT factor	3

**Resolution - Filter Image**

Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off

**Resolution - Filter Rawdata**

Raw filter	Off
Elliptical filter	Off

**Geometry - Common**

Slab group	1
Slabs	1
Position	L26.0 A0.1 F16.4 mm
Orientation	Sagittal
Phase enc. dir.	P >> A
Slab Scale	-10 %
Slices per slab	36
FoV read	186 mm
FoV phase	94.5 %
Slice thickness	0.84 mm
TR 1	68.9 ms
TR 2	6136 ms

**Geometry - AutoAlign**

Slab group	1
Position	L26.0 A0.1 F16.4 mm
Orientation	Sagittal
Phase enc. dir.	P >> A
AutoAlign	---
Initial Position	L26.0 A0.1 F16.4
L	26.0 mm
A	0.1 mm
F	16.4 mm
Initial Rotation	-180.00 deg
Initial Orientation	Sagittal

**Geometry - Saturation**

Saturation mode	Standard
Fat suppr.	Fat sat.

**Geometry - Tim Planning Suite**

Set-n-Go Protocol	Off
Table position	H

**Geometry - Tim Planning Suite**

Table position	0 mm
Inline Composing	Off

**System - Miscellaneous**

Positioning mode	REF
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
Save uncombined	Off
Matrix Optimization	Off
AutoAlign	---
Coil Select Mode	Default

**System - Adjustments**

B0 Shim mode	Brain
B1 Shim mode	TrueForm
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

**System - Adjust Volume**

! Position	L29.1 A2.4 H7.6 mm
! Orientation	Sagittal
! Rotation	0.00 deg
! A >> P	193 mm
! F >> H	124 mm
! R >> L	45 mm
Reset	Off

**System - pTx Volumes**

B1 Shim mode	TrueForm
Excitation	Slab-sel.

**System - Tx/Rx**

Frequency 1H	297.145077 MHz
Correction factor	1
Gain	Low
Img. Scale Cor.	1.000
Reset	Off
! Ref. amplitude 1H	220.000 V

**Sequence - Part 1**

Introduction	On
Dimension	3D
Reordering	Linear
Contrasts	1
Echo spacing	1.16 ms
Bandwidth	946 Hz/Px

**Sequence - Part 2**

EPI factor	52
Segmentation	3
RF pulse type	Normal
Gradient mode	Fast
Excitation	Slab-sel.
RF spoiling	On
Turbo factor	36

**Sequence - Special**

PATRef FA	3 deg
RF duration	2500 us
RF BWT product	8
Ernst T1	1200 ms
PATRef prep. shots	10
Volume dummy shots	0
Dummy Measurements	0
ETL per RTEB	1
CHECK FLIP ANGLE!	On
Invert PE	On
Min. TE if PF	On
Echo Time Shift	On
Ramp Sampling	On
NORDIC	Off
SVDPC	Off
Sym VASO	Off
Dual-pol. EPI	Off
Invert RO	On
Invert 3D	Off
Disable PF reco	Off
Disable PF reco	Off
Save sampling	Off
PE VComp	Off
Water Exc.	-none-
External PC	per Series
Saturation RF	per Shot
FIDNavs	-none-
EPI rise time factor	1.10
Mosaic DICOMs	On
Modify Ice Config	On
G-factor map	On
GRAPPA Regularization	10 10 <sup>-6</sup>
HSN RF power scale	2.50
Inversion Delay	500 ms
Relaxation Delay	0 ms
Var. FA /MAGEC	4

**Sequence - Assistant**

Mode	Off
------	-----

\\USER\FMRIF\[XT-ID:93-M-0170]\Renzo\20230327\_ptx\_used\cstfl\_wip925b\_protocolfromTobi

TA: 5:58 PM: REF Voxel size: 0.8×0.8×0.8 mmPAT: 1.0 Rel. SNR: 1.00 : WIP\_cmp

**Properties**

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	On
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	Off
Start measurements	Single measurement

**Routine**

Slab group	1
Slabs	1
Dist. factor	50 %
Position	L0.0 A11.3 F27.2 mm
Orientation	Sagittal
Phase enc. dir.	A >> P
AutoAlign	---
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	224
FoV read	230 mm
FoV phase	94.4 %
Slice thickness	0.80 mm
TR	6000.0 ms
TE	1.97 ms
Averages	1
Concatenations	1
Filter	Distortion Corr.(3D)
Coil elements	AC

**Contrast - Common**

TR	6000.0 ms
TE	1.97 ms
Magn. preparation	Non-sel. IR
T1 1	800 ms
T1 2	2700 ms
Flip angle 1	4 deg
Flip angle 2	5 deg
Fat suppr.	None
Water suppr.	None

**Contrast - Dynamic**

Averages	1
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	1
Multiple series	Each measurement

**Resolution - Common**

FoV read	230 mm
FoV phase	94.4 %
Slice thickness	0.80 mm
Base resolution	288
Phase resolution	100 %
Slice resolution	100 %
Phase partial Fourier	Off

**Resolution - Common**

Slice partial Fourier	Off
Interpolation	Off

**Resolution - iPAT**

PAT mode	None
----------	------

**Resolution - Filter Image**

Image Filter	Off
Distortion Corr.	On
Mode	3D
Unfiltered images	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off

**Resolution - Filter Rawdata**

Raw filter	Off
Elliptical filter	Off

**Geometry - Common**

Slab group	1
Slabs	1
Dist. factor	50 %
Position	L0.0 A11.3 F27.2 mm
Orientation	Sagittal
Phase enc. dir.	A >> P
Slice oversampling	0.0 %
Slices per slab	224
FoV read	230 mm
FoV phase	94.4 %
Slice thickness	0.80 mm
TR	6000.0 ms
Multi-slice mode	Single shot
Series	Interleaved
Concatenations	1

**Geometry - AutoAlign**

Slab group	1
Position	L0.0 A11.3 F27.2 mm
Orientation	Sagittal
Phase enc. dir.	A >> P
AutoAlign	---
Initial Position	L0.0 A11.3 F27.2
L	0.0 mm
A	11.3 mm
F	27.2 mm
Initial Rotation	0.00 deg
Initial Orientation	Sagittal

**Geometry - Navigator****Geometry - Tim Planning Suite**

Set-n-Go Protocol	Off
Table position	H
Table position	0 mm
Inline Composing	Off

**System - Miscellaneous**

Positioning mode	REF
------------------	-----

**System - Miscellaneous**

Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
Save uncombined	Off
Matrix Optimization	Off
AutoAlign	---
Coil Select Mode	Default

**System - Adjustments**

B0 Shim mode	Brain
B1 Shim mode	TrueForm
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

**System - Adjust Volume**

! Position	L29.1 A2.4 H7.6 mm
! Orientation	Sagittal
! Rotation	0.00 deg
! A >> P	193 mm
! F >> H	124 mm
! R >> L	45 mm
Reset	Off

**System - pTx Volumes**

B1 Shim mode	TrueForm
Excitation	Non-sel.

**System - Tx/Rx**

Frequency 1H	297.145077 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

**Physio - Signal1**

1st Signal/Mode	None
TR	6000.0 ms
Concatenations	1

**Physio - Cardiac**

Magn. preparation	Non-sel. IR
TI 1	800 ms
TI 2	2700 ms
Fat suppr.	None
Dark blood	Off
FoV read	230 mm
FoV phase	94.4 %
Phase resolution	100 %

**Physio - PACE**

Resp. control	Off
Concatenations	1

**Inline - Common**

Subtract	Off
Measurements	1

**Inline - Common**

StdDev	Off
Save original images	On

**Inline - MIP**

MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On

**Inline - Composing**

Inline Composing	Off
Distortion Corr.	On
Mode	3D
Unfiltered images	Off

**Inline - MapIt**

Save original images	On
MapIt	None
Flip angle 1	4 deg
Flip angle 2	5 deg
Measurements	1
Contrasts	1
TR	6000.0 ms
TE	1.97 ms

**Sequence - Part 1**

Introduction	On
Dimension	3D
Elliptical scanning	Off
Reordering	Linear
Asymmetric echo	Allowed
Contrasts	1
Flow comp.	No
Multi-slice mode	Single shot
Echo spacing	5.8 ms
Bandwidth	240 Hz/Px

**Sequence - Part 2**

RF pulse type	Fast
Gradient mode	Fast
Excitation	Non-sel.
RF spoiling	On
Incr. Gradient spoiling	Off
Turbo factor	224

**Sequence - pTX Pulses****Sequence - Special**

Sparse Sampling	On
US	4.0 x
Samples/TR	200
Density	0.50
Jitter Radius	1.2
Reference Scan	External
No. Ref-Lines	32
Centric	Off
Virtual Coils	Off
Shift Inv Pulse	0 Hz
No. Iterations	15
CSM RO Resolution	0
Regularisation INV1	0.00100
Regularisation INV2	0.00100

**Sequence - Special**

Uniform	On
Denoised UNI	On
FLAWS	Off
FLAWS-hc	Off
FLAWS-hc inv.	Off
Division image	Off
T1 Map	On
Synthetic T1 0	0 ms
Synthetic T1 1	0 ms
Denoise Lambda	1
Scaling	0 10^
Echo Averaging	Off
FID Monitoring	Off

**Sequence - Assistant**

Mode	Off
------	-----