

Table of contents

\USER

FMRIF

[XT-ID:93-M-0170]]Renzo

T123-EPI_DANCLA_230615

localizer_irtfl_ptx
18sl_20TRs_GRAPPA5_FA0_10
18sl_20TRs_GRAPPA3_FA0_10
18sl_20TRs_GRAPPA5_FA0_10_dualpol
18sl_20TRs_GRAPPA3_FA0_10_dualpolSeparate
18sl_20TRs_GRAPPA3_FA0_10_dualpol
18sl_20TRs_GRAPPA3_FA0_10_dualpolSeparate_PHASE
18sl_20TRs_GRAPPA3_FA6_33_forFAcomp
18sl_20TRs_GRAPPA3_FA0_10_forFAcomp

\\USER\FMRIF\XT-ID:93-M-0170\Renzo\T123-EPI_DANCLA_230615\localizer_irtfl_ptx

TA: 1:28 PM: REF Voxel size: 1.0×1.0×2.0 mmPAT: Off Rel. SNR: 1.00 : tfl

Properties

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	On
Auto close inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	Single measurement

Contrast - Dynamic

Measurements	1
Multiple series	Each measurement

Resolution - Common

FoV read	200 mm
FoV phase	100.0 %
Slice thickness	2.0 mm
Base resolution	192
Phase resolution	100 %
Phase partial Fourier	Off
Interpolation	Off

Routine

Slice group	1
Slices	6
Dist. factor	600 %
Position	L1.5 A5.0 F25.3 mm
Orientation	Sagittal
Phase enc. dir.	A >> P
Slice group	2
Slices	8
Dist. factor	250 %
Position	L1.5 P21.0 H30.0 mm
Orientation	Transversal
Phase enc. dir.	A >> P
Slice group	3
Slices	6
Dist. factor	700 %
Position	L1.5 A10.0 F25.9 mm
Orientation	Coronal
Phase enc. dir.	R >> L
AutoAlign	---
Phase oversampling	0 %
FoV read	200 mm
FoV phase	100.0 %
Slice thickness	2.0 mm
TR	4300.0 ms
TE	3.46 ms
Averages	1
Concatenations	20
Filter	None
Coil elements	AC

Resolution - iPAT

PAT mode	None
----------	------

Resolution - Filter Image

Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off

Resolution - Filter Rawdata

Raw filter	Off
Elliptical filter	Off

Geometry - Common

Slice group	1
Slices	6
Dist. factor	600 %
Position	L1.5 A5.0 F25.3 mm
Orientation	Sagittal
Phase enc. dir.	A >> P
Slice group	2
Slices	8
Dist. factor	250 %
Position	L1.5 P21.0 H30.0 mm
Orientation	Transversal
Phase enc. dir.	A >> P
Slice group	3
Slices	6
Dist. factor	700 %
Position	L1.5 A10.0 F25.9 mm
Orientation	Coronal
Phase enc. dir.	R >> L
FoV read	200 mm
FoV phase	100.0 %
Slice thickness	2.0 mm
TR	4300.0 ms
Multi-slice mode	Sequential
Series	Ascending
Concatenations	20

Contrast - Common

TR	4300.0 ms
TE	3.46 ms
TD	0 ms
Magn. preparation	Slice-sel. IR
T1 1	840 ms
T1 2	2540 ms
Flip angle 1	5.0 deg
Flip angle 2	8.0 deg
Fat suppr.	None
Water suppr.	None

Contrast - Dynamic

Averages	1
Averaging mode	Long term
Reconstruction	Magnitude

Geometry - AutoAlign

Slice group	1
Position	L1.5 A5.0 F25.3 mm
Orientation	Sagittal
Phase enc. dir.	A >> P
Slice group	2
Position	L1.5 P21.0 H30.0 mm

Geometry - AutoAlign

Orientation	Transversal
Phase enc. dir.	A >> P
Slice group	3
Position	L1.5 A10.0 F25.9 mm
Orientation	Coronal
Phase enc. dir.	R >> L
AutoAlign	---
Initial Position	L1.5 A5.0 F25.3
L	1.5 mm
A	5.0 mm
F	25.3 mm
Initial Rotation	0.00 deg
Initial Orientation	Sagittal

Geometry - Navigator**Geometry - Tim Planning Suite**

Set-n-Go Protocol	Off
Table position	H
Table position	0 mm
Inline Composing	Off

System - Miscellaneous

Positioning mode	REF
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Adaptive Combine
Save uncombined	Off
Matrix Optimization	Off
AutoAlign	---
Coil Select Mode	Default

System - Adjustments

B0 Shim mode	Tune up
B1 Shim mode	TrueForm
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

System - Adjust Volume

Position	Isocenter
Orientation	Transversal
Rotation	0.00 deg
A >> P	263 mm
R >> L	350 mm
F >> H	350 mm
Reset	Off

System - pTx Volumes

B1 Shim mode	TrueForm
Excitation	Slice-sel.

System - Tx/Rx

Frequency 1H	297.144452 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off

System - Tx/Rx

! Ref. amplitude 1H	250.000 V
---------------------	-----------

Physio - Signal1

1st Signal/Mode	None
TR	4300.0 ms
Concatenations	20

Physio - Cardiac

Magn. preparation	Slice-sel. IR
TI 1	840 ms
TI 2	2540 ms
Fat suppr.	None
Dark blood	Off
FoV read	200 mm
FoV phase	100.0 %
Phase resolution	100 %

Physio - PACE

Resp. control	Off
Concatenations	20

Inline - Common

Subtract	Off
Measurements	1
StdDev	Off
Save original images	On

Inline - MIP

MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On

Inline - Composing

Inline Composing	Off
Distortion Corr.	Off

Inline - MapIt

Save original images	On
MapIt	None
Flip angle 1	5.0 deg
Flip angle 2	8.0 deg
Measurements	1
TR	4300.0 ms
TE	3.46 ms

Sequence - Part 1

Introduction	On
Dimension	2D
Asymmetric echo	Off
Flow comp.	No
Multi-slice mode	Sequential
Echo spacing	6.7 ms
Bandwidth	240 Hz/Px

Sequence - Part 2

RF pulse type	Normal
Gradient mode	Fast
Excitation	Slice-sel.
RF spoiling	On
Incr. Gradient spoiling	Off
Turbo factor	192

Sequence - Assistant

Mode	Off
------	-----

\\USER\FMRIF\[XT-ID:93-M-0170]]Renzo\T123-EPI_DANCLA_230615\18sl_20TRs_GRAPPA5_FA0_10

TA: 8:33 PM: REF Voxel size: 0.8×0.8×0.8 mmPAT: 5 Rel. SNR: 1.00 : nih5k

Properties

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	Off
Start measurements	Single measurement

Contrast - Dynamic

Pause after meas. 13	0.0 s
Pause after meas. 14	0.0 s
Pause after meas. 15	0.0 s
Pause after meas. 16	0.0 s
Pause after meas. 17	0.0 s
Pause after meas. 18	0.0 s
Pause after meas. 19	0.0 s

Routine

Slab group	1
Slabs	1
Position	L1.6 A0.4 F25.7 mm
Orientation	Sagittal
Phase enc. dir.	A >> P
AutoAlign	---
Slices per slab	180
FoV read	210 mm
FoV phase	100.0 %
Slice thickness	0.81 mm
TR 1	35.8 ms
TR 2	25181 ms
TE 1	16.50 ms
Averages	1
Filter	None
Coil elements	AC

Resolution - Common

FoV read	210 mm
FoV phase	100.0 %
Slice thickness	0.81 mm
Base resolution	258
Phase resolution	100 %
Slice resolution	100 %
Phase partial Fourier	Off
Slice partial Fourier	Off
Interpolation	Off

Resolution - iPAT

PAT mode	CAIPIRINHA
Acc. factor PE	1
Ref. lines PE	63
Acc. factor 3D	5
Ref. lines 3D	36
CAIPI 3D Shift	2
Reference Scan Mode	GRE/separate
CAIPI Mode (tooltip)	Skipped-CAIPI
Total PAT factor	5

Contrast - Common

TR 1	35.8 ms
TR 2	25181 ms
TE 1	16.50 ms
Multi-echo spacing	29.9 ms
Magn. preparation	Non-sel. HSN IR
TI 1	854.4 ms
TI 2	2143.2 ms
Flip angle	10 deg
Fat suppr.	None
Magn. Prep. Shots	9

Resolution - Filter Image

Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off

Resolution - Filter Rawdata

Raw filter	Off
Elliptical filter	Off

Contrast - Dynamic

Averages	1
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	20
Pause after meas. 1	0.0 s
Pause after meas. 2	0.0 s
Pause after meas. 3	0.0 s
Pause after meas. 4	0.0 s
Pause after meas. 5	0.0 s
Pause after meas. 6	0.0 s
Pause after meas. 7	0.0 s
Pause after meas. 8	0.0 s
Pause after meas. 9	0.0 s
Pause after meas. 10	0.0 s
Pause after meas. 11	0.0 s
Pause after meas. 12	0.0 s

Geometry - Common

Slab group	1
Slabs	1
Position	L1.6 A0.4 F25.7 mm
Orientation	Sagittal
Phase enc. dir.	A >> P
Slices per slab	180
FoV read	210 mm
FoV phase	100.0 %
Slice thickness	0.81 mm
TR 1	35.8 ms
TR 2	25181 ms

Geometry - AutoAlign

Slab group	1
Position	L1.6 A0.4 F25.7 mm
Orientation	Sagittal
Phase enc. dir.	A >> P
AutoAlign	---

Geometry - AutoAlign

Initial Position	L1.6 A0.4 F25.7
L	1.6 mm
A	0.4 mm
F	25.7 mm
Initial Rotation	0.00 deg
Initial Orientation	Sagittal

Geometry - Saturation

Saturation mode	Standard
Fat suppr.	None

Geometry - Tim Planning Suite

Set-n-Go Protocol	Off
Table position	H
Table position	0 mm
Inline Composing	Off

System - Miscellaneous

Positioning mode	REF
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
Save uncombined	Off
Matrix Optimization	Off
AutoAlign	---
Coil Select Mode	Default

System - Adjustments

B0 Shim mode	Brain
B1 Shim mode	TrueForm
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

System - Adjust Volume

! Position	L1.0 P2.8 H21.1 mm
! Orientation	Sagittal
! Rotation	12.30 deg
! A >> P	185 mm
! F >> H	101 mm
! R >> L	129 mm
Reset	Off

System - pTx Volumes

B1 Shim mode	TrueForm
Excitation	Non-sel.

System - Tx/Rx

Frequency 1H	297.144452 MHz
Correction factor	1
Gain	Low
Img. Scale Cor.	1.000
Reset	Off
! Ref. amplitude 1H	200.000 V

Sequence - Part 1

Introduction	On
Dimension	3D

Sequence - Part 1

Reordering	Linear
Contrasts	1
Echo spacing	0.99 ms
Bandwidth	1212 Hz/Px

Sequence - Part 2

EPI factor	29
Segmentation	9
RF pulse type	Normal
Gradient mode	Fast
Excitation	Non-sel.
RF spoiling	On
Turbo factor	36

Sequence - Special

PATRef FA	3 deg
RF duration	340 us
Ernst T1	1200 ms
PATRef prep. shots	10
Volume dummy shots	0
Dummy Measurements	0
ETL per RTEB	1
Invert PE	Off
Echo Time Shift	On
Ramp Sampling	On
NORDIC	Off
SVDPC	On
Sym VASO	Off
Dual-pol. EPI	Off
Invert RO	Off
Invert 3D	Off
Disable PF reco	Off
Save sampling	Off
PE VComp	Off
Water Exc.	Binomial-11
External PC	per Series
FIDNavs	-none-
EPI rise time factor	1.10
Mosaic DICOMs	On
Modify Ice Config	On
G-factor map	Off
GRAPPA Regularization	50000 10^-6
HSN RF power scale	3.00
Inversion Delay	200 ms
Relaxation Delay	0 ms
Var. FA /MAGEC	0

Sequence - Assistant

Mode	Off
------	-----

\\USER\FMRIF\[XT-ID:93-M-0170]]Renzo\T123-EPI_DANCLA_230615\18sl_20TRs_GRAPPA3_FA0_10

TA: 13:42 PM: REF Voxel size: 0.8×0.8×0.8 mmPAT: 3 Rel. SNR: 1.00 : nih5k

Properties

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	Off
Start measurements	Single measurement

Routine

Slab group	1
Slabs	1
Position	L1.6 A0.4 F25.7 mm
Orientation	Sagittal
Phase enc. dir.	A >> P
AutoAlign	---
Slices per slab	180
FoV read	210 mm
FoV phase	100.0 %
Slice thickness	0.81 mm
TR 1	35.8 ms
TR 2	40647 ms
TE 1	16.50 ms
Averages	1
Filter	None
Coil elements	AC

Contrast - Common

TR 1	35.8 ms
TR 2	40647 ms
TE 1	16.50 ms
Multi-echo spacing	29.9 ms
Magn. preparation	Non-sel. HSN IR
TI 1	1284 ms
TI 2	3432 ms
Flip angle	10 deg
Fat suppr.	None
Magn. Prep. Shots	9

Contrast - Dynamic

Averages	1
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	20
Pause after meas. 1	0.0 s
Pause after meas. 2	0.0 s
Pause after meas. 3	0.0 s
Pause after meas. 4	0.0 s
Pause after meas. 5	0.0 s
Pause after meas. 6	0.0 s
Pause after meas. 7	0.0 s
Pause after meas. 8	0.0 s
Pause after meas. 9	0.0 s
Pause after meas. 10	0.0 s
Pause after meas. 11	0.0 s
Pause after meas. 12	0.0 s

Contrast - Dynamic

Pause after meas. 13	0.0 s
Pause after meas. 14	0.0 s
Pause after meas. 15	0.0 s
Pause after meas. 16	0.0 s
Pause after meas. 17	0.0 s
Pause after meas. 18	0.0 s
Pause after meas. 19	0.0 s

Resolution - Common

FoV read	210 mm
FoV phase	100.0 %
Slice thickness	0.81 mm
Base resolution	258
Phase resolution	100 %
Slice resolution	100 %
Phase partial Fourier	Off
Slice partial Fourier	Off
Interpolation	Off

Resolution - iPAT

PAT mode	CAIPIRINHA
Acc. factor PE	1
Ref. lines PE	63
Acc. factor 3D	3
Ref. lines 3D	36
CAIPI 3D Shift	2
Reference Scan Mode	GRE/separate
CAIPI Mode (tooltip)	Skipped-CAIPI
Total PAT factor	3

Resolution - Filter Image

Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off

Resolution - Filter Rawdata

Raw filter	Off
Elliptical filter	Off

Geometry - Common

Slab group	1
Slabs	1
Position	L1.6 A0.4 F25.7 mm
Orientation	Sagittal
Phase enc. dir.	A >> P
Slices per slab	180
FoV read	210 mm
FoV phase	100.0 %
Slice thickness	0.81 mm
TR 1	35.8 ms
TR 2	40647 ms

Geometry - AutoAlign

Slab group	1
Position	L1.6 A0.4 F25.7 mm
Orientation	Sagittal
Phase enc. dir.	A >> P
AutoAlign	---

Geometry - AutoAlign

Initial Position	L1.6 A0.4 F25.7
L	1.6 mm
A	0.4 mm
F	25.7 mm
Initial Rotation	0.00 deg
Initial Orientation	Sagittal

Geometry - Saturation

Saturation mode	Standard
Fat suppr.	None

Geometry - Tim Planning Suite

Set-n-Go Protocol	Off
Table position	H
Table position	0 mm
Inline Composing	Off

System - Miscellaneous

Positioning mode	REF
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
Save uncombined	Off
Matrix Optimization	Off
AutoAlign	---
Coil Select Mode	Default

System - Adjustments

B0 Shim mode	Brain
B1 Shim mode	TrueForm
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

System - Adjust Volume

! Position	L1.0 P2.8 H21.1 mm
! Orientation	Sagittal
! Rotation	12.30 deg
! A >> P	185 mm
! F >> H	101 mm
! R >> L	129 mm
Reset	Off

System - pTx Volumes

B1 Shim mode	TrueForm
Excitation	Non-sel.

System - Tx/Rx

Frequency 1H	297.144452 MHz
Correction factor	1
Gain	Low
Img. Scale Cor.	1.000
Reset	Off
! Ref. amplitude 1H	200.000 V

Sequence - Part 1

Introduction	On
Dimension	3D

Sequence - Part 1

Reordering	Linear
Contrasts	1
Echo spacing	0.99 ms
Bandwidth	1212 Hz/Px

Sequence - Part 2

EPI factor	29
Segmentation	9
RF pulse type	Normal
Gradient mode	Fast
Excitation	Non-sel.
RF spoiling	On
Turbo factor	60

Sequence - Special

PATRef FA	3 deg
RF duration	340 us
Ernst T1	1200 ms
PATRef prep. shots	10
Volume dummy shots	0
Dummy Measurements	0
ETL per RTEB	1
Invert PE	Off
Echo Time Shift	On
Ramp Sampling	On
NORDIC	Off
SVDPC	On
Sym VASO	Off
Dual-pol. EPI	Off
Invert RO	Off
Invert 3D	Off
Disable PF reco	Off
Save sampling	Off
PE VComp	Off
Water Exc.	Binomial-11
External PC	per Series
FIDNavs	-none-
EPI rise time factor	1.10
Mosaic DICOMs	On
Modify Ice Config	On
G-factor map	Off
GRAPPA Regularization	50000 10^-6
HSN RF power scale	3.00
Inversion Delay	200 ms
Relaxation Delay	0 ms
Var. FA /MAGEC	0

Sequence - Assistant

Mode	Off
------	-----

\\USER\FMRIF\[XT-ID:93-M-0170]]Renzo\T123-EPI_DANCLA_230615\18sl_20TRs_GRAPPA5_FA0_10
_dualpol

TA: 10:03 PM: FIX Voxel size: 0.8×0.8×0.8 mmPAT: 5 Rel. SNR: 1.00 : nih5k

Properties

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	Off
Start measurements	Single measurement

Routine

Slab group	1
Slabs	1
Position	L1.6 A0.4 F25.7 mm
Orientation	Sagittal
Phase enc. dir.	A >> P
AutoAlign	---
Slab Scale	-10 %
Slices per slab	200
FoV read	210 mm
FoV phase	100.0 %
Slice thickness	0.81 mm
TR 1	38.5 ms
TR 2	59406 ms
TE 1	16.50 ms
Averages	2
Filter	None
Coil elements	AC

Contrast - Common

TR 1	38.5 ms
TR 2	59406 ms
TE 1	16.50 ms
Multi-echo spacing	29.9 ms
Magn. preparation	Non-sel. HSN IR
TI 1	980 ms
TI 2	2520 ms
Flip angle	10 deg
Fat suppr.	None
Magn. Prep. Shots	9

Contrast - Dynamic

Averages	2
Averaging mode	Short term
Reconstruction	Magnitude
Measurements	5
Pause after meas. 1	0.0 s
Pause after meas. 2	0.0 s
Pause after meas. 3	0.0 s
Pause after meas. 4	0.0 s

Resolution - Common

FoV read	210 mm
FoV phase	100.0 %
Slice thickness	0.81 mm
Base resolution	258

Resolution - Common

Phase resolution	100 %
Slice resolution	100 %
Phase partial Fourier	Off
Slice partial Fourier	Off
Interpolation	Off

Resolution - iPAT

PAT mode	CAIPIRINHA
Acc. factor PE	1
Ref. lines PE	63
Acc. factor 3D	5
Ref. lines 3D	36
CAIPI 3D Shift	2
Reference Scan Mode	GRE/separate
CAIPI Mode (tooltip)	Skipped-CAIPI
Total PAT factor	5

Resolution - Filter Image

Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off

Resolution - Filter Rawdata

Raw filter	Off
Elliptical filter	Off

Geometry - Common

Slab group	1
Slabs	1
Position	L1.6 A0.4 F25.7 mm
Orientation	Sagittal
Phase enc. dir.	A >> P
Slab Scale	-10 %
Slices per slab	200
FoV read	210 mm
FoV phase	100.0 %
Slice thickness	0.81 mm
TR 1	38.5 ms
TR 2	59406 ms

Geometry - AutoAlign

Slab group	1
Position	L1.6 A0.4 F25.7 mm
Orientation	Sagittal
Phase enc. dir.	A >> P
AutoAlign	---
Initial Position	L1.6 A0.4 F25.7
L	1.6 mm
A	0.4 mm
F	25.7 mm
Initial Rotation	0.00 deg
Initial Orientation	Sagittal

Geometry - Saturation

Saturation mode	Standard
Fat suppr.	None

Geometry - Tim Planning Suite

Set-n-Go Protocol	Off
Table position	H
Table position	0 mm
Inline Composing	Off

System - Miscellaneous

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
Save uncombined	Off
Matrix Optimization	Off
AutoAlign	---
Coil Select Mode	Default

System - Adjustments

B0 Shim mode	Brain
B1 Shim mode	TrueForm
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

System - Adjust Volume

! Position	L1.0 P2.8 H21.1 mm
! Orientation	Sagittal
! Rotation	12.30 deg
! A >> P	185 mm
! F >> H	101 mm
! R >> L	129 mm
Reset	Off

System - pTx Volumes

B1 Shim mode	TrueForm
Excitation	Slab-sel.

System - Tx/Rx

Frequency 1H	297.144452 MHz
Correction factor	1
Gain	Low
Img. Scale Cor.	1.000
Reset	Off
! Ref. amplitude 1H	200.000 V

Sequence - Part 1

Introduction	On
Dimension	3D
Reordering	Linear
Contrasts	1
Echo spacing	0.99 ms
Bandwidth	1212 Hz/Px

Sequence - Part 2

EPI factor	29
Segmentation	9
RF pulse type	Normal
Gradient mode	Fast
Excitation	Slab-sel.
RF spoiling	On

Sequence - Part 2

Turbo factor	40
--------------	----

Sequence - Special

PATRef FA	3 deg
RF duration	340 us
RF BWT product	8
Ernst T1	1200 ms
PATRef prep. shots	10
Volume dummy shots	0
Dummy Measurements	0
ETL per RTEB	1
Invert PE	Off
Echo Time Shift	On
Ramp Sampling	On
NORDIC	Off
SVDPC	Off
Sym VASO	Off
Dual-pol. EPI	On
Invert RO	Off
Invert 3D	Off
Disable PF reco	Off
Save sampling	Off
PE VComp	Off
Water Exc.	Binomial-11
External PC	per Series
FIDNavs	-none-
EPI rise time factor	1.10
Mosaic DICOMs	On
Modify Ice Config	On
G-factor map	Off
GRAPPA Regularization	50000 10 ⁻⁶
HSN RF power scale	3.00
Inversion Delay	200 ms
Relaxation Delay	0 ms
Var. FA /MAGEC	0

Sequence - Assistant

Mode	Off
------	-----

\\USER\FMRIF\[XT-ID:93-M-0170]]Renzo\T123-EPI_DANCLA_230615\18sl_20TRs_GRAPPA3_FA0_10
_dualpolSeparate

TA: 7:24 PM: FIX Voxel size: 0.8×0.8×0.8 mmPAT: 3 Rel. SNR: 1.00 : nih5k

Properties

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	Off
Start measurements	Single measurement

Routine

Slab group	1
Slabs	1
Position	L1.6 A0.4 F25.7 mm
Orientation	Sagittal
Phase enc. dir.	A >> P
AutoAlign	---
Slab Scale	-10 %
Slices per slab	180
FoV read	210 mm
FoV phase	100.0 %
Slice thickness	0.81 mm
TR 1	38.5 ms
TR 2	43563 ms
TE 1	16.50 ms
Averages	1
Filter	None
Coil elements	AC

Contrast - Common

TR 1	38.5 ms
TR 2	43563 ms
TE 1	16.50 ms
Multi-echo spacing	29.9 ms
Magn. preparation	Non-sel. HSN IR
TI 1	1365 ms
TI 2	3675 ms
Flip angle	10 deg
Fat suppr.	None
Magn. Prep. Shots	9

Contrast - Dynamic

Averages	1
Averaging mode	Short term
Reconstruction	Magnitude
Measurements	10
Pause after meas. 1	0.0 s
Pause after meas. 2	0.0 s
Pause after meas. 3	0.0 s
Pause after meas. 4	0.0 s
Pause after meas. 5	0.0 s
Pause after meas. 6	0.0 s
Pause after meas. 7	0.0 s
Pause after meas. 8	0.0 s
Pause after meas. 9	0.0 s

Resolution - Common

FoV read	210 mm
FoV phase	100.0 %
Slice thickness	0.81 mm
Base resolution	258
Phase resolution	100 %
Slice resolution	100 %
Phase partial Fourier	Off
Slice partial Fourier	Off
Interpolation	Off

Resolution - iPAT

PAT mode	CAIPIRINHA
Acc. factor PE	1
Ref. lines PE	63
Acc. factor 3D	3
Ref. lines 3D	36
CAIPI 3D Shift	1
Reference Scan Mode	GRE/separate
CAIPI Mode (tooltip)	Skipped-CAIPI
Total PAT factor	3

Resolution - Filter Image

Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off

Resolution - Filter Rawdata

Raw filter	Off
Elliptical filter	Off

Geometry - Common

Slab group	1
Slabs	1
Position	L1.6 A0.4 F25.7 mm
Orientation	Sagittal
Phase enc. dir.	A >> P
Slab Scale	-10 %
Slices per slab	180
FoV read	210 mm
FoV phase	100.0 %
Slice thickness	0.81 mm
TR 1	38.5 ms
TR 2	43563 ms

Geometry - AutoAlign

Slab group	1
Position	L1.6 A0.4 F25.7 mm
Orientation	Sagittal
Phase enc. dir.	A >> P
AutoAlign	---
Initial Position	L1.6 A0.4 F25.7
L	1.6 mm
A	0.4 mm
F	25.7 mm
Initial Rotation	0.00 deg
Initial Orientation	Sagittal

Geometry - Saturation

Saturation mode	Standard
Fat suppr.	None

Geometry - Tim Planning Suite

Set-n-Go Protocol	Off
Table position	H
Table position	0 mm
Inline Composing	Off

System - Miscellaneous

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
Save uncombined	Off
Matrix Optimization	Off
AutoAlign	---
Coil Select Mode	Default

System - Adjustments

B0 Shim mode	Brain
B1 Shim mode	TrueForm
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

System - Adjust Volume

! Position	L1.0 P2.8 H21.1 mm
! Orientation	Sagittal
! Rotation	12.30 deg
! A >> P	185 mm
! F >> H	101 mm
! R >> L	129 mm
Reset	Off

System - pTx Volumes

B1 Shim mode	TrueForm
Excitation	Slab-sel.

System - Tx/Rx

Frequency 1H	297.144452 MHz
Correction factor	1
Gain	Low
Img. Scale Cor.	1.000
Reset	Off
! Ref. amplitude 1H	200.000 V

Sequence - Part 1

Introduction	On
Dimension	3D
Reordering	Linear
Contrasts	1
Echo spacing	0.99 ms
Bandwidth	1212 Hz/Px

Sequence - Part 2

EPI factor	29
Segmentation	9

Sequence - Part 2

RF pulse type	Normal
Gradient mode	Fast
Excitation	Slab-sel.
RF spoiling	On
Turbo factor	60

Sequence - Special

PATRef FA	3 deg
RF duration	340 us
RF BWT product	8
Ernst T1	1200 ms
PATRef prep. shots	10
Volume dummy shots	0
Dummy Measurements	0
ETL per RTEB	1
Invert PE	Off
Echo Time Shift	On
Ramp Sampling	On
NORDIC	Off
SVDPC	Off
Sym VASO	Off
Dual-pol. EPI	On
Invert RO	Off
Invert 3D	Off
Disable PF reco	Off
Save sampling	Off
PE VComp	Off
Water Exc.	Binomial-11
External PC	per Series
FIDNavs	-none-
EPI rise time factor	1.10
Mosaic DICOMs	On
Modify Ice Config	On
G-factor map	Off
GRAPPA Regularization	50000 10^-6
HSN RF power scale	3.00
Inversion Delay	200 ms
Relaxation Delay	0 ms
Var. FA /MAGEC	0

Sequence - Assistant

Mode	Off
------	-----

\\USER\FMRIF\[XT-ID:93-M-0170]]Renzo\T123-EPI_DANCLA_230615\18sl_20TRs_GRAPPA3_FA0_10
_dualpol

TA: 14:40 PM: FIX Voxel size: 0.8×0.8×0.8 mmPAT: 3 Rel. SNR: 1.00 : nih5k

Properties

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	Off
Start measurements	Single measurement

Routine

Slab group	1
Slabs	1
Position	L1.6 A0.4 F25.7 mm
Orientation	Sagittal
Phase enc. dir.	A >> P
AutoAlign	---
Slab Scale	-10 %
Slices per slab	180
FoV read	210 mm
FoV phase	100.0 %
Slice thickness	0.81 mm
TR 1	38.5 ms
TR 2	87126 ms
TE 1	16.50 ms
Averages	2
Filter	None
Coil elements	AC

Contrast - Common

TR 1	38.5 ms
TR 2	87126 ms
TE 1	16.50 ms
Multi-echo spacing	29.9 ms
Magn. preparation	Non-sel. HSN IR
TI 1	1365 ms
TI 2	3675 ms
Flip angle	10 deg
Fat suppr.	None
Magn. Prep. Shots	9

Contrast - Dynamic

Averages	2
Averaging mode	Short term
Reconstruction	Magnitude
Measurements	5
Pause after meas. 1	0.0 s
Pause after meas. 2	0.0 s
Pause after meas. 3	0.0 s
Pause after meas. 4	0.0 s

Resolution - Common

FoV read	210 mm
FoV phase	100.0 %
Slice thickness	0.81 mm
Base resolution	258

Resolution - Common

Phase resolution	100 %
Slice resolution	100 %
Phase partial Fourier	Off
Slice partial Fourier	Off
Interpolation	Off

Resolution - iPAT

PAT mode	CAIPIRINHA
Acc. factor PE	1
Ref. lines PE	63
Acc. factor 3D	3
Ref. lines 3D	36
CAIPI 3D Shift	1
Reference Scan Mode	GRE/separate
CAIPI Mode (tooltip)	Skipped-CAIPI
Total PAT factor	3

Resolution - Filter Image

Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off

Resolution - Filter Rawdata

Raw filter	Off
Elliptical filter	Off

Geometry - Common

Slab group	1
Slabs	1
Position	L1.6 A0.4 F25.7 mm
Orientation	Sagittal
Phase enc. dir.	A >> P
Slab Scale	-10 %
Slices per slab	180
FoV read	210 mm
FoV phase	100.0 %
Slice thickness	0.81 mm
TR 1	38.5 ms
TR 2	87126 ms

Geometry - AutoAlign

Slab group	1
Position	L1.6 A0.4 F25.7 mm
Orientation	Sagittal
Phase enc. dir.	A >> P
AutoAlign	---
Initial Position	L1.6 A0.4 F25.7
L	1.6 mm
A	0.4 mm
F	25.7 mm
Initial Rotation	0.00 deg
Initial Orientation	Sagittal

Geometry - Saturation

Saturation mode	Standard
Fat suppr.	None

Geometry - Tim Planning Suite

Set-n-Go Protocol	Off
Table position	H
Table position	0 mm
Inline Composing	Off

System - Miscellaneous

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
Save uncombined	Off
Matrix Optimization	Off
AutoAlign	---
Coil Select Mode	Default

System - Adjustments

B0 Shim mode	Brain
B1 Shim mode	TrueForm
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

System - Adjust Volume

! Position	L1.0 P2.8 H21.1 mm
! Orientation	Sagittal
! Rotation	12.30 deg
! A >> P	185 mm
! F >> H	101 mm
! R >> L	129 mm
Reset	Off

System - pTx Volumes

B1 Shim mode	TrueForm
Excitation	Slab-sel.

System - Tx/Rx

Frequency 1H	297.144452 MHz
Correction factor	1
Gain	Low
Img. Scale Cor.	1.000
Reset	Off
! Ref. amplitude 1H	200.000 V

Sequence - Part 1

Introduction	On
Dimension	3D
Reordering	Linear
Contrasts	1
Echo spacing	0.99 ms
Bandwidth	1212 Hz/Px

Sequence - Part 2

EPI factor	29
Segmentation	9
RF pulse type	Normal
Gradient mode	Fast
Excitation	Slab-sel.
RF spoiling	On

Sequence - Part 2

Turbo factor	60
--------------	----

Sequence - Special

PATRef FA	3 deg
RF duration	340 us
RF BWT product	8
Ernst T1	1200 ms
PATRef prep. shots	10
Volume dummy shots	0
Dummy Measurements	0
ETL per RTEB	1
Invert PE	Off
Echo Time Shift	On
Ramp Sampling	On
NORDIC	Off
SVDPC	Off
Sym VASO	Off
Dual-pol. EPI	On
Invert RO	Off
Invert 3D	Off
Disable PF reco	Off
Save sampling	Off
PE VComp	Off
Water Exc.	Binomial-11
External PC	per Series
FIDNavs	-none-
EPI rise time factor	1.10
Mosaic DICOMs	On
Modify Ice Config	On
G-factor map	Off
GRAPPA Regularization	50000 10 ⁻⁶
HSN RF power scale	3.00
Inversion Delay	200 ms
Relaxation Delay	0 ms
Var. FA /MAGEC	0

Sequence - Assistant

Mode	Off
------	-----

\\USER\FMRIF\[XT-ID:93-M-0170]]Renzo\T123-EPI_DANCLA_230615\18sl_20TRs_GRAPPA3_FA0_10
_dualpolSeparate_PHASE

TA: 7:24 PM: FIX Voxel size: 0.8×0.8×0.8 mmPAT: 3 Rel. SNR: 1.00 : nih5k

Properties

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	Off
Start measurements	Single measurement

Routine

Slab group	1
Slabs	1
Position	L1.6 A0.4 F25.7 mm
Orientation	Sagittal
Phase enc. dir.	A >> P
AutoAlign	---
Slab Scale	-10 %
Slices per slab	180
FoV read	210 mm
FoV phase	100.0 %
Slice thickness	0.81 mm
TR 1	38.5 ms
TR 2	43563 ms
TE 1	16.50 ms
Averages	1
Filter	None
Coil elements	AC

Contrast - Common

TR 1	38.5 ms
TR 2	43563 ms
TE 1	16.50 ms
Multi-echo spacing	29.9 ms
Magn. preparation	Non-sel. HSN IR
TI 1	1365 ms
TI 2	3675 ms
Flip angle	10 deg
Fat suppr.	None
Magn. Prep. Shots	9

Contrast - Dynamic

Averages	1
Averaging mode	Short term
Reconstruction	Magn./Phase
Measurements	10
Pause after meas. 1	0.0 s
Pause after meas. 2	0.0 s
Pause after meas. 3	0.0 s
Pause after meas. 4	0.0 s
Pause after meas. 5	0.0 s
Pause after meas. 6	0.0 s
Pause after meas. 7	0.0 s
Pause after meas. 8	0.0 s
Pause after meas. 9	0.0 s

Resolution - Common

FoV read	210 mm
FoV phase	100.0 %
Slice thickness	0.81 mm
Base resolution	258
Phase resolution	100 %
Slice resolution	100 %
Phase partial Fourier	Off
Slice partial Fourier	Off
Interpolation	Off

Resolution - iPAT

PAT mode	CAIPIRINHA
Acc. factor PE	1
Ref. lines PE	63
Acc. factor 3D	3
Ref. lines 3D	36
CAIPI 3D Shift	1
Reference Scan Mode	GRE/separate
CAIPI Mode (tooltip)	Skipped-CAIPI
Total PAT factor	3

Resolution - Filter Image

Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off

Resolution - Filter Rawdata

Raw filter	Off
Elliptical filter	Off

Geometry - Common

Slab group	1
Slabs	1
Position	L1.6 A0.4 F25.7 mm
Orientation	Sagittal
Phase enc. dir.	A >> P
Slab Scale	-10 %
Slices per slab	180
FoV read	210 mm
FoV phase	100.0 %
Slice thickness	0.81 mm
TR 1	38.5 ms
TR 2	43563 ms

Geometry - AutoAlign

Slab group	1
Position	L1.6 A0.4 F25.7 mm
Orientation	Sagittal
Phase enc. dir.	A >> P
AutoAlign	---
Initial Position	L1.6 A0.4 F25.7
L	1.6 mm
A	0.4 mm
F	25.7 mm
Initial Rotation	0.00 deg
Initial Orientation	Sagittal

Geometry - Saturation

Saturation mode	Standard
Fat suppr.	None

Geometry - Tim Planning Suite

Set-n-Go Protocol	Off
Table position	H
Table position	0 mm
Inline Composing	Off

System - Miscellaneous

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
Save uncombined	Off
Matrix Optimization	Off
AutoAlign	---
Coil Select Mode	Default

System - Adjustments

B0 Shim mode	Brain
B1 Shim mode	TrueForm
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

System - Adjust Volume

! Position	L1.0 P2.8 H21.1 mm
! Orientation	Sagittal
! Rotation	12.30 deg
! A >> P	185 mm
! F >> H	101 mm
! R >> L	129 mm
Reset	Off

System - pTx Volumes

B1 Shim mode	TrueForm
Excitation	Slab-sel.

System - Tx/Rx

Frequency 1H	297.144452 MHz
Correction factor	1
Gain	Low
Img. Scale Cor.	1.000
Reset	Off
! Ref. amplitude 1H	200.000 V

Sequence - Part 1

Introduction	On
Dimension	3D
Reordering	Linear
Contrasts	1
Echo spacing	0.99 ms
Bandwidth	1212 Hz/Px

Sequence - Part 2

EPI factor	29
Segmentation	9

Sequence - Part 2

RF pulse type	Normal
Gradient mode	Fast
Excitation	Slab-sel.
RF spoiling	On
Turbo factor	60

Sequence - Special

PATRef FA	3 deg
RF duration	340 us
RF BWT product	8
Ernst T1	1200 ms
PATRef prep. shots	10
Volume dummy shots	0
Dummy Measurements	0
ETL per RTEB	1
Invert PE	Off
Echo Time Shift	On
Ramp Sampling	On
NORDIC	Off
SVDPC	On
Sym VASO	Off
Dual-pol. EPI	On
Invert RO	Off
Invert 3D	Off
Disable PF reco	Off
Save sampling	Off
PE VComp	Off
Water Exc.	Binomial-11
External PC	per Series
FIDNavs	-none-
EPI rise time factor	1.10
Mosaic DICOMs	On
Modify Ice Config	On
G-factor map	Off
GRAPPA Regularization	50000 10^-6
HSN RF power scale	3.00
Inversion Delay	200 ms
Relaxation Delay	0 ms
Var. FA /MAGEC	0

Sequence - Assistant

Mode	Off
------	-----

\\USER\FMRIF\[XT-ID:93-M-0170]\Renzo\T123-EPI_DANCLA_230615\18sl_20TRs_GRAPPA3_FA6_33_forFAcomp

TA: 14:04 PM: FIX Voxel size: 0.8×0.8×0.8 mmPAT: 3 Rel. SNR: 1.00 : nih5k

Properties

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	Off
Start measurements	Single measurement

Routine

Slab group	1
Slabs	1
Position	L1.6 A0.4 F25.7 mm
Orientation	Sagittal
Phase enc. dir.	A >> P
AutoAlign	---
Slab Scale	-10 %
Slices per slab	180
FoV read	210 mm
FoV phase	100.0 %
Slice thickness	0.81 mm
TR 1	38.5 ms
TR 2	41763 ms
TE 1	16.50 ms
Averages	1
Filter	None
Coil elements	AC

Contrast - Common

TR 1	38.5 ms
TR 2	41763 ms
TE 1	16.50 ms
Multi-echo spacing	29.9 ms
Magn. preparation	Non-sel. HSN IR
TI 1	1165 ms
TI 2	3475 ms
Flip angle	33 deg
Fat suppr.	None
Magn. Prep. Shots	9

Contrast - Dynamic

Averages	1
Averaging mode	Short term
Reconstruction	Magnitude
Measurements	20
Pause after meas. 1	0.0 s
Pause after meas. 2	0.0 s
Pause after meas. 3	0.0 s
Pause after meas. 4	0.0 s
Pause after meas. 5	0.0 s
Pause after meas. 6	0.0 s
Pause after meas. 7	0.0 s
Pause after meas. 8	0.0 s
Pause after meas. 9	0.0 s
Pause after meas. 10	0.0 s

Contrast - Dynamic

Pause after meas. 11	0.0 s
Pause after meas. 12	0.0 s
Pause after meas. 13	0.0 s
Pause after meas. 14	0.0 s
Pause after meas. 15	0.0 s
Pause after meas. 16	0.0 s
Pause after meas. 17	0.0 s
Pause after meas. 18	0.0 s
Pause after meas. 19	0.0 s

Resolution - Common

FoV read	210 mm
FoV phase	100.0 %
Slice thickness	0.81 mm
Base resolution	258
Phase resolution	100 %
Slice resolution	100 %
Phase partial Fourier	Off
Slice partial Fourier	Off
Interpolation	Off

Resolution - iPAT

PAT mode	CAIPIRINHA
Acc. factor PE	1
Ref. lines PE	63
Acc. factor 3D	3
Ref. lines 3D	36
CAIPI 3D Shift	1
Reference Scan Mode	GRE/separate
CAIPI Mode (tooltip)	Skipped-CAIPI
Total PAT factor	3

Resolution - Filter Image

Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off

Resolution - Filter Rawdata

Raw filter	Off
Elliptical filter	Off

Geometry - Common

Slab group	1
Slabs	1
Position	L1.6 A0.4 F25.7 mm
Orientation	Sagittal
Phase enc. dir.	A >> P
Slab Scale	-10 %
Slices per slab	180
FoV read	210 mm
FoV phase	100.0 %
Slice thickness	0.81 mm
TR 1	38.5 ms
TR 2	41763 ms

Geometry - AutoAlign

Slab group	1
------------	---

Geometry - AutoAlign

Position	L1.6 A0.4 F25.7 mm
Orientation	Sagittal
Phase enc. dir.	A >> P
AutoAlign	---
Initial Position	L1.6 A0.4 F25.7
L	1.6 mm
A	0.4 mm
F	25.7 mm
Initial Rotation	0.00 deg
Initial Orientation	Sagittal

Geometry - Saturation

Saturation mode	Standard
Fat suppr.	None

Geometry - Tim Planning Suite

Set-n-Go Protocol	Off
Table position	H
Table position	0 mm
Inline Composing	Off

System - Miscellaneous

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
Save uncombined	Off
Matrix Optimization	Off
AutoAlign	---
Coil Select Mode	Default

System - Adjustments

B0 Shim mode	Brain
B1 Shim mode	TrueForm
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

System - Adjust Volume

! Position	L1.0 P2.8 H21.1 mm
! Orientation	Sagittal
! Rotation	12.30 deg
! A >> P	185 mm
! F >> H	101 mm
! R >> L	129 mm
Reset	Off

System - pTx Volumes

B1 Shim mode	TrueForm
Excitation	Slab-sel.

System - Tx/Rx

Frequency 1H	297.144452 MHz
Correction factor	1
Gain	Low
Img. Scale Cor.	1.000
Reset	Off
! Ref. amplitude 1H	200.000 V

Sequence - Part 1

Introduction	On
Dimension	3D
Reordering	Linear
Contrasts	1
Echo spacing	0.99 ms
Bandwidth	1212 Hz/Px

Sequence - Part 2

EPI factor	29
Segmentation	9
RF pulse type	Normal
Gradient mode	Fast
Excitation	Slab-sel.
RF spoiling	On
Turbo factor	60

Sequence - Special

PATRef FA	3 deg
RF duration	340 us
RF BWT product	8
Ernst T1	1200 ms
PATRef prep. shots	10
Volume dummy shots	0
Dummy Measurements	0
ETL per RTEB	1
CHECK FLIP ANGLE!	On
Invert PE	Off
Echo Time Shift	On
Ramp Sampling	On
NORDIC	Off
SVDPC	On
Sym VASO	Off
Dual-pol. EPI	Off
Invert RO	Off
Invert 3D	Off
Disable PF reco	Off
Save sampling	Off
PE VComp	Off
Water Exc.	Binomial-11
External PC	per Series
FIDNavs	-none-
EPI rise time factor	1.10
Mosaic DICOMs	On
Modify Ice Config	On
G-factor map	Off
GRAPPA Regularization	50000 10^-6
HSN RF power scale	3.00
Inversion Delay	0 ms
Relaxation Delay	0 ms
Var. FA /MAGEC	6

Sequence - Assistant

Mode	Off
------	-----

\\USER\FMRIF\[XT-ID:93-M-0170]]Renzo\T123-EPI_DANCLA_230615\18sl_20TRs_GRAPPA3_FA0_10_forFAcomp

TA: 14:04 PM: FIX Voxel size: 0.8×0.8×0.8 mmPAT: 3 Rel. SNR: 1.00 : nih5k

Properties

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	Off
Start measurements	Single measurement

Routine

Slab group	1
Slabs	1
Position	L1.6 A0.4 F25.7 mm
Orientation	Sagittal
Phase enc. dir.	A >> P
AutoAlign	---
Slab Scale	-10 %
Slices per slab	180
FoV read	210 mm
FoV phase	100.0 %
Slice thickness	0.81 mm
TR 1	38.5 ms
TR 2	41763 ms
TE 1	16.50 ms
Averages	1
Filter	None
Coil elements	AC

Contrast - Common

TR 1	38.5 ms
TR 2	41763 ms
TE 1	16.50 ms
Multi-echo spacing	29.9 ms
Magn. preparation	Non-sel. HSN IR
TI 1	1165 ms
TI 2	3475 ms
Flip angle	10 deg
Fat suppr.	None
Magn. Prep. Shots	9

Contrast - Dynamic

Averages	1
Averaging mode	Short term
Reconstruction	Magnitude
Measurements	20
Pause after meas. 1	0.0 s
Pause after meas. 2	0.0 s
Pause after meas. 3	0.0 s
Pause after meas. 4	0.0 s
Pause after meas. 5	0.0 s
Pause after meas. 6	0.0 s
Pause after meas. 7	0.0 s
Pause after meas. 8	0.0 s
Pause after meas. 9	0.0 s
Pause after meas. 10	0.0 s

Contrast - Dynamic

Pause after meas. 11	0.0 s
Pause after meas. 12	0.0 s
Pause after meas. 13	0.0 s
Pause after meas. 14	0.0 s
Pause after meas. 15	0.0 s
Pause after meas. 16	0.0 s
Pause after meas. 17	0.0 s
Pause after meas. 18	0.0 s
Pause after meas. 19	0.0 s

Resolution - Common

FoV read	210 mm
FoV phase	100.0 %
Slice thickness	0.81 mm
Base resolution	258
Phase resolution	100 %
Slice resolution	100 %
Phase partial Fourier	Off
Slice partial Fourier	Off
Interpolation	Off

Resolution - iPAT

PAT mode	CAIPIRINHA
Acc. factor PE	1
Ref. lines PE	63
Acc. factor 3D	3
Ref. lines 3D	36
CAIPI 3D Shift	1
Reference Scan Mode	GRE/separate
CAIPI Mode (tooltip)	Skipped-CAIPI
Total PAT factor	3

Resolution - Filter Image

Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off

Resolution - Filter Rawdata

Raw filter	Off
Elliptical filter	Off

Geometry - Common

Slab group	1
Slabs	1
Position	L1.6 A0.4 F25.7 mm
Orientation	Sagittal
Phase enc. dir.	A >> P
Slab Scale	-10 %
Slices per slab	180
FoV read	210 mm
FoV phase	100.0 %
Slice thickness	0.81 mm
TR 1	38.5 ms
TR 2	41763 ms

Geometry - AutoAlign

Slab group	1
------------	---

Geometry - AutoAlign

Position	L1.6 A0.4 F25.7 mm
Orientation	Sagittal
Phase enc. dir.	A >> P
AutoAlign	---
Initial Position	L1.6 A0.4 F25.7
L	1.6 mm
A	0.4 mm
F	25.7 mm
Initial Rotation	0.00 deg
Initial Orientation	Sagittal

Geometry - Saturation

Saturation mode	Standard
Fat suppr.	None

Geometry - Tim Planning Suite

Set-n-Go Protocol	Off
Table position	H
Table position	0 mm
Inline Composing	Off

System - Miscellaneous

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
Save uncombined	Off
Matrix Optimization	Off
AutoAlign	---
Coil Select Mode	Default

System - Adjustments

B0 Shim mode	Brain
B1 Shim mode	TrueForm
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

System - Adjust Volume

! Position	L1.0 P2.8 H21.1 mm
! Orientation	Sagittal
! Rotation	12.30 deg
! A >> P	185 mm
! F >> H	101 mm
! R >> L	129 mm
Reset	Off

System - pTx Volumes

B1 Shim mode	TrueForm
Excitation	Slab-sel.

System - Tx/Rx

Frequency 1H	297.144452 MHz
Correction factor	1
Gain	Low
Img. Scale Cor.	1.000
Reset	Off
! Ref. amplitude 1H	200.000 V

Sequence - Part 1

Introduction	On
Dimension	3D
Reordering	Linear
Contrasts	1
Echo spacing	0.99 ms
Bandwidth	1212 Hz/Px

Sequence - Part 2

EPI factor	29
Segmentation	9
RF pulse type	Normal
Gradient mode	Fast
Excitation	Slab-sel.
RF spoiling	On
Turbo factor	60

Sequence - Special

PATRef FA	3 deg
RF duration	340 us
RF BWT product	8
Ernst T1	1200 ms
PATRef prep. shots	10
Volume dummy shots	0
Dummy Measurements	0
ETL per RTEB	1
Invert PE	Off
Echo Time Shift	On
Ramp Sampling	On
NORDIC	Off
SVDPC	On
Sym VASO	Off
Dual-pol. EPI	Off
Invert RO	Off
Invert 3D	Off
Disable PF reco	Off
Save sampling	Off
PE VComp	Off
Water Exc.	Binomial-11
External PC	per Series
FIDNavs	-none-
EPI rise time factor	1.10
Mosaic DICOMs	On
Modify Ice Config	On
G-factor map	Off
GRAPPA Regularization	50000 10^-6
HSN RF power scale	3.00
Inversion Delay	0 ms
Relaxation Delay	0 ms
Var. FA /MAGEC	0

Sequence - Assistant

Mode	Off
------	-----