

\\USER\FMRIF\[XT-ID:93-M-0170]Renzo\20230404\_DPG\_T1slab\_axial\cstfl\_wip925b\_10slslab

TA: 7:22 PM: REF Voxel size: 0.4×0.4×0.4 mmPAT: 4 Rel. SNR: 1.00 : WIP\_cmp

**Properties**

|   |                    |
|---|--------------------|
| Prio recon                                    | Off                |
| Load images to viewer                         | On                 |
| Inline movie                                  | Off                |
| Auto store images                             | On                 |
| Load images to stamp segments                 | Off                |
| Load images to graphic segments               | On                 |
| Auto open inline display                      | Off                |
| Auto close inline display                     | Off                |
| Start measurement without further preparation | Off                |
| Wait for user to start                        | Off                |
| Start measurements                            | Single measurement |

**Routine**

|                    |                      |
|--------------------|----------------------|
| Slab group         | 1                    |
| Slabs              | 1                    |
| Dist. factor       | 50 %                 |
| Position           | R17.9 A6.7 F34.1 mm  |
| Orientation        | Sagittal             |
| Phase enc. dir.    | A >> P               |
| AutoAlign          | ---                  |
| Phase oversampling | 0 %                  |
| Slice oversampling | 0.0 %                |
| Slices per slab    | 88                   |
| FoV read           | 230 mm               |
| FoV phase          | 94.4 %               |
| Slice thickness    | 0.36 mm              |
| TR                 | 5000.0 ms            |
| TE                 | 4.4 ms               |
| Averages           | 1                    |
| Concatenations     | 1                    |
| Filter             | Distortion Corr.(2D) |
| Coil elements      | A32                  |

**Contrast - Common**

|                   |             |
|-------------------|-------------|
| TR                | 5000.0 ms   |
| TE                | 4.4 ms      |
| Magn. preparation | Non-sel. IR |
| T1 1              | 969 ms      |
| T1 2              | 2910 ms     |
| Flip angle 1      | 6 deg       |
| Flip angle 2      | 7 deg       |
| Fat suppr.        | None        |
| Water suppr.      | None        |

**Contrast - Dynamic**

|                 |                  |
|-----------------|------------------|
| Averages        | 1                |
| Averaging mode  | Long term        |
| Reconstruction  | Magnitude        |
| Measurements    | 1                |
| Multiple series | Each measurement |

**Resolution - Common**

|                       |         |
|-----------------------|---------|
| FoV read              | 230 mm  |
| FoV phase             | 94.4 %  |
| Slice thickness       | 0.36 mm |
| Base resolution       | 640     |
| Phase resolution      | 100 %   |
| Slice resolution      | 100 %   |
| Phase partial Fourier | Off     |

**Resolution - Common**

|                       |     |
|-----------------------|-----|
| Slice partial Fourier | Off |
| Interpolation         | Off |

**Resolution - iPAT**

|                     |            |
|---------------------|------------|
| PAT mode            | GRAPPA     |
| Accel. factor PE    | 4          |
| Ref. lines PE       | 40         |
| Accel. factor 3D    | 1          |
| Reference scan mode | Integrated |

**Resolution - Filter Image**

|                   |     |
|-------------------|-----|
| Image Filter      | Off |
| Distortion Corr.  | On  |
| Mode              | 2D  |
| Unfiltered images | Off |
| Prescan Normalize | Off |
| Normalize         | Off |
| B1 filter         | Off |

**Resolution - Filter Rawdata**

|                   |     |
|-------------------|-----|
| Raw filter        | Off |
| Elliptical filter | Off |

**Geometry - Common**

|                    |                     |
|--------------------|---------------------|
| Slab group         | 1                   |
| Slabs              | 1                   |
| Dist. factor       | 50 %                |
| Position           | R17.9 A6.7 F34.1 mm |
| Orientation        | Sagittal            |
| Phase enc. dir.    | A >> P              |
| Slice oversampling | 0.0 %               |
| Slices per slab    | 88                  |
| FoV read           | 230 mm              |
| FoV phase          | 94.4 %              |
| Slice thickness    | 0.36 mm             |
| TR                 | 5000.0 ms           |
| Multi-slice mode   | Single shot         |
| Series             | Interleaved         |
| Concatenations     | 1                   |

**Geometry - AutoAlign**

|                     |                     |
|---------------------|---------------------|
| Slab group          | 1                   |
| Position            | R17.9 A6.7 F34.1 mm |
| Orientation         | Sagittal            |
| Phase enc. dir.     | A >> P              |
| AutoAlign           | ---                 |
| Initial Position    | R17.9 A6.7 F34.1    |
| R                   | 17.9 mm             |
| A                   | 6.7 mm              |
| F                   | 34.1 mm             |
| Initial Rotation    | 0.00 deg            |
| Initial Orientation | Sagittal            |

**Geometry - Navigator****System - Miscellaneous**

|                  |           |
|------------------|-----------|
| Positioning mode | REF       |
| Table position   | H         |
| Table position   | 0 mm      |
| MSMA             | S - C - T |

**System - Miscellaneous**

|                     |                |
|---------------------|----------------|
| Sagittal            | R >> L         |
| Coronal             | A >> P         |
| Transversal         | F >> H         |
| Coil Combine Mode   | Sum of Squares |
| Save uncombined     | Off            |
| Matrix Optimization | Off            |
| AutoAlign           | ---            |
| Coil Select Mode    | Default        |

**System - Adjustments**

|                          |          |
|--------------------------|----------|
| B0 Shim mode             | Tune up  |
| B1 Shim mode             | TrueForm |
| Confirm freq. adjustment | Off      |
| Assume Dominant Fat      | Off      |
| Assume Silicone          | Off      |
| Adjustment Tolerance     | Auto     |

**System - Adjust Volume**

|               |                     |
|---------------|---------------------|
| ! Position    | R20.1 A2.7 F15.5 mm |
| ! Orientation | Transversal         |
| ! Rotation    | 0.00 deg            |
| ! A >> P      | 186 mm              |
| ! R >> L      | 41 mm               |
| ! F >> H      | 203 mm              |
| Reset         | Off                 |

**System - Tx/Rx**

|                     |                |
|---------------------|----------------|
| Frequency 1H        | 297.145675 MHz |
| Correction factor   | 1              |
| Gain                | High           |
| Img. Scale Cor.     | 1.000          |
| Reset               | Off            |
| ? Ref. amplitude 1H | 0.000 V        |

**Physio - Signal1**

|                 |           |
|-----------------|-----------|
| 1st Signal/Mode | None      |
| TR              | 5000.0 ms |
| Concatenations  | 1         |

**Physio - Cardiac**

|                   |             |
|-------------------|-------------|
| Magn. preparation | Non-sel. IR |
| TI 1              | 969 ms      |
| TI 2              | 2910 ms     |
| Fat suppr.        | None        |
| Dark blood        | Off         |
| FoV read          | 230 mm      |
| FoV phase         | 94.4 %      |
| Phase resolution  | 100 %       |

**Physio - PACE**

|                |     |
|----------------|-----|
| Resp. control  | Off |
| Concatenations | 1   |

**Inline - Common**

|                      |     |
|----------------------|-----|
| Subtract             | Off |
| Measurements         | 1   |
| StdDev               | Off |
| Save original images | On  |

**Inline - MIP**

|         |     |
|---------|-----|
| MIP-Sag | Off |
| MIP-Cor | Off |
| MIP-Tra | Off |

**Inline - MIP**

|                      |     |
|----------------------|-----|
| MIP-Time             | Off |
| Save original images | On  |

**Inline - Composing**

|                   |     |
|-------------------|-----|
| Distortion Corr.  | On  |
| Mode              | 2D  |
| Unfiltered images | Off |

**Sequence - Part 1**

|                     |             |
|---------------------|-------------|
| Introduction        | On          |
| Dimension           | 3D          |
| Elliptical scanning | Off         |
| Reordering          | Linear rot. |
| Asymmetric echo     | Allowed     |
| Contrasts           | 1           |
| Flow comp.          | No          |
| Multi-slice mode    | Single shot |
| Echo spacing        | 10.5 ms     |
| Bandwidth           | 140 Hz/Px   |

**Sequence - Part 2**

|                         |           |
|-------------------------|-----------|
| RF pulse type           | Normal    |
| Gradient mode           | Normal    |
| Excitation              | Slab-sel. |
| RF spoiling             | On        |
| Incr. Gradient spoiling | Off       |
| Turbo factor            | 603       |

**Sequence - Special**

|                 |       |
|-----------------|-------|
| Sparse Sampling | Off   |
| Shift Inv Pulse | 0 Hz  |
| Uniform         | On    |
| Denoised UNI    | On    |
| FLAWS           | Off   |
| FLAWS-hc        | Off   |
| FLAWS-hc inv.   | Off   |
| Division image  | Off   |
| T1 Map          | On    |
| Synthetic TI 0  | 0 ms  |
| Synthetic TI 1  | 0 ms  |
| Denoise Lambda  | 1     |
| Scaling         | 0 10^ |
| Echo Averaging  | Off   |
| FID Monitoring  | Off   |

**Sequence - Assistant**

|      |     |
|------|-----|
| Mode | Off |
|------|-----|