

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\RenzoHuber\Lonike\20210713_tests\Checklist_ok

TA: 1:11 PAT: Off Voxel size: 1.0x1.0x5.0 mm Rel. SNR: 1.00 SIEMENS: tfl

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	Off
Wait for user to start	Off
Start measurements	single

Routine

Slice group 1	
Slices	11
Dist. factor	150 %
Position	R1.8 A19.9 F14.3
Orientation	Sagittal
Phase enc. dir.	A >> P
Rotation	0.00 deg
Slice group 2	
Slices	5
Dist. factor	100 %
Position	R2.3 A15.9 H1.4
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	0.00 deg
Slice group 3	
Slices	7
Dist. factor	200 %
Position	R1.5 A21.4 F6.6
Orientation	Coronal
Phase enc. dir.	R >> L
Rotation	0.00 deg
Phase oversampling	0 %
FoV read	200 mm
FoV phase	100.0 %
Slice thickness	5.0 mm
TR	3000 ms
TE	2.24 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	A32

Contrast

Magn. preparation	Slice-sel. IR
TI	1100 ms
Flip angle	6 deg
Fat suppr.	None
Water suppr.	None
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	1
Multiple series	Each measurement

Resolution

Base resolution	192
Phase resolution	100 %
Phase partial Fourier	Off
Interpolation	Off

PAT mode None

Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off
Raw filter	Off
Elliptical filter	Off

Geometry

Multi-slice mode	Single shot
Series	Interleaved

Table position	H
Table position	0 mm
Inline Composing	Off

System

V32	Off
A32	On
Positioning mode	REF
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Adaptive Combine
AutoAlign	---
Auto Coil Select	Default
Shim mode	Tune up
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
! Ref. amplitude 1H	220.000 V
Adjustment Tolerance	Auto
Adjust volume	
! Position	L0.0 A23.4 F1.3
! Orientation	Transversal
! Rotation	0.00 deg
! R >> L	350 mm
! A >> P	213 mm
! F >> H	189 mm

Physio

1st Signal/Mode	None
Dark blood	Off
Resp. control	Off

Inline

Subtract	Off
Std-Dev-Sag	Off
Std-Dev-Cor	Off
Std-Dev-Tra	Off
Std-Dev-Time	Off
MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On

Sequence

Introduction	On
--------------	----

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

Dimension	2D
Asymmetric echo	Allowed
Bandwidth	240 Hz/Px
Flow comp.	No
Echo spacing	5.5 ms
<hr/>	
RF pulse type	Normal
Gradient mode	Fast
Excitation	Slice-sel.
RF spoiling	On

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\RenzoHuber\Lonike\20210713_tests\VASO_153_2D

TA: 0:48

PAT: 3

Voxel size: 0.9x0.9x0.9 mm

Rel. SNR: 1.00

USER: VASO_153

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Prescan Normalize	Off
Raw filter	Off
Elliptical filter	Off
Hamming	Off

Geometry

Multi-slice mode	Interleaved
Series	Ascending
Special sat.	Parallel F
Gap	25.0 mm
Thickness	100 mm
Table position	H
Table position	0 mm
Inline Composing	Off

Routine

Slice group 1	
Slices	12
Dist. factor	0 %
Position	L0.5 A21.9 H0.0
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
FoV read	160.0 mm
FoV phase	100.0 %
Slice thickness	0.9 mm
TR	1833.5 ms
TE	21 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	A32

System

V32	Off
A32	On
Positioning mode	REF
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
! Ref. amplitude 1H	220.000 V
Adjustment Tolerance	Auto
Adjust volume	
! Position	Isocenter
! Orientation	Transversal
! Rotation	90.00 deg
! A >> P	500 mm
! R >> L	200 mm
! F >> H	16 mm

Contrast

Perf / VASO mode	SS-SI VASO
TI2	1075 ms
TI1	50 ms
TI1s	50 ms
Flip angle	70 deg
Fat suppr.	Fat sat.
Fat sat. mode	Strong
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	23
Delay in TR	0 ms
Multiple series	Off
Perfusion mode	PICORE Q2T
Inversion time 1	50 ms
Saturation stop time	50 ms
Inversion time 2	1075 ms
Flow limit	100.0 cm/s

Physio

1st Signal/Mode	None
-----------------	------

BOLD

Sequence

Introduction	On
Contrasts	1
Bandwidth	1124 Hz/Px
Free echo spacing	Off
Echo spacing	1.01 ms
EPI factor	178
RF pulse type	Normal
Gradient mode	Normal
Ampl	100
BWDTH	300 3.1kHz
thickness	1000
Phase skip	0
Opt. TI2	1235
Volumes per TI	1
FatSat flip angle	110 deg
SMS factor	1

Resolution

Base resolution	178
Phase resolution	100 %
Phase partial Fourier	6/8
Interpolation	Off
PAT mode	GRAPPA
Accel. factor PE	3
Ref. lines PE	27
Reference scan mode	Separate

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

log physio files	Off
altern z-shim	0 uT/m
fixed z-shim	0 uT/m
EPI phase correction	local
PAT refscan mode	single-shot
RF pulse duration	2560 us
FFT scale	1.0

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\RenzoHuber\Lonike\20210713_tests\VASO_157_3D

TA: 6.4 s PAT: 3 Voxel size: 0.9x0.9x0.9 mm Rel. SNR: 1.00 USER: VASO_157

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slab group 1	
Slabs	1
Dist. factor	50 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	12
FoV read	120.0 mm
FoV phase	134.3 %
Slice thickness	0.90 mm
TR	1596.00 ms
TE	22 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	A32

Contrast

Perfusion mode	SS-SI VASO
TI2	900 ms
TI1	50 ms
TI1s	50 ms
Flip angle	4 deg
Fat suppr.	Fat sat.
Fat sat. mode	Weak
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	4
Delay in TR	0 ms
Multiple series	Off
Perfusion mode	PICORE Q2T
Inversion time 1	50 ms
Saturation stop time	50 ms
Inversion time 2	900.0 ms
Flow limit	100.0 cm/s

Resolution

Base resolution	134
Phase resolution	100 %
Slice resolution	100 %
Phase partial Fourier	6/8
Slice partial Fourier	Off
Interpolation	Off

PAT mode	GRAPPA
Accel. factor PE	3
Ref. lines PE	33
Accel. factor 3D	1
Ref. lines 3D	12
Reference scan mode	Separate
Prescan Normalize	Off
Raw filter	Off
Elliptical filter	Off
Hamming	Off

Geometry

Multi-slice mode	Interleaved
Series	Ascending
Special sat.	Parallel F
Gap	25.0 mm
Thickness	100 mm
Table position	H
Table position	0 mm
Inline Composing	Off

System

V32	Off
A32	On
Positioning mode	REF
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
! Ref. amplitude 1H	220.000 V
Adjustment Tolerance	Auto
Adjust volume	
! Position	Isocenter
! Orientation	Transversal
! Rotation	90.00 deg
! A >> P	200 mm
! R >> L	150 mm
! F >> H	11 mm

Physio

1st Signal/Mode	None
-----------------	------

BOLD

Motion correction	Off
Spatial filter	Off

Sequence

Introduction	On
Dimension	3D
Reordering	Linear
Contrasts	1
Bandwidth	1098 Hz/Px
Free echo spacing	Off
Echo spacing	1.02 ms
EPI factor	180

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

RF pulse type	Normal
Gradient mode	Fast
Excitation	Slab-sel.
RF spoiling	On
<hr/>	
Ampl	120
BWDTH	150 3.1kHz
ph.skip 4 Robert (the one)	1
use Ernst angle	Off
NORDIC	On
log physio files	Off
FFT scale	1.00
dummy prepscan time	3 s
z shim	0.00 mT/m*ms
RF duration	2560 us
RF BWTP	25.0
Renzo: Delta TI	57 ms
EFFECTIVE TR	19152 ms
PatPartitions	12
EPI phase correction	local
PAT refscan mode	Flash
FlashRef BaseRes	134
FlashRef BW	100 Hz/px
FlashRef TE	6500 us
FlashRef FA	5 deg
use CAIPI	Off

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\RenzoHuber\Lonike\20210713_tests\Checklist_ok

TA: 1:11 PAT: Off Voxel size: 1.0x1.0x5.0 mm Rel. SNR: 1.00 SIEMENS: tfl

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	Off
Wait for user to start	Off
Start measurements	single

Routine

Slice group 1	
Slices	11
Dist. factor	150 %
Position	R1.8 A19.9 F14.3
Orientation	Sagittal
Phase enc. dir.	A >> P
Rotation	0.00 deg
Slice group 2	
Slices	5
Dist. factor	100 %
Position	R2.3 A15.9 H1.4
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	0.00 deg
Slice group 3	
Slices	7
Dist. factor	200 %
Position	R1.5 A21.4 F6.6
Orientation	Coronal
Phase enc. dir.	R >> L
Rotation	0.00 deg
Phase oversampling	0 %
FoV read	200 mm
FoV phase	100.0 %
Slice thickness	5.0 mm
TR	3000 ms
TE	2.24 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	A32

Contrast

Magn. preparation	Slice-sel. IR
TI	1100 ms
Flip angle	6 deg
Fat suppr.	None
Water suppr.	None
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	1
Multiple series	Each measurement

Resolution

Base resolution	192
Phase resolution	100 %
Phase partial Fourier	Off
Interpolation	Off

PAT mode None

Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off
Raw filter	Off
Elliptical filter	Off

Geometry

Multi-slice mode	Single shot
Series	Interleaved

Table position	H
Table position	0 mm
Inline Composing	Off

System

V32	Off
A32	On
Positioning mode	REF
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Adaptive Combine
AutoAlign	---
Auto Coil Select	Default
Shim mode	Tune up
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
! Ref. amplitude 1H	220.000 V
Adjustment Tolerance	Auto
Adjust volume	
! Position	L0.0 A23.4 F1.3
! Orientation	Transversal
! Rotation	0.00 deg
! R >> L	350 mm
! A >> P	213 mm
! F >> H	189 mm

Physio

1st Signal/Mode	None
Dark blood	Off
Resp. control	Off

Inline

Subtract	Off
Std-Dev-Sag	Off
Std-Dev-Cor	Off
Std-Dev-Tra	Off
Std-Dev-Time	Off
MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On

Sequence

Introduction	On
--------------	----

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

Dimension	2D
Asymmetric echo	Allowed
Bandwidth	240 Hz/Px
Flow comp.	No
Echo spacing	5.5 ms
<hr/>	
RF pulse type	Normal
Gradient mode	Fast
Excitation	Slice-sel.
RF spoiling	On

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\RenzoHuber\Lonike\20210713_tests\VASO_153_2D_NOR_SEG_LOW

TA: 0:15 PAT: 3 Voxel size: 0.9x0.9x0.9 mm Rel. SNR: 1.00 USER: VASO_153

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Prescan Normalize	Off
Raw filter	Off
Elliptical filter	Off
Hamming	Off

Geometry

Multi-slice mode	Interleaved
Series	Ascending
Special sat.	Parallel F
Gap	25.0 mm
Thickness	100 mm
Table position	H
Table position	0 mm
Inline Composing	Off

Routine

Slice group 1	
Slices	12
Dist. factor	0 %
Position	L0.5 A33.9 F17.9
Orientation	T > C6.2
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
FoV read	160.0 mm
FoV phase	100.0 %
Slice thickness	0.9 mm
TR	1833.5 ms
TE	21 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	A32

System

V32	Off
A32	On
Positioning mode	REF
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
! Ref. amplitude 1H	220.000 V
Adjustment Tolerance	Auto
Adjust volume	
! Position	R3.3 A29.8 F12.7
! Orientation	T > C5.2
! Rotation	90.00 deg
! A >> P	225 mm
! R >> L	200 mm
! F >> H	29 mm

Contrast

Perf / VASO mode	SS-SI VASO
TI2	1075 ms
TI1	50 ms
TI1s	50 ms
Flip angle	70 deg
Fat suppr.	Fat sat.
Fat sat. mode	Strong
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	5
Delay in TR	0 ms
Multiple series	Off
Perfusion mode	PICORE Q2T
Inversion time 1	50 ms
Saturation stop time	50 ms
Inversion time 2	1075 ms
Flow limit	100.0 cm/s

Physio

1st Signal/Mode	None
-----------------	------

BOLD

Sequence

Introduction	On
Contrasts	1
Bandwidth	1124 Hz/Px
Free echo spacing	Off
Echo spacing	1.01 ms
EPI factor	178
RF pulse type	Normal
Gradient mode	Normal
Ampl	100
BWDTH	300 3.1kHz
thickness	1000
Phase skip	0
Opt. TI2	1235
Volumes per TI	1
FatSat flip angle	110 deg
SMS factor	1

Resolution

Base resolution	178
Phase resolution	100 %
Phase partial Fourier	6/8
Interpolation	Off
PAT mode	GRAPPA
Accel. factor PE	3
Ref. lines PE	27
Reference scan mode	Separate

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

log physio files	Off
altern z-shim	0 uT/m
fixed z-shim	0 uT/m
EPI phase correction	local
PAT refscan mode	segmented
RF pulse duration	2560 us
FFT scale	1.0

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\RenzoHuber\Lonike\20210713_tests\VASO_153_2D_NOR_SS_LOW

TA: 3:09 PAT: 3 Voxel size: 0.9x0.9x0.9 mm Rel. SNR: 1.00 USER: VASO_153

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Prescan Normalize	Off
Raw filter	Off
Elliptical filter	Off
Hamming	Off

Geometry

Multi-slice mode	Interleaved
Series	Ascending
Special sat.	Parallel F
Gap	25.0 mm
Thickness	100 mm
Table position	H
Table position	0 mm
Inline Composing	Off

Routine

Slice group 1	
Slices	12
Dist. factor	0 %
Position	L0.5 A33.9 F17.9
Orientation	T > C6.2
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
FoV read	160.0 mm
FoV phase	100.0 %
Slice thickness	0.9 mm
TR	1833.5 ms
TE	21 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	A32

System

V32	Off
A32	On
Positioning mode	REF
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
! Ref. amplitude 1H	220.000 V
Adjustment Tolerance	Auto
Adjust volume	
! Position	R3.3 A29.8 F12.7
! Orientation	T > C5.2
! Rotation	90.00 deg
! A >> P	225 mm
! R >> L	200 mm
! F >> H	29 mm

Contrast

Perf / VASO mode	SS-SI VASO
TI2	1075 ms
TI1	50 ms
TI1s	50 ms
Flip angle	70 deg
Fat suppr.	Fat sat.
Fat sat. mode	Strong

Averaging mode	Long term
Reconstruction	Magnitude
Measurements	100
Delay in TR	0 ms
Multiple series	Off

Perfusion mode	PICORE Q2T
Inversion time 1	50 ms
Saturation stop time	50 ms
Inversion time 2	1075 ms
Flow limit	100.0 cm/s

Resolution

Base resolution	178
Phase resolution	100 %
Phase partial Fourier	6/8
Interpolation	Off
PAT mode	GRAPPA
Accel. factor PE	3
Ref. lines PE	27
Reference scan mode	Separate

Physio

1st Signal/Mode	None
-----------------	------

BOLD

Sequence

Introduction	On
Contrasts	1
Bandwidth	1124 Hz/Px
Free echo spacing	Off
Echo spacing	1.01 ms
EPI factor	178
RF pulse type	Normal
Gradient mode	Normal
Ampl	100
BWDTH	300 3.1kHz
thickness	1000
Phase skip	0
Opt. TI2	1235
Volumes per TI	1
FatSat flip angle	110 deg
SMS factor	1

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

log physio files	Off
altern z-shim	0 uT/m
fixed z-shim	0 uT/m
EPI phase correction	local
PAT refscan mode	single-shot
RF pulse duration	2560 us
FFT scale	1.0

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\RenzoHuber\Lonike\20210713_tests\VASO_153_2D_NOR_FLEET_LOW

TA: 3:10 PAT: 3 Voxel size: 0.9x0.9x0.9 mm Rel. SNR: 1.00 USER: VASO_153

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Prescan Normalize	Off
Raw filter	Off
Elliptical filter	Off
Hamming	Off

Geometry

Multi-slice mode	Interleaved
Series	Ascending
Special sat.	Parallel F
Gap	25.0 mm
Thickness	100 mm
Table position	H
Table position	0 mm
Inline Composing	Off

Routine

Slice group 1	
Slices	12
Dist. factor	0 %
Position	L0.5 A33.9 F17.9
Orientation	T > C6.2
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
FoV read	160.0 mm
FoV phase	100.0 %
Slice thickness	0.9 mm
TR	1833.5 ms
TE	21 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	A32

System

V32	Off
A32	On
Positioning mode	REF
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
! Ref. amplitude 1H	220.000 V
Adjustment Tolerance	Auto
Adjust volume	
! Position	R3.3 A29.8 F12.7
! Orientation	T > C5.2
! Rotation	90.00 deg
! A >> P	225 mm
! R >> L	200 mm
! F >> H	29 mm

Contrast

Perf / VASO mode	SS-SI VASO
TI2	1075 ms
TI1	50 ms
TI1s	50 ms
Flip angle	70 deg
Fat suppr.	Fat sat.
Fat sat. mode	Strong

Averaging mode	Long term
Reconstruction	Magnitude
Measurements	100
Delay in TR	0 ms
Multiple series	Off

Perfusion mode	PICORE Q2T
Inversion time 1	50 ms
Saturation stop time	50 ms
Inversion time 2	1075 ms
Flow limit	100.0 cm/s

Resolution

Base resolution	178
Phase resolution	100 %
Phase partial Fourier	6/8
Interpolation	Off
PAT mode	GRAPPA
Accel. factor PE	3
Ref. lines PE	27
Reference scan mode	Separate

Physio

1st Signal/Mode	None
-----------------	------

BOLD

Sequence

Introduction	On
Contrasts	1
Bandwidth	1124 Hz/Px
Free echo spacing	Off
Echo spacing	1.01 ms
EPI factor	178
RF pulse type	Normal
Gradient mode	Normal
Ampl	120
BWDTH	300 3.1kHz
thickness	1000
Phase skip	0
Opt. TI2	1235
Volumes per TI	1
FatSat flip angle	110 deg
SMS factor	1

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

log physio files	Off
altern z-shim	0 uT/m
fixed z-shim	0 uT/m
EPI phase correction	local
PAT refscan mode	FLEET
FLEET dummies	15
FLEET flip angle	30
RF pulse duration	2560 us
FFT scale	1.0

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\RenzoHuber\Lonike\20210713_tests\VASO_153_2D_NOR_FLEET60_LOW

TA: 3:10 PAT: 3 Voxel size: 0.9x0.9x0.9 mm Rel. SNR: 1.00 USER: VASO_153

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Prescan Normalize	Off
Raw filter	Off
Elliptical filter	Off
Hamming	Off

Geometry

Multi-slice mode	Interleaved
Series	Ascending
Special sat.	Parallel F
Gap	25.0 mm
Thickness	100 mm
Table position	H
Table position	0 mm
Inline Composing	Off

Routine

Slice group 1	
Slices	12
Dist. factor	0 %
Position	L0.5 A33.9 F17.9
Orientation	T > C6.2
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
FoV read	160.0 mm
FoV phase	100.0 %
Slice thickness	0.9 mm
TR	1833.5 ms
TE	21 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	A32

System

V32	Off
A32	On
Positioning mode	REF
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
! Ref. amplitude 1H	220.000 V
Adjustment Tolerance	Auto
Adjust volume	
! Position	R3.3 A29.8 F12.7
! Orientation	T > C5.2
! Rotation	90.00 deg
! A >> P	225 mm
! R >> L	200 mm
! F >> H	29 mm

Contrast

Perf / VASO mode	SS-SI VASO
TI2	1075 ms
TI1	50 ms
TI1s	50 ms
Flip angle	70 deg
Fat suppr.	Fat sat.
Fat sat. mode	Strong
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	100
Delay in TR	0 ms
Multiple series	Off
Perfusion mode	PICORE Q2T
Inversion time 1	50 ms
Saturation stop time	50 ms
Inversion time 2	1075 ms
Flow limit	100.0 cm/s

Physio

1st Signal/Mode	None
-----------------	------

BOLD

Sequence

Introduction	On
Contrasts	1
Bandwidth	1124 Hz/Px
Free echo spacing	Off
Echo spacing	1.01 ms
EPI factor	178
RF pulse type	Normal
Gradient mode	Normal
Ampl	120
BWDTH	300 3.1kHz
thickness	1000
Phase skip	0
Opt. TI2	1235
Volumes per TI	1
FatSat flip angle	110 deg
SMS factor	1

Resolution

Base resolution	178
Phase resolution	100 %
Phase partial Fourier	6/8
Interpolation	Off
PAT mode	GRAPPA
Accel. factor PE	3
Ref. lines PE	27
Reference scan mode	Separate

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

log physio files	Off
altern z-shim	0 uT/m
fixed z-shim	0 uT/m
EPI phase correction	local
PAT refscan mode	FLEET
FLEET dummies	15
FLEET flip angle	60
RF pulse duration	2560 us
FFT scale	1.0

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\RenzoHuber\Lonike\20210713_tests\VASO_153_2D_NOR_FLEET30+_noreg_LOW

TA: 3:10 PAT: 3 Voxel size: 0.9x0.9x0.9 mm Rel. SNR: 1.00 USER: VASO_153

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Prescan Normalize	Off
Raw filter	Off
Elliptical filter	Off
Hamming	Off

Geometry

Multi-slice mode	Interleaved
Series	Ascending
Special sat.	Parallel F
Gap	25.0 mm
Thickness	100 mm
Table position	H
Table position	0 mm
Inline Composing	Off

Routine

Slice group 1	
Slices	12
Dist. factor	0 %
Position	L0.5 A33.9 F17.9
Orientation	T > C6.2
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
FoV read	160.0 mm
FoV phase	100.0 %
Slice thickness	0.9 mm
TR	1833.5 ms
TE	21 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	A32

System

V32	Off
A32	On
Positioning mode	REF
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
! Ref. amplitude 1H	220.000 V
Adjustment Tolerance	Auto
Adjust volume	
! Position	R3.3 A29.8 F12.7
! Orientation	T > C5.2
! Rotation	90.00 deg
! A >> P	225 mm
! R >> L	200 mm
! F >> H	29 mm

Contrast

Perf / VASO mode	SS-SI VASO
TI2	1075 ms
TI1	50 ms
TI1s	50 ms
Flip angle	70 deg
Fat suppr.	Fat sat.
Fat sat. mode	Strong

Averaging mode	Long term
Reconstruction	Magnitude
Measurements	100
Delay in TR	0 ms
Multiple series	Off

Perfusion mode	PICORE Q2T
Inversion time 1	50 ms
Saturation stop time	50 ms
Inversion time 2	1075 ms
Flow limit	100.0 cm/s

Resolution

Base resolution	178
Phase resolution	100 %
Phase partial Fourier	6/8
Interpolation	Off
PAT mode	GRAPPA
Accel. factor PE	3
Ref. lines PE	27
Reference scan mode	Separate

Physio

1st Signal/Mode	None
-----------------	------

BOLD

Sequence

Introduction	On
Contrasts	1
Bandwidth	1124 Hz/Px
Free echo spacing	Off
Echo spacing	1.01 ms
EPI factor	178
RF pulse type	Normal
Gradient mode	Normal
Ampl	120
BWDTH	300 3.1kHz
thickness	1000
Phase skip	0
Opt. TI2	1235
Volumes per TI	1
FatSat flip angle	110 deg
SMS factor	1

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

log physio files	Off
altern z-shim	0 uT/m
fixed z-shim	0 uT/m
EPI phase correction	local
PAT refscan mode	FLEET
FLEET dummies	15
FLEET flip angle	30
RF pulse duration	2560 us
FFT scale	1.0

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\RenzoHuber\Lonike\20210713_tests\VASO_153_2D_NOR_SEG_LOW_NOISE

TA: 9.5 s PAT: 3 Voxel size: 0.9x0.9x0.9 mm Rel. SNR: 1.00 USER: VASO_153

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slice group 1	
Slices	12
Dist. factor	0 %
Position	L0.5 A33.9 F17.9
Orientation	T > C6.2
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
FoV read	160.0 mm
FoV phase	100.0 %
Slice thickness	0.9 mm
TR	1833.5 ms
TE	21 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	A32

Contrast

Perf / VASO mode	SS-SI VASO
TI2	1075 ms
TI1	50 ms
TI1s	50 ms
Flip angle	70 deg
Fat suppr.	Fat sat.
Fat sat. mode	Strong

Averaging mode	Long term
Reconstruction	Magnitude
Measurements	2
Delay in TR	0 ms
Multiple series	Off

Perfusion mode	PICORE Q2T
Inversion time 1	50 ms
Saturation stop time	50 ms
Inversion time 2	1075 ms
Flow limit	100.0 cm/s

Resolution

Base resolution	178
Phase resolution	100 %
Phase partial Fourier	6/8
Interpolation	Off
PAT mode	GRAPPA
Accel. factor PE	3
Ref. lines PE	27
Reference scan mode	Separate

Prescan Normalize	Off
Raw filter	Off
Elliptical filter	Off
Hamming	Off

Geometry

Multi-slice mode	Interleaved
Series	Ascending
Special sat.	Parallel F
Gap	25.0 mm
Thickness	100 mm

Table position	H
Table position	0 mm
Inline Composing	Off

System

V32	Off
A32	On
Positioning mode	REF
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Default

Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
! Ref. amplitude 1H	220.000 V
Adjustment Tolerance	Auto
Adjust volume	
! Position	R3.3 A29.8 F12.7
! Orientation	T > C5.2
! Rotation	90.00 deg
! A >> P	225 mm
! R >> L	200 mm
! F >> H	29 mm

Physio

1st Signal/Mode	None
-----------------	------

BOLD

Sequence

Introduction	On
Contrasts	1
Bandwidth	1124 Hz/Px
Free echo spacing	Off
Echo spacing	1.01 ms

EPI factor	178
RF pulse type	Normal
Gradient mode	Normal

Ampl	120
BWDTH	300 3.1kHz
thickness	1000
Phase skip	0
Opt. TI2	1235
Volumes per TI	1
FatSat flip angle	110 deg
SMS factor	1

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

log physio files	Off
altern z-shim	0 uT/m
fixed z-shim	0 uT/m
EPI phase correction	local
PAT refscan mode	segmented
RF pulse duration	2560 us
FFT scale	1.0

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\RenzoHuber\Lonike\20210713_tests\VASO_157_3D_reg5000_FA26

TA: 2:47 PAT: 3 Voxel size: 0.9x0.9x0.9 mm Rel. SNR: 1.00 USER: VASO_157

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slab group 1	
Slabs	1
Dist. factor	50 %
Position	L0.5 A33.9 F17.9
Orientation	T > C6.2
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	12
FoV read	130.0 mm
FoV phase	134.7 %
Slice thickness	0.90 mm
TR	1672.70 ms
TE	24 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	A32

Contrast

Perfusion mode	SS-SI VASO
TI2	900 ms
TI1	50 ms
TI1s	50 ms
Flip angle	26 deg
Fat suppr.	Fat sat.
Fat sat. mode	Weak
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	100
Delay in TR	0 ms
Multiple series	Off
Perfusion mode	PICORE Q2T
Inversion time 1	50 ms
Saturation stop time	50 ms
Inversion time 2	900.0 ms
Flow limit	100.0 cm/s

Resolution

Base resolution	144
Phase resolution	100 %
Slice resolution	100 %
Phase partial Fourier	6/8
Slice partial Fourier	Off
Interpolation	Off

PAT mode	GRAPPA
Accel. factor PE	3
Ref. lines PE	33
Accel. factor 3D	1
Ref. lines 3D	12
Reference scan mode	Separate
Prescan Normalize	Off
Raw filter	Off
Elliptical filter	Off
Hamming	Off

Geometry

Multi-slice mode	Interleaved
Series	Ascending
Special sat.	Parallel F
Gap	25.0 mm
Thickness	100 mm
Table position	H
Table position	0 mm
Inline Composing	Off

System

V32	Off
A32	On
Positioning mode	FIX
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
! Ref. amplitude 1H	220.000 V
Adjustment Tolerance	Auto
Adjust volume	
! Position	R3.3 A29.8 F12.7
! Orientation	T > C5.2
! Rotation	90.00 deg
! A >> P	225 mm
! R >> L	200 mm
! F >> H	29 mm

Physio

1st Signal/Mode	None
-----------------	------

BOLD

Motion correction	Off
Spatial filter	Off

Sequence

Introduction	On
Dimension	3D
Reordering	Linear
Contrasts	1
Bandwidth	1052 Hz/Px
Free echo spacing	Off
Echo spacing	1.06 ms
EPI factor	194

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

RF pulse type	Normal
Gradient mode	Fast
Excitation	Slab-sel.
RF spoiling	On
<hr/>	
Ampl	120
BWDTH	150 3.1kHz
ph.skip 4 Robert (the one)	1
use Ernst angle	Off
NORDIC	On
log physio files	Off
FFT scale	1.00
dummy prepscan time	3 s
z shim	0.00 mT/m*ms
RF duration	2560 us
RF BWTP	25.0
Renzo: Delta TI	64 ms
EFFECTIVE TR	20072 ms
PatPartitions	12
EPI phase correction	local
PAT refscan mode	Flash
FlashRef BaseRes	144
FlashRef BW	100 Hz/px
FlashRef TE	6500 us
FlashRef FA	5 deg
use CAIPI	Off

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\RenzoHuber\Lonike\20210713_tests\VASO_157_3D_reg5000_FA4

TA: 2:47 PAT: 3 Voxel size: 0.9x0.9x0.9 mm Rel. SNR: 1.00 USER: VASO_157

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slab group 1	
Slabs	1
Dist. factor	50 %
Position	L0.5 A33.9 F17.9
Orientation	T > C6.2
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	12
FoV read	130.0 mm
FoV phase	134.7 %
Slice thickness	0.90 mm
TR	1672.70 ms
TE	24 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	A32

Contrast

Perfusion mode	SS-SI VASO
TI2	900 ms
TI1	50 ms
TI1s	50 ms
Flip angle	4 deg
Fat suppr.	Fat sat.
Fat sat. mode	Weak
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	100
Delay in TR	0 ms
Multiple series	Off
Perfusion mode	PICORE Q2T
Inversion time 1	50 ms
Saturation stop time	50 ms
Inversion time 2	900.0 ms
Flow limit	100.0 cm/s

Resolution

Base resolution	144
Phase resolution	100 %
Slice resolution	100 %
Phase partial Fourier	6/8
Slice partial Fourier	Off
Interpolation	Off

PAT mode	GRAPPA
Accel. factor PE	3
Ref. lines PE	33
Accel. factor 3D	1
Ref. lines 3D	12
Reference scan mode	Separate

Prescan Normalize	Off
Raw filter	Off
Elliptical filter	Off
Hamming	Off

Geometry

Multi-slice mode	Interleaved
Series	Ascending
Special sat.	Parallel F
Gap	25.0 mm
Thickness	100 mm
Table position	H
Table position	0 mm
Inline Composing	Off

System

V32	Off
A32	On
Positioning mode	FIX
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
! Ref. amplitude 1H	220.000 V
Adjustment Tolerance	Auto
Adjust volume	
! Position	R3.3 A29.8 F12.7
! Orientation	T > C5.2
! Rotation	90.00 deg
! A >> P	225 mm
! R >> L	200 mm
! F >> H	29 mm

Physio

1st Signal/Mode	None
-----------------	------

BOLD

Motion correction	Off
Spatial filter	Off

Sequence

Introduction	On
Dimension	3D
Reordering	Linear
Contrasts	1
Bandwidth	1052 Hz/Px
Free echo spacing	Off
Echo spacing	1.06 ms
EPI factor	194

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

RF pulse type	Normal
Gradient mode	Fast
Excitation	Slab-sel.
RF spoiling	On
<hr/>	
Ampl	120
BWDTH	150 3.1kHz
ph.skip 4 Robert (the one)	1
use Ernst angle	Off
NORDIC	On
log physio files	Off
FFT scale	1.00
dummy prepscan time	3 s
z shim	0.00 mT/m*ms
RF duration	2560 us
RF BWTP	25.0
Renzo: Delta TI	64 ms
EFFECTIVE TR	20072 ms
PatPartitions	12
EPI phase correction	local
PAT refscan mode	Flash
FlashRef BaseRes	144
FlashRef BW	100 Hz/px
FlashRef TE	6500 us
FlashRef FA	5 deg
use CAIPI	Off

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\RenzoHuber\Lonike\20210713_tests\VASO_157_3D_reg5000_FA4_run

TA: 15:05 PAT: 3 Voxel size: 0.9x0.9x0.9 mm Rel. SNR: 1.00 USER: VASO_157

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slab group 1	
Slabs	1
Dist. factor	50 %
Position	L0.5 A33.9 F17.9
Orientation	T > C6.2
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	12
FoV read	130.0 mm
FoV phase	134.7 %
Slice thickness	0.90 mm
TR	1672.70 ms
TE	24 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	A32

Contrast

Perfusion mode	SS-SI VASO
TI2	900 ms
TI1	50 ms
TI1s	50 ms
Flip angle	4 deg
Fat suppr.	Fat sat.
Fat sat. mode	Weak
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	541
Delay in TR	0 ms
Multiple series	Off
Perfusion mode	PICORE Q2T
Inversion time 1	50 ms
Saturation stop time	50 ms
Inversion time 2	900.0 ms
Flow limit	100.0 cm/s

Resolution

Base resolution	144
Phase resolution	100 %
Slice resolution	100 %
Phase partial Fourier	6/8
Slice partial Fourier	Off
Interpolation	Off

PAT mode	GRAPPA
Accel. factor PE	3
Ref. lines PE	33
Accel. factor 3D	1
Ref. lines 3D	12
Reference scan mode	Separate

Prescan Normalize	Off
Raw filter	Off
Elliptical filter	Off
Hamming	Off

Geometry

Multi-slice mode	Interleaved
Series	Ascending
Special sat.	Parallel F
Gap	25.0 mm
Thickness	100 mm
Table position	H
Table position	0 mm
Inline Composing	Off

System

V32	Off
A32	On
Positioning mode	FIX
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
! Ref. amplitude 1H	220.000 V
Adjustment Tolerance	Auto
Adjust volume	
! Position	R3.3 A29.8 F12.7
! Orientation	T > C5.2
! Rotation	90.00 deg
! A >> P	225 mm
! R >> L	200 mm
! F >> H	29 mm

Physio

1st Signal/Mode	None
-----------------	------

BOLD

Motion correction	Off
Spatial filter	Off

Sequence

Introduction	On
Dimension	3D
Reordering	Linear
Contrasts	1
Bandwidth	1052 Hz/Px
Free echo spacing	Off
Echo spacing	1.06 ms
EPI factor	194

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

RF pulse type	Normal
Gradient mode	Fast
Excitation	Slab-sel.
RF spoiling	On
<hr/>	
Ampl	120
BWDTH	150 3.1kHz
ph.skip 4 Robert (the one)	1
use Ernst angle	Off
NORDIC	On
log physio files	Off
FFT scale	1.00
dummy prepscan time	3 s
z shim	0.00 mT/m*ms
RF duration	2560 us
RF BWTP	25.0
Renzo: Delta TI	64 ms
EFFECTIVE TR	20072 ms
PatPartitions	12
EPI phase correction	local
PAT refscan mode	Flash
FlashRef BaseRes	144
FlashRef BW	100 Hz/px
FlashRef TE	6500 us
FlashRef FA	5 deg
use CAIPI	Off

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\RenzoHuber\Lonike\20210713_tests\VASO_153_2D_NOR_SEG_LOW_run

TA: 15:02 PAT: 3 Voxel size: 0.9x0.9x0.9 mm Rel. SNR: 1.00 USER: VASO_153

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Prescan Normalize	Off
Raw filter	Off
Elliptical filter	Off
Hamming	Off

Geometry

Multi-slice mode	Interleaved
Series	Ascending
Special sat.	Parallel F
Gap	25.0 mm
Thickness	100 mm
Table position	H
Table position	0 mm
Inline Composing	Off

Routine

Slice group 1	
Slices	12
Dist. factor	0 %
Position	L0.5 A33.9 F17.9
Orientation	T > C6.2
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
FoV read	160.0 mm
FoV phase	100.0 %
Slice thickness	0.9 mm
TR	1833.5 ms
TE	21 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	A32

System

V32	Off
A32	On
Positioning mode	REF
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
! Ref. amplitude 1H	220.000 V
Adjustment Tolerance	Auto
Adjust volume	
! Position	R3.3 A29.8 F12.7
! Orientation	T > C5.2
! Rotation	90.00 deg
! A >> P	225 mm
! R >> L	200 mm
! F >> H	29 mm

Contrast

Perf / VASO mode	SS-SI VASO
TI2	1075 ms
TI1	50 ms
TI1s	50 ms
Flip angle	70 deg
Fat suppr.	Fat sat.
Fat sat. mode	Strong

Averaging mode	Long term
Reconstruction	Magnitude
Measurements	489
Delay in TR	0 ms
Multiple series	Off

Perfusion mode	PICORE Q2T
Inversion time 1	50 ms
Saturation stop time	50 ms
Inversion time 2	1075 ms
Flow limit	100.0 cm/s

Resolution

Base resolution	178
Phase resolution	100 %
Phase partial Fourier	6/8
Interpolation	Off
PAT mode	GRAPPA
Accel. factor PE	3
Ref. lines PE	27
Reference scan mode	Separate

Physio

1st Signal/Mode	None
-----------------	------

BOLD

Sequence

Introduction	On
Contrasts	1
Bandwidth	1124 Hz/Px
Free echo spacing	Off
Echo spacing	1.01 ms
EPI factor	178
RF pulse type	Normal
Gradient mode	Normal
Ampl	120
BWDTH	300 3.1kHz
thickness	1000
Phase skip	0
Opt. TI2	1235
Volumes per TI	1
FatSat flip angle	110 deg
SMS factor	1

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

log physio files	Off
altern z-shim	0 uT/m
fixed z-shim	0 uT/m
EPI phase correction	local
PAT refscan mode	segmented
RF pulse duration	2560 us
FFT scale	1.0