

# SIEMENS MAGNETOM Investigational\_Device\_7T syngo MR B17

\\USER\RenzoHuber\KENSHU\20211011\_KEN\_movie2\Checklist\_ok

TA: 1:05 PAT: Off Voxel size: 1.0x1.0x5.0 mm Rel. SNR: 1.00 SIEMENS: tfl

## Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	Off
Wait for user to start	Off
Start measurements	single

## Routine

Slice group 1	
Slices	9
Dist. factor	80 %
Position	L0.0 A27.3 F14.3
Orientation	Sagittal
Phase enc. dir.	A >> P
Rotation	0.00 deg
Slice group 2	
Slices	5
Dist. factor	100 %
Position	R1.1 A28.9 H44.2
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	0.00 deg
Slice group 3	
Slices	7
Dist. factor	50 %
Position	R1.5 A28.2 F10.0
Orientation	Coronal
Phase enc. dir.	R >> L
Rotation	0.00 deg
Phase oversampling	0 %
FoV read	200 mm
FoV phase	100.0 %
Slice thickness	5.0 mm
TR	3000 ms
TE	2.24 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	A32

## Contrast

Magn. preparation	Slice-sel. IR
TI	1100 ms
Flip angle	6 deg
Fat suppr.	None
Water suppr.	None
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	1
Multiple series	Each measurement

## Resolution

Base resolution	192
Phase resolution	100 %
Phase partial Fourier	Off
Interpolation	Off

PAT mode	None
Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off
Raw filter	Off
Elliptical filter	Off

## Geometry

Multi-slice mode	Single shot
Series	Interleaved
Table position	H
Table position	0 mm
Inline Composing	Off

## System

V32	Off
A32	On
Positioning mode	REF
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Adaptive Combine
AutoAlign	---
Auto Coil Select	Default
Shim mode	Tune up
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
! Ref. amplitude 1H	220.000 V
Adjustment Tolerance	Auto
Adjust volume	
! Position	L0.0 A23.4 F1.3
! Orientation	Transversal
! Rotation	0.00 deg
! R >> L	350 mm
! A >> P	213 mm
! F >> H	189 mm

## Physio

1st Signal/Mode	None
Dark blood	Off
Resp. control	Off

## Inline

Subtract	Off
Std-Dev-Sag	Off
Std-Dev-Cor	Off
Std-Dev-Tra	Off
Std-Dev-Time	Off
MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On

## Sequence

Introduction	On
--------------	----

## SIEMENS MAGNETOM Investigational\_Device\_7T syngo MR B17

Dimension	2D
Asymmetric echo	Allowed
Bandwidth	240 Hz/Px
Flow comp.	No
Echo spacing	5.5 ms
<hr/>	
RF pulse type	Normal
Gradient mode	Fast
Excitation	Slice-sel.
RF spoiling	On

# SIEMENS MAGNETOM Investigational\_Device\_7T syngo MR B17

\\USER\RenzoHuber\KENSHU\20211011\_KEN\_movie2\VASO160\_G61\_nofatsat\_movie1

TA: 15:46

PAT: 6

Voxel size: 0.8x0.8x0.9 mm

Rel. SNR: 1.00

USER: VASO\_160

## Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

## Routine

Slab group 1	
Slabs	1
Dist. factor	50 %
Position	R3.9 A25.3 H16.4
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	112
FoV read	138.7 mm
FoV phase	132.9 %
Slice thickness	0.87 mm
TR	5168.90 ms
TE	18 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	A32

## Contrast

Perfusion mode	SS-SI VASO
TI2	100 ms
TI1	50 ms
TI1s	50 ms
Flip angle	40 deg
Fat suppr.	Fat sat.
Fat sat. mode	Weak
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	183
Delay in TR	0 ms
Multiple series	Off
Perfusion mode	PICORE Q2T
Inversion time 1	50 ms
Saturation stop time	50 ms
Inversion time 2	100.0 ms
Flow limit	100 cm/s

## Resolution

Base resolution	170
Phase resolution	100 %
Slice resolution	100 %
Phase partial Fourier	6/8
Slice partial Fourier	Off
Interpolation	Off

PAT mode	GRAPPA
Accel. factor PE	6
Ref. lines PE	90
Accel. factor 3D	1
Ref. lines 3D	24
Reference scan mode	Separate

Prescan Normalize	Off
Raw filter	Off
Elliptical filter	Off
Hamming	Off

## Geometry

Multi-slice mode	Interleaved
Series	Ascending
Special sat.	Parallel F
Gap	25.0 mm
Thickness	100 mm
Table position	H
Table position	0 mm
Inline Composing	Off

## System

V32	Off
A32	On
Positioning mode	FIX
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
! Ref. amplitude 1H	225.000 V
Adjustment Tolerance	Auto
Adjust volume	
! Position	R1.9 A24.9 H19.3
! Orientation	T > C-1.8
! Rotation	0.00 deg
! R >> L	137 mm
! A >> P	177 mm
! F >> H	94 mm

## Physio

1st Signal/Mode	None
-----------------	------

## BOLD

Motion correction	Off
Spatial filter	Off

## Sequence

Introduction	On
Dimension	3D
Reordering	Linear
Contrasts	1
Bandwidth	980 Hz/Px
Free echo spacing	Off
Echo spacing	1.16 ms
EPI factor	226

# SIEMENS MAGNETOM Investigational\_Device\_7T syngo MR B17

RF pulse type	Normal
Gradient mode	Fast
Excitation	Slab-sel.
RF spoiling	On
<hr/>	
Ampl	150
MAGEC FA	9 in deg
ph.skip 4 Robert (the one)	1
MAGEC SS-SI?	On
NORDIC	On
log physio files	On
FFT scale	1.00
dummy prepscan time	3 s
z shim	0.00 mT/m*ms
RF duration	2560 us
RF BWTP	8.0
Renzo: Delta TI	45 ms
EFFECTIVE TR	578916 ms
PatPartitions	112
CAIPI scale	3
EPI phase correction	local
PAT refscan mode	Flash
FlashRef BaseRes	170
FlashRef BW	103 Hz/px
FlashRef TE	6100 us
FlashRef FA	5 deg
use CAIPI	On
CAIPI shift kz	0
CAIPI shift ky	3

# SIEMENS MAGNETOM Investigational\_Device\_7T syngo MR B17

\\USER\RenzoHuber\KENSHU\20211011\_KEN\_movie2\VASO160\_G61\_nofatsat\_movie2

TA: 15:46 PAT: 6 Voxel size: 0.8x0.8x0.9 mm Rel. SNR: 1.00 USER: VASO\_160

## Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

## Routine

Slab group 1	
Slabs	1
Dist. factor	50 %
Position	R3.9 A25.3 H16.4
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	112
FoV read	138.7 mm
FoV phase	132.9 %
Slice thickness	0.87 mm
TR	5168.90 ms
TE	18 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	A32

## Contrast

Perfusion mode	SS-SI VASO
TI2	100 ms
TI1	50 ms
TI1s	50 ms
Flip angle	40 deg
Fat suppr.	Fat sat.
Fat sat. mode	Weak
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	183
Delay in TR	0 ms
Multiple series	Off
Perfusion mode	PICORE Q2T
Inversion time 1	50 ms
Saturation stop time	50 ms
Inversion time 2	100.0 ms
Flow limit	100 cm/s

## Resolution

Base resolution	170
Phase resolution	100 %
Slice resolution	100 %
Phase partial Fourier	6/8
Slice partial Fourier	Off
Interpolation	Off

PAT mode	GRAPPA
Accel. factor PE	6
Ref. lines PE	90
Accel. factor 3D	1
Ref. lines 3D	24
Reference scan mode	Separate

Prescan Normalize	Off
Raw filter	Off
Elliptical filter	Off
Hamming	Off

## Geometry

Multi-slice mode	Interleaved
Series	Ascending
Special sat.	Parallel F
Gap	25.0 mm
Thickness	100 mm
Table position	H
Table position	0 mm
Inline Composing	Off

## System

V32	Off
A32	On
Positioning mode	FIX
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
! Ref. amplitude 1H	225.000 V
Adjustment Tolerance	Auto
Adjust volume	
! Position	R1.9 A24.9 H19.3
! Orientation	T > C-1.8
! Rotation	0.00 deg
! R >> L	137 mm
! A >> P	177 mm
! F >> H	94 mm

## Physio

1st Signal/Mode	None
-----------------	------

## BOLD

Motion correction	Off
Spatial filter	Off

## Sequence

Introduction	On
Dimension	3D
Reordering	Linear
Contrasts	1
Bandwidth	980 Hz/Px
Free echo spacing	Off
Echo spacing	1.16 ms
EPI factor	226

# SIEMENS MAGNETOM Investigational\_Device\_7T syngo MR B17

RF pulse type	Normal
Gradient mode	Fast
Excitation	Slab-sel.
RF spoiling	On
<hr/>	
Ampl	150
MAGEC FA	9 in deg
ph.skip 4 Robert (the one)	1
MAGEC SS-SI?	On
NORDIC	On
log physio files	On
FFT scale	1.00
dummy prepscan time	3 s
z shim	0.00 mT/m*ms
RF duration	2560 us
RF BWTP	8.0
Renzo: Delta TI	45 ms
EFFECTIVE TR	578916 ms
PatPartitions	112
CAIPI scale	3
EPI phase correction	local
PAT refscan mode	Flash
FlashRef BaseRes	170
FlashRef BW	103 Hz/px
FlashRef TE	6100 us
FlashRef FA	5 deg
use CAIPI	On
CAIPI shift kz	0
CAIPI shift ky	3

# SIEMENS MAGNETOM Investigational\_Device\_7T syngo MR B17

\\USER\RenzoHuber\KENSU\20211011\_KEN\_movie2\VASO160\_G61\_nofatsat\_movie3

TA: 15:46 PAT: 6 Voxel size: 0.8x0.8x0.9 mm Rel. SNR: 1.00 USER: VASO\_160

## Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

## Routine

Slab group 1	
Slabs	1
Dist. factor	50 %
Position	R3.9 A25.3 H16.4
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	112
FoV read	138.7 mm
FoV phase	132.9 %
Slice thickness	0.87 mm
TR	5168.90 ms
TE	18 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	A32

## Contrast

Perfusion mode	SS-SI VASO
TI2	100 ms
TI1	50 ms
TI1s	50 ms
Flip angle	40 deg
Fat suppr.	Fat sat.
Fat sat. mode	Weak
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	183
Delay in TR	0 ms
Multiple series	Off
Perfusion mode	PICORE Q2T
Inversion time 1	50 ms
Saturation stop time	50 ms
Inversion time 2	100.0 ms
Flow limit	100 cm/s

## Resolution

Base resolution	170
Phase resolution	100 %
Slice resolution	100 %
Phase partial Fourier	6/8
Slice partial Fourier	Off
Interpolation	Off

PAT mode	GRAPPA
Accel. factor PE	6
Ref. lines PE	90
Accel. factor 3D	1
Ref. lines 3D	24
Reference scan mode	Separate

Prescan Normalize	Off
Raw filter	Off
Elliptical filter	Off
Hamming	Off

## Geometry

Multi-slice mode	Interleaved
Series	Ascending
Special sat.	Parallel F
Gap	25.0 mm
Thickness	100 mm
Table position	H
Table position	0 mm
Inline Composing	Off

## System

V32	Off
A32	On
Positioning mode	FIX
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
! Ref. amplitude 1H	225.000 V
Adjustment Tolerance	Auto
Adjust volume	
! Position	R1.9 A24.9 H19.3
! Orientation	T > C-1.8
! Rotation	0.00 deg
! R >> L	137 mm
! A >> P	177 mm
! F >> H	94 mm

## Physio

1st Signal/Mode	None
-----------------	------

## BOLD

Motion correction	Off
Spatial filter	Off

## Sequence

Introduction	On
Dimension	3D
Reordering	Linear
Contrasts	1
Bandwidth	980 Hz/Px
Free echo spacing	Off
Echo spacing	1.16 ms
EPI factor	226

# SIEMENS MAGNETOM Investigational\_Device\_7T syngo MR B17

RF pulse type	Normal
Gradient mode	Fast
Excitation	Slab-sel.
RF spoiling	On
<hr/>	
Ampl	150
MAGEC FA	9 in deg
ph.skip 4 Robert (the one)	1
MAGEC SS-SI?	On
NORDIC	On
log physio files	On
FFT scale	1.00
dummy prepscan time	3 s
z shim	0.00 mT/m*ms
RF duration	2560 us
RF BWTP	8.0
Renzo: Delta TI	45 ms
EFFECTIVE TR	578916 ms
PatPartitions	112
CAIPI scale	3
EPI phase correction	local
PAT refscan mode	Flash
FlashRef BaseRes	170
FlashRef BW	103 Hz/px
FlashRef TE	6100 us
FlashRef FA	5 deg
use CAIPI	On
CAIPI shift kz	0
CAIPI shift ky	3



# SIEMENS MAGNETOM Investigational\_Device\_7T syngo MR B17

\\USER\RenzoHuber\KENSHU\20211011\_KEN\_movie2\VASO160\_G61\_nofatsat\_movie4

TA: 15:46 PAT: 6 Voxel size: 0.8x0.8x0.9 mm Rel. SNR: 1.00 USER: VASO\_160

## Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

## Routine

Slab group 1	
Slabs	1
Dist. factor	50 %
Position	R3.9 A25.3 H16.4
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	112
FoV read	138.7 mm
FoV phase	132.9 %
Slice thickness	0.87 mm
TR	5168.90 ms
TE	18 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	A32

## Contrast

Perfusion mode	SS-SI VASO
TI2	100 ms
TI1	50 ms
TI1s	50 ms
Flip angle	40 deg
Fat suppr.	Fat sat.
Fat sat. mode	Weak
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	183
Delay in TR	0 ms
Multiple series	Off
Perfusion mode	PICORE Q2T
Inversion time 1	50 ms
Saturation stop time	50 ms
Inversion time 2	100.0 ms
Flow limit	100 cm/s

## Resolution

Base resolution	170
Phase resolution	100 %
Slice resolution	100 %
Phase partial Fourier	6/8
Slice partial Fourier	Off
Interpolation	Off

PAT mode	GRAPPA
Accel. factor PE	6
Ref. lines PE	90
Accel. factor 3D	1
Ref. lines 3D	24
Reference scan mode	Separate

Prescan Normalize	Off
Raw filter	Off
Elliptical filter	Off
Hamming	Off

## Geometry

Multi-slice mode	Interleaved
Series	Ascending
Special sat.	Parallel F
Gap	25.0 mm
Thickness	100 mm
Table position	H
Table position	0 mm
Inline Composing	Off

## System

V32	Off
A32	On
Positioning mode	FIX
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
! Ref. amplitude 1H	225.000 V
Adjustment Tolerance	Auto
Adjust volume	
! Position	R1.9 A24.9 H19.3
! Orientation	T > C-1.8
! Rotation	0.00 deg
! R >> L	137 mm
! A >> P	177 mm
! F >> H	94 mm

## Physio

1st Signal/Mode	None
-----------------	------

## BOLD

Motion correction	Off
Spatial filter	Off

## Sequence

Introduction	On
Dimension	3D
Reordering	Linear
Contrasts	1
Bandwidth	980 Hz/Px
Free echo spacing	Off
Echo spacing	1.16 ms
EPI factor	226

# SIEMENS MAGNETOM Investigational\_Device\_7T syngo MR B17

RF pulse type	Normal
Gradient mode	Fast
Excitation	Slab-sel.
RF spoiling	On
<hr/>	
Ampl	150
MAGEC FA	9 in deg
ph.skip 4 Robert (the one)	1
MAGEC SS-SI?	On
NORDIC	On
log physio files	On
FFT scale	1.00
dummy prepscan time	3 s
z shim	0.00 mT/m*ms
RF duration	2560 us
RF BWTP	8.0
Renzo: Delta TI	45 ms
EFFECTIVE TR	578916 ms
PatPartitions	112
CAIPI scale	3
EPI phase correction	local
PAT refscan mode	Flash
FlashRef BaseRes	170
FlashRef BW	103 Hz/px
FlashRef TE	6100 us
FlashRef FA	5 deg
use CAIPI	On
CAIPI shift kz	0
CAIPI shift ky	3

# SIEMENS MAGNETOM Investigational\_Device\_7T syngo MR B17

\\USER\RenzoHuber\KENSHU\20211011\_KEN\_movie2\VASO160\_G61\_nofatsat\_movie5

TA: 15:46 PAT: 6 Voxel size: 0.8x0.8x0.9 mm Rel. SNR: 1.00 USER: VASO\_160

## Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

## Routine

Slab group 1	
Slabs	1
Dist. factor	50 %
Position	R3.9 A25.3 H16.4
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	112
FoV read	138.7 mm
FoV phase	132.9 %
Slice thickness	0.87 mm
TR	5168.90 ms
TE	18 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	A32

## Contrast

Perfusion mode	SS-SI VASO
TI2	100 ms
TI1	50 ms
TI1s	50 ms
Flip angle	40 deg
Fat suppr.	Fat sat.
Fat sat. mode	Weak
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	183
Delay in TR	0 ms
Multiple series	Off
Perfusion mode	PICORE Q2T
Inversion time 1	50 ms
Saturation stop time	50 ms
Inversion time 2	100.0 ms
Flow limit	100 cm/s

## Resolution

Base resolution	170
Phase resolution	100 %
Slice resolution	100 %
Phase partial Fourier	6/8
Slice partial Fourier	Off
Interpolation	Off

PAT mode	GRAPPA
Accel. factor PE	6
Ref. lines PE	90
Accel. factor 3D	1
Ref. lines 3D	24
Reference scan mode	Separate

Prescan Normalize	Off
Raw filter	Off
Elliptical filter	Off
Hamming	Off

## Geometry

Multi-slice mode	Interleaved
Series	Ascending
Special sat.	Parallel F
Gap	25.0 mm
Thickness	100 mm
Table position	H
Table position	0 mm
Inline Composing	Off

## System

V32	Off
A32	On
Positioning mode	FIX
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
! Ref. amplitude 1H	225.000 V
Adjustment Tolerance	Auto
Adjust volume	
! Position	R1.9 A24.9 H19.3
! Orientation	T > C-1.8
! Rotation	0.00 deg
! R >> L	137 mm
! A >> P	177 mm
! F >> H	94 mm

## Physio

1st Signal/Mode	None
-----------------	------

## BOLD

Motion correction	Off
Spatial filter	Off

## Sequence

Introduction	On
Dimension	3D
Reordering	Linear
Contrasts	1
Bandwidth	980 Hz/Px
Free echo spacing	Off
Echo spacing	1.16 ms
EPI factor	226

# SIEMENS MAGNETOM Investigational\_Device\_7T syngo MR B17

RF pulse type	Normal
Gradient mode	Fast
Excitation	Slab-sel.
RF spoiling	On
<hr/>	
Ampl	150
MAGEC FA	9 in deg
ph.skip 4 Robert (the one)	1
MAGEC SS-SI?	On
NORDIC	On
log physio files	On
FFT scale	1.00
dummy prepscan time	3 s
z shim	0.00 mT/m*ms
RF duration	2560 us
RF BWTP	8.0
Renzo: Delta TI	45 ms
EFFECTIVE TR	578916 ms
PatPartitions	112
CAIPI scale	3
EPI phase correction	local
PAT refscan mode	Flash
FlashRef BaseRes	170
FlashRef BW	103 Hz/px
FlashRef TE	6100 us
FlashRef FA	5 deg
use CAIPI	On
CAIPI shift kz	0
CAIPI shift ky	3