

\\NSI\_MR\_Research\\final\_studies\\3D\_EPI\_Spielplatz\\Ausprobiert\_invivo\\RenzLay3\_functional

TA: 14:51 PM: REF Voxel size: 0.9×0.9×0.9 mmPAT: 2 Rel. SNR: 1.00 : RenzLay

**Properties**

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	Off
Start measurements	Single measurement

**Routine**

Slab group	1
Slabs	1
Dist. factor	50 %
Position	R1.2 A4.1 H50.2 mm
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
Phase oversampling	0 %
Slice oversampling	8.3 %
Slices per slab	24
FoV read	173 mm
FoV phase	100.0 %
Slice thickness	0.92 mm
TR	71 ms
TE	38.0 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	HC1-6

**Contrast - Common**

TR	71 ms
TE	38.0 ms
MTC	Off
Flip angle	19 deg
Fat suppr.	Fat sat.
SWI	Off

**Contrast - Dynamic**

Averages	1
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	303
Multiple series	Each measurement

**Resolution - Common**

FoV read	173 mm
FoV phase	100.0 %
Slice thickness	0.92 mm
Base resolution	192
Phase resolution	100 %
Slice resolution	100 %
Phase partial Fourier	6/8
Slice partial Fourier	6/8
Interpolation	Off

**Resolution - iPAT**

PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	47
Accel. factor 3D	1
Ref. lines 3D	24
Reference scan mode	GRE/separate

**Resolution - Filter Image**

Distortion Corr.	Off
Prescan Normalize	Off

**Resolution - Filter Rawdata**

Raw filter	Off
Elliptical filter	Off
Hamming	Off

**Geometry - Common**

Slab group	1
Slabs	1
Dist. factor	50 %
Position	R1.2 A4.1 H50.2 mm
Orientation	Transversal
Phase enc. dir.	A >> P
Slice oversampling	8.3 %
Slices per slab	24
FoV read	173 mm
FoV phase	100.0 %
Slice thickness	0.92 mm
TR	71 ms
Multi-slice mode	Interleaved
Series	Interleaved
Concatenations	1

**Geometry - AutoAlign**

Slab group	1
Position	R1.2 A4.1 H50.2 mm
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
Initial Position	R1.2 A4.1 H50.2
R	1.2 mm
A	4.1 mm
H	50.2 mm
Initial Rotation	0.00 deg
Initial Orientation	Transversal

**Geometry - Saturation**

Fat suppr.	Fat sat.
Special sat.	None

**Geometry - Tim Planning Suite**

Set-n-Go Protocol	Off
Table position	H
Table position	0 mm
Inline Composing	Off

**System - Miscellaneous**

Positioning mode	REF
Table position	H
Table position	0 mm

**System - Miscellaneous**

MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
Save uncombined	Off
Matrix Optimization	Off
AutoAlign	---
Coil Select Mode	Default

**System - Adjustments**

B0 Shim mode	Standard
B1 Shim mode	TrueForm
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

**System - Adjust Volume**

! Position	L0.0 A0.7 H38.4 mm
! Orientation	Transversal
! Rotation	0.00 deg
! A >> P	173 mm
! R >> L	173 mm
! F >> H	45 mm
Reset	Off

**System - pTx Volumes**

B1 Shim mode	TrueForm
--------------	----------

**System - Tx/Rx**

Frequency 1H	123.248138 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

**Physio - Signal1**

1st Signal/Mode	None
TR	71 ms
Concatenations	1

**Sequence - Part 1**

Introduction	Off
Dimension	3D
Reordering	Linear rot.
Multi-slice mode	Interleaved
Free echo spacing	On
Echo spacing	1.11 ms
Bandwidth	1302 Hz/Px

**Sequence - Part 2**

EPI factor	55
RF pulse type	Normal
Gradient mode	Fast
RF spoiling	On