

\\USER\NIMH\MARRETT\190908_HILUNA\RENZLAY4_26slices_0.5_SAMOWU

TA: 9:18 PM: FIX Voxel size: 0.6×0.6×0.6 mmPAT: 3 Rel. SNR: 1.00 : RenzLa4

Properties

| | |
|---|--------------------|
| Prio recon | Off |
| Load images to viewer | On |
| Inline movie | Off |
| Auto store images | On |
| Load images to stamp segments | Off |
| Load images to graphic segments | Off |
| Auto open inline display | Off |
| Auto close inline display | Off |
| Start measurement without further preparation | Off |
| Wait for user to start | Off |
| Start measurements | Single measurement |

Routine

| | |
|--------------------|----------------------|
| Slab group | 1 |
| Slabs | 1 |
| Position | L0.0 A5.3 H27.1 mm |
| Orientation | T > C-11.4 |
| Phase enc. dir. | A >> P |
| AutoAlign | --- |
| Slice oversampling | 0.0 % |
| Slices per slab | 144 |
| FoV read | 177 mm |
| FoV phase | 100.0 % |
| Slice thickness | 0.60 mm |
| TR 1 | 60.55 ms |
| TR 2 | 45827.28 ms |
| TE 1 | 21.0 ms |
| Averages | 1 |
| TE segmentation | 1 |
| Filter | Distortion Corr.(3D) |
| Coil elements | A32 |

Contrast - Common

| | |
|-------------------|-------------------|
| TR 1 | 60.55 ms |
| TR 2 | 45827.28 ms |
| TE 1 | 21.0 ms |
| Multi-echo dTE | 60.0 ms |
| MTC | Off |
| Magn. preparation | Non-sel. IR T1map |
| Flip angle | 33 deg |
| Fat suppr. | None |
| Number of TIs | 2 |

Contrast - Dynamic

| | |
|----------------------|-----------|
| Averages | 1 |
| Averaging mode | Long term |
| Reconstruction | Magnitude |
| Measurements | 12 |
| Pause after meas. 1 | 0.0 s |
| Pause after meas. 2 | 0.0 s |
| Pause after meas. 3 | 0.0 s |
| Pause after meas. 4 | 0.0 s |
| Pause after meas. 5 | 0.0 s |
| Pause after meas. 6 | 0.0 s |
| Pause after meas. 7 | 0.0 s |
| Pause after meas. 8 | 0.0 s |
| Pause after meas. 9 | 0.0 s |
| Pause after meas. 10 | 0.0 s |
| Pause after meas. 11 | 0.0 s |

Resolution - Common

| | |
|-----------------------|---------|
| FoV read | 177 mm |
| FoV phase | 100.0 % |
| Slice thickness | 0.60 mm |
| Base resolution | 318 |
| Phase resolution | 100 % |
| Slice resolution | 100 % |
| Phase partial Fourier | 6/8 |
| Slice partial Fourier | Off |
| Interpolation | Off |

Resolution - iPAT

| | |
|---------------------|--------------|
| PAT mode | CAIPIRINHA |
| Acc. factor PE | 3 |
| Ref. lines PE | 44 |
| Acc. factor 3D | 1 |
| Ref. lines 3D | 24 |
| CAIPI 3D Shift | 0 |
| Reference scan mode | GRE/separate |
| CAIPIRINHA mode | Free |

Resolution - Filter Image

| | |
|-------------------|-----|
| Image Filter | Off |
| Distortion Corr. | On |
| Mode | 3D |
| Unfiltered images | Off |
| Prescan Normalize | Off |
| Normalize | Off |
| B1 filter | Off |

Resolution - Filter Rawdata

| | |
|-------------------|-----|
| Raw filter | Off |
| Elliptical filter | Off |

Geometry - Common

| | |
|--------------------|--------------------|
| Slab group | 1 |
| Slabs | 1 |
| Position | L0.0 A5.3 H27.1 mm |
| Orientation | T > C-11.4 |
| Phase enc. dir. | A >> P |
| Slice oversampling | 0.0 % |
| Slices per slab | 144 |
| FoV read | 177 mm |
| FoV phase | 100.0 % |
| Slice thickness | 0.60 mm |
| TR 1 | 60.55 ms |
| TR 2 | 45827.28 ms |
| Multi-slice mode | Interleaved |
| Series | Ascending |
| TE segmentation | 1 |

Geometry - AutoAlign

| | |
|------------------|--------------------|
| Slab group | 1 |
| Position | L0.0 A5.3 H27.1 mm |
| Orientation | T > C-11.4 |
| Phase enc. dir. | A >> P |
| AutoAlign | --- |
| Initial Position | L0.0 A5.3 H27.1 |
| L | 0.0 mm |
| A | 5.3 mm |
| H | 27.1 mm |

Geometry - AutoAlign

| | |
|---------------------|----------|
| Initial Rotation | 0.00 deg |
| Initial Orientation | T > C |
| T > C | -11.4 |
| > S | 0.0 |

Geometry - Saturation

| | |
|-----------------|----------|
| Saturation mode | Standard |
| Fat suppr. | None |

System - Miscellaneous

| | |
|---------------------|----------------|
| Positioning mode | FIX |
| Table position | H |
| Table position | 0 mm |
| MSMA | S - C - T |
| Sagittal | R >> L |
| Coronal | A >> P |
| Transversal | F >> H |
| Coil Combine Mode | Sum of Squares |
| Save uncombined | Off |
| Matrix Optimization | Off |
| AutoAlign | --- |
| Coil Select Mode | Default |

System - Adjustments

| | |
|--------------------------|----------|
| B0 Shim mode | Tune up |
| B1 Shim mode | TrueForm |
| Confirm freq. adjustment | Off |
| Assume Dominant Fat | Off |
| Assume Silicone | Off |
| Adjustment Tolerance | Auto |

System - Adjust Volume

| | |
|---------------|--------------------|
| ! Position | L0.0 A8.9 H27.3 mm |
| ! Orientation | T > C-11.6 |
| ! Rotation | 0.00 deg |
| ! A >> P | 179 mm |
| ! R >> L | 184 mm |
| ! F >> H | 95 mm |
| Reset | Off |

System - Tx/Rx

| | |
|---------------------|----------------|
| Frequency 1H | 297.201506 MHz |
| Correction factor | 1 |
| Gain | Low |
| Img. Scale Cor. | 1.000 |
| Reset | Off |
| ! Ref. amplitude 1H | 220.000 V |

BOLD

| | |
|-------------------------|----------|
| GLM Statistics | Off |
| Dynamic t-maps | Off |
| Ignore meas. at start | 0 |
| Ignore after transition | 0 |
| Model transition states | On |
| Temp. highpass filter | On |
| Threshold | 4.00 |
| Paradigm size | 20 |
| Meas[1] | Baseline |
| Meas[2] | Baseline |
| Meas[3] | Baseline |
| Meas[4] | Baseline |
| Meas[5] | Baseline |
| Meas[6] | Baseline |
| Meas[7] | Baseline |

BOLD

| | |
|-------------------|----------|
| Meas[8] | Baseline |
| Meas[9] | Baseline |
| Meas[10] | Baseline |
| Meas[11] | Active |
| Meas[12] | Active |
| Meas[13] | Active |
| Meas[14] | Active |
| Meas[15] | Active |
| Meas[16] | Active |
| Meas[17] | Active |
| Meas[18] | Active |
| Meas[19] | Active |
| Meas[20] | Active |
| Motion correction | Off |
| Spatial filter | Off |
| Measurements | 12 |

Sequence - Part 1

| | |
|------------------|-------------|
| Introduction | On |
| Dimension | 3D |
| Reordering | Linear |
| Asymmetric echo | Off |
| Contrasts | 1 |
| Multi-slice mode | Interleaved |
| Echo spacing | 1.29 ms |
| Bandwidth | 874 Hz/Px |

Sequence - Part 2

| | |
|---------------|-----------|
| EPI factor | 40 |
| Segmentation | 2 |
| RF pulse type | Normal |
| Gradient mode | Fast |
| Excitation | Slab-sel. |
| RF spoiling | On |
| Turbo factor | 26 |

Sequence - Special

| | |
|--------------------|------------|
| PATRef FA | 12 deg |
| RF duration | 2500 us |
| RF BWT product | 25 |
| Slab Scale | 90 % |
| Ernst T1 | 1200 ms |
| PATRef prep. shots | 10 |
| Volume dummy shots | 0 |
| Dummy Measurements | 0 |
| CHECK FLIP ANGLE! | On |
| Invert PE | Off |
| Invert 3D | Off |
| Min. TE if PF | On |
| Are you Renzo? | Off |
| Echo Time Shift | On |
| Water Exc. | -none- |
| External PC | per Series |
| HSN RF power scale | 3.00 |
| Inversion Delay | 650000 us |
| Relaxation Delay | 0 us |

Sequence - Assistant

| | |
|------|-----|
| Mode | Off |
|------|-----|

\\USER\NIMH\MARRETT\190908_HILUNAIRENZLAY4_96 slices_MAGEC_compation

TA: 1:02 PM: FIX Voxel size: 0.8×0.8×0.8 mmPAT: 3 Rel. SNR: 1.00 : RenzLa4

Properties

| | |
|---|--------------------|
| Prio recon | Off |
| Load images to viewer | On |
| Inline movie | Off |
| Auto store images | On |
| Load images to stamp segments | Off |
| Load images to graphic segments | Off |
| Auto open inline display | Off |
| Auto close inline display | Off |
| Start measurement without further preparation | Off |
| Wait for user to start | Off |
| Start measurements | Single measurement |

Routine

| | |
|--------------------|----------------------|
| Slab group | 1 |
| Slabs | 1 |
| Position | L0.0 A5.3 H27.1 mm |
| Orientation | T > C-11.4 |
| Phase enc. dir. | A >> P |
| AutoAlign | --- |
| Slice oversampling | 0.0 % |
| Slices per slab | 96 |
| FoV read | 177 mm |
| FoV phase | 100.0 % |
| Slice thickness | 0.85 mm |
| TR 1 | 61.56 ms |
| TR 2 | 13830.54 ms |
| TE 1 | 22.0 ms |
| Averages | 1 |
| TE segmentation | 1 |
| Filter | Distortion Corr.(3D) |
| Coil elements | A32 |

Contrast - Common

| | |
|-------------------|-------------------|
| TR 1 | 61.56 ms |
| TR 2 | 13830.54 ms |
| TE 1 | 22.0 ms |
| Multi-echo dTE | 60.0 ms |
| MTC | Off |
| Magn. preparation | Non-sel. IR T1map |
| Flip angle | 33 deg |
| Fat suppr. | None |
| Number of TIs | 2 |

Contrast - Dynamic

| | |
|---------------------|-----------|
| Averages | 1 |
| Averaging mode | Long term |
| Reconstruction | Magnitude |
| Measurements | 4 |
| Pause after meas. 1 | 0.0 s |
| Pause after meas. 2 | 0.0 s |
| Pause after meas. 3 | 0.0 s |

Resolution - Common

| | |
|------------------|---------|
| FoV read | 177 mm |
| FoV phase | 100.0 % |
| Slice thickness | 0.85 mm |
| Base resolution | 216 |
| Phase resolution | 100 % |
| Slice resolution | 100 % |

Resolution - Common

| | |
|-----------------------|-----|
| Phase partial Fourier | 6/8 |
| Slice partial Fourier | Off |
| Interpolation | Off |

Resolution - iPAT

| | |
|---------------------|--------------|
| PAT mode | CAIPIRINHA |
| Acc. factor PE | 3 |
| Ref. lines PE | 44 |
| Acc. factor 3D | 1 |
| Ref. lines 3D | 24 |
| CAIPI 3D Shift | 0 |
| Reference scan mode | GRE/separate |
| CAIPIRINHA mode | Free |

Resolution - Filter Image

| | |
|-------------------|-----|
| Image Filter | Off |
| Distortion Corr. | On |
| Mode | 3D |
| Unfiltered images | Off |
| Prescan Normalize | Off |
| Normalize | Off |
| B1 filter | Off |

Resolution - Filter Rawdata

| | |
|-------------------|-----|
| Raw filter | Off |
| Elliptical filter | Off |

Geometry - Common

| | |
|--------------------|--------------------|
| Slab group | 1 |
| Slabs | 1 |
| Position | L0.0 A5.3 H27.1 mm |
| Orientation | T > C-11.4 |
| Phase enc. dir. | A >> P |
| Slice oversampling | 0.0 % |
| Slices per slab | 96 |
| FoV read | 177 mm |
| FoV phase | 100.0 % |
| Slice thickness | 0.85 mm |
| TR 1 | 61.56 ms |
| TR 2 | 13830.54 ms |
| Multi-slice mode | Interleaved |
| Series | Ascending |
| TE segmentation | 1 |

Geometry - AutoAlign

| | |
|---------------------|--------------------|
| Slab group | 1 |
| Position | L0.0 A5.3 H27.1 mm |
| Orientation | T > C-11.4 |
| Phase enc. dir. | A >> P |
| AutoAlign | --- |
| Initial Position | L0.0 A5.3 H27.1 |
| L | 0.0 mm |
| A | 5.3 mm |
| H | 27.1 mm |
| Initial Rotation | 0.00 deg |
| Initial Orientation | T > C |
| T > C | -11.4 |
| > S | 0.0 |

Geometry - Saturation

| | |
|-----------------|----------|
| Saturation mode | Standard |
| Fat suppr. | None |

System - Miscellaneous

| | |
|---------------------|----------------|
| Positioning mode | FIX |
| Table position | H |
| Table position | 0 mm |
| MSMA | S - C - T |
| Sagittal | R >> L |
| Coronal | A >> P |
| Transversal | F >> H |
| Coil Combine Mode | Sum of Squares |
| Save uncombined | Off |
| Matrix Optimization | Off |
| AutoAlign | --- |
| Coil Select Mode | Default |

System - Adjustments

| | |
|--------------------------|----------|
| B0 Shim mode | Tune up |
| B1 Shim mode | TrueForm |
| Confirm freq. adjustment | Off |
| Assume Dominant Fat | Off |
| Assume Silicone | Off |
| Adjustment Tolerance | Auto |

System - Adjust Volume

| | |
|---------------|--------------------|
| ! Position | L0.0 A8.9 H27.3 mm |
| ! Orientation | T > C-11.6 |
| ! Rotation | 0.00 deg |
| ! A >> P | 179 mm |
| ! R >> L | 184 mm |
| ! F >> H | 95 mm |
| Reset | Off |

System - Tx/Rx

| | |
|---------------------|----------------|
| Frequency 1H | 297.201506 MHz |
| Correction factor | 1 |
| Gain | Low |
| Img. Scale Cor. | 1.000 |
| Reset | Off |
| ! Ref. amplitude 1H | 220.000 V |

BOLD

| | |
|-------------------------|----------|
| GLM Statistics | Off |
| Dynamic t-maps | Off |
| Ignore meas. at start | 0 |
| Ignore after transition | 0 |
| Model transition states | On |
| Temp. highpass filter | On |
| Threshold | 4.00 |
| Paradigm size | 20 |
| Meas[1] | Baseline |
| Meas[2] | Baseline |
| Meas[3] | Baseline |
| Meas[4] | Baseline |
| Meas[5] | Baseline |
| Meas[6] | Baseline |
| Meas[7] | Baseline |
| Meas[8] | Baseline |
| Meas[9] | Baseline |
| Meas[10] | Baseline |
| Meas[11] | Active |
| Meas[12] | Active |
| Meas[13] | Active |

BOLD

| | |
|-------------------|--------|
| Meas[14] | Active |
| Meas[15] | Active |
| Meas[16] | Active |
| Meas[17] | Active |
| Meas[18] | Active |
| Meas[19] | Active |
| Meas[20] | Active |
| Motion correction | Off |
| Spatial filter | Off |
| Measurements | 4 |

Sequence - Part 1

| | |
|------------------|-------------|
| Introduction | On |
| Dimension | 3D |
| Reordering | Linear |
| Asymmetric echo | Off |
| Contrasts | 1 |
| Multi-slice mode | Interleaved |
| Echo spacing | 1.01 ms |
| Bandwidth | 1102 Hz/Px |

Sequence - Part 2

| | |
|---------------|-----------|
| EPI factor | 54 |
| Segmentation | 1 |
| RF pulse type | Normal |
| Gradient mode | Fast |
| Excitation | Slab-sel. |
| RF spoiling | On |
| Turbo factor | 32 |

Sequence - Special

| | |
|--------------------|------------|
| PATRef FA | 12 deg |
| RF duration | 2500 us |
| RF BWT product | 25 |
| Slab Scale | 90 % |
| Ernst T1 | 1200 ms |
| PATRef prep. shots | 10 |
| Volume dummy shots | 0 |
| Dummy Measurements | 0 |
| CHECK FLIP ANGLE! | On |
| Invert PE | Off |
| Invert 3D | Off |
| Min. TE if PF | On |
| Are you Renzo? | Off |
| Echo Time Shift | On |
| Water Exc. | -none- |
| External PC | per Series |
| HSN RF power scale | 3.00 |
| Inversion Delay | 650000 us |
| Relaxation Delay | 0 us |

Sequence - Assistant

| | |
|------|-----|
| Mode | Off |
|------|-----|