

\\USER\Test\Renzo\WholeBrain\_invivo20220710\WB\_0p6\_ipat8\_S2\_nolce\_movie

TA: 14:32 PM: FIX Voxel size: 0.6×0.6×0.6 mmPAT: 8 Rel. SNR: 1.00 : 8036cbcd

**Properties**

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	Off
Start measurements	Single measurement

**Routine**

Slab group	1
Slabs	1
Position	L0.0 P17.4 H29.1 mm
Orientation	T > C-12.8
Phase enc. dir.	A >> P
AutoAlign	---
Slab Scale	-10 %
Slices per slab	180
FoV read	200 mm
FoV phase	100.0 %
Slice thickness	0.63 mm
TR 1	45.9 ms
TR 2	9402 ms
TE 1	16.10 ms
Averages	1
Multi-echo Shots	1
Filter	Distortion Corr.(3D)
Coil elements	R96

**Contrast - Common**

TR 1	45.9 ms
TR 2	9402 ms
TE 1	16.10 ms
Multi-echo spacing	41.39 ms
Magn. preparation	Non-sel. HSN IR
TI 1	1569.8 ms
TI 2	3635.3 ms
Flip angle	12 deg
Fat suppr.	None
Magn. Prep. Shots	2

**Contrast - Dynamic**

Averages	1
Averaging mode	Long term
Reconstruction	Magn./Phase
Measurements	92
Pause after meas.	0.0 s

**Resolution - Common**

FoV read	200 mm
FoV phase	100.0 %
Slice thickness	0.63 mm
Base resolution	314
Phase resolution	100 %
Slice resolution	100 %
Phase partial Fourier	6/8

**Resolution - Common**

Slice partial Fourier	Off
Interpolation	Off

**Resolution - iPAT**

PAT mode	CAIPIRINHA
Acc. factor PE	2
Ref. lines PE	75
Acc. factor 3D	4
Ref. lines 3D	24
CAIPI 3D Shift	2
Reference Scan Mode	GRE/separate
CAIPIRINHA mode	Free

**Resolution - Filter Image**

Image Filter	Off
Distortion Corr.	On
Mode	3D
Unfiltered images	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off

**Resolution - Filter Rawdata**

Raw filter	Off
Elliptical filter	Off

**Geometry - Common**

Slab group	1
Slabs	1
Position	L0.0 P17.4 H29.1 mm
Orientation	T > C-12.8
Phase enc. dir.	A >> P
Slab Scale	-10 %
Slices per slab	180
FoV read	200 mm
FoV phase	100.0 %
Slice thickness	0.63 mm
TR 1	45.9 ms
TR 2	9402 ms
Multi-slice mode	Interleaved
Series	Ascending
Multi-echo Shots	1

**Geometry - AutoAlign**

Slab group	1
Position	L0.0 P17.4 H29.1 mm
Orientation	T > C-12.8
Phase enc. dir.	A >> P
AutoAlign	---
Initial Position	L0.0 P17.4 H29.1
L	0.0 mm
P	17.4 mm
H	29.1 mm
Initial Rotation	0.00 deg
Initial Orientation	T > C
T > C	-12.8
> S	0.0

**Geometry - Saturation**

Saturation mode	Standard
-----------------	----------

**Geometry - Saturation**

Fat suppr.	None
------------	------

**Geometry - Tim Planning Suite**

Set-n-Go Protocol	Off
Table position	H
Table position	0 mm
Inline Composing	Off

**System - Miscellaneous**

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
Save uncombined	Off
Matrix Optimization	Off
AutoAlign	---
Coil Select Mode	Default

**System - Adjustments**

B0 Shim mode	Tune up
B1 Shim mode	TrueForm
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

**System - Adjust Volume**

! Position	L0.0 P21.6 H33.0 mm
! Orientation	T > C-16.2
! Rotation	0.00 deg
! A >> P	199 mm
! R >> L	189 mm
! F >> H	107 mm
Reset	Off

**System - pTx Volumes**

B1 Shim mode	TrueForm
Excitation	Slab-sel.

**System - Tx/Rx**

Frequency 1H	297.192191 MHz
Correction factor	1
Gain	Low
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

**Sequence - Part 1**

Introduction	On
Dimension	3D
Reordering	Linear
Asymmetric echo	Off
Contrasts	1
Multi-slice mode	Interleaved
Echo spacing	0.69 ms
Bandwidth	1592 Hz/Px

**Sequence - Part 2**

EPI factor	59
------------	----

**Sequence - Part 2**

Segmentation	2
RF pulse type	Normal
Gradient mode	Normal
Excitation	Slab-sel.
RF spoiling	On
Turbo factor	45

**Sequence - Special**

PATRef FA	3 deg
RF duration	1400 us
RF BWT product	8
Ernst T1	1200 ms
PATRef prep. shots	10
Volume dummy shots	0
Dummy Measurements	0
Invert PE	Off
Min. TE if PF	On
Echo Time Shift	On
Ramp Sampling	On
NORDIC	On
SVDPC	On
Invert 3D	Off
Invert RO	Off
Alternate RO	On
Disable PF reco	Off
Disable PF reco	Off
Save sampling	Off
PE VComp	Off
Water Exc.	Long bino-11
External PC	per Series
EPI rise time factor	1.10
Mosaic DICOMs	On
Modify Ice Config	Off
HSN RF power scale	2.50
Inversion Delay	550 ms
Relaxation Delay	0 ms
Var. FA /MAGEC	1

**Sequence - Assistant**

Mode	Off
------	-----