

Table of contents

\\MARTINOS DEVELOPER

<table border="1"> <tr> <td>HUBER</td> </tr> </table>	HUBER
HUBER	

<table border="1"> <tr> <td>FRISGO_functional</td> </tr> </table>	FRISGO_functional
FRISGO_functional	

<table border="1"> <tr> <td>FRISGO_20251106_COD_M1_2nd</td> </tr> </table>	FRISGO_20251106_COD_M1_2nd
FRISGO_20251106_COD_M1_2nd	

<table border="1"> <tr> <td> <table border="1"> <tr> <td> AAhead_scout wip_ShimToolbox_B0map_3mm Scout_alternative AAsetter--192.168.2.1--ep2d rslh_ep3d_0.35_test rslh_ep3d_0.35_run1 Scout_alternative AAsetter--192.168.2.1--ep2d rslh_ep3d_0.35_run2_RO Scout_alternative AAsetter--192.168.2.1--ep2d rslh_ep3d_0.35_run3 Scout_alternative AAsetter--192.168.2.1--ep2d rslh_ep3d_0.35_run4_RO Scout_alternative AAsetter--192.168.2.1--ep2d rslh_ep3d_0.35_run5 mesovein_WE_RO1_use_me_princess mesovein_WE_RO1_use_me_princess </td> </tr> </table> </td> </tr> </table>	<table border="1"> <tr> <td> AAhead_scout wip_ShimToolbox_B0map_3mm Scout_alternative AAsetter--192.168.2.1--ep2d rslh_ep3d_0.35_test rslh_ep3d_0.35_run1 Scout_alternative AAsetter--192.168.2.1--ep2d rslh_ep3d_0.35_run2_RO Scout_alternative AAsetter--192.168.2.1--ep2d rslh_ep3d_0.35_run3 Scout_alternative AAsetter--192.168.2.1--ep2d rslh_ep3d_0.35_run4_RO Scout_alternative AAsetter--192.168.2.1--ep2d rslh_ep3d_0.35_run5 mesovein_WE_RO1_use_me_princess mesovein_WE_RO1_use_me_princess </td> </tr> </table>	AAhead_scout wip_ShimToolbox_B0map_3mm Scout_alternative AAsetter--192.168.2.1--ep2d rslh_ep3d_0.35_test rslh_ep3d_0.35_run1 Scout_alternative AAsetter--192.168.2.1--ep2d rslh_ep3d_0.35_run2_RO Scout_alternative AAsetter--192.168.2.1--ep2d rslh_ep3d_0.35_run3 Scout_alternative AAsetter--192.168.2.1--ep2d rslh_ep3d_0.35_run4_RO Scout_alternative AAsetter--192.168.2.1--ep2d rslh_ep3d_0.35_run5 mesovein_WE_RO1_use_me_princess mesovein_WE_RO1_use_me_princess
<table border="1"> <tr> <td> AAhead_scout wip_ShimToolbox_B0map_3mm Scout_alternative AAsetter--192.168.2.1--ep2d rslh_ep3d_0.35_test rslh_ep3d_0.35_run1 Scout_alternative AAsetter--192.168.2.1--ep2d rslh_ep3d_0.35_run2_RO Scout_alternative AAsetter--192.168.2.1--ep2d rslh_ep3d_0.35_run3 Scout_alternative AAsetter--192.168.2.1--ep2d rslh_ep3d_0.35_run4_RO Scout_alternative AAsetter--192.168.2.1--ep2d rslh_ep3d_0.35_run5 mesovein_WE_RO1_use_me_princess mesovein_WE_RO1_use_me_princess </td> </tr> </table>	AAhead_scout wip_ShimToolbox_B0map_3mm Scout_alternative AAsetter--192.168.2.1--ep2d rslh_ep3d_0.35_test rslh_ep3d_0.35_run1 Scout_alternative AAsetter--192.168.2.1--ep2d rslh_ep3d_0.35_run2_RO Scout_alternative AAsetter--192.168.2.1--ep2d rslh_ep3d_0.35_run3 Scout_alternative AAsetter--192.168.2.1--ep2d rslh_ep3d_0.35_run4_RO Scout_alternative AAsetter--192.168.2.1--ep2d rslh_ep3d_0.35_run5 mesovein_WE_RO1_use_me_princess mesovein_WE_RO1_use_me_princess	
AAhead_scout wip_ShimToolbox_B0map_3mm Scout_alternative AAsetter--192.168.2.1--ep2d rslh_ep3d_0.35_test rslh_ep3d_0.35_run1 Scout_alternative AAsetter--192.168.2.1--ep2d rslh_ep3d_0.35_run2_RO Scout_alternative AAsetter--192.168.2.1--ep2d rslh_ep3d_0.35_run3 Scout_alternative AAsetter--192.168.2.1--ep2d rslh_ep3d_0.35_run4_RO Scout_alternative AAsetter--192.168.2.1--ep2d rslh_ep3d_0.35_run5 mesovein_WE_RO1_use_me_princess mesovein_WE_RO1_use_me_princess		

\\MARTINOS DEVELOPER\HUBER\FRISGO_functional\FRISGO_20251106_COD_M1_2nd\AAhead_scout

TA: 18 sec Coil Selection: Manual Voxel Size: 1.6×1.6×1.6 mm³ Acc:: 2 Rel. SNR: 1.00**Properties**

Start measurement without further preparation	Off
Wait for User to Start	Off
Start measurements	Single Measurement
Prio Recon	Off
Auto Open Inline Display	Off
Auto Close Inline Display	Off
Load Images to MR View&GO	On
Auto Store Images	On
Disable auto transfer to PACS	Off
Load Images to Stamp Segments	Off
Load Images to Graphic Segments	On
Graphic segment	Default
Inline Movie	Off

Routine

Slab Group	1
Slabs	1
Distance Factor	20 %
Position	L0.0 A30.0 H0.0 mm
Orientation	Sagittal
Phase Encoding Dir.	A >> P
Slices per Slab	128
Phase Oversampling	0 %
Slice Oversampling	0.0 %
FOV Read	260 mm
FOV Phase	100.0 %
Slice Thickness	1.6 mm
TR	3.25 ms
TE	1.53 ms
Averages	1
Concatenations	1
AutoAlign	Head
Coil Elements	AC;VOL

Contrast - Common

TR	3.25 ms
TE	1.53 ms
Flip Angle	16 deg
Fat-Water Contrast	Standard
Contrasts	1
Reconstruction	Magnitude

Contrast - Dynamic

Dynamic Mode	Standard
Measurements	1
Time to Center	7.7 s

Resolution - Common

FOV Read	260 mm
----------	--------

Resolution - Common

FOV Phase	100.0 %
Slice Thickness	1.6 mm
Base Resolution	160
Phase Resolution	100 %
Slice Resolution	69 %
Trajectory	Cartesian

Resolution - Acceleration

Acceleration Mode	GRAPPA
Reference Scans	Integrated
Acceleration Factor PE	2
Reference Lines PE	24
Acceleration Factor 3D	1
Phase Partial Fourier	6/8
Slice Partial Fourier	6/8
Asymmetric Echo	Weak

Resolution - Filter

Raw Filter	Off
Elliptical Filter	Off
Distortion Correction	3D
Normalize	B1 Filter
Image Filter	Off

Geometry - Common

Slab Group	1
Slabs	1
Distance Factor	20 %
Position	L0.0 A30.0 H0.0 mm
Orientation	Sagittal
Phase Encoding Dir.	A >> P
Slices per Slab	128
Phase Oversampling	0 %
Slice Oversampling	0.0 %
FOV Read	260 mm
FOV Phase	100.0 %
Slice Thickness	1.6 mm
TR	3.25 ms
Multi-Slice Mode	Sequential
Series	Ascending
Concatenations	1

Geometry - AutoAlign

Slab Group	1
Position	L0.0 A30.0 H0.0 mm
Orientation	Sagittal
Phase Encoding Dir.	A >> P
AutoAlign	Head
Initial Position	Isocenter
L	0.0 mm

Geometry - AutoAlign

P	0.0 mm
H	0.0 mm
Initial Orientation	Transversal
Initial Rotation	0.00 deg

Geometry - Tim Planning Suite

Set-n-Go Protocol	Off
Table Position	0 mm
Table Position	H

System - Miscellaneous

Coil Selection	Manual
Radial Sorting	Off
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combination	Sum of Squares
Matrix Optimization	Off

System - Adjustments

Adjustment Strategy	Standard
B0 Shim	Tune up
B1 Shim	TrueForm
Adjustment Tolerance	Auto
Confirm Frequency	Never
Assume Silicone	Off

System - Adjust Volume

Position	Isocenter
Orientation	Transversal
Rotation	0.00 deg
A >> P	263 mm
R >> L	350 mm
F >> H	350 mm
Reset	Off

System - pTx

B1 Shim	TrueForm
Excitation	Non-sel.

System - Tx/Rx

Frequency 1H	297.117721 MHz
? Ref. Amplitude 1H	0.000 V
Reset	Off
Correction Factor	1.00
Image Scaling	1.000

Physio - PACE

Resp. Control	Off
Concatenations	1

Inline - Subtraction

Subtract	Off
----------	-----

Inline - Subtraction

Measurements	1
StdDev	Off
Save Original Images	On

Inline - MIP

MIP Sag	Off
MIP Cor	Off
MIP Tra	Off
MIP Time	Off
Radial MIP	Off
Save Original Images	On
MPR Sag	Off
MPR Cor	Off
MPR Tra	Off

Inline - Composing**Inline - MapIt**

MapIt	None
Flip Angle	16 deg
Measurements	1
Contrasts	1
TE	1.53 ms
TR	3.25 ms
Save Original Images	On

Sequence - Part 1

Sequence Name	fl
Dimension	3D
Excitation	Non-sel.
RF Pulse Type	Fast
Gradient Mode	Normal
Bandwidth	540 Hz/Px
Asymmetric Echo	Weak

Sequence - Part 2

Introduction	On
RF Spoiling	On
Breast Application	Off
Phase Enc. Order	Automatic

Sequence - Assistant

SAR Assistant	Off
---------------	-----

\\MARTINOS DEVELOPER\HUBER\FRISGO_functional\FRISGO_20251106_COD_M1_2nd\wip_ShimToolbox_B0map_3mm

TA: 44 sec Coil Selection: Auto Voxel Size: 3.0x3.0x3.0 mm³ Acc:: None Rel. SNR: 1.00

Properties

Start measurement without further preparation	On
Wait for User to Start	Off
Start measurements	Single Measurement
Prio Recon	Off
Auto Open Inline Display	Off
Auto Close Inline Display	Off
Load Images to MR View&GO	On
Auto Store Images	On
Disable auto transfer to PACS	Off
Load Images to Stamp Segments	Off
Load Images to Graphic Segments	Off
Graphic segment	Default
Inline Movie	Off

Routine

Slab Group	1
Slabs	1
Distance Factor	20 %
Position	Isocenter
Orientation	Sagittal
Phase Encoding Dir.	A >> P
Slices per Slab	52
Phase Oversampling	0 %
Slice Oversampling	0.0 %
FOV Read	216 mm
FOV Phase	100.0 %
Slice Thickness	3.00 mm
TR	15.0 ms
TE 1	2.20 ms
TE 2	5.40 ms
Concatenations	1
AutoAlign	---

Contrast - Common

TR	15.0 ms
TE 1	2.20 ms
TE 2	5.40 ms
MTC	Off
Magn. Preparation	None
Flip Angle	8 deg
Fat-Water Contrast	Water Excitation
Dark Blood	Off
Contrasts	2
SWI	Off
Reconstruction	Magn./Phase

Contrast - Dynamic

Dynamic Mode	Standard
--------------	----------

Contrast - Dynamic

Measurements	1
Multiple Series	Each Measurement

Resolution - Common

FOV Read	216 mm
FOV Phase	100.0 %
Slice Thickness	3.00 mm
Base Resolution	72
Phase Resolution	100 %
Slice Resolution	100 %

Resolution - Acceleration

Acceleration Mode	None
Phase Partial Fourier	Off
Slice Partial Fourier	Off
Asymmetric Echo	Off
Elliptical Scanning	On

Resolution - Filter

Raw Filter	Off
Elliptical Filter	Off
Distortion Correction	2D
Normalize	Off
Image Filter	Off

Geometry - Common

Slab Group	1
Slabs	1
Distance Factor	20 %
Position	Isocenter
Orientation	Sagittal
Phase Encoding Dir.	A >> P
Slices per Slab	52
Phase Oversampling	0 %
Slice Oversampling	0.0 %
FOV Read	216 mm
FOV Phase	100.0 %
Slice Thickness	3.00 mm
TR	15.0 ms
Multi-Slice Mode	Interleaved
Series	Interleaved
Concatenations	1

Geometry - AutoAlign

Slab Group	1
Position	Isocenter
Orientation	Sagittal
Phase Encoding Dir.	A >> P
AutoAlign	---

Geometry - AutoAlign

Initial Position	Isocenter
L	0.0 mm
P	0.0 mm
H	0.0 mm
Initial Orientation	Sagittal
Initial Rotation	0.00 deg

Geometry - Saturation

Saturation Mode	Standard
Special Saturation	None

Geometry - Tim Planning Suite

Set-n-Go Protocol	Off
Table Position	0 mm
Table Position	H

System - Miscellaneous

Coil Selection	Auto Coil Select
Radial Sorting	Off
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combination	Adaptive Combine
Matrix Optimization	Off
Coil Focus	Flat

System - Adjustments

Adjustment Strategy	Standard
B0 Shim	Standard
B1 Shim	TrueForm
Adjustment Tolerance	Auto
Adjust with Body Coil	Off
Confirm Frequency	Never
Assume Silicone	Off

System - Adjust Volume

Position	Isocenter
Orientation	Sagittal
Rotation	0.00 deg
A >> P	216 mm
F >> H	216 mm
R >> L	156 mm
Reset	Off

System - pTx

B1 Shim	TrueForm
Excitation	Non-sel.
LR Balancing	Off

System - Tx/Rx

Frequency 1H	297.117721 MHz
? Ref. Amplitude 1H	0.000 V
Reset	Off

System - Tx/Rx

Image Scaling	1.000
---------------	-------

Physio - Signal

1st Signal/Mode	None
TR	15.0 ms
Segments	1
Concatenations	1

Physio - Cardiac

Tagging	None
Fat-Water Contrast	Water Excitation
Magn. Preparation	None
Dark Blood	Off
FOV Read	216 mm
FOV Phase	100.0 %
Phase Resolution	100 %

Physio - PACE

Resp. Control	Off
Concatenations	1

Inline - Liver

Liver Registration	Off
Save Original Images	On

Inline - Subtraction

Subtract	Off
Measurements	1
StdDev	Off
Save Original Images	On

Inline - MIP

MIP Sag	Off
MIP Cor	Off
MIP Tra	Off
MIP Time	Off
Radial MIP	Off
Save Original Images	On
MPR Sag	Off
MPR Cor	Off
MPR Tra	Off

Inline - Soft Tissue

Wash-in	Off
Wash-out	Off
TTP	Off
PEI	Off
MIP Time	Off
Measurements	1

Inline - Composing**Inline - MapIt**

MapIt	None
-------	------

Inline - MapIt

Flip Angle	8 deg
Measurements	1
Contrasts	2
TE 1	2.20 ms
TE 2	5.40 ms
TR	15.0 ms
Save Original Images	On

Sequence - Part 1

Sequence Name	WIP_GRE
Dimension	3D
Excitation	Non-sel.
RF Pulse Type	Normal
Readout Mode	Monopolar
Gradient Mode	Whisper
Flow Compensation 1	None
Flow Compensation 2	None
Bandwidth 1	500 Hz/Px
Bandwidth 2	500 Hz/Px
Asymmetric Echo	Off
Segments	1

Sequence - Part 2

Introduction	Off
RF Spoiling	On
Acoustic noise reduction	Off

Sequence - Nuclei

TX/RX Nucleus	1H
TX/RX Delta Frequency	0 Hz
TX Nucleus	None
TX Delta Frequency	0 Hz

Sequence - Special

Mode	B0 Mapping
Masking	Brain Extraction
Phase Unwrapping	Spatial

Sequence - Assistant

SAR Assistant	Off
Allowed Delay	0 s

\\MARTINOS DEVELOPER\HUBER\FRISGO_functional\FRISGO_20251106_COD_M1_2nd\Scout_alternative

TA: 13 sec Coil Selection: Manual Voxel Size: 1.6×1.6×1.6 mm³ Acc:: 3 Rel. SNR: 1.00

Properties

Start measurement without further preparation	Off
Wait for User to Start	Off
Start measurements	Single Measurement
Prio Recon	Off
Auto Open Inline Display	Off
Auto Close Inline Display	Off
Load Images to MR View&GO	On
Auto Store Images	On
Disable auto transfer to PACS	Off
Load Images to Stamp Segments	Off
Load Images to Graphic Segments	On
Graphic segment	2nd Segment
Inline Movie	Off

Routine

Slab Group	1
Slabs	1
Distance Factor	20 %
Position	L0.0 A16.0 H0.0 mm
Orientation	Sagittal
Phase Encoding Dir.	A >> P
Slices per Slab	128
Phase Oversampling	0 %
Slice Oversampling	0.0 %
FOV Read	260 mm
FOV Phase	100.0 %
Slice Thickness	1.600 mm
TR	3.5 ms
TE	1.51 ms
Averages	1
Concatenations	1
AutoAlign	---
Coil Elements	AC

Contrast - Common

TR	3.5 ms
TE	1.51 ms
MTC	Off
Magn. Preparation	None
Flip Angle	15 deg
Fat-Water Contrast	Standard
Dark Blood	Off
Contrasts	1
SWI	Off
Reconstruction	Magnitude

Contrast - Dynamic

Dynamic Mode	Standard
--------------	----------

Contrast - Dynamic

Measurements	1
Multiple Series	Each Measurement

Resolution - Common

FOV Read	260 mm
FOV Phase	100.0 %
Slice Thickness	1.600 mm
Base Resolution	160
Phase Resolution	100 %
Slice Resolution	69 %
Interpolation	Off

Resolution - Acceleration

Acceleration Mode	GRAPPA
Reference Scans	Integrated
Acceleration Factor PE	3
Reference Lines PE	24
Acceleration Factor 3D	1
Phase Partial Fourier	6/8
Slice Partial Fourier	6/8
Asymmetric Echo	Allowed
Elliptical Scanning	Off

Resolution - Filter

Raw Filter	Off
Elliptical Filter	Off
Distortion Correction	3D
Normalize	B1 Filter
Image Filter	Off

Geometry - Common

Slab Group	1
Slabs	1
Distance Factor	20 %
Position	L0.0 A16.0 H0.0 mm
Orientation	Sagittal
Phase Encoding Dir.	A >> P
Slices per Slab	128
Phase Oversampling	0 %
Slice Oversampling	0.0 %
FOV Read	260 mm
FOV Phase	100.0 %
Slice Thickness	1.600 mm
TR	3.5 ms
Multi-Slice Mode	Interleaved
Series	Interleaved
Concatenations	1

Geometry - AutoAlign

Slab Group	1
Position	L0.0 A16.0 H0.0 mm
Orientation	Sagittal
Phase Encoding Dir.	A >> P
AutoAlign	---
Initial Position	L0.0 A16.0 H0.0
L	0.0 mm
A	16.0 mm
H	0.0 mm
Initial Orientation	Sagittal
Initial Rotation	0.00 deg

Geometry - Saturation

Saturation Mode	Standard
Special Saturation	None

Geometry - Tim Planning Suite

Set-n-Go Protocol	Off
Table Position	0 mm
Table Position	H

System - Miscellaneous

Coil Selection	Manual
Radial Sorting	Off
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combination	Adaptive Combine
Matrix Optimization	Off

System - Adjustments

Adjustment Strategy	Standard
B0 Shim	Tune up
B1 Shim	TrueForm
Adjustment Tolerance	Auto
Confirm Frequency	Never
Assume Silicone	Off

System - Adjust Volume

Position	Isocenter
Orientation	Transversal
Rotation	0.00 deg
A >> P	263 mm
R >> L	350 mm
F >> H	350 mm
Reset	Off

System - pTx

B1 Shim	TrueForm
Excitation	Non-sel.

System - Tx/Rx

Frequency 1H	297.117721 MHz
--------------	----------------

System - Tx/Rx

? Ref. Amplitude 1H	0.000 V
Reset	Off
Correction Factor	1.00
Image Scaling	1.000

Physio - Signal

1st Signal/Mode	None
TR	3.5 ms
Segments	1
Concatenations	1

Physio - Cardiac

Tagging	None
Fat-Water Contrast	Standard
Magn. Preparation	None
Dark Blood	Off
FOV Read	260 mm
FOV Phase	100.0 %
Phase Resolution	100 %

Physio - PACE

Resp. Control	Off
Concatenations	1

Inline - Liver

Liver Registration	Off
Save Original Images	On

Inline - Subtraction

Subtract	Off
Measurements	1
StdDev	Off
Save Original Images	On

Inline - MIP

MIP Sag	Off
MIP Cor	Off
MIP Tra	Off
MIP Time	Off
Radial MIP	Off
Save Original Images	On
MPR Sag	Off
MPR Cor	Off
MPR Tra	Off

Inline - Soft Tissue

Wash-in	Off
Wash-out	Off
TTP	Off
PEI	Off
MIP Time	Off
Measurements	1

Inline - Composing

Inline - MapIt

MapIt	None
Flip Angle	15 deg
Measurements	1
Contrasts	1
TE	1.51 ms
TR	3.5 ms
Save Original Images	On

Sequence - Part 1

Sequence Name	fl
Dimension	3D
Excitation	Non-sel.
RF Pulse Type	Fast
Gradient Mode	Normal
Flow Compensation	None
Bandwidth	540 Hz/Px
Asymmetric Echo	Allowed
Segments	1

Sequence - Part 2

Introduction	Off
RF Spoiling	On
Acoustic noise reduction	Off

Sequence - Nuclei

TX/RX Nucleus	1H
TX/RX Delta Frequency	0 Hz
TX Nucleus	None
TX Delta Frequency	0 Hz
Coil Elements	AC

Sequence - Special

Readout polarity	Positive
Image processing	Standard
Apply echo spacing	Off
Echo spacing	0 us
Delta echo spacing	0 us
Dummy scans	0 ms
RF pulse duration	100 us
Gradient spoiling	Siemens
Gradient moment factor	1.00
Receiver gain mode	Siemens
Number of segments	1
Current segment	0
Lines before/after seg	0

Sequence - Assistant

SAR Assistant	Off
Allowed Delay	0 s

\\MARTINOS DEVELOPER\HUBER\FRISGO_functional\FRISGO_20251106_COD_M1_2nd\AAsetter--192.16
8.2.1--ep2d
TA: 2 sec Coil Selection: Auto Voxel Size: 3.1×3.1×5.0 mm ³ Acc:: None Rel. SNR: 1.00

Properties

Start measurement without further preparation	On
Wait for User to Start	Off
Start measurements	Single Measurement
Prio Recon	Off
Auto Open Inline Display	Off
Auto Close Inline Display	Off
Load Images to MR View&GO	On
Auto Store Images	On
Disable auto transfer to PACS	Off
Load Images to Stamp Segments	On
Load Images to Graphic Segments	On
Graphic segment	Default
Inline Movie	Off

Resolution - Common

Base Resolution	160
Phase Resolution	91 %
Interpolation	Off

Resolution - Acceleration

Acceleration Mode	None
Phase Partial Fourier	Off

Resolution - Filter

Raw Filter	Off
Elliptical Filter	Off
Hamming	Off
Distortion Correction	Off
Normalize	Off

Routine

Slice Group	1
Slices	2
Distance Factor	0 %
Position	Isocenter
Orientation	Transversal
Phase Encoding Dir.	A >> P
Phase Oversampling	0 %
FOV Read	500 mm
FOV Phase	20.0 %
Slice Thickness	5.0 mm
TR	53.0 ms
TE	15.00 ms
Averages	1
Concatenations	1
AutoAlign	---

Geometry - Common

Slice Group	1
Slices	2
Distance Factor	0 %
Position	Isocenter
Orientation	Transversal
Phase Encoding Dir.	A >> P
Phase Oversampling	0 %
FOV Read	500 mm
FOV Phase	20.0 %
Slice Thickness	5.0 mm
TR	53.0 ms
Multi-Slice Mode	Interleaved
Series	Interleaved
Concatenations	1

Contrast - Common

TR	53.0 ms
TE	15.00 ms
MTC	Off
Flip Angle	10 deg
Fat-Water Contrast	Standard
Reconstruction	Magnitude

Contrast - Dynamic

Dynamic Mode	Standard
Measurements	1
Delay in TR	0.00 ms

Resolution - Common

FOV Read	500 mm
FOV Phase	20.0 %
Slice Thickness	5.0 mm

Geometry - AutoAlign

Slice Group	1
Position	Isocenter
Orientation	Transversal
Phase Encoding Dir.	A >> P
AutoAlign	---
Initial Position	Isocenter
L	0.0 mm
A	0.0 mm
H	0.0 mm
Initial Orientation	Transversal
Initial Rotation	0.00 deg

Geometry - Saturation

Special Saturation	None
--------------------	------

Geometry - Tim Planning Suite

Set-n-Go Protocol	Off
-------------------	-----

Geometry - Tim Planning Suite

Table Position	0 mm
Table Position	H

System - Miscellaneous

Coil Selection	Auto Coil Select
Radial Sorting	Off
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combination	Adaptive Combine
Matrix Optimization	Off
Coil Focus	Flat

System - Adjustments

Adjustment Strategy	Standard
B0 Shim	Tune up
B1 Shim	TrueForm
Adjustment Tolerance	None
Adjust with Body Coil	Off
Confirm Frequency	Never
Assume Silicone	Off

System - Adjust Volume

! Position	L2.2 P0.0 H8.8 mm
! Orientation	T > S-2.9 > C0.1
! Rotation	0.65 deg
! A >> P	263 mm
! R >> L	350 mm
! F >> H	350 mm
Reset	Off

System - pTx

B1 Shim	TrueForm
Excitation	Standard

System - Tx/Rx

Frequency 1H	297.117721 MHz
? Ref. Amplitude 1H	0.000 V
Reset	Off
Image Scaling	1.000

Physio - Signal

1st Signal/Mode	None
TR	53.0 ms
Log Signals	Off
Concatenations	1

BOLD

GLM Statistics	Off
Ignore Meas. at Start	0
Ignore After Transition	0
Model Transition States	On
Temp. Highpass Filter	On

BOLD

Threshold	4.00
Paradigm Size	3
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Active
Motion Correction	Off
Spatial Filter	Off
Measurements	1
Delay in TR	0.00 ms

Sequence - Part 1

Sequence Name	epfid
Excitation	Standard
RF Pulse Type	Fast
Gradient Mode	Fast
Bandwidth	1838 Hz/Px
Echo Spacing	0.63 ms
Free Echo Spacing	Off
EPI Factor	29

Sequence - Part 2

Introduction	On
--------------	----

Sequence - Assistant

SAR Assistant	Off
---------------	-----

\\MARTINOS DEVELOPER\HUBER\FRISGO_functional\FRISGO_20251106_COD_M1_2nd\rslh_ep3d_0.35_ test

TA: 36 sec Coil Selection: Manual Voxel Size: 0.4x0.4x0.4 mm³ Acc:: 3 Rel. SNR: 1.00

Properties

Start measurement without further preparation	On
Wait for User to Start	Off
Start measurements	Single Measurement
Prio Recon	Off
Auto Open Inline Display	Off
Auto Close Inline Display	Off
Load Images to MR View&GO	On
Auto Store Images	On
Disable auto transfer to PACS	Off
Load Images to Stamp Segments	Off
Load Images to Graphic Segments	Off
Graphic segment	Default
Inline Movie	Off

Routine

Slab Group	1
Slabs	1
Position	R1.0 P15.4 H25.4 mm
Orientation	T > C-19.5 > S-0.3
Phase Encoding Dir.	A >> P
Slices per Slab	18
Phase Oversampling	0.0 %
Slice Oversampling	0.0 %
FOV Read	130 mm
FOV Phase	100.0 %
Slice Thickness	0.35 mm
TR	54.1 ms
Vol. TR	5193.6 ms
TE 1	17.60 ms
Averages	1
Multi-echo Shots	1
AutoAlign	Head > Brain
Coil Elements	AC

Contrast - Common

TR	54.1 ms
Vol. TR	5193.6 ms
TE 1	17.60 ms
Multi-echo spacing	39.80 ms
MTC	Off
Flip Angle	18 deg
Fat-Water Contrast	Fast Fat Saturation
Magn. Prep. Shots	1
Contrasts	1
Reconstruction	Magn./Phase

Contrast - Dynamic

Dynamic Mode	Standard
--------------	----------

Contrast - Dynamic

Measurements	4
Reordering	Linear

Resolution - Common

FOV Read	130 mm
FOV Phase	100.0 %
Slice Thickness	0.35 mm
Base Resolution	368
Phase Resolution	100 %
Slice Resolution	100 %
Interpolation	Off

Resolution - Acceleration

Acceleration Mode	CAIPIRINHA
Reference Scans	GRE/Separate
Acceleration Factor PE	1
Reference Lines PE	75
Acceleration Factor 3D	3
Reference Lines 3D	18
Reordering Shift 3D	1
Phase Partial Fourier	6/8
Slice Partial Fourier	Off

Resolution - Filter

Raw Filter	Off
Elliptical Filter	Off
Distortion Correction	Off
Normalize	Off
Image Filter	Off

Geometry - Common

Slab Group	1
Slabs	1
Position	R1.0 P15.4 H25.4 mm
Orientation	T > C-19.5 > S-0.3
Phase Encoding Dir.	A >> P
Slices per Slab	18
Phase Oversampling	0.0 %
Slice Oversampling	0.0 %
FOV Read	130 mm
FOV Phase	100.0 %
Slice Thickness	0.35 mm
TR	54.1 ms
Vol. TR	5193.6 ms
Multi-echo Shots	1

Geometry - AutoAlign

Slab Group	1
Position	R1.0 P15.4 H25.4 mm

Geometry - AutoAlign

Orientation	T > C-19.5 > S-0.3
Phase Encoding Dir.	A >> P
AutoAlign	Head > Brain
Initial Position	R1.6 P9.1 H28.4
R	1.6 mm
P	9.1 mm
H	28.4 mm
Initial Orientation	T > C
T > C	-18.10
> S	0.00
Initial Rotation	1.33 deg

Geometry - Saturation

Saturation Mode	Quick
-----------------	-------

Geometry - Tim Planning Suite

Set-n-Go Protocol	Off
Table Position	0 mm
Table Position	H

System - Miscellaneous

Coil Selection	Manual
Radial Sorting	Off
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combination	Adaptive Combine
Matrix Optimization	Off

System - Adjustments

Adjustment Strategy	Standard
B0 Shim	Brain
B1 Shim	TrueForm
Adjustment Tolerance	Auto
Confirm Frequency	Never
Assume Silicone	Off

System - Adjust Volume

I Position	R1.0 P14.7 H21.7 mm
I Orientation	T > C-19.5 > S-0.3
I Rotation	6.76 deg
! A >> P	130 mm
! R >> L	130 mm
! F >> H	29 mm
Reset	Off

System - pTx

B1 Shim	TrueForm
Excitation	Slab-sel.

System - Tx/Rx

Frequency 1H	297.117721 MHz
! Ref. Amplitude 1H	370.000 V

System - Tx/Rx

Reset	Off
Correction Factor	1.00
Image Scaling	1.000

Sequence - Part 1

Sequence Name	vaso a2d6453
Dimension	3D
Excitation	Slab-sel.
RF Pulse Type	Normal
Gradient Mode	Fast
Reordering	Linear
Bandwidth	590 Hz/Px
Echo Spacing	2.12 ms
Turbo Factor	96
Segmentation	16
EPI Factor	18

Sequence - Part 2

Introduction	On
RF Spoiling	On

Sequence - Special

RF duration	1540 us
RF time x BW	8
PAT ref. FA	3 deg
Fat sat. FA	110 deg
T1 (Ernst FA)	1200 ms
Invert PE	Off
Min. TE w/ PF	On
Ramp Sampling	On
Trigger per shot	Off
Noise image	Off
Relax spoilers	Off
Round up Vol. TR	Off
MT flip angle	500 deg
MT off-res.	2000 Hz
MT RF duration	10240 us
Custom Water Exc.	-none-
Phase Correction	per Series
Saturation RF	per Shot
EPI rise time factor	1.47
G. spoil dephasing[1]	0.0 pi
G. spoil dephasing[2]	4.0 pi
G. spoil dephasing[3]	2.0 pi
Modify Ice Config	On
G-factor map	Off
GRAPPA Regularization	5000 /10^6
Slab Scale	-10 %
RF spoil scheme	Conventional
Read polarity	Inverted RO

Sequence - Assistant

SAR Assistant	Off
---------------	-----

\\MARTINOS DEVELOPER\HUBER\FRISGO_functional\FRISGO_20251106_COD_M1_2nd\rslh_ep3d_0.35_run1

TA: 12:43 min Coil Selection: Manual Voxel Size: 0.4x0.4x0.4 mm³ Acc:: 3 Rel. SNR: 1.00

Properties

Start measurement without further preparation	On
Wait for User to Start	Off
Start measurements	Single Measurement
Prio Recon	Off
Auto Open Inline Display	Off
Auto Close Inline Display	Off
Load Images to MR View&GO	On
Auto Store Images	On
Disable auto transfer to PACS	Off
Load Images to Stamp Segments	Off
Load Images to Graphic Segments	Off
Graphic segment	Default
Inline Movie	Off

Routine

Slab Group	1
Slabs	1
Position	R1.0 P15.4 H25.4 mm
Orientation	T > C-19.5 > S-0.3
Phase Encoding Dir.	A >> P
Slices per Slab	18
Phase Oversampling	0.0 %
Slice Oversampling	0.0 %
FOV Read	130 mm
FOV Phase	100.0 %
Slice Thickness	0.35 mm
TR	54.1 ms
Vol. TR	5193.6 ms
TE 1	17.60 ms
Averages	1
Multi-echo Shots	1
AutoAlign	Head > Brain
Coil Elements	AC

Contrast - Common

TR	54.1 ms
Vol. TR	5193.6 ms
TE 1	17.60 ms
Multi-echo spacing	39.80 ms
MTC	Off
Flip Angle	18 deg
Fat-Water Contrast	Fast Fat Saturation
Magn. Prep. Shots	1
Contrasts	1
Reconstruction	Magn./Phase

Contrast - Dynamic

Dynamic Mode	Standard
--------------	----------

Contrast - Dynamic

Measurements	144
Reordering	Linear

Resolution - Common

FOV Read	130 mm
FOV Phase	100.0 %
Slice Thickness	0.35 mm
Base Resolution	368
Phase Resolution	100 %
Slice Resolution	100 %
Interpolation	Off

Resolution - Acceleration

Acceleration Mode	CAIPIRINHA
Reference Scans	GRE/Separate
Acceleration Factor PE	1
Reference Lines PE	75
Acceleration Factor 3D	3
Reference Lines 3D	18
Reordering Shift 3D	1
Phase Partial Fourier	6/8
Slice Partial Fourier	Off

Resolution - Filter

Raw Filter	Off
Elliptical Filter	Off
Distortion Correction	Off
Normalize	Off
Image Filter	Off

Geometry - Common

Slab Group	1
Slabs	1
Position	R1.0 P15.4 H25.4 mm
Orientation	T > C-19.5 > S-0.3
Phase Encoding Dir.	A >> P
Slices per Slab	18
Phase Oversampling	0.0 %
Slice Oversampling	0.0 %
FOV Read	130 mm
FOV Phase	100.0 %
Slice Thickness	0.35 mm
TR	54.1 ms
Vol. TR	5193.6 ms
Multi-echo Shots	1

Geometry - AutoAlign

Slab Group	1
Position	R1.0 P15.4 H25.4 mm

Geometry - AutoAlign

Orientation	T > C-19.5 > S-0.3
Phase Encoding Dir.	A >> P
AutoAlign	Head > Brain
Initial Position	R1.6 P9.1 H28.4
R	1.6 mm
P	9.1 mm
H	28.4 mm
Initial Orientation	T > C
T > C	-18.10
> S	0.00
Initial Rotation	1.33 deg

Geometry - Saturation

Saturation Mode	Quick
-----------------	-------

Geometry - Tim Planning Suite

Set-n-Go Protocol	Off
Table Position	0 mm
Table Position	H

System - Miscellaneous

Coil Selection	Manual
Radial Sorting	Off
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combination	Adaptive Combine
Matrix Optimization	Off

System - Adjustments

Adjustment Strategy	Standard
B0 Shim	Brain
B1 Shim	TrueForm
Adjustment Tolerance	Auto
Confirm Frequency	Never
Assume Silicone	Off

System - Adjust Volume

I Position	R1.0 P14.7 H21.7 mm
I Orientation	T > C-19.5 > S-0.3
I Rotation	6.76 deg
! A >> P	130 mm
! R >> L	130 mm
! F >> H	29 mm
Reset	Off

System - pTx

B1 Shim	TrueForm
Excitation	Slab-sel.

System - Tx/Rx

Frequency 1H	297.117721 MHz
! Ref. Amplitude 1H	370.000 V

System - Tx/Rx

Reset	Off
Correction Factor	1.00
Image Scaling	1.000

Sequence - Part 1

Sequence Name	vaso a2d6453
Dimension	3D
Excitation	Slab-sel.
RF Pulse Type	Normal
Gradient Mode	Fast
Reordering	Linear
Bandwidth	590 Hz/Px
Echo Spacing	2.12 ms
Turbo Factor	96
Segmentation	16
EPI Factor	18

Sequence - Part 2

Introduction	On
RF Spoiling	On

Sequence - Special

RF duration	1540 us
RF time x BW	8
PAT ref. FA	3 deg
Fat sat. FA	110 deg
T1 (Ernst FA)	1200 ms
Invert PE	Off
Min. TE w/ PF	On
Ramp Sampling	On
Trigger per shot	Off
Noise image	Off
Relax spoilers	Off
Round up Vol. TR	Off
MT flip angle	500 deg
MT off-res.	2000 Hz
MT RF duration	10240 us
Custom Water Exc.	-none-
Phase Correction	per Series
Saturation RF	per Shot
EPI rise time factor	1.47
G. spoil dephasing[1]	0.0 pi
G. spoil dephasing[2]	4.0 pi
G. spoil dephasing[3]	2.0 pi
Modify Ice Config	On
G-factor map	Off
GRAPPA Regularization	5000 /10^6
Slab Scale	-10 %
RF spoil scheme	Conventional
Read polarity	Inverted RO

Sequence - Assistant

SAR Assistant	Off
---------------	-----

\\MARTINOS DEVELOPER\HUBER\FRISGO_functional\FRISGO_20251106_COD_M1_2nd\Scout_alternative

TA: 13 sec Coil Selection: Manual Voxel Size: 1.6×1.6×1.6 mm³ Acc:: 3 Rel. SNR: 1.00

Properties

Start measurement without further preparation	Off
Wait for User to Start	Off
Start measurements	Single Measurement
Prio Recon	Off
Auto Open Inline Display	Off
Auto Close Inline Display	Off
Load Images to MR View&GO	On
Auto Store Images	On
Disable auto transfer to PACS	Off
Load Images to Stamp Segments	Off
Load Images to Graphic Segments	On
Graphic segment	2nd Segment
Inline Movie	Off

Routine

Slab Group	1
Slabs	1
Distance Factor	20 %
Position	L0.0 A16.0 H0.0 mm
Orientation	Sagittal
Phase Encoding Dir.	A >> P
Slices per Slab	128
Phase Oversampling	0 %
Slice Oversampling	0.0 %
FOV Read	260 mm
FOV Phase	100.0 %
Slice Thickness	1.600 mm
TR	3.5 ms
TE	1.51 ms
Averages	1
Concatenations	1
AutoAlign	---
Coil Elements	AC

Contrast - Common

TR	3.5 ms
TE	1.51 ms
MTC	Off
Magn. Preparation	None
Flip Angle	15 deg
Fat-Water Contrast	Standard
Dark Blood	Off
Contrasts	1
SWI	Off
Reconstruction	Magnitude

Contrast - Dynamic

Dynamic Mode	Standard
--------------	----------

Contrast - Dynamic

Measurements	1
Multiple Series	Each Measurement

Resolution - Common

FOV Read	260 mm
FOV Phase	100.0 %
Slice Thickness	1.600 mm
Base Resolution	160
Phase Resolution	100 %
Slice Resolution	69 %
Interpolation	Off

Resolution - Acceleration

Acceleration Mode	GRAPPA
Reference Scans	Integrated
Acceleration Factor PE	3
Reference Lines PE	24
Acceleration Factor 3D	1
Phase Partial Fourier	6/8
Slice Partial Fourier	6/8
Asymmetric Echo	Allowed
Elliptical Scanning	Off

Resolution - Filter

Raw Filter	Off
Elliptical Filter	Off
Distortion Correction	3D
Normalize	B1 Filter
Image Filter	Off

Geometry - Common

Slab Group	1
Slabs	1
Distance Factor	20 %
Position	L0.0 A16.0 H0.0 mm
Orientation	Sagittal
Phase Encoding Dir.	A >> P
Slices per Slab	128
Phase Oversampling	0 %
Slice Oversampling	0.0 %
FOV Read	260 mm
FOV Phase	100.0 %
Slice Thickness	1.600 mm
TR	3.5 ms
Multi-Slice Mode	Interleaved
Series	Interleaved
Concatenations	1

Geometry - AutoAlign

Slab Group	1
Position	L0.0 A16.0 H0.0 mm
Orientation	Sagittal
Phase Encoding Dir.	A >> P
AutoAlign	---
Initial Position	L0.0 A16.0 H0.0
L	0.0 mm
A	16.0 mm
H	0.0 mm
Initial Orientation	Sagittal
Initial Rotation	0.00 deg

Geometry - Saturation

Saturation Mode	Standard
Special Saturation	None

Geometry - Tim Planning Suite

Set-n-Go Protocol	Off
Table Position	0 mm
Table Position	H

System - Miscellaneous

Coil Selection	Manual
Radial Sorting	Off
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combination	Adaptive Combine
Matrix Optimization	Off

System - Adjustments

Adjustment Strategy	Standard
B0 Shim	Tune up
B1 Shim	TrueForm
Adjustment Tolerance	Auto
Confirm Frequency	Never
Assume Silicone	Off

System - Adjust Volume

Position	Isocenter
Orientation	Transversal
Rotation	0.00 deg
A >> P	263 mm
R >> L	350 mm
F >> H	350 mm
Reset	Off

System - pTx

B1 Shim	TrueForm
Excitation	Non-sel.

System - Tx/Rx

Frequency 1H	297.117721 MHz
--------------	----------------

System - Tx/Rx

? Ref. Amplitude 1H	0.000 V
Reset	Off
Correction Factor	1.00
Image Scaling	1.000

Physio - Signal

1st Signal/Mode	None
TR	3.5 ms
Segments	1
Concatenations	1

Physio - Cardiac

Tagging	None
Fat-Water Contrast	Standard
Magn. Preparation	None
Dark Blood	Off
FOV Read	260 mm
FOV Phase	100.0 %
Phase Resolution	100 %

Physio - PACE

Resp. Control	Off
Concatenations	1

Inline - Liver

Liver Registration	Off
Save Original Images	On

Inline - Subtraction

Subtract	Off
Measurements	1
StdDev	Off
Save Original Images	On

Inline - MIP

MIP Sag	Off
MIP Cor	Off
MIP Tra	Off
MIP Time	Off
Radial MIP	Off
Save Original Images	On
MPR Sag	Off
MPR Cor	Off
MPR Tra	Off

Inline - Soft Tissue

Wash-in	Off
Wash-out	Off
TTP	Off
PEI	Off
MIP Time	Off
Measurements	1

Inline - Composing

Inline - MapIt

MapIt	None
Flip Angle	15 deg
Measurements	1
Contrasts	1
TE	1.51 ms
TR	3.5 ms
Save Original Images	On

Sequence - Part 1

Sequence Name	fl
Dimension	3D
Excitation	Non-sel.
RF Pulse Type	Fast
Gradient Mode	Normal
Flow Compensation	None
Bandwidth	540 Hz/Px
Asymmetric Echo	Allowed
Segments	1

Sequence - Part 2

Introduction	Off
RF Spoiling	On
Acoustic noise reduction	Off

Sequence - Nuclei

TX/RX Nucleus	1H
TX/RX Delta Frequency	0 Hz
TX Nucleus	None
TX Delta Frequency	0 Hz
Coil Elements	AC

Sequence - Special

Readout polarity	Positive
Image processing	Standard
Apply echo spacing	Off
Echo spacing	0 us
Delta echo spacing	0 us
Dummy scans	0 ms
RF pulse duration	100 us
Gradient spoiling	Siemens
Gradient moment factor	1.00
Receiver gain mode	Siemens
Number of segments	1
Current segment	0
Lines before/after seg	0

Sequence - Assistant

SAR Assistant	Off
Allowed Delay	0 s

\\MARTINOS DEVELOPER\HUBER\FRISGO_functional\FRISGO_20251106_COD_M1_2nd\AAsetter--192.16
8.2.1--ep2d

TA: 2 sec Coil Selection: Auto Voxel Size: 3.1×3.1×5.0 mm³ Acc:: None Rel. SNR: 1.00

Properties

Start measurement without further preparation	On
Wait for User to Start	Off
Start measurements	Single Measurement
Prio Recon	Off
Auto Open Inline Display	Off
Auto Close Inline Display	Off
Load Images to MR View&GO	On
Auto Store Images	On
Disable auto transfer to PACS	Off
Load Images to Stamp Segments	On
Load Images to Graphic Segments	On
Graphic segment	Default
Inline Movie	Off

Resolution - Common

Base Resolution	160
Phase Resolution	91 %
Interpolation	Off

Resolution - Acceleration

Acceleration Mode	None
Phase Partial Fourier	Off

Resolution - Filter

Raw Filter	Off
Elliptical Filter	Off
Hamming	Off
Distortion Correction	Off
Normalize	Off

Routine

Slice Group	1
Slices	2
Distance Factor	0 %
Position	Isocenter
Orientation	Transversal
Phase Encoding Dir.	A >> P
Phase Oversampling	0 %
FOV Read	500 mm
FOV Phase	20.0 %
Slice Thickness	5.0 mm
TR	53.0 ms
TE	15.00 ms
Averages	1
Concatenations	1
AutoAlign	---

Geometry - Common

Slice Group	1
Slices	2
Distance Factor	0 %
Position	Isocenter
Orientation	Transversal
Phase Encoding Dir.	A >> P
Phase Oversampling	0 %
FOV Read	500 mm
FOV Phase	20.0 %
Slice Thickness	5.0 mm
TR	53.0 ms
Multi-Slice Mode	Interleaved
Series	Interleaved
Concatenations	1

Contrast - Common

TR	53.0 ms
TE	15.00 ms
MTC	Off
Flip Angle	10 deg
Fat-Water Contrast	Standard
Reconstruction	Magnitude

Contrast - Dynamic

Dynamic Mode	Standard
Measurements	1
Delay in TR	0.00 ms

Resolution - Common

FOV Read	500 mm
FOV Phase	20.0 %
Slice Thickness	5.0 mm

Geometry - AutoAlign

Slice Group	1
Position	Isocenter
Orientation	Transversal
Phase Encoding Dir.	A >> P
AutoAlign	---
Initial Position	Isocenter
L	0.0 mm
A	0.0 mm
H	0.0 mm
Initial Orientation	Transversal
Initial Rotation	0.00 deg

Geometry - Saturation

Special Saturation	None
--------------------	------

Geometry - Tim Planning Suite

Set-n-Go Protocol	Off
-------------------	-----

Geometry - Tim Planning Suite

Table Position	0 mm
Table Position	H

System - Miscellaneous

Coil Selection	Auto Coil Select
Radial Sorting	Off
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combination	Adaptive Combine
Matrix Optimization	Off
Coil Focus	Flat

System - Adjustments

Adjustment Strategy	Standard
B0 Shim	Tune up
B1 Shim	TrueForm
Adjustment Tolerance	None
Adjust with Body Coil	Off
Confirm Frequency	Never
Assume Silicone	Off

System - Adjust Volume

! Position	L2.2 P0.0 H8.8 mm
! Orientation	T > S-2.9 > C0.1
! Rotation	0.65 deg
! A >> P	263 mm
! R >> L	350 mm
! F >> H	350 mm
Reset	Off

System - pTx

B1 Shim	TrueForm
Excitation	Standard

System - Tx/Rx

Frequency 1H	297.117721 MHz
? Ref. Amplitude 1H	0.000 V
Reset	Off
Image Scaling	1.000

Physio - Signal

1st Signal/Mode	None
TR	53.0 ms
Log Signals	Off
Concatenations	1

BOLD

GLM Statistics	Off
Ignore Meas. at Start	0
Ignore After Transition	0
Model Transition States	On
Temp. Highpass Filter	On

BOLD

Threshold	4.00
Paradigm Size	3
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Active
Motion Correction	Off
Spatial Filter	Off
Measurements	1
Delay in TR	0.00 ms

Sequence - Part 1

Sequence Name	epfid
Excitation	Standard
RF Pulse Type	Fast
Gradient Mode	Fast
Bandwidth	1838 Hz/Px
Echo Spacing	0.63 ms
Free Echo Spacing	Off
EPI Factor	29

Sequence - Part 2

Introduction	On
--------------	----

Sequence - Assistant

SAR Assistant	Off
---------------	-----

\\MARTINOS DEVELOPER\HUBER\FRISGO_functional\FRISGO_20251106_COD_M1_2nd\rslh_ep3d_0.35_run2_RO

TA: 12:43 min Coil Selection: Manual Voxel Size: 0.4x0.4x0.4 mm³ Acc:: 3 Rel. SNR: 1.00

Properties

Start measurement without further preparation	On
Wait for User to Start	Off
Start measurements	Single Measurement
Prio Recon	Off
Auto Open Inline Display	Off
Auto Close Inline Display	Off
Load Images to MR View&GO	On
Auto Store Images	On
Disable auto transfer to PACS	Off
Load Images to Stamp Segments	Off
Load Images to Graphic Segments	Off
Graphic segment	Default
Inline Movie	Off

Routine

Slab Group	1
Slabs	1
Position	R0.4 P15.8 H27.4 mm
Orientation	T > C-19.3 > S-0.5
Phase Encoding Dir.	A >> P
Slices per Slab	18
Phase Oversampling	0.0 %
Slice Oversampling	0.0 %
FOV Read	130 mm
FOV Phase	100.0 %
Slice Thickness	0.35 mm
TR	54.1 ms
Vol. TR	5193.6 ms
TE 1	17.60 ms
Averages	1
Multi-echo Shots	1
AutoAlign	Head > Brain
Coil Elements	AC

Contrast - Common

TR	54.1 ms
Vol. TR	5193.6 ms
TE 1	17.60 ms
Multi-echo spacing	39.80 ms
MTC	Off
Flip Angle	18 deg
Fat-Water Contrast	Fast Fat Saturation
Magn. Prep. Shots	1
Contrasts	1
Reconstruction	Magn./Phase

Contrast - Dynamic

Dynamic Mode	Standard
--------------	----------

Contrast - Dynamic

Measurements	144
Reordering	Linear

Resolution - Common

FOV Read	130 mm
FOV Phase	100.0 %
Slice Thickness	0.35 mm
Base Resolution	368
Phase Resolution	100 %
Slice Resolution	100 %
Interpolation	Off

Resolution - Acceleration

Acceleration Mode	CAIPIRINHA
Reference Scans	GRE/Separate
Acceleration Factor PE	1
Reference Lines PE	75
Acceleration Factor 3D	3
Reference Lines 3D	18
Reordering Shift 3D	1
Phase Partial Fourier	6/8
Slice Partial Fourier	Off

Resolution - Filter

Raw Filter	Off
Elliptical Filter	Off
Distortion Correction	Off
Normalize	Off
Image Filter	Off

Geometry - Common

Slab Group	1
Slabs	1
Position	R0.4 P15.8 H27.4 mm
Orientation	T > C-19.3 > S-0.5
Phase Encoding Dir.	A >> P
Slices per Slab	18
Phase Oversampling	0.0 %
Slice Oversampling	0.0 %
FOV Read	130 mm
FOV Phase	100.0 %
Slice Thickness	0.35 mm
TR	54.1 ms
Vol. TR	5193.6 ms
Multi-echo Shots	1

Geometry - AutoAlign

Slab Group	1
Position	R0.4 P15.8 H27.4 mm

Geometry - AutoAlign

Orientation	T > C-19.3 > S-0.5
Phase Encoding Dir.	A >> P
AutoAlign	Head > Brain
Initial Position	R1.6 P9.1 H28.4
R	1.6 mm
P	9.1 mm
H	28.4 mm
Initial Orientation	T > C
T > C	-18.10
> S	0.00
Initial Rotation	1.33 deg

Geometry - Saturation

Saturation Mode	Quick
-----------------	-------

Geometry - Tim Planning Suite

Set-n-Go Protocol	Off
Table Position	0 mm
Table Position	H

System - Miscellaneous

Coil Selection	Manual
Radial Sorting	Off
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combination	Adaptive Combine
Matrix Optimization	Off

System - Adjustments

Adjustment Strategy	Standard
B0 Shim	Brain
B1 Shim	TrueForm
Adjustment Tolerance	Auto
Confirm Frequency	Never
Assume Silicone	Off

System - Adjust Volume

I Position	R0.4 P15.1 H23.7 mm
I Orientation	T > C-19.3 > S-0.5
I Rotation	6.32 deg
! A >> P	130 mm
! R >> L	130 mm
! F >> H	29 mm
Reset	Off

System - pTx

B1 Shim	TrueForm
Excitation	Slab-sel.

System - Tx/Rx

Frequency 1H	297.117721 MHz
! Ref. Amplitude 1H	370.000 V

System - Tx/Rx

Reset	Off
Correction Factor	1.00
Image Scaling	1.000

Sequence - Part 1

Sequence Name	vaso a2d6453
Dimension	3D
Excitation	Slab-sel.
RF Pulse Type	Normal
Gradient Mode	Fast
Reordering	Linear
Bandwidth	590 Hz/Px
Echo Spacing	2.12 ms
Turbo Factor	96
Segmentation	16
EPI Factor	18

Sequence - Part 2

Introduction	On
RF Spoiling	On

Sequence - Special

RF duration	1540 us
RF time x BW	8
PAT ref. FA	3 deg
Fat sat. FA	110 deg
T1 (Ernst FA)	1200 ms
Invert PE	Off
Min. TE w/ PF	On
Ramp Sampling	On
Trigger per shot	Off
Noise image	Off
Relax spoilers	Off
Round up Vol. TR	Off
MT flip angle	500 deg
MT off-res.	2000 Hz
MT RF duration	10240 us
Custom Water Exc.	-none-
Phase Correction	per Series
Saturation RF	per Shot
EPI rise time factor	1.47
G. spoil dephasing[1]	0.0 pi
G. spoil dephasing[2]	4.0 pi
G. spoil dephasing[3]	2.0 pi
Modify Ice Config	On
G-factor map	Off
GRAPPA Regularization	5000 /10^6
Slab Scale	-10 %
RF spoil scheme	Conventional
Read polarity	Regular RO

Sequence - Assistant

SAR Assistant	Off
---------------	-----

\\MARTINOS DEVELOPER\HUBER\FRISGO_functional\FRISGO_20251106_COD_M1_2nd\Scout_alternativ
e

TA: 13 sec Coil Selection: Manual Voxel Size: 1.6×1.6×1.6 mm³ Acc:: 3 Rel. SNR: 1.00

Properties

Start measurement without further preparation	Off
Wait for User to Start	Off
Start measurements	Single Measurement
Prio Recon	Off
Auto Open Inline Display	Off
Auto Close Inline Display	Off
Load Images to MR View&GO	On
Auto Store Images	On
Disable auto transfer to PACS	Off
Load Images to Stamp Segments	Off
Load Images to Graphic Segments	On
Graphic segment	2nd Segment
Inline Movie	Off

Routine

Slab Group	1
Slabs	1
Distance Factor	20 %
Position	L0.0 A16.0 H0.0 mm
Orientation	Sagittal
Phase Encoding Dir.	A >> P
Slices per Slab	128
Phase Oversampling	0 %
Slice Oversampling	0.0 %
FOV Read	260 mm
FOV Phase	100.0 %
Slice Thickness	1.600 mm
TR	3.5 ms
TE	1.51 ms
Averages	1
Concatenations	1
AutoAlign	---
Coil Elements	AC

Contrast - Common

TR	3.5 ms
TE	1.51 ms
MTC	Off
Magn. Preparation	None
Flip Angle	15 deg
Fat-Water Contrast	Standard
Dark Blood	Off
Contrasts	1
SWI	Off
Reconstruction	Magnitude

Contrast - Dynamic

Dynamic Mode	Standard
--------------	----------

Contrast - Dynamic

Measurements	1
Multiple Series	Each Measurement

Resolution - Common

FOV Read	260 mm
FOV Phase	100.0 %
Slice Thickness	1.600 mm
Base Resolution	160
Phase Resolution	100 %
Slice Resolution	69 %
Interpolation	Off

Resolution - Acceleration

Acceleration Mode	GRAPPA
Reference Scans	Integrated
Acceleration Factor PE	3
Reference Lines PE	24
Acceleration Factor 3D	1
Phase Partial Fourier	6/8
Slice Partial Fourier	6/8
Asymmetric Echo	Allowed
Elliptical Scanning	Off

Resolution - Filter

Raw Filter	Off
Elliptical Filter	Off
Distortion Correction	3D
Normalize	B1 Filter
Image Filter	Off

Geometry - Common

Slab Group	1
Slabs	1
Distance Factor	20 %
Position	L0.0 A16.0 H0.0 mm
Orientation	Sagittal
Phase Encoding Dir.	A >> P
Slices per Slab	128
Phase Oversampling	0 %
Slice Oversampling	0.0 %
FOV Read	260 mm
FOV Phase	100.0 %
Slice Thickness	1.600 mm
TR	3.5 ms
Multi-Slice Mode	Interleaved
Series	Interleaved
Concatenations	1

Geometry - AutoAlign

Slab Group	1
Position	L0.0 A16.0 H0.0 mm
Orientation	Sagittal
Phase Encoding Dir.	A >> P
AutoAlign	---
Initial Position	L0.0 A16.0 H0.0
L	0.0 mm
A	16.0 mm
H	0.0 mm
Initial Orientation	Sagittal
Initial Rotation	0.00 deg

Geometry - Saturation

Saturation Mode	Standard
Special Saturation	None

Geometry - Tim Planning Suite

Set-n-Go Protocol	Off
Table Position	0 mm
Table Position	H

System - Miscellaneous

Coil Selection	Manual
Radial Sorting	Off
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combination	Adaptive Combine
Matrix Optimization	Off

System - Adjustments

Adjustment Strategy	Standard
B0 Shim	Tune up
B1 Shim	TrueForm
Adjustment Tolerance	Auto
Confirm Frequency	Never
Assume Silicone	Off

System - Adjust Volume

Position	Isocenter
Orientation	Transversal
Rotation	0.00 deg
A >> P	263 mm
R >> L	350 mm
F >> H	350 mm
Reset	Off

System - pTx

B1 Shim	TrueForm
Excitation	Non-sel.

System - Tx/Rx

Frequency 1H	297.117721 MHz
--------------	----------------

System - Tx/Rx

? Ref. Amplitude 1H	0.000 V
Reset	Off
Correction Factor	1.00
Image Scaling	1.000

Physio - Signal

1st Signal/Mode	None
TR	3.5 ms
Segments	1
Concatenations	1

Physio - Cardiac

Tagging	None
Fat-Water Contrast	Standard
Magn. Preparation	None
Dark Blood	Off
FOV Read	260 mm
FOV Phase	100.0 %
Phase Resolution	100 %

Physio - PACE

Resp. Control	Off
Concatenations	1

Inline - Liver

Liver Registration	Off
Save Original Images	On

Inline - Subtraction

Subtract	Off
Measurements	1
StdDev	Off
Save Original Images	On

Inline - MIP

MIP Sag	Off
MIP Cor	Off
MIP Tra	Off
MIP Time	Off
Radial MIP	Off
Save Original Images	On
MPR Sag	Off
MPR Cor	Off
MPR Tra	Off

Inline - Soft Tissue

Wash-in	Off
Wash-out	Off
TTP	Off
PEI	Off
MIP Time	Off
Measurements	1

Inline - Composing

Inline - MapIt

MapIt	None
Flip Angle	15 deg
Measurements	1
Contrasts	1
TE	1.51 ms
TR	3.5 ms
Save Original Images	On

Sequence - Part 1

Sequence Name	fl
Dimension	3D
Excitation	Non-sel.
RF Pulse Type	Fast
Gradient Mode	Normal
Flow Compensation	None
Bandwidth	540 Hz/Px
Asymmetric Echo	Allowed
Segments	1

Sequence - Part 2

Introduction	Off
RF Spoiling	On
Acoustic noise reduction	Off

Sequence - Nuclei

TX/RX Nucleus	1H
TX/RX Delta Frequency	0 Hz
TX Nucleus	None
TX Delta Frequency	0 Hz
Coil Elements	AC

Sequence - Special

Readout polarity	Positive
Image processing	Standard
Apply echo spacing	Off
Echo spacing	0 us
Delta echo spacing	0 us
Dummy scans	0 ms
RF pulse duration	100 us
Gradient spoiling	Siemens
Gradient moment factor	1.00
Receiver gain mode	Siemens
Number of segments	1
Current segment	0
Lines before/after seg	0

Sequence - Assistant

SAR Assistant	Off
Allowed Delay	0 s

\\MARTINOS DEVELOPER\HUBER\FRISGO_functional\FRISGO_20251106_COD_M1_2nd\AAsetter--192.16
8.2.1--ep2d
TA: 2 sec Coil Selection: Auto Voxel Size: 3.1×3.1×5.0 mm ³ Acc:: None Rel. SNR: 1.00

Properties

Start measurement without further preparation	On
Wait for User to Start	Off
Start measurements	Single Measurement
Prio Recon	Off
Auto Open Inline Display	Off
Auto Close Inline Display	Off
Load Images to MR View&GO	On
Auto Store Images	On
Disable auto transfer to PACS	Off
Load Images to Stamp Segments	On
Load Images to Graphic Segments	On
Graphic segment	Default
Inline Movie	Off

Resolution - Common

Base Resolution	160
Phase Resolution	91 %
Interpolation	Off

Resolution - Acceleration

Acceleration Mode	None
Phase Partial Fourier	Off

Resolution - Filter

Raw Filter	Off
Elliptical Filter	Off
Hamming	Off
Distortion Correction	Off
Normalize	Off

Routine

Slice Group	1
Slices	2
Distance Factor	0 %
Position	Isocenter
Orientation	Transversal
Phase Encoding Dir.	A >> P
Phase Oversampling	0 %
FOV Read	500 mm
FOV Phase	20.0 %
Slice Thickness	5.0 mm
TR	53.0 ms
TE	15.00 ms
Averages	1
Concatenations	1
AutoAlign	---

Geometry - Common

Slice Group	1
Slices	2
Distance Factor	0 %
Position	Isocenter
Orientation	Transversal
Phase Encoding Dir.	A >> P
Phase Oversampling	0 %
FOV Read	500 mm
FOV Phase	20.0 %
Slice Thickness	5.0 mm
TR	53.0 ms
Multi-Slice Mode	Interleaved
Series	Interleaved
Concatenations	1

Contrast - Common

TR	53.0 ms
TE	15.00 ms
MTC	Off
Flip Angle	10 deg
Fat-Water Contrast	Standard
Reconstruction	Magnitude

Contrast - Dynamic

Dynamic Mode	Standard
Measurements	1
Delay in TR	0.00 ms

Resolution - Common

FOV Read	500 mm
FOV Phase	20.0 %
Slice Thickness	5.0 mm

Geometry - AutoAlign

Slice Group	1
Position	Isocenter
Orientation	Transversal
Phase Encoding Dir.	A >> P
AutoAlign	---
Initial Position	Isocenter
L	0.0 mm
A	0.0 mm
H	0.0 mm
Initial Orientation	Transversal
Initial Rotation	0.00 deg

Geometry - Saturation

Special Saturation	None
--------------------	------

Geometry - Tim Planning Suite

Set-n-Go Protocol	Off
-------------------	-----

Geometry - Tim Planning Suite

Table Position	0 mm
Table Position	H

System - Miscellaneous

Coil Selection	Auto Coil Select
Radial Sorting	Off
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combination	Adaptive Combine
Matrix Optimization	Off
Coil Focus	Flat

System - Adjustments

Adjustment Strategy	Standard
B0 Shim	Tune up
B1 Shim	TrueForm
Adjustment Tolerance	None
Adjust with Body Coil	Off
Confirm Frequency	Never
Assume Silicone	Off

System - Adjust Volume

! Position	L2.2 P0.0 H8.8 mm
! Orientation	T > S-2.9 > C0.1
! Rotation	0.65 deg
! A >> P	263 mm
! R >> L	350 mm
! F >> H	350 mm
Reset	Off

System - pTx

B1 Shim	TrueForm
Excitation	Standard

System - Tx/Rx

Frequency 1H	297.117721 MHz
? Ref. Amplitude 1H	0.000 V
Reset	Off
Image Scaling	1.000

Physio - Signal

1st Signal/Mode	None
TR	53.0 ms
Log Signals	Off
Concatenations	1

BOLD

GLM Statistics	Off
Ignore Meas. at Start	0
Ignore After Transition	0
Model Transition States	On
Temp. Highpass Filter	On

BOLD

Threshold	4.00
Paradigm Size	3
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Active
Motion Correction	Off
Spatial Filter	Off
Measurements	1
Delay in TR	0.00 ms

Sequence - Part 1

Sequence Name	epfid
Excitation	Standard
RF Pulse Type	Fast
Gradient Mode	Fast
Bandwidth	1838 Hz/Px
Echo Spacing	0.63 ms
Free Echo Spacing	Off
EPI Factor	29

Sequence - Part 2

Introduction	On
--------------	----

Sequence - Assistant

SAR Assistant	Off
---------------	-----

\\MARTINOS DEVELOPER\HUBER\FRISGO_functional\FRISGO_20251106_COD_M1_2nd\rslh_ep3d_0.35_run3

TA: 12:43 min Coil Selection: Manual Voxel Size: 0.4x0.4x0.4 mm³ Acc:: 3 Rel. SNR: 1.00

Properties

Start measurement without further preparation	On
Wait for User to Start	Off
Start measurements	Single Measurement
Prio Recon	Off
Auto Open Inline Display	Off
Auto Close Inline Display	Off
Load Images to MR View&GO	On
Auto Store Images	On
Disable auto transfer to PACS	Off
Load Images to Stamp Segments	Off
Load Images to Graphic Segments	Off
Graphic segment	Default
Inline Movie	Off

Routine

Slab Group	1
Slabs	1
Position	R0.3 P16.5 H26.7 mm
Orientation	T > C-20.0 > S-0.6
Phase Encoding Dir.	A >> P
Slices per Slab	18
Phase Oversampling	0.0 %
Slice Oversampling	0.0 %
FOV Read	130 mm
FOV Phase	100.0 %
Slice Thickness	0.35 mm
TR	54.1 ms
Vol. TR	5193.6 ms
TE 1	17.60 ms
Averages	1
Multi-echo Shots	1
AutoAlign	Head > Brain
Coil Elements	AC

Contrast - Common

TR	54.1 ms
Vol. TR	5193.6 ms
TE 1	17.60 ms
Multi-echo spacing	39.80 ms
MTC	Off
Flip Angle	18 deg
Fat-Water Contrast	Fast Fat Saturation
Magn. Prep. Shots	1
Contrasts	1
Reconstruction	Magn./Phase

Contrast - Dynamic

Dynamic Mode	Standard
--------------	----------

Contrast - Dynamic

Measurements	144
Reordering	Linear

Resolution - Common

FOV Read	130 mm
FOV Phase	100.0 %
Slice Thickness	0.35 mm
Base Resolution	368
Phase Resolution	100 %
Slice Resolution	100 %
Interpolation	Off

Resolution - Acceleration

Acceleration Mode	CAIPIRINHA
Reference Scans	GRE/Separate
Acceleration Factor PE	1
Reference Lines PE	75
Acceleration Factor 3D	3
Reference Lines 3D	18
Reordering Shift 3D	1
Phase Partial Fourier	6/8
Slice Partial Fourier	Off

Resolution - Filter

Raw Filter	Off
Elliptical Filter	Off
Distortion Correction	Off
Normalize	Off
Image Filter	Off

Geometry - Common

Slab Group	1
Slabs	1
Position	R0.3 P16.5 H26.7 mm
Orientation	T > C-20.0 > S-0.6
Phase Encoding Dir.	A >> P
Slices per Slab	18
Phase Oversampling	0.0 %
Slice Oversampling	0.0 %
FOV Read	130 mm
FOV Phase	100.0 %
Slice Thickness	0.35 mm
TR	54.1 ms
Vol. TR	5193.6 ms
Multi-echo Shots	1

Geometry - AutoAlign

Slab Group	1
Position	R0.3 P16.5 H26.7 mm

Geometry - AutoAlign

Orientation	T > C-20.0 > S-0.6
Phase Encoding Dir.	A >> P
AutoAlign	Head > Brain
Initial Position	R1.6 P9.1 H28.4
R	1.6 mm
P	9.1 mm
H	28.4 mm
Initial Orientation	T > C
T > C	-18.10
> S	0.00
Initial Rotation	1.33 deg

Geometry - Saturation

Saturation Mode	Quick
-----------------	-------

Geometry - Tim Planning Suite

Set-n-Go Protocol	Off
Table Position	0 mm
Table Position	H

System - Miscellaneous

Coil Selection	Manual
Radial Sorting	Off
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combination	Adaptive Combine
Matrix Optimization	Off

System - Adjustments

Adjustment Strategy	Standard
B0 Shim	Brain
B1 Shim	TrueForm
Adjustment Tolerance	Auto
Confirm Frequency	Never
Assume Silicone	Off

System - Adjust Volume

I Position	R0.3 P15.8 H23.0 mm
I Orientation	T > C-20.0 > S-0.6
I Rotation	6.21 deg
! A >> P	130 mm
! R >> L	130 mm
! F >> H	29 mm
Reset	Off

System - pTx

B1 Shim	TrueForm
Excitation	Slab-sel.

System - Tx/Rx

Frequency 1H	297.117721 MHz
! Ref. Amplitude 1H	370.000 V

System - Tx/Rx

Reset	Off
Correction Factor	1.00
Image Scaling	1.000

Sequence - Part 1

Sequence Name	vaso a2d6453
Dimension	3D
Excitation	Slab-sel.
RF Pulse Type	Normal
Gradient Mode	Fast
Reordering	Linear
Bandwidth	590 Hz/Px
Echo Spacing	2.12 ms
Turbo Factor	96
Segmentation	16
EPI Factor	18

Sequence - Part 2

Introduction	On
RF Spoiling	On

Sequence - Special

RF duration	1540 us
RF time x BW	8
PAT ref. FA	3 deg
Fat sat. FA	110 deg
T1 (Ernst FA)	1200 ms
Invert PE	Off
Min. TE w/ PF	On
Ramp Sampling	On
Trigger per shot	Off
Noise image	Off
Relax spoilers	Off
Round up Vol. TR	Off
MT flip angle	500 deg
MT off-res.	2000 Hz
MT RF duration	10240 us
Custom Water Exc.	-none-
Phase Correction	per Series
Saturation RF	per Shot
EPI rise time factor	1.47
G. spoil dephasing[1]	0.0 pi
G. spoil dephasing[2]	4.0 pi
G. spoil dephasing[3]	2.0 pi
Modify Ice Config	On
G-factor map	Off
GRAPPA Regularization	5000 /10^6
Slab Scale	-10 %
RF spoil scheme	Conventional
Read polarity	Inverted RO

Sequence - Assistant

SAR Assistant	Off
---------------	-----

\\MARTINOS DEVELOPER\HUBER\FRISGO_functional\FRISGO_20251106_COD_M1_2nd\Scout_alternative

TA: 13 sec Coil Selection: Manual Voxel Size: 1.6×1.6×1.6 mm³ Acc:: 3 Rel. SNR: 1.00

Properties

Start measurement without further preparation	Off
Wait for User to Start	Off
Start measurements	Single Measurement
Prio Recon	Off
Auto Open Inline Display	Off
Auto Close Inline Display	Off
Load Images to MR View&GO	On
Auto Store Images	On
Disable auto transfer to PACS	Off
Load Images to Stamp Segments	Off
Load Images to Graphic Segments	On
Graphic segment	2nd Segment
Inline Movie	Off

Routine

Slab Group	1
Slabs	1
Distance Factor	20 %
Position	L0.0 A16.0 H0.0 mm
Orientation	Sagittal
Phase Encoding Dir.	A >> P
Slices per Slab	128
Phase Oversampling	0 %
Slice Oversampling	0.0 %
FOV Read	260 mm
FOV Phase	100.0 %
Slice Thickness	1.600 mm
TR	3.5 ms
TE	1.51 ms
Averages	1
Concatenations	1
AutoAlign	---
Coil Elements	AC

Contrast - Common

TR	3.5 ms
TE	1.51 ms
MTC	Off
Magn. Preparation	None
Flip Angle	15 deg
Fat-Water Contrast	Standard
Dark Blood	Off
Contrasts	1
SWI	Off
Reconstruction	Magnitude

Contrast - Dynamic

Dynamic Mode	Standard
--------------	----------

Contrast - Dynamic

Measurements	1
Multiple Series	Each Measurement

Resolution - Common

FOV Read	260 mm
FOV Phase	100.0 %
Slice Thickness	1.600 mm
Base Resolution	160
Phase Resolution	100 %
Slice Resolution	69 %
Interpolation	Off

Resolution - Acceleration

Acceleration Mode	GRAPPA
Reference Scans	Integrated
Acceleration Factor PE	3
Reference Lines PE	24
Acceleration Factor 3D	1
Phase Partial Fourier	6/8
Slice Partial Fourier	6/8
Asymmetric Echo	Allowed
Elliptical Scanning	Off

Resolution - Filter

Raw Filter	Off
Elliptical Filter	Off
Distortion Correction	3D
Normalize	B1 Filter
Image Filter	Off

Geometry - Common

Slab Group	1
Slabs	1
Distance Factor	20 %
Position	L0.0 A16.0 H0.0 mm
Orientation	Sagittal
Phase Encoding Dir.	A >> P
Slices per Slab	128
Phase Oversampling	0 %
Slice Oversampling	0.0 %
FOV Read	260 mm
FOV Phase	100.0 %
Slice Thickness	1.600 mm
TR	3.5 ms
Multi-Slice Mode	Interleaved
Series	Interleaved
Concatenations	1

Geometry - AutoAlign

Slab Group	1
Position	L0.0 A16.0 H0.0 mm
Orientation	Sagittal
Phase Encoding Dir.	A >> P
AutoAlign	---
Initial Position	L0.0 A16.0 H0.0
L	0.0 mm
A	16.0 mm
H	0.0 mm
Initial Orientation	Sagittal
Initial Rotation	0.00 deg

Geometry - Saturation

Saturation Mode	Standard
Special Saturation	None

Geometry - Tim Planning Suite

Set-n-Go Protocol	Off
Table Position	0 mm
Table Position	H

System - Miscellaneous

Coil Selection	Manual
Radial Sorting	Off
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combination	Adaptive Combine
Matrix Optimization	Off

System - Adjustments

Adjustment Strategy	Standard
B0 Shim	Tune up
B1 Shim	TrueForm
Adjustment Tolerance	Auto
Confirm Frequency	Never
Assume Silicone	Off

System - Adjust Volume

Position	Isocenter
Orientation	Transversal
Rotation	0.00 deg
A >> P	263 mm
R >> L	350 mm
F >> H	350 mm
Reset	Off

System - pTx

B1 Shim	TrueForm
Excitation	Non-sel.

System - Tx/Rx

Frequency 1H	297.117721 MHz
--------------	----------------

System - Tx/Rx

? Ref. Amplitude 1H	0.000 V
Reset	Off
Correction Factor	1.00
Image Scaling	1.000

Physio - Signal

1st Signal/Mode	None
TR	3.5 ms
Segments	1
Concatenations	1

Physio - Cardiac

Tagging	None
Fat-Water Contrast	Standard
Magn. Preparation	None
Dark Blood	Off
FOV Read	260 mm
FOV Phase	100.0 %
Phase Resolution	100 %

Physio - PACE

Resp. Control	Off
Concatenations	1

Inline - Liver

Liver Registration	Off
Save Original Images	On

Inline - Subtraction

Subtract	Off
Measurements	1
StdDev	Off
Save Original Images	On

Inline - MIP

MIP Sag	Off
MIP Cor	Off
MIP Tra	Off
MIP Time	Off
Radial MIP	Off
Save Original Images	On
MPR Sag	Off
MPR Cor	Off
MPR Tra	Off

Inline - Soft Tissue

Wash-in	Off
Wash-out	Off
TTP	Off
PEI	Off
MIP Time	Off
Measurements	1

Inline - Composing

Inline - MapIt

MapIt	None
Flip Angle	15 deg
Measurements	1
Contrasts	1
TE	1.51 ms
TR	3.5 ms
Save Original Images	On

Sequence - Part 1

Sequence Name	fl
Dimension	3D
Excitation	Non-sel.
RF Pulse Type	Fast
Gradient Mode	Normal
Flow Compensation	None
Bandwidth	540 Hz/Px
Asymmetric Echo	Allowed
Segments	1

Sequence - Part 2

Introduction	Off
RF Spoiling	On
Acoustic noise reduction	Off

Sequence - Nuclei

TX/RX Nucleus	1H
TX/RX Delta Frequency	0 Hz
TX Nucleus	None
TX Delta Frequency	0 Hz
Coil Elements	AC

Sequence - Special

Readout polarity	Positive
Image processing	Standard
Apply echo spacing	Off
Echo spacing	0 us
Delta echo spacing	0 us
Dummy scans	0 ms
RF pulse duration	100 us
Gradient spoiling	Siemens
Gradient moment factor	1.00
Receiver gain mode	Siemens
Number of segments	1
Current segment	0
Lines before/after seg	0

Sequence - Assistant

SAR Assistant	Off
Allowed Delay	0 s

\\MARTINOS DEVELOPER\HUBER\FRISGO_functional\FRISGO_20251106_COD_M1_2nd\AAsetter--192.16 8.2.1--ep2d

TA: 2 sec Coil Selection: Manual Voxel Size: 3.1×3.1×5.0 mm³ Acc:: None Rel. SNR: 1.00

Properties

Start measurement without further preparation	On
Wait for User to Start	Off
Start measurements	Single Measurement
Prio Recon	Off
Auto Open Inline Display	Off
Auto Close Inline Display	Off
Load Images to MR View&GO	On
Auto Store Images	On
Disable auto transfer to PACS	Off
Load Images to Stamp Segments	On
Load Images to Graphic Segments	On
Graphic segment	Default
Inline Movie	Off

Routine

Slice Group	1
Slices	2
Distance Factor	0 %
Position	Isocenter
Orientation	Transversal
Phase Encoding Dir.	A >> P
Phase Oversampling	0 %
FOV Read	500 mm
FOV Phase	20.0 %
Slice Thickness	5.0 mm
TR	53.0 ms
TE	15.00 ms
Averages	1
Concatenations	1
AutoAlign	---
Coil Elements	AC

Contrast - Common

TR	53.0 ms
TE	15.00 ms
MTC	Off
Flip Angle	10 deg
Fat-Water Contrast	Standard
Reconstruction	Magnitude

Contrast - Dynamic

Dynamic Mode	Standard
Measurements	1
Delay in TR	0.00 ms

Resolution - Common

FOV Read	500 mm
FOV Phase	20.0 %

Resolution - Common

Slice Thickness	5.0 mm
Base Resolution	160
Phase Resolution	91 %
Interpolation	Off

Resolution - Acceleration

Acceleration Mode	None
Phase Partial Fourier	Off

Resolution - Filter

Raw Filter	Off
Elliptical Filter	Off
Hamming	Off
Distortion Correction	Off
Normalize	Off

Geometry - Common

Slice Group	1
Slices	2
Distance Factor	0 %
Position	Isocenter
Orientation	Transversal
Phase Encoding Dir.	A >> P
Phase Oversampling	0 %
FOV Read	500 mm
FOV Phase	20.0 %
Slice Thickness	5.0 mm
TR	53.0 ms
Multi-Slice Mode	Interleaved
Series	Interleaved
Concatenations	1

Geometry - AutoAlign

Slice Group	1
Position	Isocenter
Orientation	Transversal
Phase Encoding Dir.	A >> P
AutoAlign	---
Initial Position	Isocenter
L	0.0 mm
A	0.0 mm
H	0.0 mm
Initial Orientation	Transversal
Initial Rotation	0.00 deg

Geometry - Saturation

Special Saturation	None
--------------------	------

Geometry - Tim Planning Suite

Set-n-Go Protocol	Off
Table Position	0 mm
Table Position	H

System - Miscellaneous

Coil Selection	Manual
Radial Sorting	Off
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combination	Adaptive Combine
Matrix Optimization	Off

System - Adjustments

Adjustment Strategy	Standard
B0 Shim	Tune up
B1 Shim	TrueForm
Adjustment Tolerance	None
Confirm Frequency	Never
Assume Silicone	Off

System - Adjust Volume

! Position	L2.2 P0.0 H8.8 mm
! Orientation	T > S-2.9 > C0.1
! Rotation	0.65 deg
! A >> P	263 mm
! R >> L	350 mm
! F >> H	350 mm
Reset	Off

System - pTx

B1 Shim	TrueForm
Excitation	Standard

System - Tx/Rx

Frequency 1H	297.117721 MHz
? Ref. Amplitude 1H	0.000 V
Reset	Off
Correction Factor	1.00
Image Scaling	1.000

Physio - Signal

1st Signal/Mode	None
TR	53.0 ms
Log Signals	Off
Concatenations	1

BOLD

GLM Statistics	Off
Ignore Meas. at Start	0
Ignore After Transition	0
Model Transition States	On
Temp. Highpass Filter	On

BOLD

Threshold	4.00
Paradigm Size	3
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Active
Motion Correction	Off
Spatial Filter	Off
Measurements	1
Delay in TR	0.00 ms

Sequence - Part 1

Sequence Name	epfid
Excitation	Standard
RF Pulse Type	Fast
Gradient Mode	Fast
Bandwidth	1838 Hz/Px
Echo Spacing	0.63 ms
Free Echo Spacing	Off
EPI Factor	29

Sequence - Part 2

Introduction	On
--------------	----

Sequence - Assistant

SAR Assistant	Off
---------------	-----

\\MARTINOS DEVELOPER\HUBER\FRISGO_functional\FRISGO_20251106_COD_M1_2nd\rslh_ep3d_0.35_run4_RO

TA: 12:43 min Coil Selection: Manual Voxel Size: 0.4x0.4x0.4 mm³ Acc:: 3 Rel. SNR: 1.00

Properties

Start measurement without further preparation	On
Wait for User to Start	Off
Start measurements	Single Measurement
Prio Recon	Off
Auto Open Inline Display	Off
Auto Close Inline Display	Off
Load Images to MR View&GO	On
Auto Store Images	On
Disable auto transfer to PACS	Off
Load Images to Stamp Segments	Off
Load Images to Graphic Segments	Off
Graphic segment	Default
Inline Movie	Off

Routine

Slab Group	1
Slabs	1
Position	R0.9 P16.8 H26.7 mm
Orientation	T > C-20.3
Phase Encoding Dir.	A >> P
Slices per Slab	18
Phase Oversampling	0.0 %
Slice Oversampling	0.0 %
FOV Read	130 mm
FOV Phase	100.0 %
Slice Thickness	0.35 mm
TR	54.1 ms
Vol. TR	5193.6 ms
TE 1	17.60 ms
Averages	1
Multi-echo Shots	1
AutoAlign	Head > Brain
Coil Elements	AC

Contrast - Common

TR	54.1 ms
Vol. TR	5193.6 ms
TE 1	17.60 ms
Multi-echo spacing	39.80 ms
MTC	Off
Flip Angle	18 deg
Fat-Water Contrast	Fast Fat Saturation
Magn. Prep. Shots	1
Contrasts	1
Reconstruction	Magn./Phase

Contrast - Dynamic

Dynamic Mode	Standard
--------------	----------

Contrast - Dynamic

Measurements	144
Reordering	Linear

Resolution - Common

FOV Read	130 mm
FOV Phase	100.0 %
Slice Thickness	0.35 mm
Base Resolution	368
Phase Resolution	100 %
Slice Resolution	100 %
Interpolation	Off

Resolution - Acceleration

Acceleration Mode	CAIPIRINHA
Reference Scans	GRE/Separate
Acceleration Factor PE	1
Reference Lines PE	75
Acceleration Factor 3D	3
Reference Lines 3D	18
Reordering Shift 3D	1
Phase Partial Fourier	6/8
Slice Partial Fourier	Off

Resolution - Filter

Raw Filter	Off
Elliptical Filter	Off
Distortion Correction	Off
Normalize	Off
Image Filter	Off

Geometry - Common

Slab Group	1
Slabs	1
Position	R0.9 P16.8 H26.7 mm
Orientation	T > C-20.3
Phase Encoding Dir.	A >> P
Slices per Slab	18
Phase Oversampling	0.0 %
Slice Oversampling	0.0 %
FOV Read	130 mm
FOV Phase	100.0 %
Slice Thickness	0.35 mm
TR	54.1 ms
Vol. TR	5193.6 ms
Multi-echo Shots	1

Geometry - AutoAlign

Slab Group	1
Position	R0.9 P16.8 H26.7 mm

Geometry - AutoAlign

Orientation	T > C-20.3
Phase Encoding Dir.	A >> P
AutoAlign	Head > Brain
Initial Position	R1.6 P9.1 H28.4
R	1.6 mm
P	9.1 mm
H	28.4 mm
Initial Orientation	T > C
T > C	-18.10
> S	0.00
Initial Rotation	1.33 deg

Geometry - Saturation

Saturation Mode	Quick
-----------------	-------

Geometry - Tim Planning Suite

Set-n-Go Protocol	Off
Table Position	0 mm
Table Position	H

System - Miscellaneous

Coil Selection	Manual
Radial Sorting	Off
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combination	Adaptive Combine
Matrix Optimization	Off

System - Adjustments

Adjustment Strategy	Standard
B0 Shim	Brain
B1 Shim	TrueForm
Adjustment Tolerance	Auto
Confirm Frequency	Never
Assume Silicone	Off

System - Adjust Volume

I Position	R0.8 P16.1 H23.0 mm
I Orientation	T > C-20.3
I Rotation	6.15 deg
! A >> P	130 mm
! R >> L	130 mm
! F >> H	29 mm
Reset	Off

System - pTx

B1 Shim	TrueForm
Excitation	Slab-sel.

System - Tx/Rx

Frequency 1H	297.117721 MHz
! Ref. Amplitude 1H	370.000 V

System - Tx/Rx

Reset	Off
Correction Factor	1.00
Image Scaling	1.000

Sequence - Part 1

Sequence Name	vaso a2d6453
Dimension	3D
Excitation	Slab-sel.
RF Pulse Type	Normal
Gradient Mode	Fast
Reordering	Linear
Bandwidth	590 Hz/Px
Echo Spacing	2.12 ms
Turbo Factor	96
Segmentation	16
EPI Factor	18

Sequence - Part 2

Introduction	On
RF Spoiling	On

Sequence - Special

RF duration	1540 us
RF time x BW	8
PAT ref. FA	3 deg
Fat sat. FA	110 deg
T1 (Ernst FA)	1200 ms
Invert PE	Off
Min. TE w/ PF	On
Ramp Sampling	On
Trigger per shot	Off
Noise image	Off
Relax spoilers	Off
Round up Vol. TR	Off
MT flip angle	500 deg
MT off-res.	2000 Hz
MT RF duration	10240 us
Custom Water Exc.	-none-
Phase Correction	per Series
Saturation RF	per Shot
EPI rise time factor	1.47
G. spoil dephasing[1]	0.0 pi
G. spoil dephasing[2]	4.0 pi
G. spoil dephasing[3]	2.0 pi
Modify Ice Config	On
G-factor map	Off
GRAPPA Regularization	5000 /10^6
Slab Scale	-10 %
RF spoil scheme	Conventional
Read polarity	Regular RO

Sequence - Assistant

SAR Assistant	Off
---------------	-----

\\MARTINOS DEVELOPER\HUBER\FRISGO_functional\FRISGO_20251106_COD_M1_2nd\Scout_alternative

TA: 13 sec Coil Selection: Manual Voxel Size: 1.6×1.6×1.6 mm³ Acc:: 3 Rel. SNR: 1.00

Properties

Start measurement without further preparation	Off
Wait for User to Start	Off
Start measurements	Single Measurement
Prio Recon	Off
Auto Open Inline Display	Off
Auto Close Inline Display	Off
Load Images to MR View&GO	On
Auto Store Images	On
Disable auto transfer to PACS	Off
Load Images to Stamp Segments	Off
Load Images to Graphic Segments	On
Graphic segment	2nd Segment
Inline Movie	Off

Routine

Slab Group	1
Slabs	1
Distance Factor	20 %
Position	L0.0 A16.0 H0.0 mm
Orientation	Sagittal
Phase Encoding Dir.	A >> P
Slices per Slab	128
Phase Oversampling	0 %
Slice Oversampling	0.0 %
FOV Read	260 mm
FOV Phase	100.0 %
Slice Thickness	1.600 mm
TR	3.5 ms
TE	1.51 ms
Averages	1
Concatenations	1
AutoAlign	---
Coil Elements	AC

Contrast - Common

TR	3.5 ms
TE	1.51 ms
MTC	Off
Magn. Preparation	None
Flip Angle	15 deg
Fat-Water Contrast	Standard
Dark Blood	Off
Contrasts	1
SWI	Off
Reconstruction	Magnitude

Contrast - Dynamic

Dynamic Mode	Standard
--------------	----------

Contrast - Dynamic

Measurements	1
Multiple Series	Each Measurement

Resolution - Common

FOV Read	260 mm
FOV Phase	100.0 %
Slice Thickness	1.600 mm
Base Resolution	160
Phase Resolution	100 %
Slice Resolution	69 %
Interpolation	Off

Resolution - Acceleration

Acceleration Mode	GRAPPA
Reference Scans	Integrated
Acceleration Factor PE	3
Reference Lines PE	24
Acceleration Factor 3D	1
Phase Partial Fourier	6/8
Slice Partial Fourier	6/8
Asymmetric Echo	Allowed
Elliptical Scanning	Off

Resolution - Filter

Raw Filter	Off
Elliptical Filter	Off
Distortion Correction	3D
Normalize	B1 Filter
Image Filter	Off

Geometry - Common

Slab Group	1
Slabs	1
Distance Factor	20 %
Position	L0.0 A16.0 H0.0 mm
Orientation	Sagittal
Phase Encoding Dir.	A >> P
Slices per Slab	128
Phase Oversampling	0 %
Slice Oversampling	0.0 %
FOV Read	260 mm
FOV Phase	100.0 %
Slice Thickness	1.600 mm
TR	3.5 ms
Multi-Slice Mode	Interleaved
Series	Interleaved
Concatenations	1

Geometry - AutoAlign

Slab Group	1
Position	L0.0 A16.0 H0.0 mm
Orientation	Sagittal
Phase Encoding Dir.	A >> P
AutoAlign	---
Initial Position	L0.0 A16.0 H0.0
L	0.0 mm
A	16.0 mm
H	0.0 mm
Initial Orientation	Sagittal
Initial Rotation	0.00 deg

Geometry - Saturation

Saturation Mode	Standard
Special Saturation	None

Geometry - Tim Planning Suite

Set-n-Go Protocol	Off
Table Position	0 mm
Table Position	H

System - Miscellaneous

Coil Selection	Manual
Radial Sorting	Off
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combination	Adaptive Combine
Matrix Optimization	Off

System - Adjustments

Adjustment Strategy	Standard
B0 Shim	Tune up
B1 Shim	TrueForm
Adjustment Tolerance	Auto
Confirm Frequency	Never
Assume Silicone	Off

System - Adjust Volume

Position	Isocenter
Orientation	Transversal
Rotation	0.00 deg
A >> P	263 mm
R >> L	350 mm
F >> H	350 mm
Reset	Off

System - pTx

B1 Shim	TrueForm
Excitation	Non-sel.

System - Tx/Rx

Frequency 1H	297.117721 MHz
--------------	----------------

System - Tx/Rx

? Ref. Amplitude 1H	0.000 V
Reset	Off
Correction Factor	1.00
Image Scaling	1.000

Physio - Signal

1st Signal/Mode	None
TR	3.5 ms
Segments	1
Concatenations	1

Physio - Cardiac

Tagging	None
Fat-Water Contrast	Standard
Magn. Preparation	None
Dark Blood	Off
FOV Read	260 mm
FOV Phase	100.0 %
Phase Resolution	100 %

Physio - PACE

Resp. Control	Off
Concatenations	1

Inline - Liver

Liver Registration	Off
Save Original Images	On

Inline - Subtraction

Subtract	Off
Measurements	1
StdDev	Off
Save Original Images	On

Inline - MIP

MIP Sag	Off
MIP Cor	Off
MIP Tra	Off
MIP Time	Off
Radial MIP	Off
Save Original Images	On
MPR Sag	Off
MPR Cor	Off
MPR Tra	Off

Inline - Soft Tissue

Wash-in	Off
Wash-out	Off
TTP	Off
PEI	Off
MIP Time	Off
Measurements	1

Inline - Composing

Inline - MapIt

MapIt	None
Flip Angle	15 deg
Measurements	1
Contrasts	1
TE	1.51 ms
TR	3.5 ms
Save Original Images	On

Sequence - Part 1

Sequence Name	fl
Dimension	3D
Excitation	Non-sel.
RF Pulse Type	Fast
Gradient Mode	Normal
Flow Compensation	None
Bandwidth	540 Hz/Px
Asymmetric Echo	Allowed
Segments	1

Sequence - Part 2

Introduction	Off
RF Spoiling	On
Acoustic noise reduction	Off

Sequence - Nuclei

TX/RX Nucleus	1H
TX/RX Delta Frequency	0 Hz
TX Nucleus	None
TX Delta Frequency	0 Hz
Coil Elements	AC

Sequence - Special

Readout polarity	Positive
Image processing	Standard
Apply echo spacing	Off
Echo spacing	0 us
Delta echo spacing	0 us
Dummy scans	0 ms
RF pulse duration	100 us
Gradient spoiling	Siemens
Gradient moment factor	1.00
Receiver gain mode	Siemens
Number of segments	1
Current segment	0
Lines before/after seg	0

Sequence - Assistant

SAR Assistant	Off
Allowed Delay	0 s

\\MARTINOS DEVELOPER\HUBER\FRISGO_functional\FRISGO_20251106_COD_M1_2nd\AAsetter--192.16 8.2.1--ep2d

TA: 2 sec Coil Selection: Manual Voxel Size: 3.1×3.1×5.0 mm³ Acc:: None Rel. SNR: 1.00

Properties

Start measurement without further preparation	On
Wait for User to Start	Off
Start measurements	Single Measurement
Prio Recon	Off
Auto Open Inline Display	Off
Auto Close Inline Display	Off
Load Images to MR View&GO	On
Auto Store Images	On
Disable auto transfer to PACS	Off
Load Images to Stamp Segments	On
Load Images to Graphic Segments	On
Graphic segment	Default
Inline Movie	Off

Routine

Slice Group	1
Slices	2
Distance Factor	0 %
Position	Isocenter
Orientation	Transversal
Phase Encoding Dir.	A >> P
Phase Oversampling	0 %
FOV Read	500 mm
FOV Phase	20.0 %
Slice Thickness	5.0 mm
TR	53.0 ms
TE	15.00 ms
Averages	1
Concatenations	1
AutoAlign	---
Coil Elements	AC

Contrast - Common

TR	53.0 ms
TE	15.00 ms
MTC	Off
Flip Angle	10 deg
Fat-Water Contrast	Standard
Reconstruction	Magnitude

Contrast - Dynamic

Dynamic Mode	Standard
Measurements	1
Delay in TR	0.00 ms

Resolution - Common

FOV Read	500 mm
FOV Phase	20.0 %

Resolution - Common

Slice Thickness	5.0 mm
Base Resolution	160
Phase Resolution	91 %
Interpolation	Off

Resolution - Acceleration

Acceleration Mode	None
Phase Partial Fourier	Off

Resolution - Filter

Raw Filter	Off
Elliptical Filter	Off
Hamming	Off
Distortion Correction	Off
Normalize	Off

Geometry - Common

Slice Group	1
Slices	2
Distance Factor	0 %
Position	Isocenter
Orientation	Transversal
Phase Encoding Dir.	A >> P
Phase Oversampling	0 %
FOV Read	500 mm
FOV Phase	20.0 %
Slice Thickness	5.0 mm
TR	53.0 ms
Multi-Slice Mode	Interleaved
Series	Interleaved
Concatenations	1

Geometry - AutoAlign

Slice Group	1
Position	Isocenter
Orientation	Transversal
Phase Encoding Dir.	A >> P
AutoAlign	---
Initial Position	Isocenter
L	0.0 mm
A	0.0 mm
H	0.0 mm
Initial Orientation	Transversal
Initial Rotation	0.00 deg

Geometry - Saturation

Special Saturation	None
--------------------	------

Geometry - Tim Planning Suite

Set-n-Go Protocol	Off
Table Position	0 mm
Table Position	H

System - Miscellaneous

Coil Selection	Manual
Radial Sorting	Off
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combination	Adaptive Combine
Matrix Optimization	Off

System - Adjustments

Adjustment Strategy	Standard
B0 Shim	Tune up
B1 Shim	TrueForm
Adjustment Tolerance	None
Confirm Frequency	Never
Assume Silicone	Off

System - Adjust Volume

! Position	L2.2 P0.0 H8.8 mm
! Orientation	T > S-2.9 > C0.1
! Rotation	0.65 deg
! A >> P	263 mm
! R >> L	350 mm
! F >> H	350 mm
Reset	Off

System - pTx

B1 Shim	TrueForm
Excitation	Standard

System - Tx/Rx

Frequency 1H	297.117721 MHz
? Ref. Amplitude 1H	0.000 V
Reset	Off
Correction Factor	1.00
Image Scaling	1.000

Physio - Signal

1st Signal/Mode	None
TR	53.0 ms
Log Signals	Off
Concatenations	1

BOLD

GLM Statistics	Off
Ignore Meas. at Start	0
Ignore After Transition	0
Model Transition States	On
Temp. Highpass Filter	On

BOLD

Threshold	4.00
Paradigm Size	3
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Active
Motion Correction	Off
Spatial Filter	Off
Measurements	1
Delay in TR	0.00 ms

Sequence - Part 1

Sequence Name	epfid
Excitation	Standard
RF Pulse Type	Fast
Gradient Mode	Fast
Bandwidth	1838 Hz/Px
Echo Spacing	0.63 ms
Free Echo Spacing	Off
EPI Factor	29

Sequence - Part 2

Introduction	On
--------------	----

Sequence - Assistant

SAR Assistant	Off
---------------	-----

\\MARTINOS DEVELOPER\HUBER\FRISGO_functional\FRISGO_20251106_COD_M1_2nd\rslh_ep3d_0.35_run5

TA: 12:43 min Coil Selection: Manual Voxel Size: 0.4x0.4x0.4 mm³ Acc:: 3 Rel. SNR: 1.00

Properties

Start measurement without further preparation	On
Wait for User to Start	Off
Start measurements	Single Measurement
Prio Recon	Off
Auto Open Inline Display	Off
Auto Close Inline Display	Off
Load Images to MR View&GO	On
Auto Store Images	On
Disable auto transfer to PACS	Off
Load Images to Stamp Segments	Off
Load Images to Graphic Segments	Off
Graphic segment	Default
Inline Movie	Off

Routine

Slab Group	1
Slabs	1
Position	R1.0 P17.0 H27.6 mm
Orientation	T > C-20.4
Phase Encoding Dir.	A >> P
Slices per Slab	18
Phase Oversampling	0.0 %
Slice Oversampling	0.0 %
FOV Read	130 mm
FOV Phase	100.0 %
Slice Thickness	0.35 mm
TR	54.1 ms
Vol. TR	5193.6 ms
TE 1	17.60 ms
Averages	1
Multi-echo Shots	1
AutoAlign	Head > Brain
Coil Elements	AC

Contrast - Common

TR	54.1 ms
Vol. TR	5193.6 ms
TE 1	17.60 ms
Multi-echo spacing	39.80 ms
MTC	Off
Flip Angle	18 deg
Fat-Water Contrast	Fast Fat Saturation
Magn. Prep. Shots	1
Contrasts	1
Reconstruction	Magn./Phase

Contrast - Dynamic

Dynamic Mode	Standard
--------------	----------

Contrast - Dynamic

Measurements	144
Reordering	Linear

Resolution - Common

FOV Read	130 mm
FOV Phase	100.0 %
Slice Thickness	0.35 mm
Base Resolution	368
Phase Resolution	100 %
Slice Resolution	100 %
Interpolation	Off

Resolution - Acceleration

Acceleration Mode	CAIPIRINHA
Reference Scans	GRE/Separate
Acceleration Factor PE	1
Reference Lines PE	75
Acceleration Factor 3D	3
Reference Lines 3D	18
Reordering Shift 3D	1
Phase Partial Fourier	6/8
Slice Partial Fourier	Off

Resolution - Filter

Raw Filter	Off
Elliptical Filter	Off
Distortion Correction	Off
Normalize	Off
Image Filter	Off

Geometry - Common

Slab Group	1
Slabs	1
Position	R1.0 P17.0 H27.6 mm
Orientation	T > C-20.4
Phase Encoding Dir.	A >> P
Slices per Slab	18
Phase Oversampling	0.0 %
Slice Oversampling	0.0 %
FOV Read	130 mm
FOV Phase	100.0 %
Slice Thickness	0.35 mm
TR	54.1 ms
Vol. TR	5193.6 ms
Multi-echo Shots	1

Geometry - AutoAlign

Slab Group	1
Position	R1.0 P17.0 H27.6 mm

Geometry - AutoAlign

Orientation	T > C-20.4
Phase Encoding Dir.	A >> P
AutoAlign	Head > Brain
Initial Position	R1.6 P9.1 H28.4
R	1.6 mm
P	9.1 mm
H	28.4 mm
Initial Orientation	T > C
T > C	-18.10
> S	0.00
Initial Rotation	1.32 deg

Geometry - Saturation

Saturation Mode	Quick
-----------------	-------

Geometry - Tim Planning Suite

Set-n-Go Protocol	Off
Table Position	0 mm
Table Position	H

System - Miscellaneous

Coil Selection	Manual
Radial Sorting	Off
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combination	Adaptive Combine
Matrix Optimization	Off

System - Adjustments

Adjustment Strategy	Standard
B0 Shim	Brain
B1 Shim	TrueForm
Adjustment Tolerance	Auto
Confirm Frequency	Never
Assume Silicone	Off

System - Adjust Volume

I Position	R0.9 P16.2 H23.9 mm
I Orientation	T > C-20.4
I Rotation	6.43 deg
! A >> P	130 mm
! R >> L	130 mm
! F >> H	29 mm
Reset	Off

System - pTx

B1 Shim	TrueForm
Excitation	Slab-sel.

System - Tx/Rx

Frequency 1H	297.117721 MHz
! Ref. Amplitude 1H	370.000 V

System - Tx/Rx

Reset	Off
Correction Factor	1.00
Image Scaling	1.000

Sequence - Part 1

Sequence Name	vaso a2d6453
Dimension	3D
Excitation	Slab-sel.
RF Pulse Type	Normal
Gradient Mode	Fast
Reordering	Linear
Bandwidth	590 Hz/Px
Echo Spacing	2.12 ms
Turbo Factor	96
Segmentation	16
EPI Factor	18

Sequence - Part 2

Introduction	On
RF Spoiling	On

Sequence - Special

RF duration	1540 us
RF time x BW	8
PAT ref. FA	3 deg
Fat sat. FA	110 deg
T1 (Ernst FA)	1200 ms
Invert PE	Off
Min. TE w/ PF	On
Ramp Sampling	On
Trigger per shot	Off
Noise image	Off
Relax spoilers	Off
Round up Vol. TR	Off
MT flip angle	500 deg
MT off-res.	2000 Hz
MT RF duration	10240 us
Custom Water Exc.	-none-
Phase Correction	per Series
Saturation RF	per Shot
EPI rise time factor	1.47
G. spoil dephasing[1]	0.0 pi
G. spoil dephasing[2]	4.0 pi
G. spoil dephasing[3]	2.0 pi
Modify Ice Config	On
G-factor map	Off
GRAPPA Regularization	5000 /10^6
Slab Scale	-10 %
RF spoil scheme	Conventional
Read polarity	Inverted RO

Sequence - Assistant

SAR Assistant	Off
---------------	-----

\\MARTINOS DEVELOPER\HUBER\FRISGO_functional\FRISGO_20251106_COD_M1_2nd\mesovein_WE_R
O1_use_me_princess

TA: 6:27 min Coil Selection: Manual Voxel Size: 0.4x0.4x0.4 mm³ Acc:: 6 Rel. SNR: 1.00

Properties

Start measurement without further preparation	Off
Wait for User to Start	Off
Start measurements	Single Measurement
Prio Recon	Off
Auto Open Inline Display	Off
Auto Close Inline Display	Off
Load Images to MR View&GO	On
Auto Store Images	On
Disable auto transfer to PACS	Off
Load Images to Stamp Segments	Off
Load Images to Graphic Segments	Off
Graphic segment	Default
Inline Movie	Off

Routine

Slab Group	1
Slabs	1
Position	L0.2 A2.1 F10.2 mm
Orientation	T > C-31.0
Phase Encoding Dir.	A >> P
Slices per Slab	352
Phase Oversampling	2.9 %
Slice Oversampling	0.0 %
FOV Read	165 mm
FOV Phase	112.7 %
Slice Thickness	0.35 mm
TR	60.0 ms
Vol. TR	380.16 s
TE 1	11.00 ms
TE 2	27.60 ms
TE 3	44.20 ms
Averages	1
Multi-echo Shots	1
AutoAlign	---
Coil Elements	AC

Contrast - Common

TR	60.0 ms
Vol. TR	380.16 s
TE 1	11.00 ms
TE 2	27.60 ms
TE 3	44.20 ms
Multi-echo spacing	16.60 ms
MTC	Off
Flip Angle	15 deg
Fat-Water Contrast	Water Excitation
Magn. Prep. Shots	1
Contrasts	3

Contrast - Common

Reconstruction	Magn./Phase
----------------	-------------

Contrast - Dynamic

Dynamic Mode	Standard
Measurements	1
Reordering	Linear

Resolution - Common

FOV Read	165 mm
FOV Phase	112.7 %
Slice Thickness	0.35 mm
Base Resolution	466
Phase Resolution	100 %
Slice Resolution	100 %
Interpolation	Off

Resolution - Acceleration

Acceleration Mode	CAIPIRINHA
Reference Scans	GRE/Separate
Acceleration Factor PE	3
Reference Lines PE	24
Acceleration Factor 3D	2
Reference Lines 3D	24
Reordering Shift 3D	1
Phase Partial Fourier	Off
Slice Partial Fourier	Off

Resolution - Filter

Raw Filter	Off
Elliptical Filter	Off
Distortion Correction	Off
Normalize	Off
Image Filter	Off

Geometry - Common

Slab Group	1
Slabs	1
Position	L0.2 A2.1 F10.2 mm
Orientation	T > C-31.0
Phase Encoding Dir.	A >> P
Slices per Slab	352
Phase Oversampling	2.9 %
Slice Oversampling	0.0 %
FOV Read	165 mm
FOV Phase	112.7 %
Slice Thickness	0.35 mm
TR	60.0 ms
Vol. TR	380.16 s
Multi-echo Shots	1

Geometry - AutoAlign

Slab Group	1
Position	L0.2 A2.1 F10.2 mm
Orientation	T > C-31.0
Phase Encoding Dir.	A >> P
AutoAlign	---
Initial Position	L0.2 A2.1 F10.2
L	0.2 mm
A	2.1 mm
F	10.2 mm
Initial Orientation	T > C
T > C	-31.00
> S	0.00
Initial Rotation	0.00 deg

Geometry - Saturation

Saturation Mode	Standard
-----------------	----------

Geometry - Tim Planning Suite

Set-n-Go Protocol	Off
Table Position	0 mm
Table Position	H

System - Miscellaneous

Coil Selection	Manual
Radial Sorting	Off
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combination	Adaptive Combine
Matrix Optimization	Off

System - Adjustments

Adjustment Strategy	Standard
B0 Shim	Brain
B1 Shim	TrueForm
Adjustment Tolerance	Auto
Confirm Frequency	Never
Assume Silicone	Off

System - Adjust Volume

! Position	L0.9 A1.3 F9.0 mm
! Orientation	T > C-27.1 > S0.4
! Rotation	90.14 deg
! R >> L	123 mm
! A >> P	163 mm
! F >> H	87 mm
Reset	Off

System - pTx

B1 Shim	TrueForm
Excitation	Slab-sel.

System - Tx/Rx

Frequency 1H	297.117721 MHz
? Ref. Amplitude 1H	0.000 V
Reset	Off
Correction Factor	1.00
Image Scaling	1.000

Sequence - Part 1

Sequence Name	vaso a2d6453
Dimension	3D
Excitation	Slab-sel.
RF Pulse Type	Normal
Gradient Mode	Normal
Reordering	Linear
Bandwidth	466 Hz/Px
Echo Spacing	2.83 ms
Turbo Factor	6336
Segmentation	36
EPI Factor	5

Sequence - Part 2

Introduction	On
RF Spoiling	On

Sequence - Special

RF duration	340 us
RF time x BW	12
PAT ref. FA	3 deg
Fat sat. FA	110 deg
T1 (Ernst FA)	1200 ms
Invert PE	Off
Min. TE w/ PF	On
Ramp Sampling	On
Trigger per shot	Off
Noise image	Off
Relax spoilers	Off
Round up Vol. TR	Off
MT flip angle	500 deg
MT off-res.	2000 Hz
MT RF duration	10240 us
Custom Water Exc.	-none-
Phase Correction	per Series
Saturation RF	per Shot
EPI rise time factor	1.20
G. spoil dephasing[1]	0.0 pi
G. spoil dephasing[2]	4.0 pi
G. spoil dephasing[3]	2.0 pi
Modify Ice Config	On
G-factor map	Off
GRAPPA Regularization	10 / 10^6
Slab Scale	-10 %
RF spoil scheme	Conventional
Read polarity	Regular RO

Sequence - Assistant

SAR Assistant	Off
---------------	-----

\\MARTINOS DEVELOPER\HUBER\FRISGO_functional\FRISGO_20251106_COD_M1_2nd\mesovein_WE_R
O1_use_me_princess

TA: 6:27 min Coil Selection: Manual Voxel Size: 0.4x0.4x0.4 mm³ Acc:: 6 Rel. SNR: 1.00

Properties

Start measurement without further preparation	Off
Wait for User to Start	Off
Start measurements	Single Measurement
Prio Recon	Off
Auto Open Inline Display	Off
Auto Close Inline Display	Off
Load Images to MR View&GO	On
Auto Store Images	On
Disable auto transfer to PACS	Off
Load Images to Stamp Segments	Off
Load Images to Graphic Segments	Off
Graphic segment	Default
Inline Movie	Off

Routine

Slab Group	1
Slabs	1
Position	L0.2 A2.1 F10.2 mm
Orientation	T > C-31.0
Phase Encoding Dir.	A >> P
Slices per Slab	352
Phase Oversampling	2.9 %
Slice Oversampling	0.0 %
FOV Read	165 mm
FOV Phase	112.7 %
Slice Thickness	0.35 mm
TR	60.0 ms
Vol. TR	380.16 s
TE 1	11.00 ms
TE 2	27.60 ms
TE 3	44.20 ms
Averages	1
Multi-echo Shots	1
AutoAlign	---
Coil Elements	AC

Contrast - Common

TR	60.0 ms
Vol. TR	380.16 s
TE 1	11.00 ms
TE 2	27.60 ms
TE 3	44.20 ms
Multi-echo spacing	16.60 ms
MTC	Off
Flip Angle	15 deg
Fat-Water Contrast	Water Excitation
Magn. Prep. Shots	1
Contrasts	3

Contrast - Common

Reconstruction	Magn./Phase
----------------	-------------

Contrast - Dynamic

Dynamic Mode	Standard
Measurements	1
Reordering	Linear

Resolution - Common

FOV Read	165 mm
FOV Phase	112.7 %
Slice Thickness	0.35 mm
Base Resolution	466
Phase Resolution	100 %
Slice Resolution	100 %
Interpolation	Off

Resolution - Acceleration

Acceleration Mode	CAIPIRINHA
Reference Scans	GRE/Separate
Acceleration Factor PE	3
Reference Lines PE	24
Acceleration Factor 3D	2
Reference Lines 3D	24
Reordering Shift 3D	1
Phase Partial Fourier	Off
Slice Partial Fourier	Off

Resolution - Filter

Raw Filter	Off
Elliptical Filter	Off
Distortion Correction	Off
Normalize	Off
Image Filter	Off

Geometry - Common

Slab Group	1
Slabs	1
Position	L0.2 A2.1 F10.2 mm
Orientation	T > C-31.0
Phase Encoding Dir.	A >> P
Slices per Slab	352
Phase Oversampling	2.9 %
Slice Oversampling	0.0 %
FOV Read	165 mm
FOV Phase	112.7 %
Slice Thickness	0.35 mm
TR	60.0 ms
Vol. TR	380.16 s
Multi-echo Shots	1

Geometry - AutoAlign

Slab Group	1
Position	L0.2 A2.1 F10.2 mm
Orientation	T > C-31.0
Phase Encoding Dir.	A >> P
AutoAlign	---
Initial Position	L0.2 A2.1 F10.2
L	0.2 mm
A	2.1 mm
F	10.2 mm
Initial Orientation	T > C
T > C	-31.00
> S	0.00
Initial Rotation	0.00 deg

Geometry - Saturation

Saturation Mode	Standard
-----------------	----------

Geometry - Tim Planning Suite

Set-n-Go Protocol	Off
Table Position	0 mm
Table Position	H

System - Miscellaneous

Coil Selection	Manual
Radial Sorting	Off
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combination	Adaptive Combine
Matrix Optimization	Off

System - Adjustments

Adjustment Strategy	Standard
B0 Shim	Brain
B1 Shim	TrueForm
Adjustment Tolerance	Auto
Confirm Frequency	Never
Assume Silicone	Off

System - Adjust Volume

! Position	L0.9 A1.3 F9.0 mm
! Orientation	T > C-27.1 > S0.4
! Rotation	90.14 deg
! R >> L	123 mm
! A >> P	163 mm
! F >> H	87 mm
Reset	Off

System - pTx

B1 Shim	TrueForm
Excitation	Slab-sel.

System - Tx/Rx

Frequency 1H	297.117721 MHz
? Ref. Amplitude 1H	0.000 V
Reset	Off
Correction Factor	1.00
Image Scaling	1.000

Sequence - Part 1

Sequence Name	vaso a2d6453
Dimension	3D
Excitation	Slab-sel.
RF Pulse Type	Normal
Gradient Mode	Normal
Reordering	Linear
Bandwidth	466 Hz/Px
Echo Spacing	2.83 ms
Turbo Factor	6336
Segmentation	36
EPI Factor	5

Sequence - Part 2

Introduction	On
RF Spoiling	On

Sequence - Special

RF duration	340 us
RF time x BW	12
PAT ref. FA	3 deg
Fat sat. FA	110 deg
T1 (Ernst FA)	1200 ms
Invert PE	Off
Min. TE w/ PF	On
Ramp Sampling	On
Trigger per shot	Off
Noise image	Off
Relax spoilers	Off
Round up Vol. TR	Off
MT flip angle	500 deg
MT off-res.	2000 Hz
MT RF duration	10240 us
Custom Water Exc.	-none-
Phase Correction	per Series
Saturation RF	per Shot
EPI rise time factor	1.20
G. spoil dephasing[1]	0.0 pi
G. spoil dephasing[2]	4.0 pi
G. spoil dephasing[3]	2.0 pi
Modify Ice Config	On
G-factor map	Off
GRAPPA Regularization	10 / 10^6
Slab Scale	-10 %
RF spoil scheme	Conventional
Read polarity	Inverted RO

Sequence - Assistant

SAR Assistant	Off
---------------	-----