

Table of contents

\\USER

<table border="1"><tr><td>THimpuls</td></tr></table>	THimpuls
THimpuls	

<table border="1"><tr><td><table border="1"><tr><td>v0p2</td></tr></table></td></tr></table>	<table border="1"><tr><td>v0p2</td></tr></table>	v0p2
<table border="1"><tr><td>v0p2</td></tr></table>	v0p2	
v0p2		

<table border="1"><tr><td><table border="1"><tr><td><table border="1"><tr><td>BO_B1_MP2RAGE</td></tr></table></td></tr></table></td></tr></table>	<table border="1"><tr><td><table border="1"><tr><td>BO_B1_MP2RAGE</td></tr></table></td></tr></table>	<table border="1"><tr><td>BO_B1_MP2RAGE</td></tr></table>	BO_B1_MP2RAGE
<table border="1"><tr><td><table border="1"><tr><td>BO_B1_MP2RAGE</td></tr></table></td></tr></table>	<table border="1"><tr><td>BO_B1_MP2RAGE</td></tr></table>	BO_B1_MP2RAGE	
<table border="1"><tr><td>BO_B1_MP2RAGE</td></tr></table>	BO_B1_MP2RAGE		
BO_B1_MP2RAGE			

<table border="1"><tr><td><table border="1"><tr><td><table border="1"><tr><td><table border="1"><tr><td>AAhead_scout</td></tr><tr><td>gre_b0map_4iso_sag</td></tr><tr><td>dzne_b1map_5iso_sag</td></tr><tr><td>mp2rage_0p7iso</td></tr></table></td></tr></table></td></tr></table></td></tr></table>	<table border="1"><tr><td><table border="1"><tr><td><table border="1"><tr><td>AAhead_scout</td></tr><tr><td>gre_b0map_4iso_sag</td></tr><tr><td>dzne_b1map_5iso_sag</td></tr><tr><td>mp2rage_0p7iso</td></tr></table></td></tr></table></td></tr></table>	<table border="1"><tr><td><table border="1"><tr><td>AAhead_scout</td></tr><tr><td>gre_b0map_4iso_sag</td></tr><tr><td>dzne_b1map_5iso_sag</td></tr><tr><td>mp2rage_0p7iso</td></tr></table></td></tr></table>	<table border="1"><tr><td>AAhead_scout</td></tr><tr><td>gre_b0map_4iso_sag</td></tr><tr><td>dzne_b1map_5iso_sag</td></tr><tr><td>mp2rage_0p7iso</td></tr></table>	AAhead_scout	gre_b0map_4iso_sag	dzne_b1map_5iso_sag	mp2rage_0p7iso
<table border="1"><tr><td><table border="1"><tr><td><table border="1"><tr><td>AAhead_scout</td></tr><tr><td>gre_b0map_4iso_sag</td></tr><tr><td>dzne_b1map_5iso_sag</td></tr><tr><td>mp2rage_0p7iso</td></tr></table></td></tr></table></td></tr></table>	<table border="1"><tr><td><table border="1"><tr><td>AAhead_scout</td></tr><tr><td>gre_b0map_4iso_sag</td></tr><tr><td>dzne_b1map_5iso_sag</td></tr><tr><td>mp2rage_0p7iso</td></tr></table></td></tr></table>	<table border="1"><tr><td>AAhead_scout</td></tr><tr><td>gre_b0map_4iso_sag</td></tr><tr><td>dzne_b1map_5iso_sag</td></tr><tr><td>mp2rage_0p7iso</td></tr></table>	AAhead_scout	gre_b0map_4iso_sag	dzne_b1map_5iso_sag	mp2rage_0p7iso	
<table border="1"><tr><td><table border="1"><tr><td>AAhead_scout</td></tr><tr><td>gre_b0map_4iso_sag</td></tr><tr><td>dzne_b1map_5iso_sag</td></tr><tr><td>mp2rage_0p7iso</td></tr></table></td></tr></table>	<table border="1"><tr><td>AAhead_scout</td></tr><tr><td>gre_b0map_4iso_sag</td></tr><tr><td>dzne_b1map_5iso_sag</td></tr><tr><td>mp2rage_0p7iso</td></tr></table>	AAhead_scout	gre_b0map_4iso_sag	dzne_b1map_5iso_sag	mp2rage_0p7iso		
<table border="1"><tr><td>AAhead_scout</td></tr><tr><td>gre_b0map_4iso_sag</td></tr><tr><td>dzne_b1map_5iso_sag</td></tr><tr><td>mp2rage_0p7iso</td></tr></table>	AAhead_scout	gre_b0map_4iso_sag	dzne_b1map_5iso_sag	mp2rage_0p7iso			
AAhead_scout							
gre_b0map_4iso_sag							
dzne_b1map_5iso_sag							
mp2rage_0p7iso							

\\USER\THimpuls\lv0p2\B0_B1_MP2RAGE\AAhead_scout

TA: 15 sec Coil Selection: Auto Voxel Size: 1.6×1.6×1.6 mm³ Acc:: 3 Rel. SNR: 1.00**Properties**

Start measurement without further preparation	On
Wait for User to Start	Off
Start measurements	Single Measurement
Prio Recon	Off
Auto Open Inline Display	Off
Auto Close Inline Display	Off
Load Images to MR View&GO	On
Auto Store Images	On
Disable auto transfer to PACS	Off
Load Images to Stamp Segments	Off
Load Images to Graphic Segments	On
	Default
Inline Movie	Off
Before Measurement	
After Measurement	

Routine

Slab Group	1
Slabs	1
Distance Factor	20 %
Position	L0.0 A30.0 H0.0 mm
Orientation	Sagittal
Phase Encoding Dir.	A >> P
Slices per Slab	128
Phase Oversampling	0 %
Slice Oversampling	0.0 %
FOV Read	260 mm
FOV Phase	100.0 %
Slice Thickness	1.6 mm
TR	3.53 ms
TE	1.53 ms
Averages	1
Concatenations	1
AutoAlign	Head

Contrast - Common

TR	3.53 ms
TE	1.53 ms
Flip Angle	16 deg
Fat-Water Contrast	Standard
Contrasts	1
Reconstruction	Magnitude

Contrast - Dynamic

Dynamic Mode	Standard
Measurements	1
Time to Center	6.7 s

Resolution - Common

FOV Read	260 mm
FOV Phase	100.0 %
Slice Thickness	1.6 mm
Base Resolution	160
Phase Resolution	100 %
Slice Resolution	69 %
Trajectory	Cartesian

Resolution - Acceleration

Acceleration Mode	GRAPPA
Reference Scans	Integrated
Acceleration Factor PE	3
Reference Lines PE	24
Acceleration Factor 3D	1
Phase Partial Fourier	6/8
Slice Partial Fourier	6/8
Asymmetric Echo	Weak

Resolution - Filter

Raw Filter	Off
Elliptical Filter	Off
Distortion Correction	3D
Normalize	B1 Filter
Image Filter	Off

Geometry - Common

Slab Group	1
Slabs	1
Distance Factor	20 %
Position	L0.0 A30.0 H0.0 mm
Orientation	Sagittal
Phase Encoding Dir.	A >> P
Slices per Slab	128
Phase Oversampling	0 %
Slice Oversampling	0.0 %
FOV Read	260 mm
FOV Phase	100.0 %
Slice Thickness	1.6 mm
TR	3.53 ms
Multi-Slice Mode	Sequential
Series	Ascending
Concatenations	1

Geometry - AutoAlign

Slab Group	1
Position	L0.0 A30.0 H0.0 mm
Orientation	Sagittal
Phase Encoding Dir.	A >> P
Initial Position	Isocenter
L	0.0 mm

Geometry - AutoAlign

P	0.0 mm
H	0.0 mm
Initial Orientation	Transversal
Initial Rotation	0.00 deg
AutoAlign	Head

Geometry - Tim Planning Suite

Set-n-Go Protocol	Off
Label	
Table Position	0 mm
Table Position	H

System - Miscellaneous

Coil Selection	Auto Coil Select
Radial Sorting	Off
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combination	Sum of Squares
Matrix Optimization	Off
Coil Focus	Flat

System - Adjustments

Adjustment Strategy	Standard
B0 Shim	Tune up
B1 Shim	TrueForm
Adjustment Tolerance	Auto
Adjust with Body Coil	Off
Confirm Frequency	Never
Assume Silicone	Off

System - Adjust Volume

Position	Isocenter
Orientation	Transversal
Rotation	0.00 deg
A >> P	263 mm
R >> L	350 mm
F >> H	350 mm
Reset	Off

System - pTx

B1 Shim	TrueForm
Excitation	Non-sel.

System - Tx/Rx

Frequency 1H	297.041339 MHz
? Ref. Amplitude 1H	0.000 V
Reset	Off
Image Scaling	1.000

Physio - PACE

Resp. Control	Off
Concatenations	1

Inline - Subtraction

Subtract	Off
Measurements	1
StdDev	Off
Save Original Images	On

Inline - MIP

MIP Sag	Off
MIP Cor	Off
MIP Tra	Off
MIP Time	Off
Radial MIP	Off
Save Original Images	On
MPR Sag	Off
MPR Cor	Off
MPR Tra	Off
AutoAlign MPR	On

Inline - Composing**Inline - MapIt**

MapIt	None
Flip Angle	16 deg
Measurements	1
Contrasts	1
TE	1.53 ms
TR	3.53 ms
Save Original Images	On

Inline - Open Recon

Algorithm	None
-----------	------

Sequence - Part 1

Sequence Name	fl
Dimension	3D
Excitation	Non-sel.
RF Pulse Type	Fast
Gradient Mode	Normal
Bandwidth	540 Hz/Px
Asymmetric Echo	Weak

Sequence - Part 2

Introduction	On
RF Spoiling	On
Breast Application	Off
Phase Enc. Order	Automatic

Sequence - Assistant

SAR Assistant	Off
---------------	-----

\\USER\THimpuls\lv0p2\B0_B1_MP2RAGE\gre_b0map_4iso_sag

TA: 31 sec Coil Selection: Auto Voxel Size: 4.0×4.0×4.0 mm³ Rel. SNR: 1.00**Properties**

Start measurement without further preparation	Off
Wait for User to Start	Off
Start measurements	Single Measurement
Prio Recon	Off
Auto Open Inline Display	Off
Auto Close Inline Display	Off
Load Images to MR View&GO	On
Auto Store Images	On
Disable auto transfer to PACS	Off
Load Images to Stamp Segments	On
Load Images to Graphic Segments	Off
	Default
Inline Movie	Off
Before Measurement	
After Measurement	

Routine

Slice Group	1
Slices	44
Distance Factor	0 %
Position	Isocenter
Orientation	Sagittal
Phase Encoding Dir.	A >> P
AutoAlign	---
Phase Oversampling	0 %
FOV Read	256 mm
FOV Phase	100.0 %
Slice Thickness	4.0 mm
TR	238.0 ms
TE 1	3.06 ms
TE 2	4.08 ms
Averages	1
Concatenations	1

Contrast - Common

TR	238.0 ms
TE 1	3.06 ms
TE 2	4.08 ms
MTC	Off
Flip Angle	25 deg
Fat-Water Contrast	Standard
Contrasts	2
Reconstruction	Magn./Phase

Contrast - Dynamic

Dynamic Mode	Standard
Measurements	1
Multiple Series	Off

Resolution - Common

FOV Read	256 mm
FOV Phase	100.0 %
Slice Thickness	4.0 mm
Base Resolution	64
Phase Resolution	100 %
Interpolation	Off

Resolution - Acceleration

Phase Partial Fourier	Off
Asymmetric Echo	Off

Resolution - Filter

Raw Filter	Off
Elliptical Filter	Off
Distortion Correction	3D
Normalize	Off
Image Filter	Off

Geometry - Common

Slice Group	1
Slices	44
Distance Factor	0 %
Position	Isocenter
Orientation	Sagittal
Phase Encoding Dir.	A >> P
Phase Oversampling	0 %
FOV Read	256 mm
FOV Phase	100.0 %
Slice Thickness	4.0 mm
TR	238.0 ms
Multi-Slice Mode	Interleaved
Series	Interleaved
Concatenations	1

Geometry - AutoAlign

Slice Group	1
Position	Isocenter
Orientation	Sagittal
Phase Encoding Dir.	A >> P
AutoAlign	---
Initial Position	Isocenter
L	0.0 mm
P	0.0 mm
H	0.0 mm
Initial Orientation	Sagittal
Initial Rotation	0.00 deg

Geometry - Saturation

Special Saturation	None
--------------------	------

Geometry - Tim Planning Suite

Set-n-Go Protocol	Off
Label	
Table Position	0 mm
Table Position	H

System - Miscellaneous

Coil Selection	Auto Coil Select
Radial Sorting	Off
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combination	Sum of Squares
Matrix Optimization	Off
Coil Focus	Flat

System - Adjustments

Adjustment Strategy	Standard
B0 Shim	Brain
B1 Shim	TrueForm
Adjustment Tolerance	Auto
Adjust with Body Coil	Off
Confirm Frequency	Never
Assume Silicone	Off

System - Adjust Volume

Position	Isocenter
Orientation	Sagittal
Rotation	0.00 deg
A >> P	256 mm
F >> H	256 mm
R >> L	176 mm
Reset	Off

System - pTx

B1 Shim	TrueForm
---------	----------

System - Tx/Rx

Frequency 1H	297.041339 MHz
? Ref. Amplitude 1H	0.000 V
Reset	Off
Image Scaling	1.000

Inline - Open Recon

Algorithm	None
-----------	------

Sequence - Part 1

Sequence Name	fm_r
Dimension	2D
RF Pulse Type	Fast
Gradient Mode	Fast
Flow Compensation	On
Bandwidth	908 Hz/Px
Asymmetric Echo	Off

Sequence - Part 2

Introduction	Off
RF Spoiling	On

Sequence - Assistant

SAR Assistant	Off
---------------	-----

\\USER\THimpuls\w0p2\B0_B1_MP2RAGE\ldzne_b1map_5iso_sag

TA: 33 sec Coil Selection: Manual Voxel Size: 5.0x5.0x5.0 mm³ Acc:: 4 Rel. SNR: 1.00**Properties**

Start measurement without further preparation	Off
Wait for User to Start	Off
Start measurements	Single Measurement
Prio Recon	Off
Auto Open Inline Display	Off
Auto Close Inline Display	Off
Load Images to MR View&GO	On
Auto Store Images	On
Disable auto transfer to PACS	Off
Load Images to Stamp Segments	Off
Load Images to Graphic Segments	Off
	Default
Inline Movie	Off
Before Measurement	
After Measurement	

Routine

Slab Group	1
Slabs	1
Distance Factor	20 %
Position	Isocenter
Orientation	Sagittal
Phase Encoding Dir.	A >> P
AutoAlign	---
Slices per Slab	40
Slice Oversampling	0.0 %
FOV Read	220 mm
FOV Phase	100.0 %
Slice Thickness	5.00 mm
TR	3000.0 ms
TE 1	0.74 ms
TE 2	0.86 ms
TE 3	1.99 ms
Averages	1
Concatenations	1
Coil Elements	R64

Contrast - Common

TR	3000.0 ms
TE 1	0.74 ms
TE 2	0.86 ms
TE 3	1.99 ms
Flip Angle 1	60 deg
Flip Angle 2	6 deg
Fat-Water Contrast	Standard
Reconstruction	Magnitude

Contrast - Dynamic

Dynamic Mode	Standard
Measurements	1

Contrast - Dynamic

Multiple Series	Off
-----------------	-----

Resolution - Common

FOV Read	220 mm
FOV Phase	100.0 %
Slice Thickness	5.00 mm
Base Resolution	44
Phase Resolution	100 %
Slice Resolution	100 %
Interpolation	Off

Resolution - Acceleration

Acceleration Mode	CAIPIRINHA
CAIPIRINHA Mode	Free
Reference Scans	GRE/Separate
Acceleration Factor PE	2
Reference Lines PE	24
Acceleration Factor 3D	2
Reference Lines 3D	24
Reordering Shift 3D	1
Phase Partial Fourier	Off
Slice Partial Fourier	Off
Asymmetric Echo	Off
Elliptical Scanning	On

Resolution - Filter

Raw Filter	Off
Elliptical Filter	Off
Distortion Correction	Off
Normalize	Off
Image Filter	Off

Geometry - Common

Slab Group	1
Slabs	1
Distance Factor	20 %
Position	Isocenter
Orientation	Sagittal
Phase Encoding Dir.	A >> P
Slices per Slab	40
Slice Oversampling	0.0 %
FOV Read	220 mm
FOV Phase	100.0 %
Slice Thickness	5.00 mm
TR	3000.0 ms
Multi-Slice Mode	Sequential
Series	Ascending
Concatenations	1

Geometry - AutoAlign

Slab Group	1
Position	Isocenter
Orientation	Sagittal
Phase Encoding Dir.	A >> P
AutoAlign	---
Initial Position	Isocenter
L	0.0 mm
P	0.0 mm
H	0.0 mm
Initial Orientation	Sagittal
Initial Rotation	0.00 deg

Geometry - Tim Planning Suite

Set-n-Go Protocol	Off
Label	
Table Position	0 mm
Table Position	H

System - Coils**System - Miscellaneous**

Coil Selection	Manual
Radial Sorting	Off
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combination	Adaptive Combine
Matrix Optimization	Off

System - Adjustments

Adjustment Strategy	Standard
B0 Shim	Standard
B1 Shim	TrueForm
Adjustment Tolerance	Auto
Confirm Frequency	Never
Assume Silicone	Off

System - Adjust Volume

Position	Isocenter
Orientation	Sagittal
Rotation	0.00 deg
A >> P	220 mm
F >> H	220 mm
R >> L	200 mm
Reset	Off

System - pTx

B1 Shim	TrueForm
Excitation	Non-sel.

System - Tx/Rx

Frequency 1H	297.041339 MHz
? Ref. Amplitude 1H	0.000 V

System - Tx/Rx

Reset	Off
Correction Factor	1.00
Image Scaling	1.000

Inline - Open Recon

Algorithm	None
-----------	------

Sequence - Part 1

Sequence Name	b1map 9e96902
Dimension	3D
Excitation	Non-sel.
RF Pulse Type	Normal
Gradient Mode	Fast
Bandwidth	1000 Hz/Px
Asymmetric Echo	Off
Turbo Factor	440
Echo Train Duration	1298 ms

Sequence - Part 2

Introduction	On
RF Spoiling	On

Sequence - Special

Mapping Technique	DREAM
Timing Scheme	STE*
Mixing Time	1130 us
Multichannel mode	One-Inv
Scale risetime	1.2

Sequence - Assistant

SAR Assistant	Off
---------------	-----

\\USER\THimpuls\lv0p2\B0_B1_MP2RAGE\mp2rage_0p7iso

TA: 8:27 min Coil Selection: Auto Voxel Size: 0.7×0.7×0.7 mm³ Acc:: 3 Rel. SNR: 1.00**Properties**

Start measurement without further preparation	Off
Wait for User to Start	Off
Start measurements	Single Measurement
Prio Recon	Off
Auto Open Inline Display	Off
Auto Close Inline Display	Off
Load Images to MR View&GO	On
Auto Store Images	On
Disable auto transfer to PACS	Off
Load Images to Stamp Segments	Off
Load Images to Graphic Segments	Off
	Default
Inline Movie	Off
Before Measurement	
After Measurement	

Routine

Slab Group	1
Slabs	1
Distance Factor	50 %
Position	Isocenter
Orientation	Sagittal
Phase Encoding Dir.	A >> P
AutoAlign	---
Slices per Slab	256
Phase Oversampling	0 %
Slice Oversampling	0.0 %
FOV Read	224 mm
FOV Phase	100.0 %
Slice Thickness	0.70 mm
TR	5000.0 ms
TE	3.37 ms
Averages	1
Concatenations	1

Contrast - Common

TR	5000.0 ms
TE	3.37 ms
Magn. Preparation	Non-sel. IR
TI 1	800 ms
TI 2	2700 ms
Flip Angle 1	4 deg
Flip Angle 2	5 deg
Fat-Water Contrast	Fast Water Excitation
Dark Blood	Off
Contrasts	1
Reconstruction	Magnitude

Contrast - Dynamic

Dynamic Mode	Standard
--------------	----------

Contrast - Dynamic

Measurements	1
Multiple Series	Each Measurement
Reordering	Linear

Resolution - Common

FOV Read	224 mm
FOV Phase	100.0 %
Slice Thickness	0.70 mm
Base Resolution	320
Phase Resolution	100 %
Slice Resolution	100 %
Interpolation	Off

Resolution - Acceleration

Acceleration Mode	GRAPPA
Reference Scans	Integrated
Acceleration Factor PE	3
Reference Lines PE	32
Acceleration Factor 3D	1
Deep Resolve / ID Gain	Off
Phase Partial Fourier	6/8
Slice Partial Fourier	7/8
Asymmetric Echo	Off
Elliptical Scanning	Off

Resolution - Filter

Raw Filter	Off
Elliptical Filter	Off
Distortion Correction	3D
Normalize	Off
Image Filter	Off

Geometry - Common

Slab Group	1
Slabs	1
Distance Factor	50 %
Position	Isocenter
Orientation	Sagittal
Phase Encoding Dir.	A >> P
Slices per Slab	256
Phase Oversampling	0 %
Slice Oversampling	0.0 %
FOV Read	224 mm
FOV Phase	100.0 %
Slice Thickness	0.70 mm
TR	5000.0 ms
Multi-Slice Mode	Single Shot
Series	Interleaved
Concatenations	1

Geometry - AutoAlign

Slab Group	1
Position	Isocenter
Orientation	Sagittal
Phase Encoding Dir.	A >> P
AutoAlign	---
Initial Position	Isocenter
L	0.0 mm
P	0.0 mm
H	0.0 mm
Initial Orientation	Transversal
Initial Rotation	0.00 deg

Geometry - Navigator**Geometry - Tim Planning Suite**

Set-n-Go Protocol	Off
Label	
Table Position	0 mm
Table Position	H

System - Miscellaneous

Coil Selection	Auto Coil Select
Radial Sorting	Off
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combination	Sum of Squares
Matrix Optimization	Off
Coil Focus	Flat

System - Adjustments

Adjustment Strategy	Standard
B0 Shim	Brain
B1 Shim	TrueForm
Adjustment Tolerance	Auto
Adjust with Body Coil	Off
Confirm Frequency	Never
Assume Silicone	Off

System - Adjust Volume

Position	Isocenter
Orientation	Sagittal
Rotation	0.00 deg
A >> P	224 mm
F >> H	224 mm
R >> L	180 mm
Reset	Off

System - pTx

B1 Shim	TrueForm
Excitation	Non-sel.

System - Tx/Rx

Frequency 1H	297.041339 MHz
? Ref. Amplitude 1H	0.000 V
Reset	Off
Image Scaling	1.000

Physio - Signal

1st Signal/Mode	None
TR	5000.0 ms
Concatenations	1

Physio - Cardiac

Fat-Water Contrast	Fast Water Excitation
Magn. Preparation	Non-sel. IR
TI 1	800 ms
TI 2	2700 ms
Dark Blood	Off
FOV Read	224 mm
FOV Phase	100.0 %
Phase Resolution	100 %
Motion Correction	None

Physio - PACE

Resp. Control	Off
Concatenations	1

Inline - Subtraction

Subtract	Off
Measurements	1
StdDev	Off
Save Original Images	On

Inline - MIP

MIP Sag	Off
MIP Cor	Off
MIP Tra	Off
MIP Time	Off
Radial MIP	Off
Save Original Images	On
MPR Sag	Off
MPR Cor	Off
MPR Tra	Off

Inline - Composing**Inline - MapIt**

MapIt	T1 Map
Flip Angle 1	4 deg
Flip Angle 2	5 deg
Measurements	1
Contrasts	1
TE	3.37 ms
TR	5000.0 ms
Save Original Images	On

Inline - Open Recon

Algorithm	None
-----------	------

Sequence - Part 1

Sequence Name	tfl
Dimension	3D
Excitation	Non-sel.
RF Pulse Type	Fast
Gradient Mode	Fast
Flow Compensation	None
Reordering	Linear
Bandwidth	200 Hz/Px
Echo Spacing	7.68 ms
Asymmetric Echo	Off
Turbo Factor	224

Sequence - Part 2

Introduction	On
RF Spoiling	On
Incr. Gradient Spoiling	Off
Motion Correction	None

Sequence - Assistant

SAR Assistant	Off
---------------	-----