

\\MARTINOS DEVELOPER\HUBER\Physics_playground\Executed_20260603_Skopelwip_3DEPI_seg_fid_E
S1_4_Nicolas

TA: 7 sec Coil Selection: Manual Voxel Size: 0.9x0.9x0.9 mm³ Acc:: 3 Rel. SNR: 1.00

Properties

Start measurement without further preparation	On
Wait for User to Start	Off
Start measurements	Single Measurement
Prio Recon	Off
Auto Open Inline Display	Off
Auto Close Inline Display	Off
Load Images to MR View&GO	On
Auto Store Images	On
Disable auto transfer to PACS	Off
Load Images to Stamp Segments	Off
Load Images to Graphic Segments	Off
Graphic segment	Default
Inline Movie	Off

Routine

Slab Group	1
Slabs	1
Distance Factor	-50 %
Position	Isocenter
Orientation	Transversal
Phase Encoding Dir.	A >> P
Slices per Slab	30
Phase Oversampling	0 %
Slice Oversampling	0.0 %
FOV Read	190 mm
FOV Phase	100.0 %
Slice Thickness	0.88 mm
TR	50.0 ms
TE	19.00 ms
Averages	1
Concatenations	1
AutoAlign	---
Coil Elements	AC

Contrast - Common

TR	50.0 ms
TE	19.00 ms
MTC	Off
Flip Angle	10 deg
Fat-Water Contrast	Fat Saturation
SWI	Off
Reconstruction	Magnitude

Contrast - Dynamic

Dynamic Mode	Standard
Measurements	1
Multiple Series	Off
3D Reordering	Standard

Contrast - Dynamic

Reordering	Linear
------------	--------

Resolution - Common

FOV Read	190 mm
FOV Phase	100.0 %
Slice Thickness	0.88 mm
Base Resolution	208
Phase Resolution	98 %
Slice Resolution	100 %
Interpolation	Off

Resolution - Acceleration

Acceleration Mode	CAIPIRINHA
CAIPIRINHA Mode	Free
Reference Scans	GRE/Separate
Acceleration Factor PE	1
Reference Lines PE	24
Acceleration Factor 3D	3
Reference Lines 3D	24
Reordering Shift 3D	1
Deep Resolve	Off
Phase Partial Fourier	Off
Slice Partial Fourier	6/8

Resolution - Filter

Raw Filter	Off
Elliptical Filter	Off
Hamming	Off
Distortion Correction	2D
Static Field Correction	Off
Normalize	Off
Image Filter	Off

Geometry - Common

Slab Group	1
Slabs	1
Distance Factor	-50 %
Position	Isocenter
Orientation	Transversal
Phase Encoding Dir.	A >> P
Slices per Slab	30
Phase Oversampling	0 %
Slice Oversampling	0.0 %
FOV Read	190 mm
FOV Phase	100.0 %
Slice Thickness	0.88 mm
TR	50.0 ms
Multi-Slice Mode	Sequential
Series	Ascending

Geometry - Common

Concatenations	1
----------------	---

Geometry - AutoAlign

Slab Group	1
Position	Isocenter
Orientation	Transversal
Phase Encoding Dir.	A >> P
AutoAlign	---
Initial Position	Isocenter
L	0.0 mm
P	0.0 mm
H	0.0 mm
Initial Orientation	Transversal
Initial Rotation	0.00 deg

Geometry - Saturation

Special Saturation	None
--------------------	------

Geometry - Tim Planning Suite

Set-n-Go Protocol	Off
Table Position	0 mm
Table Position	H

System - Miscellaneous

Coil Selection	Manual
Radial Sorting	Off
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combination	Adaptive Combine
Matrix Optimization	Off

System - Adjustments

Adjustment Strategy	Standard
B0 Shim	Tune up
B1 Shim	TrueForm
Adjustment Tolerance	Auto
Confirm Frequency	Never
Assume Silicone	Off

System - Adjust Volume

Position	Isocenter
Orientation	Transversal
Rotation	0.00 deg
A >> P	263 mm
R >> L	350 mm
F >> H	350 mm
Reset	Off

System - pTx

B1 Shim	TrueForm
Excitation	Slab-sel.

System - Tx/Rx

Frequency 1H	297.135464 MHz
? Ref. Amplitude 1H	0.000 V
Reset	Off
Correction Factor	1.00
Image Scaling	1.000

Physio - Signal

1st Signal/Mode	None
TR	50.0 ms
Segments	3
Concatenations	1

Sequence - Part 1

Sequence Name	WIP_epi
Dimension	3D
Excitation	Slab-sel.
RF Pulse Type	Normal
Gradient Mode	Fast*
Flow Compensation	None
Reordering	Linear
Bandwidth	1202 Hz/Px
Echo Spacing	1.00 ms
Free Echo Spacing	Off
Segments	3
EPI Factor	25

Sequence - Part 2

Introduction	Off
RF Spoiling	On
Phase Correction	Internal

Sequence - Special

Prep/dummy scan duration	2 sec
Selective RF Duration	1024 us
Desired RF TBW	8.0
Flyback	No
QuickSat	1
Excitation Orientation	SS - Slice Selective
Ernst Angle T1	1220 ms
Reduce Slab Thickness?	Weak
skipped-CAIPI?	Off
Log Signals	Off

Sequence - Assistant

SAR Assistant	Off
---------------	-----