

Table of contents

\\BTM5 TERRA 7T

<table border="1"> <tr> <td>BRAIN BWH</td> </tr> </table>	BRAIN BWH
BRAIN BWH	

<table border="1"> <tr> <td>Camilo_Peds</td> </tr> </table>	Camilo_Peds
Camilo_Peds	

<table border="1"> <tr> <td>Epilepsy</td> </tr> </table>	Epilepsy
Epilepsy	

<table border="1"> <tr> <td> <table border="1"> <tr> <td>AAhead_scout</td> </tr> </table> </td> </tr> </table>	<table border="1"> <tr> <td>AAhead_scout</td> </tr> </table>	AAhead_scout
<table border="1"> <tr> <td>AAhead_scout</td> </tr> </table>	AAhead_scout	
AAhead_scout		

<table border="1"> <tr> <td> <table border="1"> <tr> <td>t2_flair_DRB_pTx_sag_p10_0.8mm_MT</td> </tr> </table> </td> </tr> </table>	<table border="1"> <tr> <td>t2_flair_DRB_pTx_sag_p10_0.8mm_MT</td> </tr> </table>	t2_flair_DRB_pTx_sag_p10_0.8mm_MT
<table border="1"> <tr> <td>t2_flair_DRB_pTx_sag_p10_0.8mm_MT</td> </tr> </table>	t2_flair_DRB_pTx_sag_p10_0.8mm_MT	
t2_flair_DRB_pTx_sag_p10_0.8mm_MT		

<table border="1"> <tr> <td> <table border="1"> <tr> <td>t1_mp2rage_sag_p3_0.75mm</td> </tr> </table> </td> </tr> </table>	<table border="1"> <tr> <td>t1_mp2rage_sag_p3_0.75mm</td> </tr> </table>	t1_mp2rage_sag_p3_0.75mm
<table border="1"> <tr> <td>t1_mp2rage_sag_p3_0.75mm</td> </tr> </table>	t1_mp2rage_sag_p3_0.75mm	
t1_mp2rage_sag_p3_0.75mm		

<table border="1"> <tr> <td> <table border="1"> <tr> <td>Cor</td> <td>T2</td> <td>2mm</td> </tr> </table> </td> </tr> </table>	<table border="1"> <tr> <td>Cor</td> <td>T2</td> <td>2mm</td> </tr> </table>	Cor	T2	2mm
<table border="1"> <tr> <td>Cor</td> <td>T2</td> <td>2mm</td> </tr> </table>	Cor	T2	2mm	
Cor	T2	2mm		

<table border="1"> <tr> <td> <table border="1"> <tr> <td>Post</td> <td>Sag</td> <td>t1fs</td> <td>space</td> </tr> </table> </td> </tr> </table>	<table border="1"> <tr> <td>Post</td> <td>Sag</td> <td>t1fs</td> <td>space</td> </tr> </table>	Post	Sag	t1fs	space
<table border="1"> <tr> <td>Post</td> <td>Sag</td> <td>t1fs</td> <td>space</td> </tr> </table>	Post	Sag	t1fs	space	
Post	Sag	t1fs	space		

<table border="1"> <tr> <td> <table border="1"> <tr> <td>t1_space_DRB_pTx_sag__matchproduct</td> </tr> </table> </td> </tr> </table>	<table border="1"> <tr> <td>t1_space_DRB_pTx_sag__matchproduct</td> </tr> </table>	t1_space_DRB_pTx_sag__matchproduct
<table border="1"> <tr> <td>t1_space_DRB_pTx_sag__matchproduct</td> </tr> </table>	t1_space_DRB_pTx_sag__matchproduct	
t1_space_DRB_pTx_sag__matchproduct		

<table border="1"> <tr> <td> <table border="1"> <tr> <td>t1_space_DRB_pTx_sag_p14_0.55</td> </tr> </table> </td> </tr> </table>	<table border="1"> <tr> <td>t1_space_DRB_pTx_sag_p14_0.55</td> </tr> </table>	t1_space_DRB_pTx_sag_p14_0.55
<table border="1"> <tr> <td>t1_space_DRB_pTx_sag_p14_0.55</td> </tr> </table>	t1_space_DRB_pTx_sag_p14_0.55	
t1_space_DRB_pTx_sag_p14_0.55		

<table border="1"> <tr> <td> <table border="1"> <tr> <td>t1_space_DRB_NopTx_sag_p14_0.6</td> </tr> </table> </td> </tr> </table>	<table border="1"> <tr> <td>t1_space_DRB_NopTx_sag_p14_0.6</td> </tr> </table>	t1_space_DRB_NopTx_sag_p14_0.6
<table border="1"> <tr> <td>t1_space_DRB_NopTx_sag_p14_0.6</td> </tr> </table>	t1_space_DRB_NopTx_sag_p14_0.6	
t1_space_DRB_NopTx_sag_p14_0.6		

\\BTM5 TERRA 7T\BRAIN BWH\Camilo_Peds\Epilepsy\AAhead_scout

TA: 14 sec Coil Selection: Auto Voxel Size: 1.6×1.6×1.6 mm³ Acc:: 3 Rel. SNR: 1.00**Properties**

Start measurement without further preparation	Off
Wait for User to Start	Off
Start measurements	Single Measurement
Prio Recon	Off
Auto Open Inline Display	Off
Auto Close Inline Display	Off
Load Images to MR View&GO	On
Auto Store Images	On
Disable auto transfer to PACS	Off
Load Images to Stamp Segments	Off
Load Images to Graphic Segments	On
Graphic segment	Default
Inline Movie	Off

Routine

Slab Group	1
Slabs	1
Distance Factor	20 %
Position	L0.0 A30.0 H0.0 mm
Orientation	Sagittal
Phase Encoding Dir.	A >> P
Slices per Slab	128
Phase Oversampling	0 %
Slice Oversampling	0.0 %
FOV Read	260 mm
FOV Phase	100.0 %
Slice Thickness	1.6 mm
TR	3.25 ms
TE	1.53 ms
Averages	1
Concatenations	1
AutoAlign	Head

Contrast - Common

TR	3.25 ms
TE	1.53 ms
Flip Angle	16 deg
Fat-Water Contrast	Standard
Contrasts	1
Reconstruction	Magnitude

Contrast - Dynamic

Dynamic Mode	Standard
Measurements	1
Time to Center	6.3 s

Resolution - Common

FOV Read	260 mm
FOV Phase	100.0 %

Resolution - Common

Slice Thickness	1.6 mm
Base Resolution	160
Phase Resolution	100 %
Slice Resolution	69 %
Trajectory	Cartesian

Resolution - Acceleration

Acceleration Mode	GRAPPA
Reference Scans	Integrated
Acceleration Factor PE	3
Reference Lines PE	24
Acceleration Factor 3D	1
Phase Partial Fourier	6/8
Slice Partial Fourier	6/8
Asymmetric Echo	Weak

Resolution - Filter

Raw Filter	Off
Elliptical Filter	Off
Distortion Correction	3D
Normalize	B1 Filter
Image Filter	Off

Geometry - Common

Slab Group	1
Slabs	1
Distance Factor	20 %
Position	L0.0 A30.0 H0.0 mm
Orientation	Sagittal
Phase Encoding Dir.	A >> P
Slices per Slab	128
Phase Oversampling	0 %
Slice Oversampling	0.0 %
FOV Read	260 mm
FOV Phase	100.0 %
Slice Thickness	1.6 mm
TR	3.25 ms
Multi-Slice Mode	Sequential
Series	Ascending
Concatenations	1

Geometry - AutoAlign

Slab Group	1
Position	L0.0 A30.0 H0.0 mm
Orientation	Sagittal
Phase Encoding Dir.	A >> P
AutoAlign	Head
Initial Position	Isocenter
L	0.0 mm
P	0.0 mm

Geometry - AutoAlign

H	0.0 mm
Initial Orientation	Transversal
Initial Rotation	0.00 deg

Geometry - Tim Planning Suite

Set-n-Go Protocol	Off
Table Position	0 mm
Table Position	H

System - Miscellaneous

Coil Selection	Auto Coil Select
Radial Sorting	Off
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combination	Sum of Squares
Matrix Optimization	Off
Coil Focus	Flat

System - Adjustments

Adjustment Strategy	Standard
B0 Shim	Tune up
B1 Shim	TrueForm
Adjustment Tolerance	Auto
Adjust with Body Coil	Off
Confirm Frequency	Never
Assume Silicone	Off

System - Adjust Volume

Position	Isocenter
Orientation	Transversal
Rotation	0.00 deg
A >> P	263 mm
R >> L	350 mm
F >> H	350 mm
Reset	Off

System - pTx

B1 Shim	TrueForm
Excitation	Non-sel.

System - Tx/Rx

Frequency 1H	297.208927 MHz
? Ref. Amplitude 1H	0.000 V
Reset	Off
Image Scaling	1.000

Physio - PACE

Resp. Control	Off
Concatenations	1

Inline - Subtraction

Subtract	Off
----------	-----

Inline - Subtraction

Measurements	1
StdDev	Off
Save Original Images	On

Inline - MIP

MIP Sag	Off
MIP Cor	Off
MIP Tra	Off
MIP Time	Off
Radial MIP	Off
Save Original Images	On
MPR Sag	Off
MPR Cor	Off
MPR Tra	Off

Inline - Composing**Inline - MapIt**

MapIt	None
Flip Angle	16 deg
Measurements	1
Contrasts	1
TE	1.53 ms
TR	3.25 ms
Save Original Images	On

Sequence - Part 1

Sequence Name	fl
Dimension	3D
Excitation	Non-sel.
RF Pulse Type	Fast
Gradient Mode	Normal
Bandwidth	540 Hz/Px
Asymmetric Echo	Weak

Sequence - Part 2

Introduction	On
RF Spoiling	On
Breast Application	Off
Phase Enc. Order	Automatic

Sequence - Assistant

SAR Assistant	Off
---------------	-----

\\BTM5 TERRA 7T\BRAIN BWH\Camilo_Peds\Epilepsy\t2_flair_DRB_pTx_sag_p10_0.8mm_MT

TA: 4:26 min Coil Selection: Auto Voxel Size: 0.8×0.8×0.8 mm³ Acc:: 10.0 Rel. SNR: 1.00**Properties**

Start measurement without further preparation	On
Wait for User to Start	Off
Start measurements	Single Measurement
Prio Recon	Off
Auto Open Inline Display	Off
Auto Close Inline Display	Off
Load Images to MR View&GO	On
Auto Store Images	On
Disable auto transfer to PACS	Off
Load Images to Stamp Segments	Off
Load Images to Graphic Segments	Off
Graphic segment	Default
Inline Movie	Off

Routine

Slab Group	1
Slabs	1
Position	R7.3 A7.6 F15.9 mm
Orientation	Sagittal
Phase Encoding Dir.	A >> P
Slices per Slab	224
Phase Oversampling	15 %
Slice Oversampling	0.0 %
FOV Read	220 mm
FOV Phase	100.0 %
Slice Thickness	0.80 mm
TR	8000.0 ms
TE	302.00 ms
Averages	1.0
Concatenations	1
AutoAlign	---

Contrast - Common

TR	8000.0 ms
TE	302.00 ms
Magn. Preparation	Non-sel. T2 prep. IR
T1	2250 ms
T2 Prep. Duration	125 ms
Flip Angle Mode	T2 Var
Fat-Water Contrast	Fat Saturation
Fat Saturation	Strong
Dark Blood	Off
Blood Suppression	Off
Excitation offset frequency	0 Hz
Wrap-up Magn.	None
Refoc. offset frequency	0 Hz
Reconstruction	Magnitude

Contrast - Dynamic

goKART	Off
--------	-----

Contrast - Dynamic

Measurements	1
Multiple Series	Each Measurement
Reordering	Linear

Contrast - ASL

MTC	On
MTC frequency sign	Positive
MTC # RF	7
MTC pulse duration	7680.0 µs
MTC frequency offset	1200 Hz
MTC pause between RF	0 ms

Resolution - Common

FOV Read	220 mm
FOV Phase	100.0 %
Slice Thickness	0.80 mm
Base Resolution	272
Phase Resolution	100 %
Slice Resolution	100 %
Interpolation	Off

Resolution - Acceleration

Acceleration Mode	CS
Total Factor	10.0
Reference Scans	Integrated
Acceleration Factor PE	3
Reference Lines PE	24
Acceleration Factor 3D	2
Reference Lines 3D	24
Reordering Shift 3D	1
Deep Resolve	Off
Phase Partial Fourier	Allowed
Slice Partial Fourier	Off
Elliptical Scanning	On

Resolution - Filter

Raw Filter	On
Elliptical Filter	Off
Distortion Correction	3D
Normalize	B1 Filter
Image Filter	Off

Geometry - Common

Slab Group	1
Slabs	1
Position	R7.3 A7.6 F15.9 mm
Orientation	Sagittal
Phase Encoding Dir.	A >> P
Slices per Slab	224
Phase Oversampling	15 %

Geometry - Common

Slice Oversampling	0.0 %
FOV Read	220 mm
FOV Phase	100.0 %
Slice Thickness	0.80 mm
TR	8000.0 ms
Concatenations	1

Geometry - AutoAlign

Slab Group	1
Position	R7.3 A7.6 F15.9 mm
Orientation	Sagittal
Phase Encoding Dir.	A >> P
AutoAlign	---
Initial Position	R7.3 A7.6 F15.9
R	7.3 mm
A	7.6 mm
F	15.9 mm
Initial Orientation	Sagittal
Initial Rotation	0.00 deg

Geometry - Navigator**Geometry - Saturation**

Special Saturation	None
--------------------	------

Geometry - Tim Planning Suite

Set-n-Go Protocol	Off
Table Position	0 mm
Table Position	H

System - Miscellaneous

Coil Selection	Auto Coil Select
Radial Sorting	Off
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combination	Adaptive Combine
Matrix Optimization	Off
Coil Focus	Flat

System - Adjustments

Adjustment Strategy	Standard
B0 Shim	Brain
B1 Shim	TrueForm
Adjustment Tolerance	Auto
Adjust with Body Coil	Off
Confirm Frequency	Never
Assume Silicone	Off
Adj Freq. protocol related	On

System - Adjust Volume

Position	R7.3 A7.6 F15.9 mm
Orientation	Sagittal

System - Adjust Volume

Rotation	0.00 deg
A >> P	220 mm
F >> H	220 mm
R >> L	180 mm
Reset	Off

System - pTx

B1 Shim	TrueForm
Excitation	Non-sel.
pTx Pulse	1

System - Tx/Rx

Frequency 1H	297.208927 MHz
? Ref. Amplitude 1H	0.000 V
Reset	Off
Image Scaling	1.000
Gain	High

Physio - Signal

1st Signal/Mode	None
Trigger Delay	0 ms
TR	8000.0 ms
Rd crusher scale	893
Concatenations	1

Physio - Cardiac

Fat-Water Contrast	Fat Saturation
Magn. Preparation	Non-sel. T2 prep. IR
TI	2250 ms
T2 Prep. Duration	125 ms
Dark Blood	Off
FOV Read	220 mm
FOV Phase	100.0 %
Phase Resolution	100 %

Physio - PACE

Resp. Control	Off
Concatenations	1

Inline - Subtraction

Subtract	Off
Measurements	1
StdDev	Off
Save Original Images	On

Inline - MIP

MIP Sag	Off
MIP Cor	Off
MIP Tra	Off
MIP Time	Off
Radial MIP	Off
Save Original Images	On
MPR Sag	Off
MPR Cor	Off

Inline - MIP

MPR Tra	Off
---------	-----

Inline - Composing**Sequence - Part 1**

Sequence Name	wip-spc-t2p+ir
Dimension	3D
Excitation	Non-sel.
RF Pulse Type	Low SAR
Refoc. duration	1000
Gradient Mode	Normal
Flow Compensation	None
Reordering	Linear
Bandwidth	735 Hz/Px
Echo Spacing	4.08 ms
Concomitant field corr. R	Product
Turbo Factor	170
Echo Train Duration	645 ms
Rd crusher scale	893

Sequence - Part 2

Introduction	On
Rel. crusher mom. P	0 %
Rel crusher mom. S	0 %

Sequence - Special

IR refoc pulse	IR12800H180.IR180_HS3
IR RF duration	12800 us
IR HS mu	5.00
IR HS BWTP	10.99
IR HS B1 peak	22.5 μ T
T2-Prep. excit/tip-up pulse	RECT
T2-Prep. refoc pulse	IR10240H180.IR180_36B1_1
T2-Prep. # refoc. RF	4
T2-Prep. excit. duration	800 us
T2-Prep. refoc. duration	10240 us
T2-Prep. spoil mom. $10^3 \times$	-71.20 usmT/m
T2-Prep. B1 scale	1.00
T2-Prep. crusher mom. rel.	0.00
T2-Prep. Refoc shift	0.00 %
Excit. flip angle	90 deg
Blood Suppression	Off
Trajectory Optimization	8

Sequence - Assistant

SAR Assistant	Off
Allowed Delay	0 s

\\BTM5 TERRA 7T\\BRAIN BWH\\Camilo_Peds\\Epilepsy\\t1_mp2rage_sag_p3_0.75mm

TA: 8:50 min Coil Selection: Auto Voxel Size: 0.8×0.8×0.8 mm³ Acc:: 3 Rel. SNR: 1.00**Properties**

Start measurement without further preparation	On
Wait for User to Start	Off
Start measurements	Single Measurement
Prio Recon	Off
Auto Open Inline Display	Off
Auto Close Inline Display	Off
Load Images to MR View&GO	On
Auto Store Images	On
Disable auto transfer to PACS	Off
Load Images to Stamp Segments	Off
Load Images to Graphic Segments	Off
Graphic segment	Default
Inline Movie	Off

Routine

Slab Group	1
Slabs	1
Distance Factor	50 %
Position	R7.3 A7.6 F15.9 mm
Orientation	Sagittal
Phase Encoding Dir.	A >> P
Slices per Slab	224
Phase Oversampling	0 %
Slice Oversampling	7.1 %
FOV Read	240 mm
FOV Phase	93.8 %
Slice Thickness	0.75 mm
TR	4300.0 ms
TE	1.99 ms
Averages	1
Concatenations	1
AutoAlign	---

Contrast - Common

TR	4300.0 ms
TE	1.99 ms
Magn. Preparation	Non-sel. IR
TI 1	840 ms
TI 2	2370 ms
Flip Angle 1	5 deg
Flip Angle 2	6 deg
Fat-Water Contrast	Fast Water Excitation
Dark Blood	Off
Contrasts	1
Reconstruction	Magnitude

Contrast - Dynamic

Dynamic Mode	Standard
Measurements	1
Multiple Series	Each Measurement

Contrast - Dynamic

Reordering	Linear
------------	--------

Resolution - Common

FOV Read	240 mm
FOV Phase	93.8 %
Slice Thickness	0.75 mm
Base Resolution	320
Phase Resolution	100 %
Slice Resolution	100 %
Interpolation	Off

Resolution - Acceleration

Acceleration Mode	GRAPPA
Reference Scans	Integrated
Acceleration Factor PE	3
Reference Lines PE	37
Acceleration Factor 3D	1
Phase Partial Fourier	Off
Slice Partial Fourier	6/8
Asymmetric Echo	Allowed
Elliptical Scanning	Off

Resolution - Filter

Raw Filter	Off
Elliptical Filter	Off
Distortion Correction	2D
Normalize	Off
Image Filter	Off

Geometry - Common

Slab Group	1
Slabs	1
Distance Factor	50 %
Position	R7.3 A7.6 F15.9 mm
Orientation	Sagittal
Phase Encoding Dir.	A >> P
Slices per Slab	224
Phase Oversampling	0 %
Slice Oversampling	7.1 %
FOV Read	240 mm
FOV Phase	93.8 %
Slice Thickness	0.75 mm
TR	4300.0 ms
Multi-Slice Mode	Single Shot
Series	Interleaved
Concatenations	1

Geometry - AutoAlign

Slab Group	1
Position	R7.3 A7.6 F15.9 mm

Geometry - AutoAlign

Orientation	Sagittal
Phase Encoding Dir.	A >> P
AutoAlign	---
Initial Position	R7.3 A7.6 F15.9
R	7.3 mm
A	7.6 mm
F	15.9 mm
Initial Orientation	Sagittal
Initial Rotation	0.00 deg

Geometry - Navigator**Geometry - Tim Planning Suite**

Set-n-Go Protocol	Off
Table Position	0 mm
Table Position	H

System - Miscellaneous

Coil Selection	Auto Coil Select
Radial Sorting	Off
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combination	Sum of Squares
Matrix Optimization	Off
Coil Focus	Flat

System - Adjustments

Adjustment Strategy	Standard
B0 Shim	Brain
B1 Shim	TrueForm
Adjustment Tolerance	Auto
Adjust with Body Coil	Off
Confirm Frequency	Never
Assume Silicone	Off

System - Adjust Volume

Position	R7.3 A7.6 F15.9 mm
Orientation	Sagittal
Rotation	0.00 deg
A >> P	225 mm
F >> H	240 mm
R >> L	168 mm
Reset	Off

System - pTx

B1 Shim	TrueForm
Excitation	Non-sel.

System - Tx/Rx

Frequency 1H	297.208927 MHz
? Ref. Amplitude 1H	0.000 V
Reset	Off

System - Tx/Rx

Image Scaling	1.000
---------------	-------

Physio - Signal

1st Signal/Mode	None
TR	4300.0 ms
Concatenations	1

Physio - Cardiac

Fat-Water Contrast	Fast Water Excitation
Magn. Preparation	Non-sel. IR
TI 1	840 ms
TI 2	2370 ms
Dark Blood	Off
FOV Read	240 mm
FOV Phase	93.8 %
Phase Resolution	100 %

Physio - PACE

Resp. Control	Off
Concatenations	1

Inline - Subtraction

Subtract	Off
Measurements	1
StdDev	Off
Save Original Images	On

Inline - MIP

MIP Sag	Off
MIP Cor	Off
MIP Tra	Off
MIP Time	Off
Radial MIP	Off
Save Original Images	On
MPR Sag	Off
MPR Cor	Off
MPR Tra	Off

Inline - Composing**Inline - MapIt**

MapIt	None
Flip Angle 1	5 deg
Flip Angle 2	6 deg
Measurements	1
Contrasts	1
TE	1.99 ms
TR	4300.0 ms
Save Original Images	On

Sequence - Part 1

Sequence Name	tfl
Dimension	3D
Excitation	Non-sel.

Sequence - Part 1

RF Pulse Type	Fast
Gradient Mode	Fast*
Flow Compensation	None
Reordering	Linear
Bandwidth	250 Hz/Px
Echo Spacing	6.44 ms
Asymmetric Echo	Allowed
Turbo Factor	180

Sequence - Part 2

Introduction	On
RF Spoiling	On
Incr. Gradient Spoiling	Off

Sequence - Assistant

SAR Assistant	Off
---------------	-----

\\BTM5 TERRA 7T\BRAIN BWH\Camilo_Peds\Epilepsy\Cor T2 2mmTA: 4:13 min Coil Selection: Auto Voxel Size: 0.1×0.1×2.0 mm³ Acc:: 3 Rel. SNR: 1.00**Properties**

Start measurement without further preparation	On
Wait for User to Start	On
Start measurements	Single Measurement
Prio Recon	Off
Auto Open Inline Display	Off
Auto Close Inline Display	Off
Load Images to MR View&GO	On
Auto Store Images	On
Disable auto transfer to PACS	Off
Load Images to Stamp Segments	Off
Load Images to Graphic Segments	Off
Graphic segment	Default
Inline Movie	Off

Resolution - Common

FOV Read	160 mm
FOV Phase	100.0 %
Slice Thickness	2.0 mm
Base Resolution	560
Phase Resolution	75 %
Interpolation	On

Resolution - Acceleration

Acceleration Mode	GRAPPA
Reference Scans	TSE/Separate
Acceleration Factor PE	3
Reference Lines PE	32
Deep Resolve	On
Phase Partial Fourier	Off

Routine

Slice Group	1
Slices	52
Distance Factor	14 %
Position	R7.7 A12.6 F7.6 mm
Orientation	C > T-19.5 > S-0.2
Phase Encoding Dir.	R >> L
Phase Oversampling	50 %
FOV Read	160 mm
FOV Phase	100.0 %
Slice Thickness	2.0 mm
TR	4890.0 ms
TE	44.00 ms
Averages	1
Concatenations	2
AutoAlign	Head > Basis

Resolution - Filter

Raw Filter	Off
Elliptical Filter	Off
Distortion Correction	2D
Normalize	Off
Image Filter	Off

Geometry - Common

Slice Group	1
Slices	52
Distance Factor	14 %
Position	R7.7 A12.6 F7.6 mm
Orientation	C > T-19.5 > S-0.2
Phase Encoding Dir.	R >> L
Phase Oversampling	50 %
FOV Read	160 mm
FOV Phase	100.0 %
Slice Thickness	2.0 mm
TR	4890.0 ms
Multi-Slice Mode	Interleaved
Series	Interleaved
Concatenations	2

Contrast - Common

TR	4890.0 ms
TE	44.00 ms
TD	0.00 ms
MTC	Off
Magn. Preparation	None
Flip Angle Mode	Hyperecho
Flip Angle	120 deg
Fat-Water Contrast	Standard
Dark Blood	Off
Contrasts	1
Wrap-up Magn.	None
Reconstruction	Magnitude

Contrast - Dynamic

Dynamic Mode	Standard
Measurements	1
Multiple Series	Each Measurement

Geometry - AutoAlign

Slice Group	1
Position	R7.7 A12.6 F7.6 mm
Orientation	C > T-19.5 > S-0.2
Phase Encoding Dir.	R >> L
AutoAlign	Head > Basis
Initial Position	R1.8 P1.9 H27.6
R	1.8 mm
P	1.9 mm
H	27.6 mm
Initial Orientation	C > T
C > T	-14.90

Geometry - AutoAlign

> S	-0.40
Initial Rotation	0.03 deg

Geometry - Navigator**Geometry - Saturation**

Special Saturation	None
--------------------	------

Geometry - Tim Planning Suite

Set-n-Go Protocol	Off
Table Position	0 mm
Table Position	H

System - Miscellaneous

Coil Selection	Auto Coil Select
Radial Sorting	Off
MSMA	S - C - T
Sagittal	L >> R
Coronal	P >> A
Transversal	F >> H
Coil Combination	Adaptive Combine
Matrix Optimization	Off
Coil Focus	Flat

System - Adjustments

Adjustment Strategy	Standard
B0 Shim	Brain
B1 Shim	TrueForm
Adjustment Tolerance	Auto
Adjust with Body Coil	Off
Confirm Frequency	Never
Assume Silicone	Off

System - Adjust Volume

Position	R7.7 A12.6 F7.6 mm
Orientation	C > T-19.5 > S-0.2
Rotation	2.14 deg
R >> L	160 mm
F >> H	160 mm
A >> P	119 mm
Reset	Off

System - pTx

B1 Shim	TrueForm
LR Balancing	Off

System - Tx/Rx

Frequency 1H	297.208927 MHz
? Ref. Amplitude 1H	0.000 V
Reset	Off
Image Scaling	1.000

Physio - Signal

1st Signal/Mode	None
-----------------	------

Physio - Signal

TR	4890.0 ms
Concatenations	2

Physio - Cardiac

Fat-Water Contrast	Standard
Magn. Preparation	None
Dark Blood	Off
FOV Read	160 mm
FOV Phase	100.0 %
Phase Resolution	75 %
Motion Correction	None

Physio - PACE

Resp. Control	Off
Concatenations	2

Inline - Subtraction

Subtract	Off
Measurements	1
StdDev	Off
Motion Correction	None
Save Original Images	On

Inline - MIP

MIP Sag	Off
MIP Cor	Off
MIP Tra	Off
MIP Time	Off
Radial MIP	Off
Save Original Images	On
MPR Sag	Off
MPR Cor	Off
MPR Tra	Off

Inline - Composing**Sequence - Part 1**

Sequence Name	tse_rr
Dimension	2D
RF Pulse Type	Low SAR
Gradient Mode	Fast
Flow Compensation	Read
Bandwidth	255 Hz/Px
Echo Spacing	11.1 ms
Free Echo Spacing	Off
Define	Turbo Factor
Turbo Factor	9
Echo Trains per Slice	24

Sequence - Part 2

Introduction	On
Incr. Gradient Spoiling	Off
Phase Correction	Automatic
Compensate T2 Decay	Off

Sequence - Part 2

Fast Mode	Off
WARP	Off
Acoustic noise reduction	Off
Reduce Motion Sens.	On
Motion Correction	None

Sequence - Assistant

SAR Assistant	Off
Allowed Delay	60 s

\\BTM5 TERRA 7T\BRAIN BWH\Camilo_Peds\Epilepsy\Post Sag t1fs spaceTA: 6:28 min Coil Selection: Auto Voxel Size: 0.3×0.3×0.5 mm³ Acc:: 8 Rel. SNR: 1.00**Properties**

Start measurement without further preparation	On
Wait for User to Start	Off
Start measurements	Single Measurement
Prio Recon	Off
Auto Open Inline Display	Off
Auto Close Inline Display	Off
Load Images to MR View&GO	On
Auto Store Images	On
Disable auto transfer to PACS	Off
Load Images to Stamp Segments	Off
Load Images to Graphic Segments	Off
Graphic segment	Default
Inline Movie	Off

Routine

Slab Group	1
Slabs	1
Position	R5.3 P1.6 F4.1 mm
Orientation	S > T-2.0 > C0.5
Phase Encoding Dir.	A >> P
Slices per Slab	352
Phase Oversampling	10 %
Slice Oversampling	0.0 %
FOV Read	136 mm
FOV Phase	132.8 %
Slice Thickness	0.50 mm
TR	1160.0 ms
TE	23.00 ms
Averages	1.0
Concatenations	1
AutoAlign	Head > Brain

Contrast - Common

TR	1160.0 ms
TE	23.00 ms
MTC	Off
Magn. Preparation	None
Flip Angle Mode	T1 Var
Fat-Water Contrast	SPAIR
Fat Saturation	Strong
Dark Blood	Off
Blood Suppression	Off
Wrap-up Magn.	Restore
Reconstruction	Magnitude

Contrast - Dynamic

Dynamic Mode	Standard
Measurements	1
Multiple Series	Each Measurement
Reordering	Linear

Resolution - Common

FOV Read	136 mm
FOV Phase	132.8 %
Slice Thickness	0.50 mm
Base Resolution	256
Phase Resolution	100 %
Slice Resolution	80 %
Interpolation	On

Resolution - Acceleration

Acceleration Mode	CAIPIRINHA
Total Factor	8
Reference Scans	Integrated
Acceleration Factor PE	2
Reference Lines PE	64
Acceleration Factor 3D	4
Reference Lines 3D	68
Reordering Shift 3D	2
Deep Resolve	Off
Phase Partial Fourier	Off
Slice Partial Fourier	Off
Elliptical Scanning	On

Resolution - Filter

Raw Filter	On
Elliptical Filter	Off
Distortion Correction	2D
Normalize	B1 Filter
Image Filter	Off

Geometry - Common

Slab Group	1
Slabs	1
Position	R5.3 P1.6 F4.1 mm
Orientation	S > T-2.0 > C0.5
Phase Encoding Dir.	A >> P
Slices per Slab	352
Phase Oversampling	10 %
Slice Oversampling	0.0 %
FOV Read	136 mm
FOV Phase	132.8 %
Slice Thickness	0.50 mm
TR	1160.0 ms
Concatenations	1

Geometry - AutoAlign

Slab Group	1
Position	R5.3 P1.6 F4.1 mm
Orientation	S > T-2.0 > C0.5
Phase Encoding Dir.	A >> P
AutoAlign	Head > Brain

Geometry - AutoAlign

Initial Position	L0.0 A2.5 F0.4
L	0.0 mm
A	2.5 mm
F	0.4 mm
Initial Orientation	Sagittal
Initial Rotation	0.00 deg

Geometry - Navigator**Geometry - Saturation**

Special Saturation	None
--------------------	------

Geometry - Tim Planning Suite

Set-n-Go Protocol	Off
Table Position	0 mm
Table Position	H

System - Miscellaneous

Coil Selection	ACS All but spine
Radial Sorting	Off
MSMA	S - C - T
Sagittal	L >> R
Coronal	P >> A
Transversal	F >> H
Coil Combination	Adaptive Combine
Matrix Optimization	Off
Coil Focus	Flat

System - Adjustments

Adjustment Strategy	Standard
B0 Shim	Standard
B1 Shim	TrueForm
Adjustment Tolerance	Auto
Adjust with Body Coil	Off
Confirm Frequency	Never
Assume Silicone	Off

System - Adjust Volume

Position	R5.3 P1.6 F4.1 mm
Orientation	S > T-2.0 > C0.5
Rotation	95.63 deg
F >> H	136 mm
A >> P	181 mm
R >> L	176 mm
Reset	Off

System - pTx

B1 Shim	TrueForm
Excitation	Non-sel.

System - Tx/Rx

Frequency 1H	297.208927 MHz
? Ref. Amplitude 1H	0.000 V
Reset	Off

System - Tx/Rx

Image Scaling	1.000
Gain	High

Physio - Signal

1st Signal/Mode	None
Trigger Delay	0 ms
TR	1160.0 ms
Concatenations	1

Physio - Cardiac

Fat-Water Contrast	SPAIR
Magn. Preparation	None
Dark Blood	Off
FOV Read	136 mm
FOV Phase	132.8 %
Phase Resolution	100 %

Physio - PACE

Resp. Control	Off
Concatenations	1

Inline - Subtraction

Subtract	Off
Measurements	1
StdDev	Off
Save Original Images	On

Inline - MIP

MIP Sag	Off
MIP Cor	Off
MIP Tra	Off
MIP Time	Off
Radial MIP	Off
Save Original Images	On
MPR Sag	Off
MPR Cor	Off
MPR Tra	Off

Inline - Composing**Sequence - Part 1**

Sequence Name	spcR
Dimension	3D
Excitation	Non-sel.
RF Pulse Type	Normal
Gradient Mode	Fast
Flow Compensation	None
Reordering	Linear
Bandwidth	610 Hz/Px
Echo Spacing	3.84 ms
Turbo Factor	43
Echo Train Duration	188 ms

Sequence - Part 2

Introduction	On
--------------	----

Sequence - Assistant

SAR Assistant	Off
Allowed Delay	30 s

\\BTM5 TERRA 7T\\BRAIN BWH\\Camilo_Peds\\Epilepsy\\t1_space_DRB_pTx_sag_matchproduct

TA: 6:20 min Coil Selection: Auto Voxel Size: 0.5×0.5×0.6 mm³ Acc.: 14.0 Rel. SNR: 1.00**Properties**

Start measurement without further preparation	On
Wait for User to Start	Off
Start measurements	Single Measurement
Prio Recon	Off
Auto Open Inline Display	Off
Auto Close Inline Display	Off
Load Images to MR View&GO	On
Auto Store Images	On
Disable auto transfer to PACS	Off
Load Images to Stamp Segments	Off
Load Images to Graphic Segments	Off
Graphic segment	Default
Inline Movie	Off

Routine

Slab Group	1
Slabs	1
Position	R4.7 A13.6 F19.4 mm
Orientation	Sagittal
Phase Encoding Dir.	A >> P
Slices per Slab	288
Phase Oversampling	15 %
Slice Oversampling	0.0 %
FOV Read	136 mm
FOV Phase	132.8 %
Slice Thickness	0.55 mm
TR	900.0 ms
TE	25.00 ms
Averages	1.0
Concatenations	1
AutoAlign	---

Contrast - Common

TR	900.0 ms
TE	25.00 ms
Magn. Preparation	None
Flip Angle Mode	T1 Var
Fat-Water Contrast	Standard
Dark Blood	Off
Blood Suppression	Off
Excitation offset frequency	0 Hz
Wrap-up Magn.	Restore
Refoc. offset frequency	0 Hz
Reconstruction	Magnitude

Contrast - Dynamic

goKART	Off
Measurements	1
Multiple Series	Each Measurement
Reordering	Linear

Contrast - ASL

MTC	Off
MTC frequency sign	Positive

Resolution - Common

FOV Read	136 mm
FOV Phase	132.8 %
Slice Thickness	0.55 mm
Base Resolution	256
Phase Resolution	100 %
Slice Resolution	80 %
Interpolation	Off

Resolution - Acceleration

Acceleration Mode	CS
Total Factor	14.0
Reference Scans	Integrated
Acceleration Factor PE	3
Reference Lines PE	24
Acceleration Factor 3D	2
Reference Lines 3D	24
Reordering Shift 3D	1
Deep Resolve	Off
Phase Partial Fourier	Off
Slice Partial Fourier	Off
Elliptical Scanning	On

Resolution - Filter

Raw Filter	On
Elliptical Filter	Off
Distortion Correction	3D
Normalize	B1 Filter
Image Filter	Off

Geometry - Common

Slab Group	1
Slabs	1
Position	R4.7 A13.6 F19.4 mm
Orientation	Sagittal
Phase Encoding Dir.	A >> P
Slices per Slab	288
Phase Oversampling	15 %
Slice Oversampling	0.0 %
FOV Read	136 mm
FOV Phase	132.8 %
Slice Thickness	0.55 mm
TR	900.0 ms
Concatenations	1

Geometry - AutoAlign

Slab Group	1
------------	---

Geometry - AutoAlign

Position	R4.7 A13.6 F19.4 mm
Orientation	Sagittal
Phase Encoding Dir.	A >> P
AutoAlign	---
Initial Position	R4.7 A13.6 F19.4
R	4.7 mm
A	13.6 mm
F	19.4 mm
Initial Orientation	Sagittal
Initial Rotation	0.00 deg

Geometry - Navigator**Geometry - Saturation**

Special Saturation	None
--------------------	------

Geometry - Tim Planning Suite

Set-n-Go Protocol	Off
Table Position	0 mm
Table Position	H

System - Miscellaneous

Coil Selection	Auto Coil Select
Radial Sorting	Off
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combination	Adaptive Combine
Matrix Optimization	Off
Coil Focus	Flat

System - Adjustments

Adjustment Strategy	Standard
B0 Shim	Brain
B1 Shim	TrueForm
Adjustment Tolerance	Auto
Adjust with Body Coil	Off
Confirm Frequency	Never
Assume Silicone	Off
Adj Freq. protocol related	Off

System - Adjust Volume

Position	R4.7 A13.6 F19.4 mm
Orientation	Sagittal
Rotation	90.00 deg
F >> H	136 mm
A >> P	181 mm
R >> L	159 mm
Reset	Off

System - pTx

B1 Shim	TrueForm
Excitation	Non-sel.

System - pTx

pTx Pulse	1
-----------	---

System - Tx/Rx

Frequency 1H	297.208927 MHz
? Ref. Amplitude 1H	0.000 V
Reset	Off
Image Scaling	1.000
Gain	High

Physio - Signal

1st Signal/Mode	None
Trigger Delay	0 ms
TR	900.0 ms
Rd crusher scale	893
Concatenations	1

Physio - Cardiac

Fat-Water Contrast	Standard
Magn. Preparation	None
Dark Blood	Off
FOV Read	136 mm
FOV Phase	132.8 %
Phase Resolution	100 %

Physio - PACE

Resp. Control	Off
Concatenations	1

Inline - Subtraction

Subtract	Off
Measurements	1
StdDev	Off
Save Original Images	On

Inline - MIP

MIP Sag	Off
MIP Cor	Off
MIP Tra	Off
MIP Time	Off
Radial MIP	Off
Save Original Images	On
MPR Sag	Off
MPR Cor	Off
MPR Tra	Off

Inline - Composing**Sequence - Part 1**

Sequence Name	wip-spc-R
Dimension	3D
Excitation	Non-sel.
RF Pulse Type	Low SAR
Refoc. duration	1000
Gradient Mode	Fast

Sequence - Part 1

Flow Compensation	None
Reordering	Linear
Bandwidth	444 Hz/Px
Echo Spacing	5.06 ms
Concomitant field corr. R	Product
Turbo Factor	16
Echo Train Duration	86 ms
Rd crusher scale	893

Sequence - Part 2

Introduction	On
Rel. crusher mom. P	0 %
Rel crusher mom. S	0 %

Sequence - Special

Excit. flip angle	90 deg
Blood Suppression	Off
Trajectory Optimization	8

Sequence - Assistant

SAR Assistant	Off
Allowed Delay	0 s

\\BTM5 TERRA 7T\BRAIN BWH\Camilo_Peds\Epilepsy\t1_space_DRB_pTx_sag_p14_0.55TA: 8:44 min Coil Selection: Auto Voxel Size: 0.6×0.6×0.6 mm³ Acc:: 14.0 Rel. SNR: 1.00**Properties**

Start measurement without further preparation	On
Wait for User to Start	Off
Start measurements	Single Measurement
Prio Recon	Off
Auto Open Inline Display	Off
Auto Close Inline Display	Off
Load Images to MR View&GO	On
Auto Store Images	On
Disable auto transfer to PACS	Off
Load Images to Stamp Segments	Off
Load Images to Graphic Segments	Off
Graphic segment	Default
Inline Movie	Off

Routine

Slab Group	1
Slabs	1
Position	R4.7 A13.6 F19.4 mm
Orientation	Sagittal
Phase Encoding Dir.	A >> P
Slices per Slab	288
Phase Oversampling	0 %
Slice Oversampling	0.0 %
FOV Read	220 mm
FOV Phase	100.0 %
Slice Thickness	0.55 mm
TR	900.0 ms
TE	15.00 ms
Averages	1.0
Concatenations	1
AutoAlign	---

Contrast - Common

TR	900.0 ms
TE	15.00 ms
Magn. Preparation	None
Flip Angle Mode	T1 Var
Fat-Water Contrast	Standard
Dark Blood	Off
Blood Suppression	Off
Excitation offset frequency	0 Hz
Wrap-up Magn.	Restore
Refoc. offset frequency	0 Hz
Reconstruction	Magnitude

Contrast - Dynamic

goKART	Off
Measurements	1
Multiple Series	Each Measurement
Reordering	Linear

Contrast - ASL

MTC	Off
MTC frequency sign	Positive

Resolution - Common

FOV Read	220 mm
FOV Phase	100.0 %
Slice Thickness	0.55 mm
Base Resolution	384
Phase Resolution	100 %
Slice Resolution	100 %
Interpolation	Off

Resolution - Acceleration

Acceleration Mode	CS
Total Factor	14.0
Reference Scans	Integrated
Acceleration Factor PE	3
Reference Lines PE	24
Acceleration Factor 3D	2
Reference Lines 3D	24
Reordering Shift 3D	1
Deep Resolve	Off
Phase Partial Fourier	Off
Slice Partial Fourier	Off
Elliptical Scanning	On

Resolution - Filter

Raw Filter	On
Elliptical Filter	Off
Distortion Correction	3D
Normalize	B1 Filter
Image Filter	Off

Geometry - Common

Slab Group	1
Slabs	1
Position	R4.7 A13.6 F19.4 mm
Orientation	Sagittal
Phase Encoding Dir.	A >> P
Slices per Slab	288
Phase Oversampling	0 %
Slice Oversampling	0.0 %
FOV Read	220 mm
FOV Phase	100.0 %
Slice Thickness	0.55 mm
TR	900.0 ms
Concatenations	1

Geometry - AutoAlign

Slab Group	1
------------	---

Geometry - AutoAlign

Position	R4.7 A13.6 F19.4 mm
Orientation	Sagittal
Phase Encoding Dir.	A >> P
AutoAlign	---
Initial Position	R4.7 A13.6 F19.4
R	4.7 mm
A	13.6 mm
F	19.4 mm
Initial Orientation	Sagittal
Initial Rotation	0.00 deg

Geometry - Navigator**Geometry - Saturation**

Special Saturation	None
--------------------	------

Geometry - Tim Planning Suite

Set-n-Go Protocol	Off
Table Position	0 mm
Table Position	H

System - Miscellaneous

Coil Selection	Auto Coil Select
Radial Sorting	Off
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combination	Adaptive Combine
Matrix Optimization	Off
Coil Focus	Flat

System - Adjustments

Adjustment Strategy	Standard
B0 Shim	Brain
B1 Shim	TrueForm
Adjustment Tolerance	Auto
Adjust with Body Coil	Off
Confirm Frequency	Never
Assume Silicone	Off
Adj Freq. protocol related	Off

System - Adjust Volume

Position	R4.7 A13.6 F19.4 mm
Orientation	Sagittal
Rotation	0.00 deg
A >> P	220 mm
F >> H	220 mm
R >> L	159 mm
Reset	Off

System - pTx

B1 Shim	TrueForm
Excitation	Non-sel.

System - pTx

pTx Pulse	1
-----------	---

System - Tx/Rx

Frequency 1H	297.208927 MHz
? Ref. Amplitude 1H	0.000 V
Reset	Off
Image Scaling	1.000
Gain	High

Physio - Signal

1st Signal/Mode	None
Trigger Delay	0 ms
TR	900.0 ms
Rd crusher scale	893
Concatenations	1

Physio - Cardiac

Fat-Water Contrast	Standard
Magn. Preparation	None
Dark Blood	Off
FOV Read	220 mm
FOV Phase	100.0 %
Phase Resolution	100 %

Physio - PACE

Resp. Control	Off
Concatenations	1

Inline - Subtraction

Subtract	Off
Measurements	1
StdDev	Off
Save Original Images	On

Inline - MIP

MIP Sag	Off
MIP Cor	Off
MIP Tra	Off
MIP Time	Off
Radial MIP	Off
Save Original Images	On
MPR Sag	Off
MPR Cor	Off
MPR Tra	Off

Inline - Composing**Sequence - Part 1**

Sequence Name	wip-spc-R
Dimension	3D
Excitation	Non-sel.
RF Pulse Type	Low SAR
Refoc. duration	1000
Gradient Mode	Fast

Sequence - Part 1

Flow Compensation	None
Reordering	Linear
Bandwidth	449 Hz/Px
Echo Spacing	5.08 ms
Concomitant field corr. R	Product
Turbo Factor	14
Echo Train Duration	76 ms
Rd crusher scale	893

Sequence - Part 2

Introduction	On
Rel. crusher mom. P	0 %
Rel crusher mom. S	0 %

Sequence - Special

Excit. flip angle	90 deg
Blood Suppression	Off
Trajectory Optimization	8

Sequence - Assistant

SAR Assistant	Off
Allowed Delay	0 s

\\BTM5 TERRA 7T\\BRAIN BWH\\Camilo_Peds\\Epilepsy\\t1_space_DRB_NopTx_sag_p14_0.6

TA: 5:01 min Coil Selection: Auto Voxel Size: 0.6×0.6×0.6 mm³ Acc:: 14.0 Rel. SNR: 1.00**Properties**

Start measurement without further preparation	On
Wait for User to Start	Off
Start measurements	Single Measurement
Prio Recon	Off
Auto Open Inline Display	Off
Auto Close Inline Display	Off
Load Images to MR View&GO	On
Auto Store Images	On
Disable auto transfer to PACS	Off
Load Images to Stamp Segments	Off
Load Images to Graphic Segments	Off
Graphic segment	Default
Inline Movie	Off

Routine

Slab Group	1
Slabs	1
Position	R4.7 A13.6 F19.4 mm
Orientation	Sagittal
Phase Encoding Dir.	A >> P
Slices per Slab	288
Phase Oversampling	0 %
Slice Oversampling	0.0 %
FOV Read	220 mm
FOV Phase	100.0 %
Slice Thickness	0.55 mm
TR	700.0 ms
TE	15.00 ms
Averages	1.0
Concatenations	1
AutoAlign	---

Contrast - Common

TR	700.0 ms
TE	15.00 ms
Magn. Preparation	None
Flip Angle Mode	T1 Var
Fat-Water Contrast	Standard
Dark Blood	Off
Blood Suppression	Off
Excitation offset frequency	0 Hz
Wrap-up Magn.	Restore
Refoc. offset frequency	0 Hz
Reconstruction	Magnitude

Contrast - Dynamic

goKART	Off
Measurements	1
Multiple Series	Each Measurement
Reordering	Linear

Contrast - ASL

MTC	Off
MTC frequency sign	Positive

Resolution - Common

FOV Read	220 mm
FOV Phase	100.0 %
Slice Thickness	0.55 mm
Base Resolution	384
Phase Resolution	100 %
Slice Resolution	100 %
Interpolation	Off

Resolution - Acceleration

Acceleration Mode	CS
Total Factor	14.0
Reference Scans	Integrated
Acceleration Factor PE	3
Reference Lines PE	24
Acceleration Factor 3D	2
Reference Lines 3D	24
Reordering Shift 3D	1
Deep Resolve	Off
Phase Partial Fourier	Off
Slice Partial Fourier	Off
Elliptical Scanning	On

Resolution - Filter

Raw Filter	On
Elliptical Filter	Off
Distortion Correction	3D
Normalize	B1 Filter
Image Filter	Off

Geometry - Common

Slab Group	1
Slabs	1
Position	R4.7 A13.6 F19.4 mm
Orientation	Sagittal
Phase Encoding Dir.	A >> P
Slices per Slab	288
Phase Oversampling	0 %
Slice Oversampling	0.0 %
FOV Read	220 mm
FOV Phase	100.0 %
Slice Thickness	0.55 mm
TR	700.0 ms
Concatenations	1

Geometry - AutoAlign

Slab Group	1
------------	---

Geometry - AutoAlign

Position	R4.7 A13.6 F19.4 mm
Orientation	Sagittal
Phase Encoding Dir.	A >> P
AutoAlign	---
Initial Position	R4.7 A13.6 F19.4
R	4.7 mm
A	13.6 mm
F	19.4 mm
Initial Orientation	Sagittal
Initial Rotation	0.00 deg

Geometry - Navigator**Geometry - Saturation**

Special Saturation	None
--------------------	------

Geometry - Tim Planning Suite

Set-n-Go Protocol	Off
Table Position	0 mm
Table Position	H

System - Miscellaneous

Coil Selection	Auto Coil Select
Radial Sorting	Off
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combination	Adaptive Combine
Matrix Optimization	Off
Coil Focus	Flat

System - Adjustments

Adjustment Strategy	Standard
B0 Shim	Brain
B1 Shim	TrueForm
Adjustment Tolerance	Auto
Adjust with Body Coil	Off
Confirm Frequency	Never
Assume Silicone	Off
Adj Freq. protocol related	Off

System - Adjust Volume

Position	R4.7 A13.6 F19.4 mm
Orientation	Sagittal
Rotation	0.00 deg
A >> P	220 mm
F >> H	220 mm
R >> L	159 mm
Reset	Off

System - pTx

B1 Shim	TrueForm
Excitation	Non-sel.

System - Tx/Rx

Frequency 1H	297.208927 MHz
? Ref. Amplitude 1H	0.000 V
Reset	Off
Image Scaling	1.000
Gain	High

Physio - Signal

1st Signal/Mode	None
Trigger Delay	0 ms
TR	700.0 ms
Rd crusher scale	893
Concatenations	1

Physio - Cardiac

Fat-Water Contrast	Standard
Magn. Preparation	None
Dark Blood	Off
FOV Read	220 mm
FOV Phase	100.0 %
Phase Resolution	100 %

Physio - PACE

Resp. Control	Off
Concatenations	1

Inline - Subtraction

Subtract	Off
Measurements	1
StdDev	Off
Save Original Images	On

Inline - MIP

MIP Sag	Off
MIP Cor	Off
MIP Tra	Off
MIP Time	Off
Radial MIP	Off
Save Original Images	On
MPR Sag	Off
MPR Cor	Off
MPR Tra	Off

Inline - Composing**Sequence - Part 1**

Sequence Name	wip-spc-R
Dimension	3D
Excitation	Non-sel.
RF Pulse Type	Low SAR
Refoc. duration	1000
Gradient Mode	Fast
Flow Compensation	None
Reordering	Linear
Bandwidth	449 Hz/Px

Sequence - Part 1

Echo Spacing	4.88 ms
Concomitant field corr. R	Product
Turbo Factor	19
Echo Train Duration	98 ms
Rd crusher scale	893

Sequence - Part 2

Introduction	On
Rel. crusher mom. P	0 %
Rel crusher mom. S	0 %

Sequence - Special

Excit. flip angle	90 deg
Blood Suppression	Off
Trajectory Optimization	8

Sequence - Assistant

SAR Assistant	Off
Allowed Delay	0 s