

\\USER\NEURO_JAX\EPILEPSY NEW\Test\PTX_DL_SAG_MP2RAGE_EDGE_CS_DL_0.63MM_Q9WS105

TA: 6:51 min Coil Selection: Auto Voxel Size: 0.6x0.6x0.6 mm³ Acc:: 7.3 Rel. SNR: 1.00**Properties**

Start measurement without further preparation	On
Wait for User to Start	Off
Start measurements	Single Measurement
Prio Recon	Off
Auto Open Inline Display	On
Auto Close Inline Display	Off
Load Images to MR View&GO	On
Auto Store Images	On
Disable auto transfer to PACS	Off
Load Images to Stamp Segments	On
Load Images to Graphic Segments	On
Graphic segment	1st Segment
Inline Movie	Off

Routine

Slab Group	1
Slabs	1
Distance Factor	50 %
Position	L0.0 A12.2 F22.8 mm
Orientation	Sagittal
Phase Encoding Dir.	A >> P
Slices per Slab	256
Phase Oversampling	0 %
Slice Oversampling	12.5 %
FOV Read	240 mm
FOV Phase	93.8 %
Slice Thickness	0.63 mm
TR	6000.0 ms
TE	3.02 ms
Averages	1
Concatenations	1
AutoAlign	Head > Brain

Contrast - Common

TR	6000.0 ms
TE	3.02 ms
Magn. Preparation	Non-sel. IR
TI 1	875 ms
TI 2	2700 ms
Flip Angle 1	4 deg
Flip Angle 2	5 deg
Fat-Water Contrast	Standard
Dark Blood	Off
Contrasts	1
Reconstruction	Magnitude

Contrast - Dynamic

Dynamic Mode	Standard
Measurements	1
Multiple Series	Each Measurement

Contrast - Dynamic

Reordering	Linear
------------	--------

Resolution - Common

FOV Read	240 mm
FOV Phase	93.8 %
Slice Thickness	0.63 mm
Base Resolution	384
Phase Resolution	100 %
Slice Resolution	100 %
Interpolation	Off

Resolution - Acceleration

Acceleration Mode	CS
Total Factor	7.3
Reference Scans	GRE/Separate
Acceleration Factor PE	1
Reference Lines PE	32
Acceleration Factor 3D	1
Reordering Shift 3D	0
Phase Partial Fourier	Off
Slice Partial Fourier	Off
Asymmetric Echo	Off
Elliptical Scanning	Off

Resolution - Filter

Raw Filter	Off
Elliptical Filter	Off
Distortion Correction	3D
Normalize	B1 Filter
Image Filter	Off

Geometry - Common

Slab Group	1
Slabs	1
Distance Factor	50 %
Position	L0.0 A12.2 F22.8 mm
Orientation	Sagittal
Phase Encoding Dir.	A >> P
Slices per Slab	256
Phase Oversampling	0 %
Slice Oversampling	12.5 %
FOV Read	240 mm
FOV Phase	93.8 %
Slice Thickness	0.63 mm
TR	6000.0 ms
Multi-Slice Mode	Single Shot
Series	Interleaved
Concatenations	1

Geometry - AutoAlign

Slab Group	1
Position	L0.0 A12.2 F22.8 mm
Orientation	Sagittal
Phase Encoding Dir.	A >> P
AutoAlign	Head > Brain
Initial Position	L0.0 A12.2 F22.8
R	0.0 mm
A	12.2 mm
F	22.8 mm
Initial Orientation	Sagittal
Initial Rotation	0.00 deg

Geometry - Navigator**Geometry - Tim Planning Suite**

Set-n-Go Protocol	Off
Table Position	23 mm
Table Position	F

System - Miscellaneous

Coil Selection	Auto Coil Select
Radial Sorting	Off
MSMA	S - C - T
Sagittal	L >> R
Coronal	A >> P
Transversal	F >> H
Coil Combination	Adaptive Combine
Matrix Optimization	Off
Coil Focus	Flat

System - Adjustments

Adjustment Strategy	Standard
B0 Shim	Brain
B1 Shim	Patient-specific
Adjustment Tolerance	Auto
Adjust with Body Coil	Off
Confirm Frequency	Never
Assume Silicone	Off

System - Adjust Volume

Position	L0.0 A12.2 F22.8 mm
Orientation	Sagittal
Rotation	0.00 deg
A >> P	225 mm
F >> H	240 mm
R >> L	162 mm
Reset	Off

System - pTx

B1 Shim	Patient-specific
Excitation	Non-sel.
pTx Pulse	1

System - Tx/Rx

Frequency 1H	297.208927 MHz
? Ref. Amplitude 1H	0.000 V
Reset	Off
Image Scaling	1.000

Physio - Signal

1st Signal/Mode	None
TR	6000.0 ms
Concatenations	1

Physio - Cardiac

Fat-Water Contrast	Standard
Magn. Preparation	Non-sel. IR
TI 1	875 ms
TI 2	2700 ms
Dark Blood	Off
FOV Read	240 mm
FOV Phase	93.8 %
Phase Resolution	100 %

Physio - PACE

Resp. Control	Off
Concatenations	1

Inline - Subtraction

Subtract	Off
Measurements	1
StdDev	Off
Save Original Images	On

Inline - MIP

MIP Sag	Off
MIP Cor	Off
MIP Tra	Off
MIP Time	Off
Radial MIP	Off
Save Original Images	On
MPR Sag	Off
MPR Cor	Off
MPR Tra	Off

Inline - Composing**Inline - MapIt**

MapIt	T1 Map
Flip Angle 1	4 deg
Flip Angle 2	5 deg
Measurements	1
Contrasts	1
TE	3.02 ms
TR	6000.0 ms
Save Original Images	On

Sequence - Part 1

Sequence Name	WIP_tfl
Dimension	3D
Excitation	Non-sel.
RF Pulse Type	Fast
Gradient Mode	Fast
Flow Compensation	None
Reordering	Linear
Bandwidth	240 Hz/Px
Echo Spacing	7.12 ms
Asymmetric Echo	Off
Turbo Factor	205

Sequence - Part 2

Introduction	On
RF Spoiling	On
Incr. Gradient Spoiling	Off

Sequence - Special

Interpolation Mode	Post Volume
DL Recon	DL-Finetuned_CS-7
Run in Background	Off
Data Fidelity	3.0

Sequence - Assistant

SAR Assistant	Off
---------------	-----