

\\USER\NEURO_JAX\EPILEPSY NEW\Test\PTX_SAG_T2prep_FLAIR_SPC_Q9WS104-0.78mm

TA: 7:06 min Coil Selection: Auto Voxel Size: 0.8×0.8×0.8 mm³ Acc:: 6 Rel. SNR: 1.00**Properties**

Start measurement without further preparation	On
Wait for User to Start	Off
Start measurements	Single Measurement
Prio Recon	Off
Auto Open Inline Display	Off
Auto Close Inline Display	Off
Load Images to MR View&GO	On
Auto Store Images	On
Disable auto transfer to PACS	Off
Load Images to Stamp Segments	On
Load Images to Graphic Segments	On
Graphic segment	1st Segment
Inline Movie	Off

Routine

Slab Group	1
Slabs	1
Position	L0.0 P3.9 F29.8 mm
Orientation	S > T-0.3 > C-0.1
Phase Encoding Dir.	A >> P
Slices per Slab	208
Phase Oversampling	0 %
Slice Oversampling	7.7 %
FOV Read	225 mm
FOV Phase	100.0 %
Slice Thickness	0.78 mm
TR	8000.0 ms
TE	302.00 ms
Averages	1.0
Concatenations	1
AutoAlign	Head > Brain

Contrast - Common

TR	8000.0 ms
TE	302.00 ms
Magn. Preparation	Non-sel. T2 prep. IR
T1	2250 ms
T2 Prep. Duration	125 ms
Flip Angle Mode	T2 Var
Fat-Water Contrast	Fat Saturation
Fat Saturation	Strong
Dark Blood	Off
Blood Suppression	Off
Excitation offset frequency	0 Hz
Wrap-up Magn.	None
Refoc. offset frequency	0 Hz
Reconstruction	Magnitude

Contrast - Dynamic

goKART	Off
--------	-----

Contrast - Dynamic

Measurements	1
Multiple Series	Each Measurement
Reordering	Linear

Contrast - ASL

MTC	Off
MTC frequency sign	Positive

Resolution - Common

FOV Read	225 mm
FOV Phase	100.0 %
Slice Thickness	0.78 mm
Base Resolution	288
Phase Resolution	100 %
Slice Resolution	100 %
Interpolation	Off

Resolution - Acceleration

Acceleration Mode	CAIPIRINHA
Total Factor	6
Reference Scans	Integrated
Acceleration Factor PE	3
Reference Lines PE	32
Acceleration Factor 3D	2
Reference Lines 3D	32
Reordering Shift 3D	1
Deep Resolve	Boost
Phase Partial Fourier	Off
Slice Partial Fourier	Off
Elliptical Scanning	On

Resolution - Filter

Raw Filter	On
Elliptical Filter	Off
Distortion Correction	3D
Normalize	B1 Filter
Image Filter	Off

Geometry - Common

Slab Group	1
Slabs	1
Position	L0.0 P3.9 F29.8 mm
Orientation	S > T-0.3 > C-0.1
Phase Encoding Dir.	A >> P
Slices per Slab	208
Phase Oversampling	0 %
Slice Oversampling	7.7 %
FOV Read	225 mm
FOV Phase	100.0 %
Slice Thickness	0.78 mm

Geometry - Common

TR	8000.0 ms
Concatenations	1

Geometry - AutoAlign

Slab Group	1
Position	L0.0 P3.9 F29.8 mm
Orientation	S > T-0.3 > C-0.1
Phase Encoding Dir.	A >> P
AutoAlign	Head > Brain
Initial Position	L0.0 A2.0 F23.4
R	0.0 mm
A	2.0 mm
F	23.4 mm
Initial Orientation	Sagittal
Initial Rotation	0.00 deg

Geometry - Navigator**Geometry - Saturation**

Special Saturation	None
--------------------	------

Geometry - Tim Planning Suite

Set-n-Go Protocol	Off
Table Position	0 mm
Table Position	H

System - Miscellaneous

Coil Selection	Auto Coil Select
Radial Sorting	Off
MSMA	S - C - T
Sagittal	L >> R
Coronal	A >> P
Transversal	F >> H
Coil Combination	Adaptive Combine
Matrix Optimization	Off
Coil Focus	Flat

System - Adjustments

Adjustment Strategy	Standard
B0 Shim	Brain
B1 Shim	TrueForm
Adjustment Tolerance	Auto
Adjust with Body Coil	Off
Confirm Frequency	Always
Assume Silicone	Off
Adj Freq. protocol related	On

System - Adjust Volume

Position	L0.0 P3.9 F29.8 mm
Orientation	S > T-0.3 > C-0.1
Rotation	10.43 deg
A >> P	225 mm
F >> H	225 mm
R >> L	163 mm

System - Adjust Volume

Reset	Off
-------	-----

System - pTx

B1 Shim	TrueForm
Excitation	Non-sel.
pTx Pulse	1

System - Tx/Rx

Frequency 1H	297.208927 MHz
? Ref. Amplitude 1H	0.000 V
Reset	Off
Image Scaling	1.000
Gain	High

Physio - Signal

1st Signal/Mode	None
Trigger Delay	0 ms
TR	8000.0 ms
Rd crusher scale	893
Concatenations	1

Physio - Cardiac

Fat-Water Contrast	Fat Saturation
Magn. Preparation	Non-sel. T2 prep. IR
TI	2250 ms
T2 Prep. Duration	125 ms
Dark Blood	Off
FOV Read	225 mm
FOV Phase	100.0 %
Phase Resolution	100 %

Physio - PACE

Resp. Control	Off
Concatenations	1

Inline - Subtraction

Subtract	Off
Measurements	1
StdDev	Off
Save Original Images	On

Inline - MIP

MIP Sag	Off
MIP Cor	Off
MIP Tra	Off
MIP Time	Off
Radial MIP	Off
Save Original Images	On
MPR Sag	Off
MPR Cor	Off
MPR Tra	Off

Inline - Composing

Sequence - Part 1

Sequence Name	wip-spc-t2p+ir
Dimension	3D
Excitation	Non-sel.
RF Pulse Type	Low SAR
Refoc. duration	1000
Gradient Mode	Normal
Flow Compensation	None
Reordering	Linear
Bandwidth	755 Hz/Px
Echo Spacing	4.14 ms
Concomitant field corr. R	Product
Turbo Factor	180
Echo Train Duration	741 ms
Rd crusher scale	893

Sequence - Part 2

Introduction	On
Rel. crusher mom. P	0 %
Rel crusher mom. S	0 %

Sequence - Special

IR refoc pulse	IR12800H180.IR180_HS3
IR RF duration	12800 us
IR HS mu	5.00
IR HS BWTP	10.99
IR HS B1 peak	22.5 μ T
T2-Prep. excit/tip-up pulse	RECT
T2-Prep. refoc pulse	BIREF-1
T2-Prep. # refoc. RF	4
T2-Prep. excit. duration	800 us
T2-Prep. refoc. duration	11520 us
T2-Prep. spoil mom. 10^3 x	-71.20 usmT/m
T2-Prep. B1 scale	1.00
T2-Prep. crusher mom. rel.	0.00
T2-Prep. Refoc shift	0.00 %
Excit. flip angle	90 deg
Blood Suppression	Off
Trajectory Optimization	8

Sequence - Assistant

SAR Assistant	Off
Allowed Delay	0 s