

\\MARTINOS DEVELOPER\HUBER\resting_veins\20260818_JEN\t2_space_tra_p4_iso

TA: 6:20 min Coil Selection: Auto Voxel Size: 0.3×0.3×0.4 mm³ Acc:: 2 Rel. SNR: 1.00**Properties**

Start measurement without further preparation	On
Wait for User to Start	Off
Start measurements	Single Measurement
Prio Recon	Off
Auto Open Inline Display	Off
Auto Close Inline Display	Off
Load Images to MR View&GO	On
Auto Store Images	On
Disable auto transfer to PACS	Off
Load Images to Stamp Segments	Off
Load Images to Graphic Segments	Off
Graphic segment	Default
Inline Movie	Off

Routine

Slab Group	1
Slabs	1
Position	L0.0 A83.0 F32.7 mm
Orientation	C > T2.2
Phase Encoding Dir.	R >> L
Slices per Slab	112
Phase Oversampling	0 %
Slice Oversampling	14.3 %
FOV Read	120 mm
FOV Phase	100.0 %
Slice Thickness	0.40 mm
TR	3100.0 ms
TE	55.00 ms
Averages	1.5
Concatenations	1
AutoAlign	---

Contrast - Common

TR	3100.0 ms
TE	55.00 ms
MTC	Off
Magn. Preparation	None
Flip Angle Mode	Constant
Flip Angle	100 deg
Fat-Water Contrast	Standard
Dark Blood	Off
Blood Suppression	Off
Wrap-up Magn.	Restore
Reconstruction	Magnitude

Contrast - Dynamic

Dynamic Mode	Standard
Measurements	1
Multiple Series	Each Measurement
Reordering	Linear

Resolution - Common

FOV Read	120 mm
FOV Phase	100.0 %
Slice Thickness	0.40 mm
Base Resolution	224
Phase Resolution	100 %
Slice Resolution	100 %
Interpolation	On

Resolution - Acceleration

Acceleration Mode	CAIPIRINHA
Total Factor	2
Reference Scans	Integrated
Acceleration Factor PE	2
Reference Lines PE	24
Acceleration Factor 3D	1
Reordering Shift 3D	0
Deep Resolve	Off
Phase Partial Fourier	Off
Slice Partial Fourier	Off
Elliptical Scanning	On

Resolution - Filter

Raw Filter	On
Elliptical Filter	Off
Distortion Correction	2D
Normalize	B1 Filter
Image Filter	On

Geometry - Common

Slab Group	1
Slabs	1
Position	L0.0 A83.0 F32.7 mm
Orientation	C > T2.2
Phase Encoding Dir.	R >> L
Slices per Slab	112
Phase Oversampling	0 %
Slice Oversampling	14.3 %
FOV Read	120 mm
FOV Phase	100.0 %
Slice Thickness	0.40 mm
TR	3100.0 ms
Concatenations	1

Geometry - AutoAlign

Slab Group	1
Position	L0.0 A83.0 F32.7 mm
Orientation	C > T2.2
Phase Encoding Dir.	R >> L
AutoAlign	---
Initial Position	L0.0 A83.0 F32.7

Geometry - AutoAlign

L	0.0 mm
A	83.0 mm
F	32.7 mm
Initial Orientation	C > T
C > T	2.20
> S	0.00
Initial Rotation	0.00 deg

Geometry - Navigator**Geometry - Saturation**

Special Saturation	None
--------------------	------

Geometry - Tim Planning Suite

Set-n-Go Protocol	Off
Table Position	0 mm
Table Position	H

System - Miscellaneous

Coil Selection	Auto Coil Select
Radial Sorting	Off
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combination	Adaptive Combine
Matrix Optimization	Performance
Coil Focus	Flat

System - Adjustments

Adjustment Strategy	Standard
B0 Shim	Standard
B1 Shim	TrueForm
Adjustment Tolerance	Auto
Adjust with Body Coil	Off
Confirm Frequency	Never
Assume Silicone	Off

System - Adjust Volume

! Position	L0.0 A83.0 F32.7 mm
! Orientation	C > T2.2
! Rotation	0.00 deg
! R >> L	180 mm
! F >> H	180 mm
! A >> P	38 mm
Reset	Off

System - pTx

B1 Shim	TrueForm
Excitation	Slab-sel.

System - Tx/Rx

Frequency 1H	297.136418 MHz
! Ref. Amplitude 1H	300.000 V

System - Tx/Rx

Reset	Off
Image Scaling	1.000
Gain	High

Physio - Signal

1st Signal/Mode	None
Trigger Delay	0 ms
TR	3100.0 ms
Concatenations	1

Physio - Cardiac

Fat-Water Contrast	Standard
Magn. Preparation	None
Dark Blood	Off
FOV Read	120 mm
FOV Phase	100.0 %
Phase Resolution	100 %

Physio - PACE

Resp. Control	Off
Concatenations	1

Inline - Subtraction

Subtract	Off
Measurements	1
StdDev	Off
Save Original Images	On

Inline - MIP

MIP Sag	Off
MIP Cor	Off
MIP Tra	Off
MIP Time	Off
Radial MIP	Off
Save Original Images	On
MPR Sag	Off
MPR Cor	Off
MPR Tra	Off

Inline - Composing**Sequence - Part 1**

Sequence Name	spcR
Dimension	3D
Excitation	Slab-sel.
RF Pulse Type	Low SAR
Gradient Mode	Fast
Flow Compensation	None
Reordering	Linear
Bandwidth	378 Hz/Px
Echo Spacing	5.56 ms
Turbo Factor	159
Echo Train Duration	923 ms

Sequence - Part 2

Introduction	On
--------------	----

Sequence - Assistant

SAR Assistant	Off
Allowed Delay	30 s