

Table of contents

\\MARTINOS DEVELOPER				
	Renzo			
		Stefano_20260424		
			20260424	
				AAhead_scout cmrr_mbep2d_bold_2mm_AP cmrr_mbep2d_2mm_AP cmrr_mbep2d_2mm_PA mesovein_1_RO

\\MARTINOS DEVELOPER\Renzo\Stefano_20260424\20260424\AAhead_scout

TA: 14 sec Coil Selection: Auto Voxel Size: 1.6×1.6×1.6 mm³ Acc:: 3 Rel. SNR: 1.00**Properties**

Start measurement without further preparation	Off
Wait for User to Start	Off
Start measurements	Single Measurement
Prio Recon	Off
Auto Open Inline Display	Off
Auto Close Inline Display	Off
Load Images to MR View&GO	On
Auto Store Images	On
Disable auto transfer to PACS	Off
Load Images to Stamp Segments	Off
Load Images to Graphic Segments	On
Graphic segment	Default
Inline Movie	Off

Routine

Slab Group	1
Slabs	1
Distance Factor	20 %
Position	Isocenter
Orientation	Sagittal
Phase Encoding Dir.	A >> P
Slices per Slab	128
Phase Oversampling	0 %
Slice Oversampling	0.0 %
FOV Read	260 mm
FOV Phase	100.0 %
Slice Thickness	1.6 mm
TR	3.15 ms
TE	1.37 ms
Averages	1
Concatenations	1
AutoAlign	Head

Contrast - Common

TR	3.15 ms
TE	1.37 ms
Flip Angle	8 deg
Fat-Water Contrast	Standard
Contrasts	1
Reconstruction	Magnitude

Contrast - Dynamic

Dynamic Mode	Standard
Measurements	1
Time to Center	6.2 s

Resolution - Common

FOV Read	260 mm
FOV Phase	100.0 %

Resolution - Common

Slice Thickness	1.6 mm
Base Resolution	160
Phase Resolution	100 %
Slice Resolution	69 %
Trajectory	Cartesian

Resolution - Acceleration

Acceleration Mode	GRAPPA
Reference Scans	Integrated
Acceleration Factor PE	3
Reference Lines PE	24
Acceleration Factor 3D	1
Phase Partial Fourier	6/8
Slice Partial Fourier	6/8
Asymmetric Echo	Weak

Resolution - Filter

Raw Filter	Off
Elliptical Filter	Off
Distortion Correction	3D
Normalize	Prescan
Noise Masking	Off
Image Filter	Off

Geometry - Common

Slab Group	1
Slabs	1
Distance Factor	20 %
Position	Isocenter
Orientation	Sagittal
Phase Encoding Dir.	A >> P
Slices per Slab	128
Phase Oversampling	0 %
Slice Oversampling	0.0 %
FOV Read	260 mm
FOV Phase	100.0 %
Slice Thickness	1.6 mm
TR	3.15 ms
Multi-Slice Mode	Sequential
Series	Ascending
Concatenations	1

Geometry - AutoAlign

Slab Group	1
Position	Isocenter
Orientation	Sagittal
Phase Encoding Dir.	A >> P
AutoAlign	Head
Initial Position	Isocenter
L	0.0 mm

Geometry - AutoAlign

P	0.0 mm
H	0.0 mm
Initial Orientation	Transversal
Initial Rotation	0.00 deg

Geometry - Tim Planning Suite

Set-n-Go Protocol	Off
Table Position	0 mm
Table Position	H
Inline Composing	Off

System - Miscellaneous

Coil Selection	Auto Coil Select
Radial Sorting	Off
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combination	Adaptive Combine
Matrix Optimization	Off
Coil Focus	Flat

System - Adjustments

Adjustment Strategy	Standard
B0 Shim	Tune up
B1 Shim	TrueForm
CoilShim	Off
Adjustment Tolerance	Auto
Adjust with Body Coil	Off
Confirm Frequency	Never
Assume Silicone	Off

System - Adjust Volume

Position	Isocenter
Orientation	Transversal
Rotation	0.00 deg
A >> P	263 mm
R >> L	350 mm
F >> H	350 mm
Reset	Off

System - pTx

B1 Shim	TrueForm
Excitation	Non-sel.

System - Tx/Rx

Frequency 1H	123.221495 MHz
? Ref. Amplitude 1H	0.000 V
Reset	Off
Image Scaling	1.000

Physio - PACE

Resp. Control	Off
Concatenations	1

Inline - Subtraction

Subtract	Off
Measurements	1
StdDev	Off
Save Original Images	On

Inline - MIP

MIP Sag	Off
MIP Cor	Off
MIP Tra	Off
MIP Time	Off
Radial MIP	Off
Save Original Images	On
MPR Sag	Off
MPR Cor	Off
MPR Tra	Off

Inline - Composing

Inline Composing	Off
------------------	-----

Inline - MapIt

MapIt	None
Flip Angle	8 deg
Measurements	1
Contrasts	1
TE	1.37 ms
TR	3.15 ms
Save Original Images	On

Inline - Open Recon

Algorithm	None
-----------	------

Sequence - Part 1

Sequence Name	fl
Dimension	3D
Excitation	Non-sel.
RF Pulse Type	Fast
Gradient Mode	Normal
Bandwidth	540 Hz/Px
Asymmetric Echo	Weak

Sequence - Part 2

Introduction	On
RF Spoiling	On
Breast Application	Off
Phase Enc. Order	Automatic

Sequence - Assistant

SAR Assistant	Off
---------------	-----

\\MARTINOS DEVELOPER\Renzo\Stefano_20260424\20260424\cmrr_mbep2d_bold_2mm_AP

TA: 5:18 min Coil Selection: Manual Voxel Size: 2.0x2.0x2.0 mm³ Acc:: 3 Rel. SNR: 1.00**Properties**

Start measurement without further preparation	On
Wait for User to Start	Off
Start measurements	Single Measurement
Prio Recon	Off
Auto Open Inline Display	Off
Auto Close Inline Display	Off
Load Images to MR View&GO	On
Auto Store Images	On
Disable auto transfer to PACS	Off
Load Images to Stamp Segments	Off
Load Images to Graphic Segments	Off
Graphic segment	Default
Inline Movie	Off

Routine

Slice Group	1
Slices	1
Distance Factor	50 %
Position	L9.2 A25.6 F0.3 mm
Orientation	S > C-0.4
Phase Encoding Dir.	A >> P
Phase Oversampling	0 %
FOV Read	200 mm
FOV Phase	100.0 %
Slice Thickness	2.0 mm
TR	35.0 ms
TE	9.62 ms
Averages	1
Multi-band accel. factor	1
AutoAlign	---
Coil Elements	HC1-7;NC1,2

Contrast - Common

TR	35.0 ms
TE	9.62 ms
MTC	Off
Magn. Preparation	None
Flip Angle	13 deg
Fat-Water Contrast	Fat Saturation
Contrasts	1
Reconstruction	Magnitude

Contrast - Dynamic

Dynamic Mode	Standard
Measurements	9000
Delay in TR	0.00 ms

Resolution - Common

FOV Read	200 mm
----------	--------

Resolution - Common

FOV Phase	100.0 %
Slice Thickness	2.0 mm
Base Resolution	100
Phase Resolution	100 %
Interpolation	Off

Resolution - Acceleration

Acceleration Mode	GRAPPA
Reference scan mode	Segmented
Acceleration Factor PE	3
Reference Lines PE	48
Phase Partial Fourier	6/8

Resolution - Filter

Raw Filter	Off
Elliptical Filter	Off
Hamming	Off
Distortion Correction	Off
Static Field Correction	Off
Normalize	Off

Geometry - Common

Slice Group	1
Slices	1
Distance Factor	50 %
Position	L9.2 A25.6 F0.3 mm
Orientation	S > C-0.4
Phase Encoding Dir.	A >> P
Phase Oversampling	0 %
FOV Read	200 mm
FOV Phase	100.0 %
Slice Thickness	2.0 mm
TR	35.0 ms
Multi-Slice Mode	Interleaved
Series	Interleaved
Multi-band accel. factor	1

Geometry - AutoAlign

Slice Group	1
Position	L9.2 A25.6 F0.3 mm
Orientation	S > C-0.4
Phase Encoding Dir.	A >> P
AutoAlign	---
Initial Position	L9.2 A25.6 F0.3
L	9.2 mm
A	25.6 mm
F	0.3 mm
Initial Orientation	S > C
S > C	-0.40
> T	0.00

Geometry - AutoAlign

Initial Rotation	0.00 deg
------------------	----------

Geometry - Saturation

Special Saturation	None
--------------------	------

Geometry - Tim Planning Suite

Set-n-Go Protocol	Off
Table Position	0 mm
Table Position	H
Inline Composing	Off

System - Miscellaneous

Coil Selection	Manual
Radial Sorting	Off
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combination	Sum of Squares
Matrix Optimization	Off

System - Adjustments

Adjustment Strategy	Standard
B0 Shim	Standard
B1 Shim	TrueForm
CoilShim	Off
Adjustment Tolerance	Auto
Adjust with Body Coil	Off
Confirm Frequency	Never
Assume Silicone	Off

System - Adjust Volume

! Position	L7.8 A19.8 H29.2 mm
! Orientation	S > C0.7 > T-0.7
! Rotation	99.44 deg
! F >> H	113 mm
! A >> P	190 mm
! R >> L	41 mm
Reset	Off

System - pTx

B1 Shim	TrueForm
Excitation	Standard

System - Tx/Rx

Frequency 1H	123.221495 MHz
? Ref. Amplitude 1H	0.000 V
Reset	Off
Correction Factor	1.00
Image Scaling	1.000

Physio - Signal

1st Signal/Mode	None
TR	35.0 ms

Physio - Signal

Log Signals	Off
Multi-band accel. factor	1

BOLD

GLM Statistics	Off
Ignore Meas. at Start	0
Ignore After Transition	0
Model Transition States	On
Temp. Highpass Filter	On
Threshold	4.00
Paradigm Size	20
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active
Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Meas[16]	Active
Meas[17]	Active
Meas[18]	Active
Meas[19]	Active
Meas[20]	Active
Motion Correction	Off
Spatial Filter	Off
Measurements	9000
Delay in TR	0.00 ms

Inline - Open Recon

Algorithm	None
-----------	------

Sequence - Part 1

Sequence Name	epfid
Dimension	2D
Excitation	Standard
Gradient Mode	Performance*
Flow Compensation	None
Bandwidth	2084 Hz/Px
Echo Spacing	0.58 ms
Free Echo Spacing	Off
EPI Factor	100

Sequence - Part 2

Introduction	Off
RF Spoiling	Off
Ghost Reduction	Off

Sequence - Special

Excite pulse duration	2560 us
Min. prep scans	0
Delay before PC scans	0 us
EPI noise scans	2
SENSE1 coil combine	Off
Invert RO/PE polarity	Off
Limit grad. capability	Off
PF omits higher k-space	Off
Disable B1 control loop	Off
Disable freq. update	On
Suppress 16-bit DICOM	On
Force equal slice timing	Off
FFT scale factor	1.00
Fat saturation FA	110.00 deg
Fat sat. offset	0.00 Hz
Sinc exc. pulse BWTP	5.20
Physio recording	File
Triggering scheme	Standard

Sequence - Assistant

SAR Assistant	Off
---------------	-----

\\MARTINOS DEVELOPER\Renzo\Stefano_20260424\20260424\cmrr_mbep2d_2mm_AP

TA: 27 sec Coil Selection: Auto Voxel Size: 2.0x2.0x2.0 mm³ Acc:: 3 Rel. SNR: 1.00**Properties**

Start measurement without further preparation	On
Wait for User to Start	Off
Start measurements	Single Measurement
Prio Recon	Off
Auto Open Inline Display	Off
Auto Close Inline Display	Off
Load Images to MR View&GO	On
Auto Store Images	On
Disable auto transfer to PACS	Off
Load Images to Stamp Segments	Off
Load Images to Graphic Segments	Off
Graphic segment	Default
Inline Movie	Off

Resolution - Common

FOV Read	200 mm
FOV Phase	100.0 %
Slice Thickness	2.0 mm
Base Resolution	100
Phase Resolution	100 %
Interpolation	Off

Resolution - Acceleration

Acceleration Mode	GRAPPA
Reference scan mode	Segmented
Acceleration Factor PE	3
Reference Lines PE	48
Phase Partial Fourier	6/8

Routine

Slice Group	1
Slices	1
Distance Factor	50 %
Position	L9.2 A25.6 F0.3 mm
Orientation	S > C-0.4
Phase Encoding Dir.	A >> P
Phase Oversampling	0 %
FOV Read	200 mm
FOV Phase	100.0 %
Slice Thickness	2.0 mm
TR	1000.0 ms
TE	14.40 ms
Averages	1
Multi-band accel. factor	1
AutoAlign	---

Resolution - Filter

Raw Filter	Off
Elliptical Filter	Off
Hamming	Off
Distortion Correction	2D
Static Field Correction	Off
Normalize	Off

Geometry - Common

Slice Group	1
Slices	1
Distance Factor	50 %
Position	L9.2 A25.6 F0.3 mm
Orientation	S > C-0.4
Phase Encoding Dir.	A >> P
Phase Oversampling	0 %
FOV Read	200 mm
FOV Phase	100.0 %
Slice Thickness	2.0 mm
TR	1000.0 ms
Multi-Slice Mode	Interleaved
Series	Interleaved
Multi-band accel. factor	1

Contrast - Common

TR	1000.0 ms
TE	14.40 ms
MTC	Off
Magn. Preparation	None
Flip Angle	90 deg
Refocus flip angle	180 deg
Fat-Water Contrast	Fat Saturation
Grad. rev. fat suppr.	Enabled
Contrasts	1
Reconstruction	Magnitude

Contrast - Dynamic

Dynamic Mode	Standard
Measurements	20
Delay in TR	0.00 ms

Geometry - AutoAlign

Slice Group	1
Position	L9.2 A25.6 F0.3 mm
Orientation	S > C-0.4
Phase Encoding Dir.	A >> P
AutoAlign	---
Initial Position	L9.2 A25.6 F0.3
L	9.2 mm
A	25.6 mm
F	0.3 mm
Initial Orientation	S > C
S > C	-0.40

Geometry - AutoAlign

> T	0.00
Initial Rotation	0.00 deg

Geometry - Saturation

Special Saturation	None
--------------------	------

Geometry - Tim Planning Suite

Set-n-Go Protocol	Off
Table Position	0 mm
Table Position	H
Inline Composing	Off

System - Miscellaneous

Coil Selection	Auto Coil Select
Radial Sorting	Off
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combination	Sum of Squares
Matrix Optimization	Off
Coil Focus	Flat

System - Adjustments

Adjustment Strategy	Standard
B0 Shim	Standard
B1 Shim	TrueForm
CoilShim	Off
Adjustment Tolerance	Auto
Adjust with Body Coil	Off
Confirm Frequency	Never
Assume Silicone	Off

System - Adjust Volume

I Position	L7.8 A19.8 H29.2 mm
I Orientation	S > C0.7 > T-0.7
I Rotation	99.44 deg
I F >> H	113 mm
I A >> P	190 mm
I R >> L	41 mm
Reset	Off

System - pTx

B1 Shim	TrueForm
Excitation	Standard

System - Tx/Rx

Frequency 1H	123.221495 MHz
? Ref. Amplitude 1H	0.000 V
Reset	Off
Image Scaling	1.000

Physio - Signal

1st Signal/Mode	None
-----------------	------

Physio - Signal

TR	1000.0 ms
Log Signals	Off
Multi-band accel. factor	1

BOLD

GLM Statistics	On
Ignore Meas. at Start	0
Ignore After Transition	0
Model Transition States	On
Temp. Highpass Filter	On
Threshold	4.00
Paradigm Size	20
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active
Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Meas[16]	Active
Meas[17]	Active
Meas[18]	Active
Meas[19]	Active
Meas[20]	Active
Motion Correction	On
Spatial Filter	Off
Measurements	20
Delay in TR	0.00 ms

Inline - Open Recon

Algorithm	None
-----------	------

Sequence - Part 1

Sequence Name	epse
Dimension	2D
Excitation	Standard
RF Pulse Type	Normal
Gradient Mode	Performance
Bandwidth	2084 Hz/Px
Echo Spacing	0.58 ms
Free Echo Spacing	Off
EPI Factor	100

Sequence - Part 2

Introduction	Off
Ghost Reduction	Off

Sequence - Special

Min. prep scans	0
Delay before PC scans	0 us
EPI noise scans	0
SENSE1 coil combine	Off
Invert RO/PE polarity	Off
Limit grad. capability	Off
PF omits higher k-space	Off
Force GPA balance	Off
Disable B1 control loop	Off
Disable freq. update	Off
Suppress 16-bit DICOM	Off
Force equal slice timing	Off
FFT scale factor	1.00
Fat saturation FA	110.00 deg
Fat sat. offset	0.00 Hz
Physio recording	Off
Triggering scheme	Standard

Sequence - Assistant

SAR Assistant	Off
---------------	-----

\\MARTINOS DEVELOPER\Renzo\Stefano_20260424\20260424\cmrr_mbep2d_2mm_PA

TA: 27 sec Coil Selection: Auto Voxel Size: 2.0x2.0x2.0 mm³ Acc:: 3 Rel. SNR: 1.00**Properties**

Start measurement without further preparation	On
Wait for User to Start	Off
Start measurements	Single Measurement
Prio Recon	Off
Auto Open Inline Display	Off
Auto Close Inline Display	Off
Load Images to MR View&GO	On
Auto Store Images	On
Disable auto transfer to PACS	Off
Load Images to Stamp Segments	Off
Load Images to Graphic Segments	Off
Graphic segment	Default
Inline Movie	Off

Resolution - Common

FOV Read	200 mm
FOV Phase	100.0 %
Slice Thickness	2.0 mm
Base Resolution	100
Phase Resolution	100 %
Interpolation	Off

Resolution - Acceleration

Acceleration Mode	GRAPPA
Reference scan mode	Segmented
Acceleration Factor PE	3
Reference Lines PE	48
Phase Partial Fourier	6/8

Routine

Slice Group	1
Slices	1
Distance Factor	50 %
Position	L9.2 A25.6 F0.3 mm
Orientation	S > C-0.4
Phase Encoding Dir.	A >> P
Phase Oversampling	0 %
FOV Read	200 mm
FOV Phase	100.0 %
Slice Thickness	2.0 mm
TR	1000.0 ms
TE	14.40 ms
Averages	1
Multi-band accel. factor	1
AutoAlign	---

Resolution - Filter

Raw Filter	Off
Elliptical Filter	Off
Hamming	Off
Distortion Correction	2D
Static Field Correction	Off
Normalize	Off

Geometry - Common

Slice Group	1
Slices	1
Distance Factor	50 %
Position	L9.2 A25.6 F0.3 mm
Orientation	S > C-0.4
Phase Encoding Dir.	A >> P
Phase Oversampling	0 %
FOV Read	200 mm
FOV Phase	100.0 %
Slice Thickness	2.0 mm
TR	1000.0 ms
Multi-Slice Mode	Interleaved
Series	Interleaved
Multi-band accel. factor	1

Contrast - Common

TR	1000.0 ms
TE	14.40 ms
MTC	Off
Magn. Preparation	None
Flip Angle	90 deg
Refocus flip angle	180 deg
Fat-Water Contrast	Fat Saturation
Grad. rev. fat suppr.	Enabled
Contrasts	1
Reconstruction	Magnitude

Contrast - Dynamic

Dynamic Mode	Standard
Measurements	20
Delay in TR	0.00 ms

Geometry - AutoAlign

Slice Group	1
Position	L9.2 A25.6 F0.3 mm
Orientation	S > C-0.4
Phase Encoding Dir.	A >> P
AutoAlign	---
Initial Position	L9.2 A25.6 F0.3
L	9.2 mm
A	25.6 mm
F	0.3 mm
Initial Orientation	S > C
S > C	-0.40

Geometry - AutoAlign

> T	0.00
Initial Rotation	0.00 deg

Geometry - Saturation

Special Saturation	None
--------------------	------

Geometry - Tim Planning Suite

Set-n-Go Protocol	Off
Table Position	0 mm
Table Position	H
Inline Composing	Off

System - Miscellaneous

Coil Selection	Auto Coil Select
Radial Sorting	Off
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combination	Sum of Squares
Matrix Optimization	Off
Coil Focus	Flat

System - Adjustments

Adjustment Strategy	Standard
B0 Shim	Standard
B1 Shim	TrueForm
CoilShim	Off
Adjustment Tolerance	Auto
Adjust with Body Coil	Off
Confirm Frequency	Never
Assume Silicone	Off

System - Adjust Volume

I Position	L7.8 A19.8 H29.2 mm
I Orientation	S > C0.7 > T-0.7
I Rotation	99.44 deg
I F >> H	113 mm
I A >> P	190 mm
I R >> L	41 mm
Reset	Off

System - pTx

B1 Shim	TrueForm
Excitation	Standard

System - Tx/Rx

Frequency 1H	123.221495 MHz
? Ref. Amplitude 1H	0.000 V
Reset	Off
Image Scaling	1.000

Physio - Signal

1st Signal/Mode	None
-----------------	------

Physio - Signal

TR	1000.0 ms
Log Signals	Off
Multi-band accel. factor	1

BOLD

GLM Statistics	On
Ignore Meas. at Start	0
Ignore After Transition	0
Model Transition States	On
Temp. Highpass Filter	On
Threshold	4.00
Paradigm Size	20
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active
Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Meas[16]	Active
Meas[17]	Active
Meas[18]	Active
Meas[19]	Active
Meas[20]	Active
Motion Correction	On
Spatial Filter	Off
Measurements	20
Delay in TR	0.00 ms

Inline - Open Recon

Algorithm	None
-----------	------

Sequence - Part 1

Sequence Name	epse
Dimension	2D
Excitation	Standard
RF Pulse Type	Normal
Gradient Mode	Performance
Bandwidth	2084 Hz/Px
Echo Spacing	0.58 ms
Free Echo Spacing	Off
EPI Factor	100

Sequence - Part 2

Introduction	Off
Ghost Reduction	Off

Sequence - Special

Min. prep scans	0
Delay before PC scans	0 us
EPI noise scans	0
SENSE1 coil combine	Off
Invert RO/PE polarity	On
Limit grad. capability	Off
PF omits higher k-space	Off
Force GPA balance	Off
Disable B1 control loop	Off
Disable freq. update	Off
Suppress 16-bit DICOM	Off
Force equal slice timing	Off
FFT scale factor	1.00
Fat saturation FA	110.00 deg
Fat sat. offset	0.00 Hz
Physio recording	Off
Triggering scheme	Standard

Sequence - Assistant

SAR Assistant	Off
---------------	-----

\\MARTINOS DEVELOPER\Renzo\Stefano_20260424\20260424\mesovein_1_RO

TA: 6:19 min Coil Selection: Auto Voxel Size: 0.4x0.4x0.4 mm³ Acc:: 6 Rel. SNR: 1.00**Properties**

Start measurement without further preparation	Off
Wait for User to Start	Off
Start measurements	Single Measurement
Prio Recon	Off
Auto Open Inline Display	Off
Auto Close Inline Display	Off
Load Images to MR View&GO	On
Auto Store Images	On
Disable auto transfer to PACS	Off
Load Images to Stamp Segments	Off
Load Images to Graphic Segments	Off
Graphic segment	Default
Inline Movie	Off

Routine

Slab Group	1
Slabs	1
Position	L0.1 A14.3 H21.8 mm
Orientation	T > C-17.1
Phase Encoding Dir.	A >> P
Slices per Slab	416
Phase Oversampling	0.0 %
Slice Oversampling	3.8 %
FOV Read	175 mm
FOV Phase	120.0 %
Slice Thickness	0.35 mm
TR	149.6 ms
Vol. TR	367.56720 s
TE 1	22.60 ms
TE 2	67.80 ms
TE 3	113.00 ms
Averages	1
Multi-echo Shots	1
AutoAlign	---

Contrast - Common

TR	149.6 ms
Vol. TR	367.56720 s
TE 1	22.60 ms
TE 2	67.80 ms
TE 3	113.00 ms
Multi-echo spacing	45.20 ms
MTC	Off
Flip Angle	23 deg
Fat-Water Contrast	Fast Fat Saturation
Contrasts	3
Reconstruction	Magnitude

Contrast - Dynamic

Dynamic Mode	Standard
--------------	----------

Contrast - Dynamic

Measurements	1
Reordering	Linear

Resolution - Common

FOV Read	175 mm
FOV Phase	120.0 %
Slice Thickness	0.35 mm
Base Resolution	500
Phase Resolution	100 %
Slice Resolution	100 %
Interpolation	Off

Resolution - Acceleration

Acceleration Mode	CAIPIRINHA
Reference Scans	GRE/Separate
Acceleration Factor PE	3
Reference Lines PE	36
Acceleration Factor 3D	2
Reference Lines 3D	36
Reordering Shift 3D	1
Phase Partial Fourier	7/8
Slice Partial Fourier	7/8

Resolution - Filter

Raw Filter	Off
Elliptical Filter	Off
Distortion Correction	Off
Normalize	Off
Image Filter	Off

Geometry - Common

Slab Group	1
Slabs	1
Position	L0.1 A14.3 H21.8 mm
Orientation	T > C-17.1
Phase Encoding Dir.	A >> P
Slices per Slab	416
Phase Oversampling	0.0 %
Slice Oversampling	3.8 %
FOV Read	175 mm
FOV Phase	120.0 %
Slice Thickness	0.35 mm
TR	149.6 ms
Vol. TR	367.56720 s
Multi-echo Shots	1

Geometry - AutoAlign

Slab Group	1
Position	L0.1 A14.3 H21.8 mm
Orientation	T > C-17.1

Geometry - AutoAlign

Phase Encoding Dir.	A >> P
AutoAlign	---
Initial Position	L0.1 A14.3 H21.8
L	0.1 mm
A	14.3 mm
H	21.8 mm
Initial Orientation	T > C
T > C	-17.10
> S	0.00
Initial Rotation	0.00 deg

Geometry - Saturation

Saturation Mode	Quick
-----------------	-------

Geometry - Tim Planning Suite

Set-n-Go Protocol	Off
Table Position	0 mm
Table Position	H
Inline Composing	Off

System - Miscellaneous

Coil Selection	Auto Coil Select
Radial Sorting	Off
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combination	Adaptive Combine
Matrix Optimization	Off
Coil Focus	Flat

System - Adjustments

Adjustment Strategy	Standard
B0 Shim	Brain
B1 Shim	TrueForm
CoilShim	Off
Adjustment Tolerance	Auto
Adjust with Body Coil	Off
Confirm Frequency	Never
Assume Silicone	Off

System - Adjust Volume

! Position	L2.6 A16.3 H26.2 mm
! Orientation	T > C-11.9 > S0.2
! Rotation	89.83 deg
! R >> L	137 mm
! A >> P	157 mm
! F >> H	87 mm
Reset	Off

System - pTx

B1 Shim	TrueForm
Excitation	Slab-sel.

System - Tx/Rx

Frequency 1H	123.221495 MHz
? Ref. Amplitude 1H	0.000 V
Reset	Off
Image Scaling	1.000

Inline - Open Recon

Algorithm	None
-----------	------

Sequence - Part 1

Sequence Name	ep v3.1
Dimension	3D
Excitation	Slab-sel.
RF Pulse Type	Normal
Gradient Mode	Performance
Reordering	Linear
Bandwidth	370 Hz/Px
Echo Spacing	3.08 ms
Segmentation	13
EPI Factor	14

Sequence - Part 2

Introduction	On
RF Spoiling	On

Sequence - Special

RF duration	1000 us
RF time x BW	18
PAT ref. FA	7 deg
Fat sat. FA	110 deg
T1 (Ernst FA)	1200 ms
Invert PE	On
Min. TE w/ PF	On
Ramp Sampling	On
Trigger per shot	Off
Noise image	Off
Relax spoilers	On
Round up Vol. TR	Off
MT flip angle	500 deg
MT off-res.	2000 Hz
MT RF duration	10240 us
Custom Water Exc.	-none-
Phase Correction	per Series
Saturation RF	per Shot
EPI rise time factor	1.50
G. spoil dephasing[1]	0.0 pi
G. spoil dephasing[2]	2.0 pi
G. spoil dephasing[3]	0.0 pi
Modify Ice Config	On
G-factor map	Off
GRAPPA Regularization	200 /10^6
Slab Scale	0 %
RF spoil scheme	Conventional
Read polarity	Regular RO

Sequence - Assistant

SAR Assistant	Off
---------------	-----