

\\MARTINOS DEVELOPER\Renzo\mesovein\20260407\mesovein

TA: 7:18 min Coil Selection: Auto Voxel Size: 0.4x0.4x0.4 mm³ Acc:: 6 Rel. SNR: 1.00**Properties**

Start measurement without further preparation	Off
Wait for User to Start	Off
Start measurements	Single Measurement
Prio Recon	Off
Auto Open Inline Display	Off
Auto Close Inline Display	Off
Load Images to MR View&GO	On
Auto Store Images	On
Disable auto transfer to PACS	Off
Load Images to Stamp Segments	Off
Load Images to Graphic Segments	Off
Graphic segment	Default
Inline Movie	Off

Routine

Slab Group	1
Slabs	1
Position	L0.0 A16.1 H1.3 mm
Orientation	T > C-18.0
Phase Encoding Dir.	A >> P
Slices per Slab	384
Phase Oversampling	0.0 %
Slice Oversampling	4.2 %
FOV Read	175 mm
FOV Phase	120.0 %
Slice Thickness	0.35 mm
TR	152.9 ms
Vol. TR	428.12 s
TE 1	13.90 ms
TE 2	41.90 ms
TE 3	69.90 ms
TE 4	97.90 ms
TE 5	125.90 ms
Averages	1
Multi-echo Shots	1
AutoAlign	---

Contrast - Common

TR	152.9 ms
Vol. TR	428.12 s
TE 1	13.90 ms
TE 2	41.90 ms
TE 3	69.90 ms
TE 4	97.90 ms
TE 5	125.90 ms
Multi-echo spacing	28.00 ms
MTC	Off
Flip Angle	23 deg
Fat-Water Contrast	Fast Fat Saturation
Contrasts	5

Contrast - Common

Reconstruction	Magn./Phase
----------------	-------------

Contrast - Dynamic

Dynamic Mode	Standard
Measurements	1
Reordering	Linear

Resolution - Common

FOV Read	175 mm
FOV Phase	120.0 %
Slice Thickness	0.35 mm
Base Resolution	500
Phase Resolution	100 %
Slice Resolution	100 %
Interpolation	Off

Resolution - Acceleration

Acceleration Mode	CAIPIRINHA
Reference Scans	GRE/Separate
Acceleration Factor PE	3
Reference Lines PE	36
Acceleration Factor 3D	2
Reference Lines 3D	36
Reordering Shift 3D	1
Phase Partial Fourier	7/8
Slice Partial Fourier	7/8

Resolution - Filter

Raw Filter	Off
Elliptical Filter	Off
Distortion Correction	Off
Normalize	Off
Image Filter	Off

Geometry - Common

Slab Group	1
Slabs	1
Position	L0.0 A16.1 H1.3 mm
Orientation	T > C-18.0
Phase Encoding Dir.	A >> P
Slices per Slab	384
Phase Oversampling	0.0 %
Slice Oversampling	4.2 %
FOV Read	175 mm
FOV Phase	120.0 %
Slice Thickness	0.35 mm
TR	152.9 ms
Vol. TR	428.12 s
Multi-echo Shots	1

Geometry - AutoAlign

Slab Group	1
Position	L0.0 A16.1 H1.3 mm
Orientation	T > C-18.0
Phase Encoding Dir.	A >> P
AutoAlign	---
Initial Position	L0.0 A16.1 H1.3
L	0.0 mm
A	16.1 mm
H	1.3 mm
Initial Orientation	T > C
T > C	-18.00
> S	0.00
Initial Rotation	0.00 deg

Geometry - Saturation

Saturation Mode	Quick
-----------------	-------

Geometry - Tim Planning Suite

Set-n-Go Protocol	Off
Table Position	0 mm
Table Position	H
Inline Composing	Off

System - Miscellaneous

Coil Selection	Auto Coil Select
Radial Sorting	Off
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combination	Adaptive Combine
Matrix Optimization	Off
Coil Focus	Flat

System - Adjustments

Adjustment Strategy	Standard
B0 Shim	Brain
B1 Shim	TrueForm
CoilShim	Off
Adjustment Tolerance	Auto
Adjust with Body Coil	Off
Confirm Frequency	Never
Assume Silicone	Off

System - Adjust Volume

Position	L0.0 A16.1 H1.3 mm
Orientation	T > C-18.0
Rotation	90.00 deg
R >> L	175 mm
A >> P	210 mm
F >> H	135 mm
Reset	Off

System - pTx

B1 Shim	TrueForm
Excitation	Slab-sel.

System - Tx/Rx

Frequency 1H	123.221487 MHz
? Ref. Amplitude 1H	0.000 V
Reset	Off
Image Scaling	1.000

Inline - Open Recon

Algorithm	None
-----------	------

Sequence - Part 1

Sequence Name	ep v3.1
Dimension	3D
Excitation	Slab-sel.
RF Pulse Type	Normal
Gradient Mode	Performance
Reordering	Linear
Bandwidth	526 Hz/Px
Echo Spacing	2.32 ms
Segmentation	16
EPI Factor	11

Sequence - Part 2

Introduction	On
RF Spoiling	On

Sequence - Special

RF duration	1000 us
RF time x BW	18
PAT ref. FA	7 deg
Fat sat. FA	110 deg
T1 (Ernst FA)	1200 ms
Invert PE	On
Min. TE w/ PF	On
Ramp Sampling	On
Trigger per shot	Off
Noise image	Off
Relax spoilers	On
Round up Vol. TR	Off
MT flip angle	500 deg
MT off-res.	2000 Hz
MT RF duration	10240 us
Custom Water Exc.	-none-
Phase Correction	per Series
Saturation RF	per Shot
EPI rise time factor	1.50
G. spoil dephasing[1]	0.0 pi
G. spoil dephasing[2]	2.0 pi
G. spoil dephasing[3]	0.0 pi
Modify Ice Config	On
G-factor map	Off
GRAPPA Regularization	200 /10^6

Sequence - Special

Slab Scale	0 %
RF spoil scheme	Conventional
Read polarity	Regular RO

Sequence - Assistant

SAR Assistant	Off
---------------	-----