

\\USER\jon\Saskia\THS\ldzne_ep3d_vasoNIH

TA: 35 sec Coil Selection: Auto Voxel Size: 0.8x0.8x0.9 mm³ Acc.: 16 Rel. SNR: 1.00**Properties**

Start measurement without further preparation	Off
Wait for User to Start	Off
Start measurements	Single Measurement
Prio Recon	Off
Auto Open Inline Display	Off
Auto Close Inline Display	Off
Load Images to MR View&GO	On
Auto Store Images	On
Disable auto transfer to PACS	Off
Load Images to Stamp Segments	Off
Load Images to Graphic Segments	Off
	Default
Inline Movie	Off
Before Measurement	
After Measurement	

Routine

Slab Group	1
Slabs	1
Position	L0.0 P11.4 F12.6 mm
Orientation	T > C-14.7
Phase Encoding Dir.	A >> P
AutoAlign	---
Slices per Slab	208
Phase Oversampling	0.0 %
Slice Oversampling	0.0 %
FOV Read	200 mm
FOV Phase	96.6 %
Slice Thickness	0.85 mm
TR	39.5 ms
Vol. TR	2054 ms
TE 1	17.50 ms
Averages	1
Multi-echo Shots	1

Contrast - Common

TR	39.5 ms
Vol. TR	2054 ms
TE 1	17.50 ms
Multi-echo spacing	32.50 ms
MTC	Off
Flip Angle	18 deg
Fat-Water Contrast	Fat Saturation
Contrasts	1
Reconstruction	Magnitude

Contrast - Dynamic

Dynamic Mode	Standard
Measurements	10
Reordering	Linear

Resolution - Common

FOV Read	200 mm
FOV Phase	96.6 %
Slice Thickness	0.85 mm
Base Resolution	236
Phase Resolution	100 %
Slice Resolution	100 %
Interpolation	Off

Resolution - Acceleration

Acceleration Mode	CAIPIRINHA
Reference Scans	GRE/Separate
Acceleration Factor PE	4
Reference Lines PE	80
Acceleration Factor 3D	4
Reference Lines 3D	36
Reordering Shift 3D	2
Phase Partial Fourier	7/8
Slice Partial Fourier	Off
Asymmetric Echo	Off

Resolution - Filter

Raw Filter	Off
Elliptical Filter	Off
Distortion Correction	3D
Normalize	Off
Image Filter	Off

Geometry - Common

Slab Group	1
Slabs	1
Position	L0.0 P11.4 F12.6 mm
Orientation	T > C-14.7
Phase Encoding Dir.	A >> P
Slices per Slab	208
Phase Oversampling	0.0 %
Slice Oversampling	0.0 %
FOV Read	200 mm
FOV Phase	96.6 %
Slice Thickness	0.85 mm
TR	39.5 ms
Vol. TR	2054 ms
Multi-echo Shots	1

Geometry - AutoAlign

Slab Group	1
Position	L0.0 P11.4 F12.6 mm
Orientation	T > C-14.7
Phase Encoding Dir.	A >> P
AutoAlign	---
Initial Position	Isocenter

Geometry - AutoAlign

L	0.0 mm
P	0.0 mm
H	0.0 mm
Initial Orientation	Transversal
Initial Rotation	0.00 deg

Geometry - Saturation

Saturation Mode	Standard
-----------------	----------

Geometry - Tim Planning Suite

Set-n-Go Protocol	Off
Label	
Table Position	13 mm
Table Position	F

System - Miscellaneous

Coil Selection	Auto Coil Select
Radial Sorting	Off
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combination	Adaptive Combine
Matrix Optimization	Off
Coil Focus	Flat

System - Adjustments

Adjustment Strategy	Standard
B0 Shim	Brain
B1 Shim	TrueForm
Adjustment Tolerance	Auto
Adjust with Body Coil	Off
Confirm Frequency	Never
Assume Silicone	Off

System - Adjust Volume

Position	L0.0 P11.4 F12.6 mm
Orientation	T > C-14.7
Rotation	0.00 deg
A >> P	194 mm
R >> L	200 mm
F >> H	177 mm
Reset	Off

System - pTx

B1 Shim	TrueForm
Excitation	Slab-sel.

System - Tx/Rx

Frequency 1H	297.168938 MHz
? Ref. Amplitude 1H	0.000 V
Reset	Off
Image Scaling	1.000

Inline - Open Recon

Algorithm	None
-----------	------

Sequence - Part 1

Sequence Name	epi v3.1.10-dev
Dimension	3D
Excitation	Slab-sel.
RF Pulse Type	Normal
Gradient Mode	Normal
Reordering	Linear
Bandwidth	1766 Hz/Px
Echo Spacing	0.64 ms
Asymmetric Echo	Off
Turbo Factor	52
Segmentation	1
EPI Factor	50

Sequence - Part 2

Introduction	On
RF Spoiling	On

Sequence - Special

Custom Water Exc.	-none-
Phase Correction	per Series
Saturation RF	per Shot
Spoiling scheme	Conventional
Read polarity	Dual-polarity
Local Scrambling	-none-

Sequence - Assistant

SAR Assistant	Off
---------------	-----