

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\UserProtocols\Yuhui\VAPER_largeFA\A VAPER_0.8i0.9sl_Pat3_pf78_ref300v_26sl_IPS_MT

TA: 30:45 PAT: 3 Voxel size: 0.8x0.8x0.9 mm Rel. SNR: 1.00 UNKNOWN:

Properties

| | |
|---|--------|
| Prio Recon | Off |
| Before measurement | |
| After measurement | |
| Load to viewer | On |
| Inline movie | Off |
| Auto store images | On |
| Load to stamp segments | Off |
| Load images to graphic segments | Off |
| Auto open inline display | Off |
| Start measurement without further preparation | On |
| Wait for user to start | Off |
| Start measurements | single |

Routine

| | |
|--------------------|------------------|
| Slab group 1 | |
| Slabs | 1 |
| Dist. factor | 50 % |
| Position | R2.4 P23.3 H34.7 |
| Orientation | C > T-7.7 |
| Phase enc. dir. | F >> H |
| Rotation | 90.00 deg |
| Phase oversampling | 9 % |
| Slice oversampling | 7.7 % |
| Slices per slab | 26 |
| FoV read | 130.0 mm |
| FoV phase | 100.0 % |
| Slice thickness | 0.90 mm |
| TR | 3024.00 ms |
| TE | 28 ms |
| Averages | 1 |
| Concatenations | 1 |
| Filter | None |
| Coil elements | A32 |

Contrast

| | |
|----------------------|---------------|
| Perfusion mode | Picore Q2TIPS |
| TI2 | 1100 ms |
| TI1 | 50 ms |
| TI1s | 50 ms |
| Flip angle | 27.5 deg |
| Fat suppr. | None |
| Averaging mode | Long term |
| Reconstruction | Magnitude |
| Measurements | 610 |
| Delay in TR | 0 ms |
| Multiple series | Off |
| Perfusion mode | PICORE Q2T |
| Inversion time 1 | 50 ms |
| Saturation stop time | 50 ms |
| Inversion time 2 | 1100.0 ms |
| Flow limit | 100 cm/s |

Resolution

| | |
|-----------------------|--------|
| Base resolution | 162 |
| Phase resolution | 100 % |
| Slice resolution | 100 % |
| Phase partial Fourier | 7/8 |
| Slice partial Fourier | Off |
| Interpolation | Off |
| PAT mode | GRAPPA |

| | |
|---------------------|----------|
| Accel. factor PE | 3 |
| Ref. lines PE | 48 |
| Accel. factor 3D | 1 |
| Ref. lines 3D | 24 |
| Reference scan mode | Separate |
| Prescan Normalize | Off |
| Raw filter | Off |
| Elliptical filter | Off |
| Hamming | Off |

Geometry

| | |
|------------------|-------------|
| Multi-slice mode | Interleaved |
| Series | Ascending |
| Special sat. | Parallel A |
| Gap | 25.0 mm |
| Thickness | 100 mm |
| Table position | H |
| Table position | 0 mm |
| Inline Composing | Off |

System

| | |
|--------------------------|-------------------|
| V32 | Off |
| A32 | On |
| Positioning mode | FIX |
| MSMA | S - C - T |
| Sagittal | R >> L |
| Coronal | A >> P |
| Transversal | F >> H |
| Save uncombined | Off |
| Coil Combine Mode | Sum of Squares |
| AutoAlign | --- |
| Auto Coil Select | Default |
| Shim mode | Standard |
| Adjust with body coil | Off |
| Confirm freq. adjustment | Off |
| Assume Silicone | Off |
| ! Ref. amplitude 1H | 290.000 V |
| Adjustment Tolerance | Auto |
| Adjust volume | |
| ! Position | R2.8 P25.7 H42.4 |
| ! Orientation | S > T-0.4 > C-0.1 |
| ! Rotation | 0.00 deg |
| ! F >> H | 82 mm |
| ! A >> P | 61 mm |
| ! R >> L | 127 mm |

Physio

| | |
|-----------------|------|
| 1st Signal/Mode | None |
|-----------------|------|

BOLD

| | |
|-------------------|-----|
| Motion correction | Off |
| Spatial filter | Off |

Sequence

| | |
|-------------------|------------|
| Introduction | On |
| Dimension | 3D |
| Reordering | Linear |
| Contrasts | 1 |
| Bandwidth | 1102 Hz/Px |
| Free echo spacing | Off |
| Echo spacing | 1.03 ms |
| EPI factor | 162 |
| RF pulse type | Normal |

| | | |
|---------------------------|--------------|-------------------------------------|
| Gradient mode | Normal | |
| Excitation | Slab-sel. | |
| RF spoiling | On | |
| <hr/> | | |
| Read Diff Amp | 21.5 mT/m | set diff amp to 0 in all directions |
| Phase Diff Amp | 0.0 mT/m | |
| Slice Diff Amp | 21.5 mT/m | |
| Dante puls # in 1st par | 132 | largest FA allowed by SAR limit |
| Dante puls # in 2nd par | 32 | |
| Pulses FA in DANTE | 10.6 degree | |
| TAU in DANTE | 1100 us | |
| Vari readFA | 0 | |
| Blank bef/aft DANTE-RF | 50 us | |
| Grad # bef DANTE | 0 | |
| DANTE-RF dur | 120 us | |
| use Ernst angle | Off | |
| Maxwell Correction | Off | |
| log physio files | Off | |
| FFT scale | 2.00 | |
| dummy prepscan time | 3 s | |
| z shim | 0.00 mT/m*ms | |
| RF duration | 2470 us | |
| RF BWTP | 25.0 | |
| EFFECTIVE TR | 97 ms | |
| PatPartitions | 28 | |
| EPI phase correction | local | |
| PAT refscan mode | Flash | |
| FlashRef BaseRes | 162 | |
| FlashRef BW | 100 Hz/px | |
| FlashRef TE | 10000 us | |
| FlashRef FA | 5 deg | |
| use CAIPI | Off | |