

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Renzo\FAIR\neg_BOLD\MAFI			
TA: 2:17	Voxel size: 6.0x6.0x6.0 mm	Rel. SNR: 1.00	USER: Renzo\MAFI

Properties		Geometry	
Prio Recon	Off	Multi-slice mode	Sequential
Before measurement		Series	Ascending
After measurement		Special sat.	None
Load to viewer	On	Table position	H
Inline movie	Off	Table position	0 mm
Auto store images	On	Inline Composing	Off
Load to stamp segments	Off		
Load images to graphic segments	Off		
Auto open inline display	Off	System	
Start measurement without further preparation	On	V24	Off
Wait for user to start	Off	A24	On
Start measurements	single		
Routine		Positioning mode	REF
Slab group 1		MSMA	S - C - T
Slabs	1	Sagittal	R >> L
Dist. factor	0 %	Coronal	A >> P
Position	L0.0 A28.5 F1.2	Transversal	F >> H
Orientation	Sagittal	Save uncombined	Off
Phase enc. dir.	A >> P	Coil Combine Mode	Adaptive Combine
Rotation	0.00 deg	AutoAlign	---
Phase oversampling	0 %	Auto Coil Select	Default
Slice oversampling	0.0 %		
Slices per slab	28	Shim mode	Standard
FOV read	192 mm	Adjust with body coil	Off
FOV phase	100.0 %	Confirm freq. adjustment	Off
Slice thickness	6.00 mm	Assume Silicone	Off
TR	150 ms	! Ref. amplitude 1H	270.000 V
TE 1	1.080 ms	Adjustment Tolerance	Auto
TE 2	1.080 ms	Adjust volume	
TE 3	2.09 ms	Position	L0.0 A28.5 F1.2
TE 4	3.100 ms	Orientation	Sagittal
Averages	1	Rotation	0.00 deg
Concatenations	1	F >> H	192 mm
Filter	None	A >> P	192 mm
Coil elements	A24	R >> L	168 mm
Contrast		Physio	
MTC	Off	1st Signal/Mode	None
Flip angle	90 deg	Inline	
Fat suppr.	None		
Water suppr.	None	Subtract	Off
		Std-Dev-Sag	Off
Averaging mode	Short term	Std-Dev-Cor	Off
Reconstruction	Magn./Phase	Std-Dev-Tra	Off
Measurements	1	Std-Dev-Time	Off
Multiple series	Off	MIP-Sag	Off
		MIP-Cor	Off
		MIP-Tra	Off
		MIP-Time	Off
		Save original images	On
Resolution		Sequence	
Base resolution	32	Introduction	Off
Phase resolution	100 %	Dimension	3D
Slice resolution	100 %	Contrasts	4
Phase partial Fourier	Off	Bandwidth	1560 Hz/Px
Interpolation	Off		
Image Filter	Off	Gradient mode	Fast
Distortion Corr.	Off	RF spoiling	On
Prescan Normalize	Off	Online ICE	Off
Normalize	Off	RF pulse type	square
B1 filter	Off	Pulse duration	500 us
Raw filter	Off	Spoil me!	On
Elliptical filter	Off	TR2/TR1	5

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

N dummy TRs	20
Sample T1	1800 ms
Diffusion damping d= bD	0.6000
Diffusion coefficient D	2.2000 $\mu\text{m}^2/\text{ms}$
RF spoil phase increment	129.3 deg
Number of pulse shapes	1
<hr/>	
TX/RX Nucleus	1H
TX/RX delta frequency	0 Hz
TX Nucleus	None
TX delta frequency	0 Hz