

<b>Table of contents</b>
--------------------------

\\USER
--------

<table border="1"><tr><td>FMRIF</td></tr></table>	FMRIF
FMRIF	

<table border="1"><tr><td><table border="1"><tr><td>[XT-ID:93-M-0170] Renzo</td></tr></table></td></tr></table>	<table border="1"><tr><td>[XT-ID:93-M-0170] Renzo</td></tr></table>	[XT-ID:93-M-0170] Renzo
<table border="1"><tr><td>[XT-ID:93-M-0170] Renzo</td></tr></table>	[XT-ID:93-M-0170] Renzo	
[XT-ID:93-M-0170] Renzo		

<table border="1"><tr><td><table border="1"><tr><td>coil_comp_sharing</td></tr></table></td></tr></table>	<table border="1"><tr><td>coil_comp_sharing</td></tr></table>	coil_comp_sharing
<table border="1"><tr><td>coil_comp_sharing</td></tr></table>	coil_comp_sharing	
coil_comp_sharing		

<table border="1"><tr><td><table border="1"><tr><td><a href="#">coil_utils</a></td></tr></table></td></tr></table>	<table border="1"><tr><td><a href="#">coil_utils</a></td></tr></table>	<a href="#">coil_utils</a>
<table border="1"><tr><td><a href="#">coil_utils</a></td></tr></table>	<a href="#">coil_utils</a>	
<a href="#">coil_utils</a>		

\\USER\FMRIF\[XT-ID:93-M-0170]\Renzo\coil\_comp\_sharing\coil\_utils

TA: 0:13 PM: REF Voxel size: 1.2×1.2×5.0 mmRel. SNR: 1.00 : fl

**Properties**

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	Off
Start measurements	Single measurement

**Routine**

Slice group	1
Slices	1
Dist. factor	20 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
Phase oversampling	0 %
FoV read	300 mm
FoV phase	100.0 %
Slice thickness	5.0 mm
TR	10.0 ms
TE	4.8 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	AC

**Contrast - Common**

TR	10.0 ms
TE	4.8 ms
MTC	Off
Flip angle	15 deg
Fat suppr.	None
Water suppr.	None

**Contrast - Dynamic**

Averages	1
Averaging mode	Short term
Reconstruction	Magnitude
Measurements	1
Multiple series	Off

**Resolution - Common**

FoV read	300 mm
FoV phase	100.0 %
Slice thickness	5.0 mm
Base resolution	256
Phase resolution	100 %
Phase partial Fourier	Off
Interpolation	Off

**Resolution - Filter Image**

Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	Off

**Resolution - Filter Image**

Normalize	Off
B1 filter	Off

**Resolution - Filter Rawdata**

Raw filter	Off
Elliptical filter	Off

**Geometry - Common**

Slice group	1
Slices	1
Dist. factor	20 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
FoV read	300 mm
FoV phase	100.0 %
Slice thickness	5.0 mm
TR	10.0 ms
Multi-slice mode	Sequential
Series	Ascending
Concatenations	1

**Geometry - AutoAlign**

Slice group	1
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
Initial Position	Isocenter
L	0.0 mm
P	0.0 mm
H	0.0 mm
Initial Rotation	0.00 deg
Initial Orientation	Transversal

**Geometry - Saturation**

Fat suppr.	None
Water suppr.	None
Special sat.	None

**Geometry - Tim Planning Suite**

Set-n-Go Protocol	Off
Table position	H
Table position	0 mm
Inline Composing	Off

**System - Miscellaneous**

Positioning mode	REF
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
Save uncombined	Off
Matrix Optimization	Off
AutoAlign	---
Coil Select Mode	Default

**System - Adjustments**

B0 Shim mode	Tune up
B1 Shim mode	TrueForm
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

**System - Adjust Volume**

Position	Isocenter
Orientation	Transversal
Rotation	0.00 deg
A >> P	263 mm
R >> L	350 mm
F >> H	350 mm
Reset	Off

**System - pTx Volumes**

B1 Shim mode	TrueForm
--------------	----------

**System - Tx/Rx**

Frequency 1H	297.184941 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

**Physio - Signal1**

1st Signal/Mode	None
TR	10.0 ms
Concatenations	1

**Sequence - Part 1**

Introduction	Off
Dimension	2D
Contrasts	1
Multi-slice mode	Sequential
Bandwidth	390 Hz/Px

**Sequence - Part 2**

Gradient mode	Fast
RF spoiling	On

**Sequence - Special**

ICE program	IceProgramCoilUtils
Prep. scans duration	10000 ms
Optimal SNR	On
GFactor	On
Rx coil diode switching	On

**Sequence - Assistant**

Mode	Off
------	-----