

Table of contents

\USER

Test

Joseph

20210920_VASOremote

[rslh_Gz3_S6_0p39_glob_PC](#)
[202111101_executed](#)
[rslh_Gz3_S6_C1_0p39_PCglob_run4](#)
[rslh_Gz3_S6_C1_0p39_PCshot_run2](#)

\\USER\Test\Joseph\20210920_VASOremote\rslh_Gz3_S6_0p39_glob_PC

TA: 14:21 PM: FIX Voxel size: 0.4×0.4×0.4 mmPAT: 3 Rel. SNR: 1.00 : d238999c

Properties

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	Off
Start measurements	Single measurement

Routine

Slab group	1
Slabs	1
Position	L0.0 P51.0 H14.6 mm
Orientation	C > T-36.1
Phase enc. dir.	R >> L
AutoAlign	---
Slab Scale	-10 %
Slices per slab	18
FoV read	182 mm
FoV phase	100.0 %
Slice thickness	0.39 mm
TR 1	67.7 ms
TR 2	6870 ms
TE 1	23.50 ms
Averages	1
Multi-echo Shots	1
Filter	Distortion Corr.(3D)
Coil elements	R96

Contrast - Common

TR 1	67.7 ms
TR 2	6870 ms
TE 1	23.50 ms
Multi-echo spacing	63.37 ms
Magn. preparation	Non-sel. HSN IR
TI 1	1269.3 ms
TI 2	2487.9 ms
Flip angle	40 deg
Fat suppr.	Fat sat.
Magn. Prep. Shots	2

Contrast - Dynamic

Averages	1
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	124
Pause after meas.	0.0 s

Resolution - Common

FoV read	182 mm
FoV phase	100.0 %
Slice thickness	0.39 mm
Base resolution	462
Phase resolution	100 %
Slice resolution	100 %
Phase partial Fourier	6/8

Resolution - Common

Slice partial Fourier	Off
Interpolation	Off

Resolution - iPAT

PAT mode	CAIPIRINHA
Acc. factor PE	1
Ref. lines PE	75
Acc. factor 3D	3
Ref. lines 3D	18
CAIPI 3D Shift	0
Reference Scan Mode	GRE/separate
CAIPIRINHA mode	Free

Resolution - Filter Image

Image Filter	Off
Distortion Corr.	On
Mode	3D
Unfiltered images	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off

Resolution - Filter Rawdata

Raw filter	Off
Elliptical filter	Off

Geometry - Common

Slab group	1
Slabs	1
Position	L0.0 P51.0 H14.6 mm
Orientation	C > T-36.1
Phase enc. dir.	R >> L
Slab Scale	-10 %
Slices per slab	18
FoV read	182 mm
FoV phase	100.0 %
Slice thickness	0.39 mm
TR 1	67.7 ms
TR 2	6870 ms
Multi-slice mode	Interleaved
Series	Ascending
Multi-echo Shots	1

Geometry - AutoAlign

Slab group	1
Position	L0.0 P51.0 H14.6 mm
Orientation	C > T-36.1
Phase enc. dir.	R >> L
AutoAlign	---
Initial Position	L0.0 P51.0 H14.6
L	0.0 mm
P	51.0 mm
H	14.6 mm
Initial Rotation	0.00 deg
Initial Orientation	C > T
C > T	-36.1
> S	0.0

Geometry - Saturation

Saturation mode	Standard
-----------------	----------

Geometry - Saturation

Fat suppr.	Fat sat.
------------	----------

Geometry - Tim Planning Suite

Set-n-Go Protocol	Off
Table position	H
Table position	0 mm
Inline Composing	Off

System - Miscellaneous

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
Save uncombined	Off
Matrix Optimization	Off
AutoAlign	---
Coil Select Mode	Default

System - Adjustments

B0 Shim mode	Tune up
B1 Shim mode	TrueForm
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

System - Adjust Volume

! Position	L0.0 P38.8 H22.5 mm
! Orientation	C > T-36.7
! Rotation	-90.00 deg
! F >> H	180 mm
! R >> L	350 mm
! A >> P	40 mm
Reset	Off

System - pTx Volumes

B1 Shim mode	TrueForm
Excitation	Slab-sel.

System - Tx/Rx

Frequency 1H	297.210751 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	0.200
Reset	Off
! Ref. amplitude 1H	300.000 V

Sequence - Part 1

Introduction	On
Dimension	3D
Reordering	Linear
Asymmetric echo	Off
Contrasts	1
Multi-slice mode	Interleaved
Echo spacing	1.08 ms
Bandwidth	984 Hz/Px

Sequence - Part 2

EPI factor	58
------------	----

Sequence - Part 2

Segmentation	6
RF pulse type	Normal
Gradient mode	Fast
Excitation	Slab-sel.
RF spoiling	On
Turbo factor	18

Sequence - Special

PATRef FA	3 deg
RF duration	2000 us
RF BWT product	15
Ernst T1	1200 ms
PATRef prep. shots	10
Volume dummy shots	0
Dummy Measurements	0
Invert PE	Off
Min. TE if PF	On
Echo Time Shift	On
Ramp Sampling	On
NORDIC	Off
Water Exc.	-none-
External PC	per Series
Saturation RF	per Shot
EPI rise time factor	1.00
Mosaic DICOMs	On
Modify Ice Config	On
GRAPPA Regularization	50000 10 ⁻⁶
HSN RF power scale	2.50
Inversion Delay	650 ms
Relaxation Delay	0 ms
Var. FA /MAGEC	4

Sequence - Assistant

Mode	Off
------	-----

\\USER\Test\Joseph\20210920_VASOremote\202111101_executed

TA: 0:27 PM: FIX Voxel size: 0.5×0.5×3.0 mmPAT: Off Rel. SNR: 1.00 : qfl

Properties

Prio recon	On
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	On
Load images to graphic segments	On
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	Off
Start measurements	Single measurement

Routine

Slice group	1
Slices	5
Dist. factor	200 %
Position	Isocenter
Orientation	Sagittal
Phase enc. dir.	A >> P
Slice group	2
Slices	5
Dist. factor	200 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
Slice group	3
Slices	5
Dist. factor	200 %
Position	Isocenter
Orientation	Coronal
Phase enc. dir.	R >> L
AutoAlign	---
Phase oversampling	0 %
FoV read	280 mm
FoV phase	100.0 %
Slice thickness	3.0 mm
TR	10.0 ms
TE	3.00 ms
Averages	1
Concatenations	15
Filter	Elliptical filter
Coil elements	R96

Contrast - Common

TR	10.0 ms
TE	3.00 ms
TD	0 ms
MTC	Off
Magn. preparation	None
Flip angle	10 deg
Fat suppr.	None
Water suppr.	None
SWI	Off

Contrast - Dynamic

Averages	1
Averaging mode	Short term
Reconstruction	Magnitude
Measurements	1

Contrast - Dynamic

Multiple series	Each measurement
-----------------	------------------

Resolution - Common

FoV read	280 mm
FoV phase	100.0 %
Slice thickness	3.0 mm
Base resolution	256
Phase resolution	90 %
Phase partial Fourier	6/8
Interpolation	On

Resolution - iPAT

PAT mode	None
----------	------

Resolution - Filter Image

Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off

Resolution - Filter Rawdata

Raw filter	Off
Elliptical filter	On

Geometry - Common

Slice group	1
Slices	5
Dist. factor	200 %
Position	Isocenter
Orientation	Sagittal
Phase enc. dir.	A >> P
Slice group	2
Slices	5
Dist. factor	200 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
Slice group	3
Slices	5
Dist. factor	200 %
Position	Isocenter
Orientation	Coronal
Phase enc. dir.	R >> L
FoV read	280 mm
FoV phase	100.0 %
Slice thickness	3.0 mm
TR	10.0 ms
Multi-slice mode	Sequential
Series	Interleaved
Concatenations	15

Geometry - AutoAlign

Slice group	1
Position	Isocenter
Orientation	Sagittal
Phase enc. dir.	A >> P
Slice group	2
Position	Isocenter
Orientation	Transversal

Geometry - AutoAlign

Phase enc. dir.	A >> P
Slice group	3
Position	Isocenter
Orientation	Coronal
Phase enc. dir.	R >> L
AutoAlign	---
Initial Position	Isocenter
L	0.0 mm
P	0.0 mm
H	0.0 mm
Initial Rotation	0.00 deg
Initial Orientation	Sagittal

Geometry - Saturation

Saturation mode	Standard
Fat suppr.	None
Water suppr.	None
Special sat.	None

Geometry - Tim Planning Suite

Set-n-Go Protocol	Off
Table position	H
Table position	0 mm
Inline Composing	Off

Geometry - Tim CT

Tim CT mode	Off
Slices	5
Slice thickness	3.0 mm
Dist. factor	200 %
FoV read	280 mm
FoV phase	100.0 %
Segments	1

System - Miscellaneous

Positioning mode	FIX
Table position	F
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
Save uncombined	Off
Matrix Optimization	Off
AutoAlign	---
Coil Select Mode	Off - AutoCoilSelect

System - Adjustments

B0 Shim mode	Tune up
B1 Shim mode	TrueForm
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

System - Adjust Volume

Position	Isocenter
Orientation	Transversal
Rotation	0.00 deg
A >> P	263 mm
R >> L	350 mm
F >> H	350 mm
Reset	Off

System - pTx Volumes

B1 Shim mode	TrueForm
Excitation	Slice-sel.

System - Tx/Rx

Frequency 1H	297.210751 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

Physio - Signal1

1st Signal/Mode	None
TR	10.0 ms
Concatenations	15
Segments	1

Physio - Cardiac

Tagging	None
Magn. preparation	None
Fat suppr.	None
Dark blood	Off
FoV read	280 mm
FoV phase	100.0 %
Phase resolution	90 %

Physio - PACE

Resp. control	Off
Concatenations	15

Inline - Common

Subtract	Off
Measurements	1
StdDev	Off
Liver registration	Off
Save original images	On

Inline - MIP

MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On

Inline - Soft Tissue

Wash - In	Off
Wash - Out	Off
TTP	Off
PEI	Off
MIP - time	Off
Measurements	1

Inline - Composing

Inline Composing	Off
Distortion Corr.	Off

Inline - MapIt

Save original images	On
MapIt	None
Flip angle	10 deg
Measurements	1
Contrasts	1
TR	10.0 ms

Inline - MapIt

TE	3.00 ms
----	---------

Sequence - Part 1

Introduction	On
Dimension	2D
Phase stabilisation	Off
Asymmetric echo	Allowed
Contrasts	1
Flow comp.	No
Multi-slice mode	Sequential
Bandwidth	320 Hz/Px

Sequence - Part 2

Segments	1
Acoustic noise reduction	Active
RF pulse type	Normal
Gradient mode	Normal
Excitation	Slice-sel.
RF spoiling	On

Sequence - Assistant

Mode	Off
------	-----

\\USER\Test\Joseph\20210920_VASOremote\rslh_Gz3_S6_C1_0p39_PCglob_run4

TA: 14:21 PM: FIX Voxel size: 0.4×0.4×0.4 mmPAT: 3 Rel. SNR: 1.00 : d238999c

Properties

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	Off
Start measurements	Single measurement

Routine

Slab group	1
Slabs	1
Position	L0.0 P53.3 H21.0 mm
Orientation	C > T-37.4
Phase enc. dir.	R >> L
AutoAlign	---
Slab Scale	-10 %
Slices per slab	18
FoV read	182 mm
FoV phase	100.0 %
Slice thickness	0.39 mm
TR 1	67.7 ms
TR 2	6870 ms
TE 1	23.50 ms
Averages	1
Multi-echo Shots	1
Filter	Distortion Corr.(3D)
Coil elements	R96

Contrast - Common

TR 1	67.7 ms
TR 2	6870 ms
TE 1	23.50 ms
Multi-echo spacing	63.37 ms
Magn. preparation	Non-sel. HSN IR
TI 1	1269.3 ms
TI 2	2487.9 ms
Flip angle	40 deg
Fat suppr.	Fat sat.
Magn. Prep. Shots	2

Contrast - Dynamic

Averages	1
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	124
Pause after meas.	0.0 s

Resolution - Common

FoV read	182 mm
FoV phase	100.0 %
Slice thickness	0.39 mm
Base resolution	462
Phase resolution	100 %
Slice resolution	100 %
Phase partial Fourier	6/8

Resolution - Common

Slice partial Fourier	Off
Interpolation	Off

Resolution - iPAT

PAT mode	CAIPIRINHA
Acc. factor PE	1
Ref. lines PE	75
Acc. factor 3D	3
Ref. lines 3D	18
CAIPI 3D Shift	1
Reference Scan Mode	GRE/separate
CAIPIRINHA mode	Free

Resolution - Filter Image

Image Filter	Off
Distortion Corr.	On
Mode	3D
Unfiltered images	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off

Resolution - Filter Rawdata

Raw filter	Off
Elliptical filter	Off

Geometry - Common

Slab group	1
Slabs	1
Position	L0.0 P53.3 H21.0 mm
Orientation	C > T-37.4
Phase enc. dir.	R >> L
Slab Scale	-10 %
Slices per slab	18
FoV read	182 mm
FoV phase	100.0 %
Slice thickness	0.39 mm
TR 1	67.7 ms
TR 2	6870 ms
Multi-slice mode	Interleaved
Series	Ascending
Multi-echo Shots	1

Geometry - AutoAlign

Slab group	1
Position	L0.0 P53.3 H21.0 mm
Orientation	C > T-37.4
Phase enc. dir.	R >> L
AutoAlign	---
Initial Position	L0.0 P53.3 H21.0
L	0.0 mm
P	53.3 mm
H	21.0 mm
Initial Rotation	0.00 deg
Initial Orientation	C > T
C > T	-37.4
> S	0.0

Geometry - Saturation

Saturation mode	Standard
-----------------	----------

Geometry - Saturation

Fat suppr.	Fat sat.
------------	----------

Geometry - Tim Planning Suite

Set-n-Go Protocol	Off
Table position	H
Table position	0 mm
Inline Composing	Off

System - Miscellaneous

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
Save uncombined	Off
Matrix Optimization	Off
AutoAlign	---
Coil Select Mode	Default

System - Adjustments

B0 Shim mode	Tune up
B1 Shim mode	TrueForm
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

System - Adjust Volume

! Position	L0.0 P48.1 H17.8 mm
! Orientation	C > T-36.7
! Rotation	-90.00 deg
! F >> H	180 mm
! R >> L	350 mm
! A >> P	40 mm
Reset	Off

System - pTx Volumes

B1 Shim mode	TrueForm
Excitation	Slab-sel.

System - Tx/Rx

Frequency 1H	297.210751 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	0.200
Reset	Off
! Ref. amplitude 1H	320.000 V

Sequence - Part 1

Introduction	On
Dimension	3D
Reordering	Linear
Asymmetric echo	Off
Contrasts	1
Multi-slice mode	Interleaved
Echo spacing	1.08 ms
Bandwidth	984 Hz/Px

Sequence - Part 2

EPI factor	58
------------	----

Sequence - Part 2

Segmentation	6
RF pulse type	Normal
Gradient mode	Fast
Excitation	Slab-sel.
RF spoiling	On
Turbo factor	18

Sequence - Special

PATRef FA	3 deg
RF duration	2000 us
RF BWT product	15
Ernst T1	1200 ms
PATRef prep. shots	10
Volume dummy shots	0
Dummy Measurements	0
Invert PE	Off
Min. TE if PF	On
Echo Time Shift	On
Ramp Sampling	On
NORDIC	On
Water Exc.	-none-
External PC	per Series
Saturation RF	per Shot
EPI rise time factor	1.00
Mosaic DICOMs	On
Modify Ice Config	On
GRAPPA Regularization	50000 10 ⁻⁶
HSN RF power scale	2.50
Inversion Delay	650 ms
Relaxation Delay	0 ms
Var. FA /MAGEC	4

Sequence - Assistant

Mode	Off
------	-----

\\USER\Test\Joseph\20210920_VASOremote\rslh_Gz3_S6_C1_0p39_PCshot_run2

TA: 14:55 PM: FIX Voxel size: 0.4×0.4×0.4 mmPAT: 3 Rel. SNR: 1.00 : d238999c

Properties

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	Off
Start measurements	Single measurement

Routine

Slab group	1
Slabs	1
Position	L0.0 P53.3 H21.0 mm
Orientation	C > T-37.4
Phase enc. dir.	R >> L
AutoAlign	---
Slab Scale	-10 %
Slices per slab	18
FoV read	182 mm
FoV phase	100.0 %
Slice thickness	0.39 mm
TR 1	71.4 ms
TR 2	7136 ms
TE 1	26.80 ms
Averages	1
Multi-echo Shots	1
Filter	Distortion Corr.(3D)
Coil elements	R96

Contrast - Common

TR 1	71.4 ms
TR 2	7136 ms
TE 1	26.80 ms
Multi-echo spacing	63.37 ms
Magn. preparation	Non-sel. HSN IR
TI 1	1302.6 ms
TI 2	2587.8 ms
Flip angle	40 deg
Fat suppr.	Fat sat.
Magn. Prep. Shots	2

Contrast - Dynamic

Averages	1
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	124
Pause after meas.	0.0 s

Resolution - Common

FoV read	182 mm
FoV phase	100.0 %
Slice thickness	0.39 mm
Base resolution	462
Phase resolution	100 %
Slice resolution	100 %
Phase partial Fourier	6/8

Resolution - Common

Slice partial Fourier	Off
Interpolation	Off

Resolution - iPAT

PAT mode	CAIPIRINHA
Acc. factor PE	1
Ref. lines PE	75
Acc. factor 3D	3
Ref. lines 3D	18
CAIPI 3D Shift	1
Reference Scan Mode	GRE/separate
CAIPIRINHA mode	Free

Resolution - Filter Image

Image Filter	Off
Distortion Corr.	On
Mode	3D
Unfiltered images	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off

Resolution - Filter Rawdata

Raw filter	Off
Elliptical filter	Off

Geometry - Common

Slab group	1
Slabs	1
Position	L0.0 P53.3 H21.0 mm
Orientation	C > T-37.4
Phase enc. dir.	R >> L
Slab Scale	-10 %
Slices per slab	18
FoV read	182 mm
FoV phase	100.0 %
Slice thickness	0.39 mm
TR 1	71.4 ms
TR 2	7136 ms
Multi-slice mode	Interleaved
Series	Ascending
Multi-echo Shots	1

Geometry - AutoAlign

Slab group	1
Position	L0.0 P53.3 H21.0 mm
Orientation	C > T-37.4
Phase enc. dir.	R >> L
AutoAlign	---
Initial Position	L0.0 P53.3 H21.0
L	0.0 mm
P	53.3 mm
H	21.0 mm
Initial Rotation	0.00 deg
Initial Orientation	C > T
C > T	-37.4
> S	0.0

Geometry - Saturation

Saturation mode	Standard
-----------------	----------

Geometry - Saturation

Fat suppr.	Fat sat.
------------	----------

Geometry - Tim Planning Suite

Set-n-Go Protocol	Off
Table position	H
Table position	0 mm
Inline Composing	Off

System - Miscellaneous

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
Save uncombined	Off
Matrix Optimization	Off
AutoAlign	---
Coil Select Mode	Default

System - Adjustments

B0 Shim mode	Tune up
B1 Shim mode	TrueForm
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

System - Adjust Volume

! Position	L0.0 P48.1 H17.8 mm
! Orientation	C > T-36.7
! Rotation	-90.00 deg
! F >> H	180 mm
! R >> L	350 mm
! A >> P	40 mm
Reset	Off

System - pTx Volumes

B1 Shim mode	TrueForm
Excitation	Slab-sel.

System - Tx/Rx

Frequency 1H	297.210751 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	0.200
Reset	Off
! Ref. amplitude 1H	300.000 V

Sequence - Part 1

Introduction	On
Dimension	3D
Reordering	Linear
Asymmetric echo	Off
Contrasts	1
Multi-slice mode	Interleaved
Echo spacing	1.08 ms
Bandwidth	984 Hz/Px

Sequence - Part 2

EPI factor	58
------------	----

Sequence - Part 2

Segmentation	6
RF pulse type	Normal
Gradient mode	Fast
Excitation	Slab-sel.
RF spoiling	On
Turbo factor	18

Sequence - Special

PATRef FA	3 deg
RF duration	2000 us
RF BWT product	15
Ernst T1	1200 ms
PATRef prep. shots	10
Volume dummy shots	0
Dummy Measurements	0
Invert PE	Off
Min. TE if PF	On
Echo Time Shift	On
Ramp Sampling	On
NORDIC	On
Water Exc.	-none-
External PC	per Series
Saturation RF	per Shot
EPI rise time factor	1.00
Mosaic DICOMs	On
Modify Ice Config	On
GRAPPA Regularization	50000 10 ⁻⁶
HSN RF power scale	2.50
Inversion Delay	650 ms
Relaxation Delay	0 ms
Var. FA /MAGEC	4

Sequence - Assistant

Mode	Off
------	-----